

AD-A098 258

SOUTHWEST RESEARCH INST SAN ANTONIO TX ARMY FUELS AN--ETC F/G 21/4  
INFRARED ANALYSIS OF GASOLINE/ALCOHOL BLENDS.(U)

FEB 81 F M NEWMAN, S M BIELA


DAAK70-80-C-0001

UNCLASSIFIED


AFLRL-134

NL

1 1  
A  
2




END  
DATE  
FORM  
81  
DTIC



MICROCOPY RESOLUTION TEST CHART  
NATIONAL BUREAU OF STANDARDS 1963 A

AD A 098258

**LEVEL**



12

**INFRARED ANALYSIS OF  
GASOLINE/ALCOHOL BLENDS**

**INTERIM REPORT**

**AFLRL No. 134**

**DTIC**  
APR 27 1981

by

**Frank M. Newman**

**Susan M. Biela**

**U. S. Army Fuels and Lubricants Research Laboratory**

**Southwest Research Institute**

**San Antonio, Texas**

Under contract to

**U. S. Army Mobility Equipment Research**

**and Development Command**

**Fort Belvoir, Virginia**

**Contract No. DAAK70-80-C-0001**

Approved for public release; distribution unlimited

January 1981

81 4 27 050

DTIC FILE COPY

### **Disclaimers**

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorized documents.

Trade names cited in this report do not constitute an official endorsement or approval of the use of such commercial hardware or software.

### **DDC Availability Notice**

Qualified requestors may obtain copies of this report from Defense Documentation Center, Cameron Station, Alexandria, Virginia 22314.

### **Disposition Instructions**

Destroy this report when no longer needed. Do not return it to the originator.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER AFLRL <del>134</del> ✓	2. GOVT ACCESSION NO. AD-A098 258	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) INFRARED ANALYSIS OF GASOLINE/ALCOHOL BLENDS.		5. TYPE OF REPORT & PERIOD COVERED Interim Report.
7. AUTHOR(s) Frank M. Newman Susan M. Biela		6. PERFORMING ORG. REPORT NUMBER AFLRL <del>134</del>
9. PERFORMING ORGANIZATION NAME AND ADDRESSES U.S. Army Fuels & Lubricants Research Laboratory Southwest Research Institute San Antonio, TX 78284		8. CONTRACT OR GRANT NUMBER(s) DAAK70-80-C-0001
11. CONTROLLING OFFICE NAME AND ADDRESS U.S. Army Mobility Equipment Research & Development Command, Energy & Water Resources Lab. Fort Belvoir, VA 22060		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS 12 30
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE 11 Feb 1981
		13. NUMBER OF PAGES 21
		15. SECURITY CLASS. (of this report) Unclassified
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Infrared Analysis                      Gasoline Gasohol                                      Oxygenates Alcohols		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) An infrared method for qualitative and quantitative determination of alcohols in gasoline/alcohol blends has been developed. The method has high potential for speed, low cost, and specificity.  Gasoline containing methanol, ethanol, iso-propanol, t-butanol, methyl t-butyl ether, and methyl iso-butyl ketone (MIBK) can be analyzed using a computing infrared spectrophotometer. → next page		

DTIC  
APR 27 1981

DD FORM 1 JAN 73 1473

EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED

387 339

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

mt

**UNCLASSIFIED**

**SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)**

20. ABSTRACT (Cont'd)

Due to variable composition of different gasolines, each base gasoline analyzed required a new quantitative program to be established, using standards prepared with that particular gasoline.

**UNCLASSIFIED**

**SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)**

FOREWORD

The work reported herein was conducted at the U.S. Army Fuels and Lubricants Research Laboratory (USAFRLRL) located at Southwest Research Institute, San Antonio, Texas under Contract No. DAAK70-80-C-0001. The contracting officer's representative was Mr. F.W. Schaekel, Energy and Water Resources Laboratory, USAMERADCOM, DRDME-GL, Fort Belvoir, Virginia.

Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By _____	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A	

TABLE OF CONTENTS

	<u>Page</u>
I. INTRODUCTION.....	5
II. APPROACH.....	6
III. EXPERIMENTAL.....	10
A. Standards.....	10
B. Instrumental.....	10
C. Sample Preparation and Handling.....	14
IV. RESULTS.....	19
V. CONCLUSIONS.....	22

## LIST OF ILLUSTRATIONS

<u>Figure</u>		<u>Page</u>
1	Typical Infrared Spectrum of a Fuel Blend.....	7
2	Instrument Printout of Peakpick Mode.....	8
3	Infrared Spectrum of a Fuel Blend.....	9
4	Infrared Spectra of Low and High Standards for Methanol.....	11
5	Infrared Response for Methanol in Gasoline.....	12
6	Point Program Printout.....	13
7	Normalize Printouts for Two Standards.....	16
8	Array Build Mode Printout.....	17
9	Analysis Printout.....	18

## LIST OF TABLES

<u>Table</u>		<u>Page</u>
1	Alcohol Analyte Band Numbers.....	6
2	Alcohol-Gasoline Blends: IR Analysis.....	20
3	Multicomponent Analysis.....	21

## I. INTRODUCTION

The use of alcohols as fuel extenders is being considered to reduce U.S. dependence on foreign oil and to prolong the existing supply of domestic petroleum. The use of ethanol in particular has been intensively researched for possible economic and environmental advantages over neat gasoline. For example, a blend of 10 vol% ethanol in gasoline increases the antiknock index 2-5 octane numbers over that of the base gasoline and reduces the hydrocarbon and carbon monoxide exhaust emissions. However, deleterious effects can arise from incompatibilities of ethanol with different gasolines and elastomers which may be accentuated in storage, routine handling and distribution. As a result, other oxygenates such as methanol, iso-propanol, t-butanol, methyl-t-butyl ether, and methyl-iso-butyl ketone (MIBK) blended with gasoline alone or in various combinations are being considered to determine any advantages over ethanol/gasoline blends.

Several analytical methods dealing with gasohol as well as synfuels and conventional fossil fuels have been studied. Various qualitative and quantitative methods for the determination of alcohols in gasohols are being evaluated. Infrared spectrophotometry has been found to be a useful method for the quantitative determination of the volume percent of oxygenates in a gasoline blend for specification conformation.

## II. APPROACH

An infrared spectroscopic method for oxygenates has been emphasized because of its high potential for speed, low cost, and specificity. Figure 1 shows a typical IR spectrum of a fuel blend. Using the PEAKPICK MODE of a Beckman Microlab 620MX computing infrared spectrophotometer, analytical frequencies for each component were chosen where there was little or no interference from other components. Figure 2 shows an instrument printout, and Table 1 lists

TABLE 1. ALCOHOL ANALYTE BAND NUMBERS

<u>Component</u>	<u>Analytical Frequency, cm<sup>-1</sup></u>
Gasoline	967
Methanol	1030
Ethanol	882
<u>iso-propanol</u>	952
<u>t-butanol</u>	914
<u>Methyl-t-butyl ether</u>	1086

these chosen frequencies. Three background frequencies, 1255 cm<sup>-1</sup>, 1230 cm<sup>-1</sup>, and 860 cm<sup>-1</sup>, were also selected which bracket the analytical frequencies and are used for baseline correction. Figure 3 shows the IR spectrum of a fuel blend with the analytical and baseline frequencies marked. Using these background points, the net percent transmittance at each frequency is obtained, converted to absorbance values, and normalized with respect to unit path length and concentration.

Each set of values was stored in array rows of a net normalized absorbance matrix and was later used to calculate the percentages of alcohols in unknown samples.




FIGURE 1. TYPICAL INFRARED SPECTRUM OF A FUEL BLEND

SCAN ROUTINES: 1-SCAN 2-REPCAN  
3-PEAKPICK 4-TIMEDRIVE

\*\*\*\*\* PEAKPICK \*\*\*\*\*

SET %T DISCRIMINATOR

50

ENTER SCANNING PARAMETERS THEN  
PUSH START

SCAN PARAMETERS: 1-RANGE 2-SPEED

3-CHART FORMAT 4-SLIT 5-GAIN

ENTER RANGE

2000 TO 0600

1604CM-1	46.7%T
1512CM-1	49.9%T
1494CM-1	22.2%T
1459CM-1	9.5%T
1378CM-1	26.0%T
1080CM-1	49.7%T
1030CM-1	6.2%T
0967CM-1	48.7%T
0888CM-1	50.0%T
0803CM-1	44.4%T
0793CM-1	37.7%T
0767CM-1	24.1%T
0727CM-1	1.1%T
0689CM-1	1.5%T
0667CM-1	0.2%T
0640CM-1	0.5%T
0633CM-1	0.4%T
0618CM-1	0.4%T
0612CM-1	0.5%T
0608CM-1	0.5%T

FIGURE 2. INSTRUMENT PRINTOUT OF PEAKPICK MODE




FIGURE 3. INFRARED SPECTRUM OF A FUEL BLEND

### III. EXPERIMENTAL

All calculations were performed using a matrix set up for each gasoline by the Beckman Microlab 620MX computing infrared spectrophotometer. An unknown sample was analyzed using a precisely calibrated (0.025 cm nominal) sealed cell to prevent evaporation. The cell was filled completely, including the syringe fittings to ensure that absolutely no air was trapped.

A NaCl cell was used with no difficulty since no analytical frequencies were needed between 700 and 600  $\text{cm}^{-1}$ . However, the attack on the cell by methanol was too severe and required the use of Irtran cells. After the analysis was completed, the cell was cleaned with heptane and dried.


#### A. Standards

To establish a procedure for preparing standards, spectra of varying concentrations of each alcohol in a gasoline were obtained, with Figure 4 showing a low and high standard for methanol. The net peak absorbance was plotted as a function of concentration for each alcohol, and a linear response with a non-zero intercept was obtained in each case. Figure 5 shows the plot for methanol. As a result, a low and high calibration standard is used for each component in a non-zero intercept method.

#### B. Instrumental

Using the quantitative routines on the computing infrared spectrophotometer, a POINT PROGRAM MODE and NORMALIZE MODE are run for each component. A POINT PROGRAM operation obtains transmittance measurement at prescribed wavelengths and stores these data for further use. A POINT PROGRAM printout is shown in Figure 6. The NORMALIZE MODE operation retrieves the data stored by the POINT PROGRAM, and with cell pathlength and component concentration values supplied, converts all data to a "normalized" standard format for future use. These normalized data take into account variations in cell pathlengths in subsequent analyses and the concentration of components in the standardization runs. Although used as the solvent for the alcohols, the gasoline was run as an undiluted component due to some absorbance at the analytical frequencies and

HIGH



LOW




FIGURE 4. INFRARED SPECTRA OF LOW AND HIGH STANDARDS FOR METHANOL




FIGURE 5. INFRARED RESPONSE FOR METHANOL IN GASOLINE

QUANT ROUTINES: 1-PT PROGRAM 2-NORMALIZE  
3-ARRAYBUILD 4-ANALYSIS 5-ENTER ARRAY

##### POINT PROGRAM #####

ENTER NUMBER OF POINTS

5

POINTS MUST BE ENTERED IN  
DECREASING ORDER

ENTER A WAVENUMBER

1255 CM-1

ENTER A WAVENUMBER

1230 CM-1

ENTER A WAVENUMBER

1030 CM-1

ENTER A WAVENUMBER

0967 CM-1

ENTER A WAVENUMBER

0860 CM-1

ENTER NUMBER OF READINGS TO AVERAGE

09

PUSH START

#1 1255 CM-1 77.3 %T

#2 1230 CM-1 78.6 %T

#3 1030 CM-1 8.2 %T

#4 0967 CM-1 61.3 %T

#5 0860 CM-1 78.3 %T

QUANT ROUTINES: 1-PT PROGRAM 2-NORMALIZE  
3-ARRAYBUILD 4-ANALYSIS 5-ENTER ARRAY

FIGURE 6. POINT PROGRAM PRINTOUT

possible variability in concentration. The absorbance values obtained for the base gasoline were stored in Array Row 1 of the net normalized absorbance matrix and were used to correct the absorbance values at the other analytical frequencies.

A POINT PROGRAM and NORMALIZE MODE were run for each low and high standard of the alcohol components. The low standard, usually 1 percent alcohol, was run first, and the high standard, usually 10 percent alcohol, was then run as an intercept correction. The original concentration of the alcohol (1 percent) and the analyte band number of that particular alcohol were identified as required by the program. The absorbance values obtained for the high standard were entered in the same array row as the values for the low standard. NORMALIZE printouts for the two standards are shown in Figure 7. After all the components had been entered, the normalized absorbance matrix was inverted using the ARRAY BUILD MODE which formed the analysis matrix shown in Figure 8.

A sample was scanned using the ANALYSIS MODE of the quantitative routines which recalled the matrix set up by the computer, measured the data at the specified frequencies, and computed the concentrations of the component present. The sample was normalized to 100 percent since the sum of all components was 100 percent. Figure 9 shows an ANALYSIS printout.

### C. Sample Preparation and Handling

Blends of alcohols with different base fuels were analyzed to study the gasoline matrix effect. Since different gasolines have variable compositions, a new quantitative program was established for each new gasoline to be analyzed. To identify which alcohols were present as analytes, a scan of the blend was obtained and examined for peaks at the respective analyte wavelengths. If the base fuel was not readily available, it was recovered by removing the alcohols by water extraction.

A 100-ml sample was extracted using two 50-ml portions of water. After the washings, the extracted gasoline was filtered over  $\text{Na}_2\text{SO}_4$ . The recovered fuel was then used to prepare the low and high standards needed to set up the new quantitative program.

Because methanol severely attacks the NaCl cell, methanol was run using an Irtran cell. After being placed in the infrared spectrophotometer, the cell was allowed to warm up for approximately 3-5 minutes to reach a constant pathlength before analysis.

#####	NORMALIZE	#####	#####	NORMALIZE	#####	
01	1255 CH-1	81.3 XT		01	1255 CH-1	77.3 XT
02	1230 CH-1	82.8 XT		02	1230 CH-1	78.6 XT
03	1030 CH-1	51.1 XT		03	1030 CH-1	8.2 XT
04	0967 CH-1	63.9 XT		04	0967 CH-1	61.3 XT
05	0060 CH-1	81.6 XT		05	0060 CH-1	78.3 XT

ENTER BACKGROUND NUMBERS	ENTER BACKGROUND NUMBERS
01	01
ENTER BACKGROUND NUMBERS	ENTER BACKGROUND NUMBERS
02	02
ENTER BACKGROUND NUMBERS	ENTER BACKGROUND NUMBERS
05	05
ENTER PATHLENGTH	ENTER PATHLENGTH
.020	.020
ENTER CONCENTRATION	ENTER CONCENTRATION
0001	0010
ENTER SOLVENT CONCENTRATION	ENTER SOLVENT CONCENTRATION
0099	0090
ENTER ARRAY ROW NUMBER	ENTER ARRAY ROW NUMBER
2	2
IS THIS AN INTERCEPT CORRECTION?	IS THIS AN INTERCEPT CORRECTION?
1 FOR YES 0 FOR NO	1 FOR YES 0 FOR NO
NO	YES
PUSH START	IS THIS THE FIRST CORRECTION?
ROW = 0 2	1 FOR YES 0 FOR NO
PATHLENGTH = .020	YES
CONCENTRATION = 0001	ENTER ORIGINAL CONCENTRATION
SOLVENT CONCENTRATION = 0099	0001
1030 CH-1 1.0005 A	ENTER ANALYTE BAND #
0967 CH-1 0.0101 A	3
	PUSH START
	ROW = 0 2
	PATHLENGTH = .020
	CONCENTRATION = 0010
	SOLVENT CONCENTRATION = 0090
	1030 CH-1 3.0007 A
	0967 CH-1 0.0305 A
	INTERCEPT = - 0.035
	CORRECTED ABSORBTIVITY = 3.124

FIGURE 7. NORMALIZE PRINTOUTS FOR TWO STANDARDS

```

##### ARRAY BUILD #####
IS THIS AN INTERCEPT CORRECTION?
1 FOR YES 0 FOR NO
YES
ENTER ANALYSIS #
1
ENTER SLIT TYPE:
0-RESOLUTION 1-NORMAL
NORMAL SLIT
ROW = # 1
1030 CM-1 0.0553 A
0967 CM-1 0.0389 A
ROW = # 2
1030 CM-1 3.1243 A
0967 CM-1 0.0305 A
INTERCEPT ROW
1030 CM-1 - 0.0346 A
0967 CM-1 0.0000 A
MODIFY? 1-YES 0-NO
PUSH START
ROW = # 1
1030 CM-1 - 0.2543
0967 CM-1 0.3246
ROW = # 2
1030 CM-1 26.0477
0967 CM-1 - 0.4610
INVERTED INTERCEPTS
1030 CM-1 - 0.0000
0967 CM-1 0.0112

##### STANDBY #####

```

FIGURE 8. ARRAY BUILD MODE PRINTOUT

```

##### ANALYSIS #####
ENTER ANALYSIS #
1
ENTER NUMBER OF READINGS TO AVERAGE
09
ENTER PATHLENGTH
.045
NORMALIZE ANSWERS TO 100% ?
1 FOR YES      0 FOR NO
YES
PUSH START
ANALYSIS
NORMAL SLIT
GAIN SETTING = 040
#1  1255 CH-1    49.0 XT
#2  1230 CH-1    45.8 XT
#3  1030 CH-1     0.1 XT
#4  0967 CH-1    29.5 XT
#5  0860 CH-1    39.8 XT

1030 CH-1        2.6209 A
0967 CH-1        0.1404 A
COMPONENT 1      00.0206
COMPONENT 2      19.9794

```

FIGURE 9. ANALYSIS PRINTOUT

#### IV. RESULTS

Single-component standards of known concentrations were run in tests of repeatability and reproducibility. Low standard and average deviations for all alcohols are shown in Table 2. In tests of repeatability, the relative deviations remained low for the 1- and 10-vol% blends. The 20-vol% blends have higher standard deviations, but this concentration is not normally encountered. The relative deviations, which remained fairly constant for all blends, were larger for the reproducibility tests due to the decreased transmittance of the Irtran cells.

Several multicomponent mixtures of methanol, ethanol, and *t*-butanol were prepared and analyzed. Table 3 shows the results of these analyses. Compared to the single-component analysis shown in Table 2, a multicomponent analysis was not quite as accurate due to the functional group similarity of the alcohols. However, this is a good method to identify the relative concentrations of any alcohols present. A single-component analysis has been recommended to be set up for MIBK at  $1725\text{ cm}^{-1}$ , since MIBK does not use the same background points of  $1255\text{ cm}^{-1}$ ,  $1230\text{ cm}^{-1}$ , and  $860\text{ cm}^{-1}$  as used for alcohols.

Actual analysis of an individual sample required only five minutes. Approximately thirty minutes is needed to set up a single-component quantitative program. An extra twenty minutes is needed for each additional alcohol, and if necessary, a water extraction can be completed in thirty minutes.

TABLE 2. ALCOHOL-GASOLINE BLENDS: IR ANALYSIS

		Repeatability		
Alcohol	Label	Value, %	5 Runs	
			Stan. Dev.	Avg. Dev.
MeOH*	1%	1.04	0.01	0.01
	10%	9.49	0.09	0.07
	20%	19.39	0.14	0.12
EtOH	1%	1.18	0.02	0.01
	10%	9.73	0.02	0.02
	20%	19.45	0.04	0.03
IPA	1%	1.16	0.01	0.004
	10%	9.58	0.05	0.04
	20%	18.90	0.24	0.18
TBA	1%	1.01	0.01	0.01
	10%	9.73	0.03	0.02
	20%	19.87	0.06	0.05
MTBE	1%	1.01	0.05	0.04
	10%	10.37	0.06	0.04
	20%	17.97	0.21	0.17

		Reproducibility		
Alcohol	Label	Value, %	5 Runs	
			Stan. Dev.	Avg. Dev.
MeOH*	1%	0.99	0.05	0.03
	10%	8.81	0.11	0.08
	20%	20.91	0.15	0.11
EtOH	1%	1.10	0.02	0.01
	10%	9.88	0.06	0.05
	20%	19.13	0.05	0.04
IPA	1%	1.16	0.01	0.01
	10%	9.86	0.04	0.03
	20%	18.24	0.08	0.05
TBA	1%	1.00	0.01	0.004
	10%	10.06	0.05	0.04
	20%	19.55	0.04	0.05
MTBE	1%	0.96	0.003	0.003
	10%	10.13	0.04	0.04
	20%	20.11	0.08	0.06

\*Intran cell used with 4 runs.

TABLE 3. MULTICOMPONENT ANALYSIS

<u>Sample</u>	<u>Component</u>	<u>Known vol%</u>	<u>Determined vol%</u>
A	MeOH	1	1.15
	EtOH	9	9.58
	TBA	5	5.35
B	MeOH	2.75	3.16
	EtOH	10.0	12.05
	TBA	2.75	2.95
C	MeOH	2.75	3.08
	EtOH	---	0.06
	TBA	2.75	3.01

## V. CONCLUSIONS

The infrared method for oxygenates is a quick, economical way to quantitatively determine the presence of alcohols in gasoline/alcohol blends. The Beckman Microlab 620MX computing infrared spectrophotometer is easily programmed for this analysis.

Both single- and multicomponent analyses were run using several different alcohols. However, due to the functional group similarity of the alcohols, the multicomponent analysis was not as accurate. The use of Irtran cells for methanol blends increases the relative deviations due to lower transmittance. When not readily available to prepare the standards, the base fuel is recovered by water extraction of the alcohols.

Once the quantitative program is set up for a particular gasoline, an actual analysis requires only five minutes, and tests of repeatability and reproducibility have shown low standard deviations for each alcohol. Multicomponent standards are also easily set up for methanol, ethanol, and t-butanol mixtures, but a single-component program using different background points is needed for methyl-t-butyl ether.

DISTRIBUTION LIST

DEPARTMENT OF DEFENSE		CDR
DEFENSE DOCUMENTATION CTR		US ARMY MATERIAL DEVEL&READINESS
CAMERON STATION	12	COMMAND
ALEXANDRIA VA 22314		ATTN: DRCLDC (MR BENDER) 1
		DRCMM-SP (LTC O'CONNER) 1
DEPT OF DEFENSE		DRCQA-E (MR SMART) 1
OFC OF SECRETARY OF DEFENSE		DRCDE-DG (MR MCGOWAN) 1
ATTN USD(R&E) DR ADICOFF	1	DRCIS-S (MR SPRAGUE) 1
WASHINGTON DC 20301		DRCIS-C (LTC CROW) 1
		5001 EISENHOWER AVE
		ALEXANDRIA VA 22333
DEPT OF DEFENSE		CDR
ATTN: DASA(MRA&L)-ES(MR DYCKMAN)	1	US ARMY TANK-AUTOMOTIVE CMD
WASHINGTON DC 20301		ATTN DRSDA-NW (TWVMO) 1
COMMANDER		DRSTA-RG (MR HAMPARIAN) 1
DEFENSE LOGISTICS AGY		DRSTA-NS (DR PETRICK) 1
ATTN DLA-SME (MRS P MCLAIN)	1	DRSTA-J 1
CAMERON STATION		DRSTA-G (COL MILLS) 1
ALEXANDRIA VA 22314		DRSTA-M 1
		DRSTA-GBP (MR MCCARTNEY) 1
COMMANDER		WARREN MI 48090
DEFENSE FUEL SUPPLY CTR		DIRECTOR
ATTN: DFSC-T	1	US ARMY MATERIAL SYSTEMS
CAMERON STA		ANALYSIS AGENCY
ALEXANDRIA VA 22314		ATTN DRXSY-CM 1
COMMANDER		DRXSY-S 1
DEFENSE GENERAL SUPPLY CTR		DRXSY-L 1
ATTN: DGSC-SSA	1	ABERDEEN PROVING GROUND MD 21005
RICHMOND VA 23297		CDR
DEPARTMENT OF THE ARMY		US ARMY APPLIED TECH LAB
HQ, DEPT OF ARMY		ATTN DAVDL-ATL-DD 1
ATTN: DALO-TSE	1	DAVDL-ATL-ATP (MR MORROW) 1
DALO-AU	1	DAVDL-ATL-ASV 1
DAMA-CSS-P (DR BRYANT)	1	FORT EUSTIS VA 23604
DAMA-ARZ (DR CHURCH)	1	HQ, 172D INFANTRY BRIGADE (ALASKA)
DAMA-SM/E	1	ATTN AFZT-DI-L 1
WASHINGTON DC 20310		AFZT-DI-M 1
CDR		DIRECTORATE OF INDUSTRIAL
U.S. ARMY MOBILITY EQUIPMENT		OPERATIONS
R&D COMMAND		FT RICHARDSON AK 99505
Attn: DRDME-GL	10	CDR
FORT BELVOIR VA 22060		US ARMY GENERAL MATERIAL &
		PETROLEUM ACTIVITY
		ATTN STSGP-FT (MS GEORGE) 1
		STSGP-PE 1
		STSGP (COL HILL) 1
		NEW CUMBERLAND ARMY DEPOT
		NEW CUMBERLAND PA 17070

CDR  
US ARMY ARRCOM, LOG ENGR DIR  
ATTN DRSAR-LEM (MR MENKE) 1  
ROCK ISLAND ARSENAL IL 61299

CDR  
US ARMY COLD REGION TEST CENTER  
ATTN STECR-TA (MR HASLEM) 1  
APO SEATTLE 98733

CDR  
US ARMY RES & STDZN GROUP  
(EUROPE)  
ATTN DRXSN-E-RA 1  
BOX 65  
FPO NEW YORK 09510

HQ, US ARMY AVIATION R&D CMD  
ATTN DRDAV-D (MR CRAWFORD) 1  
DRDAV-N (MR BORGMAN) 1  
DRDAV-E (MR LONG) 1  
P O BOX 209  
ST LOUIS MO 63166

CDR  
US ARMY FORCES COMMAND  
ATTN AFLG-REG (MR HAMMERSTROM) 1  
AFLG-POP (MR COOK) 1  
FORT MCPHERSON GA 30330

CDR  
US ARMY ABERDEEN PROVING GROUND  
ATTN STEAP-MT 1  
STEAP-MT-U (MR DEEVER) 1  
ABERDEEN PROVING GROUND MD 21005

CDR  
US ARMY YUMA PROVING GROUND  
ATTN STEYP-MT (MR DOEBBLER) 1  
YUMA AR 85364

MICHIGAN ARMY MISSILE PLANT  
OFC OF PROJ MGR, XM-1 TANK SYS  
ATTN DRCPM-GCM-S 1  
WARREN MI 48090

MICHIGAN ARMY MISSILE PLANT  
PROG MGR, FIGHTING VEHICLE SYS  
ATTN DRCPM-FVS-SE 1  
WARREN MI 48090

PROJ MGR, M60 TANK DEVELOPMENT  
ATTN DRCPM-M60-E (MR WESALA) 1  
WARREN MI 48090

PROG MGR, M113/M113A1 FAMILY  
OF VEHICLES  
ATTN DRCPM-M113 1  
WARREN MI 48090

PROJ MGR, MOBILE ELECTRIC POWER  
ATTN DRCPM-MEP-TM 1  
7500 BACKLICK ROAD  
SPRINGFIELD VA 22150

OFC OF PROJ MGR, IMPROVED TOW  
VEHICLE  
US ARMY TANK-AUTOMOTIVE R&D CMD  
ATTN DRCPM-ITV-T 1  
WARREN MI 48090

CDR  
US ARMY EUROPE & SEVENTH ARMY  
ATTN AEAGC-FMD 1  
APO NY 09403

PROJ MGR, PATRIOT PROJ OFC  
ATTN DRCPM-MD-T-G 1  
US ARMY DARCOM  
REDSTONE ARSENAL AL 35809

CDR  
THEATER ARMY MATERIAL MGMT  
CENTER (200TH)  
DIRECTORATE FOR PETROL MGMT  
ATTN AEAGD-MM-PT-Q (MR PINZOLA) 1  
ZWEIBRUCKEN  
APO NY 09052

CDR  
US ARMY RESEARCH OFC  
ATTN DRXRO-EG 1  
DRXRO-CB (DR GHIRARDELLI) 1  
DRXRO-TT (DR SCHMIEORSHOFF) 1  
P O BOX 12211  
RSCH TRIANGLE PARK NC 27709

DIR  
US ARMY R&T LAB  
ADVANCED SYSTEMS RSCH OFC  
ATTN MR D WILSTED 1  
AMES RSCH CTR  
MOFFITT FIELD CA 94035

CDR TOBYHANNA ARMY DEPOT ATTN SDSTO-TP-S TOBYHANNA PA 18466	1	HQ, US ARMY ARMAMENT R&D CMD ATTN DRDAR-SCM-OO (MR MUFFLEY) DRDAR-TST-S DOVER NJ 07801	1 1
DIR US ARMY MATERIALS & MECHANICS RSCH CTR ATTN DRXMR-E DRXMR-T DRXMR-R WATERTOWN MA 02172	1 1 1	HQ, US ARMY TROOP SUPPORT & AVIATION MATERIAL READINESS COMMAND ATTN DRSTS-MEG (2) DRCPO-PDE (LTC FOSTER) 4300 GOODFELLOW BLVD ST LOUIS MO 63120	1 1
CDR US ARMY DEPOT SYSTEMS CMD ATTN DRSDS CHAMBERSBURG PA 17201	1	DEPARTMENT OF THE ARMY CONSTRUCTION ENG RSCH LAB ATTN CERL-EM CERL-ZT CERL-EH P O BOX 4005 CHAMPAIGN IL 61820	1 1 1
CDR US ARMY WATERVLIET ARSENAL ATTN SARWY-RDD WATERVLIET NY 12189	1	HQ US ARMY TRAINING & DOCTRINE CMD ATTN ATCD-SL (MR RAFFERTY) ATCD-TA ATCD-D FORT MONROE VA 23651	1 1 1
CDR US ARMY LEA ATTN DALO-LEP NEW CUMBERLAND ARMY DEPOT NEW CUMBERLAND PA 17070	1	DIRECTOR US ARMY RSCH & TECH LAB (AVRADCOM) PROPULSION LABORATORY ATTN DAVDL-PL-D (MR ACURIO) 21000 BROOKPARK ROAD CLEVELAND OH 44135	1 1 1
CDR US ARMY GENERAL MATERIAL & PETROLEUM ACTIVITY ATTN STSGP-PW (MR PRICE) SHARPE ARMY DEPOT LATHROP CA 95330	1	CDR US ARMY NATICK RES & DEV CMD ATTN DRDNA-YEP (DR KAPLAN) NATICK MA 01760	1
CDR US ARMY FOREIGN SCIENCE & TECH CENTER ATTN DRXST-MT1 FEDERAL BLDG CHARLOTTESVILLE VA 22901	1	CDR US ARMY TRANSPORTATION SCHOOL ATTN ATSP-CD-MS FORT EUSTIS VA 23604	1
CDR DARCOM MATERIAL READINESS SUPPORT ACTIVITY (MRSA) ATTN DRXMD-MD LEXINGTON KY 40511	1	CDR US ARMY QUARTERMASTER SCHOOL ATTN ATSM-CTD-MS ATSM-TNG-PT (COL VOLPE) FORT LEE VA 23801	1 1
HQ, US ARMY T&E COMMAND ATTN DRSTE-TO-O ABERDEEN PROVING GROUND, MD 21005	1	HQ, US ARMY ARMOR SCHOOL ATTN ATSB-TD FORT KNOX KY 40121	1

CDR  
US ARMY LOGISTICS CTR  
ATTN ATCL-MS (MR A MARSHALL) 1  
FORT LEE VA 23801

CDR  
US ARMY FIELD ARTILLERY SCHOOL  
ATTN ATSF-CD 1  
FORT SILL OK 73503

CDR  
US ARMY ORDNANCE CTR & SCHOOL  
ATTN ATSL-CTD-MS 1  
ABERDEEN PROVING GROUND MD 21005

CDR  
US ARMY ENGINEER SCHOOL  
ATTN ATSE-CDM 1  
FORT BELVOIR VA 22060

CDR  
US ARMY INFANTRY SCHOOL  
ATTN ATSH-CD-MS-M 1  
FORT BENNING GA 31905

CDR  
US ARMY AVIATION CTR & FT RUCKER  
ATTN ATZQ-D 1  
FORT RUCKER AL 36362

DEPARTMENT OF THE NAVY

CDR  
NAVAL AIR PROPULSION CENTER  
ATTN PE-71 1  
PE-72 (MR D'ORAZIO) 1  
P O BOX 7176  
TRENTON NJ 06828

CDR  
NAVAL SEA SYSTEMS CMD  
CODE 6101F (MR R LAYNE) 1  
WASHINGTON DC 20362

CDR  
DAVID TAYLOR NAVAL SHIP R&D CTR  
CODE 2830 (MR G BOSMAJIAN) 1  
CODE 2831 1  
ANNAPOLIS MD 2140

JOINT OIL ANALYSIS PROGRAM -  
TECHNICAL SUPPORT CTR 1  
BLDG 780  
NAVAL AIR STATION  
PENSACOLA FL 32508

DEPARTMENT OF THE NAVY  
HQ, US MARINE CORPS  
ATTN LPP (MAJ SANBERG) 1  
LMM (MAJ GRIGGS) 1  
WASHINGTON DC 20380

CDR  
NAVAL AIR SYSTEMS CMD  
ATTN CODE 52032E (MR WEINBURG) 1  
CODE 53645 1  
WASHINGTON DC 20361

CDR  
NAVAL AIR DEVELOPMENT CTR  
ATTN CODE 60612 (MR L STALLINGS) 1  
WARMINSTER PA 18974

CDR  
NAVAL RESEARCH LABORATORY  
ATTN CODE 6170 (MR H RAVNER) 1  
CODE 6180 1  
CODE 6110 (DR HARVEY) 1  
WASHINGTON DC 20375

CDR  
NAVAL FACILITIES ENGR CTR  
ATTN CODE 1202B (MR R BURRIS) 1  
CODE 120B (MR BUSCHELMAN) 1  
200 STOVWALL ST  
ALEXANDRIA VA 22322

CHIEF OF NAVAL RESEARCH  
ATTN CODE 473 (DR R MILLER) 1  
ARLINGTON VA 22217

CDR  
NAVAL AIR ENGR CENTER  
ATTN CODE 92727 1  
LAKEHURST NJ 08733

CDR  
NAVY FACILITIES ENGRG CMD  
CIVIL ENGR SUPPORT OFC  
CODE 15312A (ATTN EOC COOK) 1  
NAVAL CONSTRUCTION BATTALION CTR  
PORT HUENEME CA 93043

CDR, NAVAL MATERIAL COMMAND  
ATTN MAT-08E (DR A ROBERTS) 1  
CP6, RM 606 MAT-08E (MR ZIEM) 1  
WASHINGTON DC 20360

CDR  
NAVY PETROLEUM OPC  
ATTN CODE 40 1  
CAMERON STATION  
ALEXANDRIA VA 22314

CDR  
MARINE CORPS LOGISTICS SUPPORT  
BASE ATLANTIC  
ATTN CODE P841 1  
ALBANY GA 31704

DEPARTMENT OF THE AIR FORCE

HQ, USAF  
ATTN RDPT 1  
WASHINGTON DC 20330

HQ AIR FORCE SYSTEMS CMD  
ATTN AFSC/DLF (LTC RADLOF) 1  
ANDREWS AFB MD 20334

CDR  
US AIR FORCE WRIGHT AERONAUTICAL  
LAB  
ATTN AFWAL/POSF (MR CHURCHILL) 1  
AFWAL/POSL (MR JONES) 1  
WRIGHT-PATTERSON AFB OH 45433

CDR  
USAF SAN ANTONIO AIR LOGISTICS  
CTR  
ATTN SAALC/SFQ (MR MAKRIS) 1  
SAALC/MMPRR (MR ELLIOT) 1  
KELLY AIR FORCE BASE, TX 78241

CDR  
US AIR FORCE WRIGHT AERONAUTICAL  
LAB  
ATTN AFWAL/MLSE (MR MORRIS) 1  
AFWAL/MLBT 1  
WRIGHT-PATTERSON AFB OH 45433

CDR  
USAF WARNER ROBINS AIR LOGISTIC  
CTR  
ATTN WR-ALC/MMIRAB-1 (MR GRAHAM) 1  
ROBINS AFB GA 31098

OTHER GOVERNMENT AGENCIES

US DEPARTMENT OF TRANSPORTATION  
ATTN AIRCRAFT DESIGN CRITERIA  
BRANCH 2  
FEDERAL AVIATION ADMIN  
2100 2ND ST SW  
WASHINGTON DC 20590

US DEPARTMENT OF ENERGY  
DIV OF TRANS ENERGY CONSERV 2  
ALTERNATIVE FUELS UTILIZATION  
BRANCH  
20 MASSACHUSETTS AVENUE  
WASHINGTON DC 20545

DIRECTOR  
NATL MAINTENANCE TECH SUPPORT  
CTR 2  
US POSTAL SERVICE  
NORMAN OK 73069

US DEPARTMENT OF ENERGY  
BARTLESVILLE ENERGY RSCH CTR  
DIV OF PROCESSING & THERMO RES 1  
DIV OF UTILIZATION RES 1  
BOX 1398  
BARTLESVILLE OK 74003

SCI & TECH INFO FACILITY  
ATTN NASA REP (SAK/DL) 1  
P O BOX 8757  
BALTIMORE/WASH INT AIRPORT MD 21240

**DA  
FILM**