

R

USAARL Report No. 93-18

AD-A267 100



**Test and Evaluation Report
of the
Medical Data Electronics Escort Patient Monitor
Model E300A**

By

**Martin D. Quattlebaum (Project Officer)
James E. Bruckart (Project Officer)
Brendan E. Squire (Project Officer)
Leslie W. Johnston (Project Officer)**

**S DTIC
ELECTE
JUL 22 1993
E D**

Impact, Tolerance, and Protection Division

March 1993

93-16401



WPT

93 7 20 005

Approved for public release; distribution unlimited.

**United States Army Aeromedical Research Laboratory
Fort Rucker, Alabama 36362-5292**

Notice

Qualified requesters

Qualified requesters may obtain copies from the Defense Technical Information Center (DTIC), Cameron Station, Alexandria, Virginia 22314. Orders will be expedited if placed through the librarian or other person designated to request documents from DTIC.

Change of address

Organizations receiving reports from the U.S. Army Aeromedical Research Laboratory on automatic mailing lists should confirm correct address when corresponding about laboratory reports.

Disposition

Destroy this document when it is no longer needed. Do not return it to the originator.

Disclaimer

The views, opinions, and/or findings contained in this report are those of the author(s) and should not be construed as an official Department of the Army position, policy, or decision, unless so designated by other official documentation. Citation of trade names in this report does not constitute an official Department of the Army endorsement or approval of the use of such commercial items.

Reviewed:



JOHN V. BARSON
LTC, MC, SFS
Director, Impact, Tolerance,
and Protection Division

Released for publication:



ROGER W. WILEY, O.D., Ph.D.
Chairman, Scientific
Review Committee



DAVID H. KARNEY
Colonel, MC, SFS
Commanding

REPORT DOCUMENTATION PAGE

1a. REPORT SECURITY CLASSIFICATION Unclassified		1b. RESTRICTIVE MARKINGS	
2a. SECURITY CLASSIFICATION AUTHORITY		3. DISTRIBUTION/AVAILABILITY OF REPORT Approved for public release, distribution unlimited	
2b. DECLASSIFICATION/DOWNGRADING SCHEDULE		5. MONITORING ORGANIZATION REPORT NUMBER(S)	
4. PERFORMING ORGANIZATION REPORT NUMBER(S) USAARL Report No. 93-18		7a. NAME OF MONITORING ORGANIZATION U.S. Army Medical Research and Development Command	
6a. NAME OF PERFORMING ORGANIZATION U.S. Army Aeromedical Research Laboratory	6b. OFFICE SYMBOL (if applicable) SGRD-UAD-IE	7b. ADDRESS (City, State, and ZIP Code) Fort Detrick Frederick, MD 21702-5012	
6c. ADDRESS (City, State, and ZIP Code) P.O. Box 620577 Fort Rucker, AL 36362-0577		9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER	
8a. NAME OF FUNDING/SPONSORING ORGANIZATION	8b. OFFICE SYMBOL (if applicable)	10. SOURCE OF FUNDING NUMBERS	
8c. ADDRESS (City, State, and ZIP Code)		PROGRAM ELEMENT NO. 0603807A	PROJECT NO. 3M463807D836 LC
		TASK NO.	WORK UNIT ACCESSION NO. 201
11. TITLE (Include Security Classification) Test and Evaluation Report of the Medical Data Electronics Escort Patient Monitor, Model E300A			
12. PERSONAL AUTHOR(S) Martin D. Quattlebaum, James E. Bruckart, Brendan E. Squire, and Leslie W. Johnston			
13a. TYPE OF REPORT Final	13b. TIME COVERED FROM _____ TO _____	14. DATE OF REPORT (Year, Month, Day) 93 March	15. PAGE COUNT 52
16. SUPPLEMENTARY NOTATION			
17. COSATI CODES		18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number)	
FIELD	GROUP	Electromagnetic compatibility, test and evaluation, aeromedical equipment	
01	03		
14	02		
19. ABSTRACT (Continue on reverse if necessary and identify by block number) The Medical Data Electronics escort patient monitor, model E300A, was tested for environmental and electromagnetic interference/compatibility in the UH-60A helicopter under the U.S. Army Program for Testing and Evaluation of Equipment for Aeromedical Operations. The tests were conducted using current military and industrial standards and procedures for environmental tests and electromagnetic interference/compatibility and human factors. The Medical Data Electronics escort patient monitor, model E300A, did not perform properly in the vibration test, high temperature operation test, and low temperature operation test. Excess electromagnetic emissions and electromagnetic susceptibility make it ineligible for an airworthiness release. The unit is not considered compatible with the U.S. Army Medical Evacuation UH-60A Black Hawk.			
20. DISTRIBUTION/AVAILABILITY OF ABSTRACT <input checked="" type="checkbox"/> UNCLASSIFIED/UNLIMITED <input type="checkbox"/> SAME AS RPT <input type="checkbox"/> DTIC USERS		21. ABSTRACT SECURITY CLASSIFICATION Unclassified	
22a. NAME OF RESPONSIBLE INDIVIDUAL Chief, Scientific Information Center		22b. TELEPHONE (include Area Code) (205) 255-6907	22c. OFFICE SYMBOL SGRD-UAX-SI

Table of contents

SECTION	PAGE
1. <u>EXECUTIVE DIGEST</u>	
1.1 Test objectives	1-1
1.2 Testing authority	1-2
1.3 Scope	1-2
1.4 Material description	1-2
1.5 Summary	1-3
1.6 Conclusion	1-4
2. <u>SUBTESTS</u>	
2.1 Initial inspection	2-1
2.2 Battery life evaluation (laboratory)	2-1
2.3 Electrical safety evaluation	2-2
2.4 Human factors evaluation (laboratory)	2-3
2.5 Altitude (low pressure) test	2-3
2.6 Vibration test	2-4
2.7 High temperature test	2-5
2.8 Low temperature test	2-7
2.9 Humidity test	2-8
2.10 Electromagnetic characteristics test	2-9
2.11 In-flight human factors evaluation	2-12

3. SUPPORTING DOCUMENTATION

3.1 Detailed test information 3-1

3.2 Test data 3-2

3.3 Criteria, significant problems, and suggested improvements 3-28

3.4 References 3-31

3.5 Abbreviations 3-32

3.6 List of manufacturers. 3-34

3.7 Distribution list 3-35

Accession For	
NTIS CRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A-1	

DTIC QUALITY INSPECTED 8

Section 1. Executive digest

The Army program for Test and Evaluation of Aeromedical Equipment uses existing military standards (MIL-STD) and collective professional expertise to test and evaluate selected medical equipment proposed for use aboard Army aircraft. Equipment meeting these standards ensures the safety of the crew, patients, and aircraft by eliminating risks due to: (1) Interference by the medical equipment with aircraft systems/subsystems operation, (2) the aircraft systems' interference with the operation of the medical equipment, (3) the medical equipment's susceptibility to environmental exposure, or (4) physical and/or functional incompatibility while in use on board selected rotary-wing aircraft. This program tests both developmental and nondevelopmental (off the shelf) medical equipment destined for use aboard Army medical evacuation aircraft.

1.1 TEST OBJECTIVES

- 1.1.1 To determine if the medical equipment is complete and operational per the manufacturer's operating instructions.
- 1.1.2 To ensure the electrical safety of the medical equipment.
- 1.1.3 To ensure the equipment will function as designed throughout the rated battery operation time.
- 1.1.4 To ensure the safety of the operator, the patient, and the aircrew.
- 1.1.5 To assess design considerations which potentially could contribute to an operator error.
- 1.1.6 To determine if the medical equipment can function as designed in a low pressure environment.
- 1.1.7 To determine the ability of the medical equipment to withstand the vibrational stresses expected in a rotary-wing flight environment without degradation or malfunction.
- 1.1.8 To determine the ability of the medical equipment to be stored and operated in a high-temperature environment.
- 1.1.9 To determine the ability of the medical equipment to be stored and operated in a low-temperature environment.
- 1.1.10 To determine the ability of the medical equipment to operate satisfactorily for short periods during exposure to high humidity conditions.

1.1.11 To assess the levels of electromagnetic emissions produced by the medical equipment within selected frequency ranges.

1.1.12 To assess the minimum electromagnetic susceptibility levels of the medical equipment within selected frequency ranges.

1.1.13 To assess the physical and/or functional compatibility of the medical equipment while in use on board the aircraft.

1.1.14 To assess the electromagnetic interference (EMI) and electromagnetic compatibility (EMC) characteristics of the medical equipment with the host aircraft and its installed systems.

1.2 TESTING AUTHORITY

Research and Technology Work Unit Summary, dated 5 October 1989. Project number 3M463807D836, titled, Army Program for Testing and Evaluation of Equipment for Aeromedical Operations.

1.3 SCOPE

1.3.1 This test was conducted at the United States Army Aeromedical Research Laboratory (USAARL) at Fort Rucker, Alabama.

1.3.2 The USAARL UH-60A aircraft, serial number 88-26069, with subsystems delineated in paragraph 3.2.2, was available for flight test of the MDE® patient monitor*, model E300A. The unit was not approved for flight testing.

1.3.3 The following tests were accomplished: Acceptance inspection, equipment training, electromagnetic compatibility, human factors and safety, and environmental compatibility.

1.4 MATERIAL DESCRIPTION

The MDE® model E300A patient monitor is electrically driven and designed to measure and monitor a patient's electrocardiograph (ECG), noninvasive blood pressure (NIBP), respiration (RESP) and temperature (TEMP). Six soft keys (membrane switches) on the bottom front panel of the patient monitor allow the operator to select the desired function. Each function's membrane switch has menu driven screens for adjustments and selections which are displayed on the bottom of the backlight liquid crystal display (LCD). Each menu and submenu screen allows the operator to select parameters and other menus ("NEXT PAGE") or

* See list of manufacturers

return to the original menu ("PAGE HOME"). Parameter measurements and corresponding wave forms, if selected, are displayed on the LCD. Each measured parameter has a separate alarm indicator and warning setting which are selected via the front panel membrane switches. If activated, an "ALARM SUSPEND" button on the front panel will suspend all alarms for 180 seconds in the adult mode and 90 seconds in the neonatal mode. The side panel of the E300A houses the patient input connectors for the ECG leads, NIBP, SAO2, BP, and TEMP. The rear of the patient monitor incorporates the ac receptacle, on/off ac line switch, ac voltage select switch, fuse holder, auxiliary output, two rechargeable batteries, and a defibrillator interface.

1.5 SUMMARY

1.5.1 Laboratory testing

1.5.1.1 Battery life evaluation: The MDE model E300A was operated on a fully-charged battery and connected to a Valmedix ECG simulator, at room temperature (75° F and 51% RH). The fully-charged MDE® model E300A averaged 2 hours and 39 minutes of operation. The battery is rated for up to a maximum of 3 hours of operation.

1.5.1.2 Electrical safety evaluation: All measurements were within acceptable limits. No unsafe qualities were found in the MDE® model E300A. The limits for currents and resistances were in accordance with (IAW) the limits specified in TB-38-750-2, April 1987.

1.5.1.3 Human factors evaluation: The MDE® model E300A was found to be satisfactory in all major categories of the evaluation criteria with one exception.

1.5.1.4 Environmental tests: The MDE® model E300A may be unsatisfactory for use in certain environmental conditions. Its performance was found to be unsatisfactory in stages of the environmental testing. The MDE® E300A patient monitor performance was found to be unsatisfactory during the vibration, high temperature operation, and low temperature operation tests. During vibration on the Z-axis, the patient monitor malfunctioned internally and would not operate. During the high temperature operation test, the unit's backlight LCD screen failed. The unit was sent back to the manufacturer for repairs on both occasions. At the manufacturer's request, the low temperature operation test was not performed. The unit will not operate in an environment below 5°C. The requirements for environmental tests are established in MIL-STD-810E, methods 500.2 (altitude), 514.3 (vibration), 501.2 (high temperature), 502.2 (low temperature), and 507.2 (humidity).

1.5.1.5 Radiated emissions tests (RE02): The MDE® model E300A may be unsatisfactory for use in certain EMI sensitive environments. Narrowband (NB) and broadband (BB) radiated emissions were detected in the test frequency ranges. Some narrowband and broadband emissions exceeded the test limits. Emission limits are set forth in MIL-STD-461C, Notice 5.

1.5.1.6 Radiated susceptibility test (RS03): The MDE® model E300A was found to be susceptible to radio frequency interference in the testing range and magnitude.

1.5.1.7 Conducted emissions test (CE01, CE02, and CE04): Narrowband and broadband emissions exceeding specification limits were detected.

1.5.1.8 Conducted susceptibility test (CS02 and CS06): The MDE model E300A was not found to be susceptible to the test-generated signals and spikes on its power lines.

1.5.2 In-flight testing

In-flight testing of the MDE® model E300A patient monitor could not be accomplished. As a result of excessive failures in laboratory environmental (vibration, high temperature, and low temperature) and electromagnetic compatibility (RE02 and RS03) a request for airworthiness release from the U.S. Army Aviation Systems Command will not be pursued.

1.6 CONCLUSION

The MDE® model E300A did not perform properly in the vibration test, high temperature operation test, and low temperature operation test. This, combined with excess electromagnetic emissions and electromagnetic susceptibility, make it ineligible for an airworthiness release. The unit is not considered compatible with the U.S. Army medical evacuation UH-60A Black Hawk with the subsystems listed in paragraph 3.2.2.

Section 2. Subtests

2.1 INITIAL INSPECTION

2.1.1 Objective

To determine if the MDE® model E300A is complete and operational for testing per the manufacturer's operating instructions.

2.1.2 Criteria

2.1.2.1 The physical inventory is conducted solely for investigation and documentation.

2.1.2.2 The MDE® model E300A will display consistent and accurate measurements as an acceptable performance test.

2.1.3 Test procedure

2.1.3.1 A complete physical inventory of the MDE® model E300A was completed per the manufacturer's equipment list.

2.1.3.2 An operational validation test of the MDE® model E300A was conducted per the manufacturer's operating instructions by USAARL's medical maintenance personnel.

2.1.4 Test findings

2.1.4.1 The MDE® model E300A was inventoried and found to be complete.

2.1.4.2 The MDE® model E300A operated as prescribed in the manufacturer's operating manual. Criteria met.

2.2 BATTERY LIFE EVALUATION (laboratory)

2.2.1 Objective

To ensure the equipment will function as designed throughout the rated battery operation time.

2.2.2 Criterion

Verify manufacturer's specified full power battery life expectancy of up to 3 hours operation.

2.2.3 Test procedure

2.2.3.1 Charging and operation cycles were conducted in ambient room conditions.

2.2.3.2 The MDE® model E300A was operated continuously using its fully-charged battery until a low battery indication occurred. The depletion time was noted and the battery was recharged. This procedure was repeated three times.

2.2.4 Test findings

The test was conducted using the fully-charged battery. The average operating time in testing was 2 hours and 39 minutes at room temperature. Criterion met.

2.3 ELECTRICAL SAFETY EVALUATION

2.3.1 Objective

To ensure the electrical safety of the MDE® model E300A by evaluation of case-to-ground resistance and case-to-ground current leakage.

2.3.2 Criterion

The MDE® model E300A shall meet the standards established in Technical Bulletin Number 38-750-2, April 1987 for electrical safety of medical equipment.

2.3.3 Test procedure

Measurements in the electrical safety evaluation were made with a Neurodyne-Dempsey model 431F electrical safety analyzer*, IAW the procedures described in Technical Bulletin (TB) Number 38-750-2. Case-to-ground resistance and various case-to-ground leakage currents were measured. Checks were made for safety concerns such as case integrity, breaks in power cord insulation, and connectors.

2.3.4 Test findings

Grounding conductor resistance was 80.2 milliohms and maximum case leakage current was 28.5 microamperes. These measurements are below the limits specified in TB 38-750-2 and NAFP 99. Criterion met.

2.4 HUMAN FACTORS EVALUATION (laboratory)

2.4.1 Objectives

2.4.1.1 To assure the safety of the operator, the potential patient, and the aircrew.

2.4.1.2 To assess the design considerations which could potentially contribute to an operator error.

2.4.2 Criterion

The MDE® model E300A must be rated satisfactory in all major categories of the evaluation. These include visual displays, controls, maintainability, conductors, fasteners, test points, test equipment, fuses and circuit breakers, labels and coding, and safety.

2.4.3 Test procedure

2.4.3.1 The evaluation was conducted in a laboratory under fluorescent lighting and ambient room conditions.

2.4.3.2 The MDE® model E300A was operated according to prescribed instructions through its full range of functions.

2.4.4 Test finding

The MDE® model E300A was found to be satisfactory in all of the evaluation criteria Criterion met.

2.5 ALTITUDE (LOW PRESSURE) TEST [IAW MIL-STD-810E, METHOD 500.2]

2.5.1 Objective

To determine if the MDE® model E300A can function as designed in a low-pressure environment.

2.5.2 Criterion

The MDE® model E300A will display consistent and accurate measurements while exposed to an altitude equivalency of 15,000 feet above sea level.

2.5.3 Test procedure

2.5.3.1 A pretest performance check was conducted to ensure proper operation of the MDE® model E300A.

2.5.3.2 The altitude test was performed in a Tenney Engineering model 64S altitude chamber*. This test is based on MIL-STD-810E, Method 500.2. The MDE® model E300A was operating and placed on the floor of the chamber. Chamber pressure was decreased to 420 mmHg (15,000 ft equivalent altitude) over a 15-minute period, held constant for 60 minutes, then raised, at 1500 fpm, to ambient conditions (760 mmHg) over a 10-minute period. There were no provisions for the control of temperature or humidity inside this chamber.

2.5.3.3 A posttest performance check was conducted to ensure proper operation of the MDE® model E300A after the exposure to low pressure.

2.5.4 Test findings

2.5.4.1 The pretest performance check met criterion 2.1.2.2.

2.5.4.2 The posttest performance check met criterion 2.1.2.2.

2.6 VIBRATION TEST [IAW MIL-STD-810E, METHOD 514.3]

2.6.1 Objective

To determine the ability of the MDE® model E300A to withstand the vibrational stresses expected in a rotary-wing environment without degradation or malfunction.

2.6.2 Criterion

The MDE® model E300A will remain operational and be able to display consistent and accurate measurements while exposed to vibrational stresses.

2.6.3 Test procedure

2.6.3.1 A pretest performance check was conducted to ensure proper operation of the MDE® model E300A.

2.6.3.2 The vibration test was performed using an Unholtz-Dickey model TA115-40/CSTA vibration test system*. It is a single-axis system with an electromagnetic driver unit. The test consisted of sinusoidal vibrations superimposed on random vibrations over a frequency range of 500 Hz, as shown below. These vibrations are derived from measurements taken on the floor under the copilot's seat in a UH-1 helicopter traveling at 120 knots. The reference spectrum breakpoints are from MIL-STD-810E, Method 514.3; Independent tests were conducted in the X-, Y-, and Z-axes.

X-, Y-, and Z-axes

duration: 60 minutes
broadband intensity: 2.116 G_{rms}
random vibration: initial slope : 99.00 dB/oct
5 Hz level: 0.00007 $G_{sqr/Hz}$
100 Hz level: 0.0007 $G_{sqr/Hz}$
300 Hz level: 0.0007 $G_{sqr/Hz}$
500 Hz level: 0.00007 $G_{sqr/Hz}$
final slope: -99.00 dB/oct

sinusoidal vibration:

0.9500 G_{pk} at 17.00 Hz
1.4000 G_{pk} at 34.00 Hz
0.9000 G_{pk} at 51.00 Hz
0.9000 G_{pk} at 68.01 Hz

The MDE® model E300A was strapped to the vibration table fixture, and its performance was evaluated before, during, and after exposure to vibration.

2.6.3.3 A posttest performance check was conducted to ensure proper operation of the MDE® model E300A.

2.6.4 Test findings

2.6.4.1 The pretest performance check met criterion 2.1.2.2.

2.6.4.2 Failures in the performance of the MDE® model E300A occurred during exposure to vibration. The E300A internally malfunctioned and would not operate. The unit was sent back to the manufacturer for repairs. Criterion not met.

2.6.4.3 The posttest performance check met criterion 2.1.2.2.

2.7 HIGH TEMPERATURE TEST [IAW MIL-STD-810E, METHOD 501.2]

2.7.1 Objective

To determine the ability of the MDE® model E300A to be stored and operated in a high-temperature environment.

2.7.2 Criteria

2.7.2.1 The MDE® model E300A will display consistent and accurate measurements during the high-temperature operation check.

2.7.2.2 The MDE® model E300A will display consistent and accurate measurements after the high-temperature storage cycle.

2.7.3 Test procedure

2.7.3.1 A pretest performance check was conducted to ensure proper operation of the MDE® model E300A.

2.7.3.2 The high-temperature test was conducted in a Tenney Engineering model ZWUL-10107D walk-in controlled environment chamber*. This test is based on MIL-STD-810E, Method 501.2. For the high-temperature operation test, the MDE® model E300A was operating and connected to a Valmedix ECG simulator and placed on the floor of the environmental chamber. The chamber temperature was raised to 49°C and the humidity was stabilized at a maximum of 20 percent RH within 15 minutes. The environmental control system is capable of regulating temperature within $\pm 2^\circ\text{C}$ and humidity within ± 5 percent RH. Temperature and humidity were held constant for 2 hours. At 30-minute intervals, the chamber door was opened briefly to minimize the change in chamber conditions during performance checks. After the operational test, the MDE® model E300A was allowed to return to ambient conditions over a 30-minute period.

2.7.3.3 A posttest performance check was conducted to ensure proper operation of the MDE® model E300A.

2.7.3.4 The MDE® model E300A was stored (not operated) at temperatures of 63°C for 1 hour, 71°C for 4 hours, then again at 63°C for 1 hour. The chamber and MDE® model E300A then were returned to ambient conditions over a 30-minute period.

2.7.3.5 A poststorage performance check was conducted to ensure proper performance of the MDE® model E300A.

2.7.4 Test findings

2.7.4.1 The pretest performance check met criterion 2.1.2.2.

2.7.4.2 The MDE® model E300A failed to function correctly during the high temperature operation test. The backlight LCD failed and the measured parameters could not be monitored. The unit was sent back to the manufacturer for repairs. Criterion not met.

2.7.4.3 The posttest performance check did not meet criterion 2.1.2.2.

2.7.4.4 The MDE® model E300A functioned properly after the high-temperature storage test. Criterion met.

2.8 LOW TEMPERATURE TEST [IAW MIL-STD-810E, METHOD 502.2]

2.8.1 Objective

To determine the ability of the MDE® model E300A to be stored and operated in a low-temperature environment.

2.8.2 Criteria

2.8.2.1 The MDE® model E300A will display consistent and accurate measurements during the low-temperature operation check.

2.8.2.2 The MDE® model E300A will display consistent and accurate measurements after the low-temperature storage cycle.

2.8.3 Test procedure

2.8.3.1 A pretest performance check was conducted to ensure proper operation of the MDE® model E300A.

2.8.3.2 The MDE® model E300A was not evaluated in the low temperature operation test at request from the manufacture. The manufacturer indicates the unit will not operate in an environment below 5°C.

2.8.3.3 A posttest performance check was conducted to ensure proper operation of the MDE® model E300A.

2.8.3.4 The MDE® model E300A was "stored" in a nonoperational mode. The MDE® model E300A was placed on the floor of the environmental test chamber and the temperature was lowered to -46°C for 6 hours. The chamber then was raised to ambient temperature over a 30-minute period.

2.8.3.5 A poststorage performance check was conducted to ensure proper operation of the MDE® model E300A.

2.8.4 Test findings

2.8.4.1 The pretest performance check met criterion 2.1.2.2.

2.8.4.2 The unit was not evaluated in the low temperature operation test at request from the manufacture. Criterion not met.

2.8.4.3 The posttest performance check did not meet criterion 2.1.2.2.

2.8.4.4 The MDE® model E300A functioned properly after the low-temperature storage test. Criterion met.

2.9 HUMIDITY TEST [IAW MIL-STD-810E, METHOD 507.2]

2.9.1 Objective

To determine the ability of the MDE® model E300A to operate satisfactorily for short periods of time during exposure to highly humid conditions.

2.9.2 Criterion

The MDE® model E300A will display consistent and accurate measurements while exposed to a high humidity environment.

2.9.3 Test procedure

2.9.3.1 A pretest performance check was conducted to ensure the proper operation of the MDE® model E300A.

2.9.3.2 The humidity test was conducted in a Tenney Engineering model ZWUL-10107D walk-in controlled environment chamber*. This test is based on MIL-STD-810D, Method 507.2. For the humidity test, the MDE® model E300A was placed in operation on the floor of the environmental chamber. The chamber temperature was raised to a temperature of 30°C and a relative humidity of 95 percent within 25 minutes. Temperature and relative humidity were maintained for 4 hours. The environmental control system is capable of regulating temperature within $\pm 2^\circ\text{C}$ and humidity within ± 5 percent RH. At 45-minute intervals, the performance of the BP monitor was checked. The chamber door was opened briefly to minimize the change in chamber conditions. The chamber and the MDE® model E300A were returned to ambient conditions before the posttest performance validation check was conducted.

2.9.3.3 A posttest performance check was conducted to ensure the proper operation of the MDE® model E300A.

2.9.4 Test findings

2.9.4.1 The pretest performance check met criterion 2.1.2.2.

2.9.4.2 The MDE® model E300A was not adversely affected by the high humid environment. Criterion met.

2.9.4.3 The posttest performance check met criterion 2.1.2.2.

2.10 ELECTROMAGNETIC CHARACTERISTICS TEST [IAW MIL-STD-461C, Notice 5, and MIL-STD-462, Notice 3]

2.10.1 Objectives

2.10.1.1 To assess the maximum levels of radiated electromagnetic emissions produced by the MDE® model E300A in the 14 kHz to 12.4 GHz frequency range.

2.10.1.2 To assess the tolerances of radiated electromagnetic susceptibility of the MDE® model E300A within the 10 kHz to 10 GHz electric field.

2.10.1.3 To assess the maximum levels of conducted electromagnetic emissions produced by the MDE® model E300A in the 10 kHz to 50 MHz frequency range.

2.10.1.4 To assess the tolerances of conducted electromagnetic susceptibility of the MDE® model E300A within the range of 50 kHz to 400 MHz and power spikes.

2.10.2 Criteria

2.10.2.1 The MDE® model E300A will not produce emissions in excess of the limits set forth in MIL-STD-461C, Notice 5, paragraph 6.13.

2.10.2.2 The MDE® model E300A will not malfunction when it is subjected to radiated emissions as specified in MIL-STD-461C, Notice 5, paragraph 6.20.

2.10.2.3 The MDE® model E300A shall not conduct emissions in excess of the limits set forth in MIL-STD-461C, Notice 5, paragraphs 6.1 and 6.2.

2.10.2.4 The MDE® model E300A shall not malfunction when it is subjected to conducted emissions as specified in MIL-STD-461C, Notice 5, paragraphs 6.7 and 6.10.

2.10.3 Test procedure

2.10.3.1 The radiated emissions test was performed according to MIL-STD-462, Notice 3, Method RE02. The MDE® model E300A was positioned on a wooden test stand inside the EMI chamber, 1 meter away from the receiving antennas. The antennas were mounted for both vertical and horizontal polarities and connected to EMI receivers. While the MDE® model E300A was operating, the frequency spectrum (14 kHz to 12.4 GHz) was scanned for emissions. The MDE® model E300A was operated with both ac and battery power.

2.10.3.2 The radiated susceptibility test was performed according to MIL-STD-462, Notice 3, Method RS03. The MDE® model E300A was positioned on a wooden test stand inside the EMI chamber 1 meter away from the transmitting antennas. The antennas were mounted for both vertical and horizontal polarities and connected to radio frequency (RF) transmitters. While the MDE® model E300A was operating, it was monitored for faulty operation during exposures to fields of 20 V/m from 10 kHz to 10 GHz. The MDE® model E300A was operated with ac and battery power.

2.10.3.3 The conducted emissions tests were performed according to MIL-STD-462, Notice 3, Methods CE02 and CE04. The MDE® model E300A was placed on a grounded, copper-covered workbench. The top of the workbench was 1 meter from floor level, 1.37 meters long and 0.81 meters wide. Power was supplied via a pair of line impedance stabilization networks (LISN) and a test jig. The test jig is a wooden tray with two power receptacles and two slots to hold current probes in place around power supply conductors. While the MDE® model E300A was operating, the frequency range (10 kHz to 50 MHz) was scanned for emissions conducted in the power cable.

2.10.3.4 The conducted susceptibility spike test was performed according to MIL-STD-462, Notice 3, Method CS06, on a chemical resistant counter top. Power was supplied via a customized metal connection box. The connection box has two power receptacles and four banana jacks on its front panel. Connections to the individual power lines are made in series through the banana jacks. Transient spikes of 200 volts, 10 microseconds were generated with a Solar Electronics model 8282-1 transient pulse generator* and induced onto the power leads at the connection box banana jacks. The spikes were monitored with a Tektronix 2235 oscilloscope* connected to a power receptacle on the connection box. The MDE® model E300A was plugged into the other receptacle on the connection box and placed in operation. It was observed for correct operation and visual displays while it was subjected to the power line spikes.

2.10.3.5 The conducted susceptibility test was performed according to MIL-STD-462, Notice 3, Method CS02. The MDE® model E300A was placed on a grounded, copper-covered workbench. Radio frequency interference was induced on the power leads and measured at the power cable. The frequency of the interference was incremented over the 50 kHz to 400 MHz range while the MDE® model E300A was operated. It was observed for correct operation and visual displays while it was subjected to the radio interference on the power leads. Each frequency was held for 15 seconds.

2.10.4 Test findings

2.10.4.1 During the radiated emissions test, emissions which exceeded specification limits of MIL-STD-461C, Notice 5, were detected. These included:

ac operating failure data:

<u>Frequency (MHz)</u>	<u>Amount of failure (dB)</u>
NB: 0.024 - 363.037	0.5 - 34.6
BB: 0.016 - 3.131	0.2 - 20.9
BB: 25.609	9.3
BB: 30.000	11.8

Battery operating failure data:

<u>Frequency (MHz)</u>	<u>Amount of failure (dB)</u>
NB: 0.126 - 425	0.1 - 35.0
BB: 0.016 - 1.688	0.1 - 21.5
BB: 30.0	10.1

Criterion not met.

2.10.4.2 The MDE® model E300A was found to be susceptible to radio frequency interference in the testing range and magnitude.

ac operating failure data:

<u>Frequency (MHz)</u>	<u>Threshold of failure (V/m)</u>
10.0 - 14.2	2.12 - 10.0
29.6 - 38.5	1.89 - 7.10
69.1	6.32
9.7 - 145.6	1.5 - 15.0

Battery operating failure data:

<u>Frequency (MHz)</u>	<u>Threshold of failure (V/m)</u>
10.0 - 14.0	4.75 - 12.6
16.2 - 25.6	6.32 - 15.0
30.0 - 145.6	1.12 - 11.9
174.5	1.68

Criterion partially met.

2.10.4.3 The MDE® model E300A did not produce conducted emissions in excess of the limits. Criterion met.

2.10.4.4 The MDE® model E300A was not susceptible to the test-generated signals or test-generated spikes produced on the power lines. Criterion met.

2.11 IN-FLIGHT HUMAN FACTORS EVALUATION

2.11.1 Objective

To assess the physical and/or functional compatibility of the MDE® model E300A while in use on board the aircraft.

2.11.2 Criterion

The flight surgeon will be able to operate the MDE® model E300A without physical or functional restrictions aboard the aircraft. Major areas of concern include: proper operation, visual displays, controls, maintainability, conductors, fasteners, test points, test equipment, fuses and circuit breakers, labels and coding, and safety.

2.11.3 Test procedure

Not applicable since unit not available for in-flight testing.

2.11.4 Test findings

Not available since unit not available for in-flight testing. Criterion not evaluated.

Section 3. Supporting documentation

3.1 DETAILED TEST INFORMATION

3.1.1 General information

3.1.1.1 MDE® model E300A testing is not considered a major action significantly affecting the quality of the human environment and, therefore, qualifies for categorical exclusion A-28, appendix A, AR 200-1.

3.1.1.2 A safety pilot will be designated for each flight. Flight operations will be conducted IAW the aircraft operator's manual, appropriate aircrew training manuals, and test item technical data.

3.1.2 Material description

3.1.2.1 The MDE® model E300A patient monitor is electrically driven and designed to measure and monitor a patient's electrocardiograph (ECG) noninvasive blood pressure (NIBP), respiration (RESP) and temperature (TEMP). Six soft keys (membrane switches) on the bottom front panel of the patient monitor allow the operator to select the desired function. Each function's membrane switch has menu-driven screens for adjustments and selections which are displayed on the bottom of the backlight liquid crystal display (LCD). Each menu and submenu screen allows the operator to select parameters and other menus ("NEXT PAGE") or return to the original menu ("PAGE HOME"). Parameter measurements and corresponding wave forms, if selected, are displayed on the LCD. Each measured parameter has a separate alarm indicator and warning setting which are selected via the front panel membrane switches. If activated, an "ALARM SUSPEND" button on the front panel will suspend all alarms for 180 seconds in the adult mode and 90 seconds in the neonatal mode. The side panel of the E300A houses the patient input connectors for the ECG leads, NIBP, SAO2, BP, and TEMP. The rear of the patient monitor incorporates the ac receptacle, on/off ac line switch, ac voltage select switch, fuse holder, auxiliary output, two rechargeable batteries, and a defibrillator interface.

3.1.2.2 Dimensions: 7 x 7 x 8.5 in. (17.78 x 17.78 x 21.59 cm.)

3.1.2.3 Weight: 14 lb (6.35 kg).

3.1.2.4 Power requirements: 110 Vac, 60 Hz, 0.5 amp, 3-wire grounded system. External rechargeable batteries provide up to 3 hours of rated operation. Charge time for a low battery to a full charge approximately 14 hours with the unit off. Power cord is 18-3, type SJT, PE CO, of approximately 119 inches in length.

3.2 TEST DATA

3.2.1 Photographic description

Photo not available

3.2.2 Aircraft equipment list

Item No.	Nomenclature
1	Receiver radio -- R-1496A/ARN-89 (automatic direction finder)
2	Displacement gyro -- CN-1314/A
3	Gyro directional -- CN-998/ASN-43
4	Signal data converter -- CV-3338/ASN-128
5	Receiver -- R-2139/ARN-123 (VOR/LOC/MB/GS)
6	Command instrument system processor -- 70600-01038-101
7	SAS amplifier -- 70901-02908-104 (flight control stability augmentation system)
8	Rate gyro -- TRU-2A/A
9	Amplifier, impedance -- AM-4859A/ARN-89
10	Cargo hook -- FE-7590-145
11	Receiver, radar -- RT-1193/ASN-128 (doppler navigation receiver)
12	Barometric altimeter -- AAU-31/A-1
13	Barometric altimeter -- AAU-32A
14	Receiver/transmitter -- RT-1300/ARC-186 (VHF-AM and/or FM radio)
15	UHF-AM radio set -- RT-1518/ARC-164
16	Interphone control -- C6533/ARC (aircraft intercom control)
17	Receiver/transmitter -- RT-1115D/APN-209 (radar altimeter)
18	Indicator altimeter -- ID-1917C/APN-209 (radar altimeter)
19	Control radio set -- C-7392A/ARN-89 (automatic direction finder)
20	Comparator signal data -- CM-482/ARC-186 (comparator for ARC-186)
21	Receiver/transmitter -- RT-1296A/APX-100 (transponder with IFF)
22	Computer display unit -- CP-1252/ASN-128 (doppler navigation system)
23	Compass set infusion pump -- C-8021E/ASN75
24	Magnetic compass - standby -- MS-17983-4

3.2.3 Battery life evaluation

Battery life evaluation
report form

Nomenclature: Patient monitor
Manufacturer: MDE[®] Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Manufacturer battery life specification: Up to 3 hours

Specified battery recharge time: 14 hours (full charge).

Specified mode of operation under battery power: ECG connected to simulator. SA02 and NIBP connected to the test engineer.

Overall performance: Pass

Measurements:

Dates of first test: 9/29/92
Temperature: 23°C
Humidity: 51% RH
Start time: 1226
End time: 1507
Operating time: 2 hours 41 minutes

TOTAL OPERATING TIME: 2 hours 41 minutes
PERFORMANCE: Pass

Dates of second test: 9/30/92
Temperature: 24°C
Humidity: 46% RH
Start time: 0730
End time: 1007
Operating time: 2 hour 37 minutes

TOTAL OPERATING TIME: 2 hour 37 minutes
PERFORMANCE: Pass

Dates of third test: 10/1/92
Temperature: 22°C
Humidity: 46% RH
Start time: 0725
End time: 1005
Operating time: 2 hours 40 minutes

TOTAL OPERATING TIME: 2 hours 40 minutes
PERFORMANCE: Pass

Comments: The unit averaged 2 hours and 39 minutes of operation.

3.2.4 Electrical safety test

Electrical safety test report form

Nomenclature: Patient monitor
Manufacturer: MDE® Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: 9 March 1992

Measurements:

Grounding conductor resistance (milliohms): 80.5

Leakage current - case to ground (microamperes):

unit off, grounded, normal polarity	0.1
unit off, ungrounded, normal polarity	0.1
unit off, ungrounded, reverse polarity	28.5
unit on, grounded, normal polarity	0.1
unit on, ungrounded, normal polarity	21.2
unit on, ungrounded, reverse polarity	28.5

MAXIMUM LIMITS:

ground resistance (milliohms):	150
current (grounded, type A unit):	10
current (ungrounded, type A unit):	100
current (grounded, type B unit):	50
current (ungrounded, type B unit):	500

Comments on item setup or checks: None

Comments on test run (including interruptions): None

Comments on other data: None

Additional measurements for equipment with patient leads

EUT always ON, set for Lead I

Current (μ A): All leads to ground
(G, P)0.3 (UG, NP)6.8 (UG, RP)9.1

Current (μ A): Right Arm to ground
(G, NP)2.8 (UG, NP)7.6 (UG, RP)9.4

Current (μ A): Left Arm to ground
(G, NP)2.8 (UG, NP)7.8 (UG, RP)9.8

Current (μ A): Right Leg to ground
(G, NP) NA (UG, NP) NA (UG, RP)NA

Current (μ A): Left Leg to ground
(G, NP)0.3 (UG, NP)6.4 (UG, RP)8.6

Comments on item setup or checks: None

Comments on test run (including interruptions): None

Comments on other data: None

3.2.5 Human factors evaluation

Human factors evaluation report form

Nomenclature: Patient monitor
Manufacturer: HDE® Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: 23 November 1992

Item configuration during test: Item prepared for operation,
sitting on a counter top and connected to simulators.

Checklist for HFE

RESULTS

VISUAL DISPLAYS:

Satisfactory

display type, format, content
location of displays
indicator lights
scalar displays
color coding
legends and labels
cathode ray tubes
counters
flags, go-no-go, center-null indicators

Comments:

CONTROLS:

Satisfactory

location
characteristics of controls
labeling
control - display relationships

Comments: None

TIME REQUIRED TO PREPARE FOR OPERATION (list in comment)

Comments: 2.5 to 3 minutes required to prepare for operation.

MAINTAINABILITY:

Satisfactory

component location
component characteristics
rests and stands
covers, cases, access doors
handles
lubrication
component mounting
cord storage provisions
external accessibility
internal accessibility
list special tools required
list realistic inspection requirements
list realistic inspection intervals

Comments: There is an internal self-test with each power-up. There is no lubrication or need for internal accessibility.

CONDUCTORS:

Satisfactory

binding and securing
length
protection
routing
conductor coding
fabrication
connectors

Comments: None

FASTENERS:

Satisfactory

access through inspection panel covers
enclosure fasteners
device mounting bolts and fasteners

Comments: None

TEST POINTS:

NA

general
location and mounting
test point labeling and coding

Comments: No test points.

TEST EQUIPMENT:

Satisfactory

general
equipment self-test
indicators (list in comments)
controls
positive indication of proper operation

Comments: An internal self-test outputs codes for normal status, leads off, and for failures.

FUSES AND CIRCUIT BREAKERS:

Satisfactory

external accessibility
easy replacement or reset by operator

Comments: None

LABELS and CODING:

Satisfactory

placed above controls and displays
near or on the items they identify
not obscured by other equipment components
describe the function of the items they identify
readable from normal operating distance
conspicuous placards adjacent to hazardous items

Comments: None

SAFETY:

Satisfactory

manual
materials
fire and explosive protection
operator protection from mechanical hazards
patient protection from mechanical hazards
electrical safety (operator and patient)

Comments: None

3.2.6 Altitude test

Altitude test report form

Nomenclature: Patient monitor
Manufacturer: MDE[®] Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: 31 March 1992

Item configuration during test: Item operating, sitting on chamber floor connected to ECG simulator.

Performance test criteria: Consistent and accurate displays and measurements.

Ambient conditions outside chamber:

Temperature	76°F
Humidity	65% RH
Barometric pressure	1 atm

PRETEST DATA

Pretest performance check:

Item functional (based on performance test criteria): Yes

Installation of item in test facility:

list connections to power	None (battery)
list connections to simulators	Valmedix ECG
list connections to dummy loads	None
list unconnected terminals	ac

IN-TEST DATA

Time of test start: 0935

POSTTEST DATA

Posttest performance check (complete check of item and accessories):

Time of test end: 1100

Item functional (based on performance test criteria): Yes

Deviation from pretest : None

Comments on item setup or checks: None

Comments on test run (including interruptions):

Comments on other data:

3.2.7 Vibration test

Vibration test
report form

Nomenclature: Patient monitor
Manufacturer: MDE® Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: 30 March 1992

Item configuration during test: Item strapped down on vibration table fixture. Unit operating on ac and battery power.

Performance test criteria: Consistent and accurate displays and measurements.

PRETEST DATA

Pretest performance check:

Item functional (based on performance test criteria): Yes

Installation of item in test facility:

list connections to power	120 Vac
list connections to simulators	Valmedix ECG
list connections to dummy loads	None
list unconnected terminals	None

Ambient conditions

Temperature	80°F
Humidity	63% RH
Barometric pressure	1 atm

IN-TEST DATA

Data and performance checks during test:

Time at first check:

X: 1320 (3/30/92) Y: 1432 (3/30/92) Z: 1245 (3/30/92)

Item functional (based on performance test criteria): Yes

Deviation from pretest: None

Time at second check:

X: 1410

Y: 1522

Z: NA

Item functional (based on performance test criteria): No

Deviation from pretest: Unit malfunctioned during the z-axis shake. Unit would not operate and was sent back to the manufacturer for repairs.

POSTTEST DATA

Time at test end:

X: 1420

Y: 1532

Z: NA

Posttest performance check (complete check of item and accessories):

Item functional (based on performance test criteria): No

Item intact: Yes

Deviation from pretest: None

Comments on item setup or checks: The unit would not operate during or after the z-axis shake in each of the three test runs.

Comments on test run (including interruptions): The unit was retested in the z-axis on 23 November and 1 December 1993. During each retest, the unit failed and would not operate correctly. The unit was again sent back to the manufacturer for repairs.

Comments on other data: The unit was not tested after 1 December 1992 for vibration compatibility.

3.2.8 High temperature test

High temperature test (equipment operating) report form

Nomenclature: Patient monitor
Manufacturer: MDE® Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: 9 March 1992

Item configuration during test: Unit was sitting on the wire test stand, operated on ac and battery power.

Performance test criteria: Consistent and accurate displays and measurements.

Ambient conditions outside chamber:
Temperature 24°C
Humidity 53% RH
Barometric pressure 1 atm

PRETEST DATA

Pretest performance check :
Item functional (based on performance test criteria): Yes

Installation of item in test facility:

list connections to power	120 Vac
list connections to simulators	Valmedix ECG
list connections to dummy loads	None
list unconnected terminals	Temperature probe
distance from north wall (meters)	0.6223
distance from south wall (meters)	0.5334
distance from east wall (meters)	1.4986
distance from west wall (meters)	1.2446
distance from ceiling (meters)	1.3208
distance from floor (meters)	0.4826

IN-TEST DATA

Time of test start: 1050
Performance checks during test:

First check:

Time: 1120
Temperature: 49°C
Humidity: 15% RH
Barometric pressure: 1 atm
Item functional (based on performance test criteria):
all OK: Pass
Deviation from pretest: None

Second check:

Time: 1150
Temperature: 49°C
Humidity: 15% RH
Barometric pressure: 1 atm
Item functional (based on performance test criteria):
all OK: Pass
Deviation from pretest: None

Third check:

Time: 1220
Temperature: 49°C
Humidity: 15% RH
Barometric pressure: 1 atm
Item functional (based on performance test criteria):
all OK: Pass
Deviation from pretest: None

POSTTEST DATA

Posttest performance check:

(complete check of item and accessories)
Time of test end: 1330
Item functional (based on performance test criteria):
Fail
Deviation from pretest: None

Comments on item setup or checks: The illumination of the LED failed prior to the end of the test. The screen was unreadable.

Comments on test run (including interruptions): None

Comments on other data: None

3.2.9 High temperature storage test

High temperature test
(equipment in storage)
report form

Nomenclature: Patient monitor
Manufacturer: MDE® Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: 20 March 1992

Item configuration during test: Sitting on wire test stand in storage, not operating.

Performance test criteria: Consistent and accurate displays and measurements.

Ambient conditions outside chamber:

Temperature	22°C
Humidity	48% RH
Barometric pressure	1 atm

PRETEST DATA

Pretest performance check:

Item functional (based on performance test criteria): Yes

Installation of item in test facility:

list connections to power	None
list connections to simulators	None
list connections to dummy loads	None
list unconnected terminals	All
distance from north wall (meters)	0.6223
distance from south wall (meters)	0.5334
distance from east wall (meters)	1.4986
distance from west wall (meters)	1.2446
distance from ceiling (meters)	1.3208
distance from floor (meters)	0.4826

Time of test start: 0830

POSTTEST DATA

Posttest performance check:
(complete check of item and accessories)

Time of test end: 1430

Item functional (based on performance test criteria): Yes

Deviation from pretest: None

Comments on item setup or checks:

The unit was allowed to cool overnight, before the posttest performance check was conducted.

Comments on test run (including interruptions): None

Comments on other data: None

3.2.10 Low temperature test

Low temperature test
(equipment operating)
report form

Nomenclature: Patient monitor
Manufacturer: MDE[®] Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: NA

Item configuration during test: NA

Performance test criteria: NA

Ambient conditions outside chamber:

Temperature	NA
Humidity	NA
Barometric pressure	NA

PRETEST DATA

Pretest performance check:

Item functional (based on performance test criteria): The E300A patient monitor was not evaluated at the MIL-STD-810E temperature of -25°C at the manufacturer's request. The unit will not function correctly at a temperature below +5°C. This represents a failure of the low temperature operational test.

Installation of item in test facility:

- list connections to power
- list connections to simulators
- list connections to dummy loads
- list unconnected terminals
- distance from north wall (meters)
- distance from south wall (meters)
- distance from east wall (meters)
- distance from west wall (meters)
- distance from ceiling (meters)
- distance from floor (meters)

Time of test start: NA

Performance checks during test:

First check:

Time: NA
Temperature: NA
Humidity: NA
Barometric pressure: NA
Item functional (based on performance test criteria): NA
Deviation from pretest: NA

Second check:

Time: NA
Temperature: NA
Humidity: NA
Barometric pressure: NA
Item functional (based on performance test criteria): NA
Deviation from pretest: NA

Third check:

Time: NA
Temperature: NA
Humidity: NA
Barometric pressure: NA
Item functional (based on performance test criteria): NA
Deviation from pretest: NA

POSTTEST DATA

Posttest performance check:
(complete check of item and accessories)

Time of test end: NA
Item functional (based on performance test criteria): NA
Deviation from pretest:

Comments on item setup or checks: NA

Comments on test run (including interruptions): NA

Comments on other data: NA

3.2.11 Low temperature storage test

Low temperature test
(equipment in storage)
report form

Nomenclature: Patient monitor
Manufacturer: MDE® Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: 19 March 1992

Item configuration during test: Sitting on the wire test stand,
not operating, unit is in storage.

Performance test criteria: Consistent and accurate displays and
measurements

Ambient conditions outside chamber:

Temperature	23°C
Humidity	52% RH
Barometric pressure	1 atm

PRETEST DATA

Pretest performance check:

Item functional (based on performance test criteria): Yes

Installation of item in test facility:

list connections to power	None
list connections to simulators	None
list connections to dummy loads	None
list unconnected terminals	All
distance from north wall (meters)	0.6223
distance from south wall (meters)	0.5334
distance from east wall (meters)	1.4986
distance from west wall (meters)	1.2446
distance from ceiling (meters)	1.3208
distance from floor (meters)	0.4826

Time of test start:	0830
Midtest time:	1130
Midtest temperature:	-46°C

POSTTEST DATA

Posttest performance check:
(complete check of item and accessories)

Time of test end: 1430

Item functional (based on performance test criteria): Yes

Deviation from pretest: None

Comments on item setup or checks: None

Comments on test run (including interruptions): None

Comments on other data: The unit was allowed to return to
ambient conditions overnight before final performance check.

3.2.12 Humidity test

Humidity test
report form

Nomenclature: Patient monitor
Manufacturer: MDE® Inc.
Model number: E300A
Serial number: 1427
Military item number: None

Options installed: None

Date of test: 10 March 1992

Item configuration during test: Sitting on the wire test
stand, operating on ac and battery power.

Performance test criteria: Consistent and accurate displays and
measurements.

Ambient conditions outside chamber:

Temperature	22°C
Humidity	50% RH
Barometric pressure	1 atm

PRETEST DATA

Pretest performance check:

Item functional (based on performance test criteria): Yes

Installation of item in test facility:

list connections to power	120 Vac
list connections to simulators	Valmedix ECG
list connections to dummy loads	None
list unconnected terminals	Temperature probe
distance from north wall (meters)	0.6223
distance from south wall (meters)	0.5334
distance from east wall (meters)	1.4986
distance from west wall (meters)	1.2446
distance from ceiling (meters)	1.3208
distance from floor (meters)	0.4826

IN-TEST DATA

Time of test start: 1000

Performance checks during test:

First check:

Time: 1045
Temperature: 29.5°C
Humidity: 95% RH
Barometric pressure: 1 atm
Item functional (based on performance test criteria): Yes
Deviation from pretest: None

Second check:

Time: 1130
Temperature: 29.5°C
Humidity: 95% RH
Barometric pressure: 1 atm
Item functional (based on performance test criteria): Yes
Deviation from pretest: None

Third check:

Time: 1215
Temperature: 29.5°C
Humidity: 95% RH
Barometric pressure: 1 atm
Item functional (based on performance test criteria): Yes
Deviation from pretest: None

Fourth check:

Time: 1300
Temperature: 29.5°C
Humidity: 95% RH
Barometric pressure: 1 atm
Item functional (based on performance test criteria): Yes
Deviation from pretest: None

Fifth check:

Time: 1345
Temperature: 29.5°C
Humidity: 95% RH
Barometric pressure: 1 atm
Item functional (based on performance test criteria): Yes
Deviation from pretest: None

POSTTEST DATA

Posttest performance check:

(complete check of item and accessories)

Time of test end: 1430

Item functional (based on performance test criteria): Yes

Deviation from pretest: None

Comments on item setup or checks: None

Comments on test run (including interruptions): None

Comments on other data: None

3.2.13 Electromagnetic characteristics test

Electromagnetic characteristics testing
evaluation of performance

T & E item number: 14

Date: 23 March 1992

Nomenclature: Patient monitor

Manufacturer: MDE[®] Inc.

Model number: E300A

Serial number: 1427

Military item number: NA

Conducted emissions tests

CE01 Testing configuration(s): NA
Performance (pass/fail): NA

Comments: NA

CE02 Testing configuration(s): Operating on test bench,
connected to ECG simulator.
Performance (pass/fail): Pass

Comments: No emissions in excess of limits were
detected.

CE04 Testing configuration(s): Operating on test bench,
connected to ECG simulator.
Performance (pass/fail): Pass

Comments: No emissions in excess of limits were
detected.

Conducted susceptibility tests

CS02 Testing configuration(s): Operating on test bench,
connected to ECG simulator and the test jig.
Performance (pass/fail): Pass

Comments: Not susceptible to test generated sig-
nals.

CS06 Testing configuration(s): Operating on counter
top, connected to ECG simulator and to test jig.
Performance (pass/fail): Pass

Comments: Not susceptible to test spikes

Radiated emissions tests

RE02 Testing configuration(s): Operating on the wooden test stand in the EMC chamber, ac and battery power. Performance (pass/fail): Fail

Comments:

ac operating failure data:

<u>Frequency (MHz)</u>	<u>Amount of failure (dB)</u>
NB: 0.024 - 363.037	0.5 - 34.6
BB: 0.016 - 3.131	0.2 - 20.9
BB: 25.609	9.3
BB: 30.000	11.8

Battery operating failure data:

<u>Frequency (MHz)</u>	<u>Amount of failure (dB)</u>
NB: 0.126 - 425	0.1 - 35.0
BB: 0.016 - 1.688	0.1 - 21.5
BB: 30.0	10.1

Radiated susceptibility tests

RS03 Testing configuration(s): Operating on the wooden test stand in the EMC chamber, ac and battery power. Performance (pass/fail): Fail

Comments:

ac operating failure data:

<u>Frequency (MHz)</u>	<u>Threshold of failure (V/m)</u>
10.0 - 14.2	2.12 - 10.0
29.6 - 38.5	1.89 - 7.10
69.1	6.32
99.7 - 145.6	1.5 - 15.0

Battery operating failure data:

<u>Frequency (MHz)</u>	<u>Threshold of failure (V/m)</u>
10.0 - 14.0	4.75 - 12.6
16.2 - 25.6	6.32 - 15.0
30.0 - 145.6	1.12 - 11.9
174.5	1.68

3.3 CRITERIA, SIGNIFICANT PROBLEMS, AND SUGGESTED IMPROVEMENTS

3.3.1 Criteria

<u>Item</u>			<u>Applicable</u>
<u>No.</u>	<u>Criteria (source)</u>	<u>Remarks</u>	<u>subparagraph</u>
1	The physical inventory is conducted solely for investigation and documentation.	NA	2.1.2.1
2	The E300A will display consistent and accurate displays and measurements.	met	2.1.2.2
3	Verify manufacturer's specified full power internal battery life expectancy of 6 hours.	met	2.2.2
4	The E300A will meet the limits established in TB 38-750-2 for electrical safety of medical equipment.	met	2.3.2
5	The E300A will be rated satisfactory in all major categories of the evaluation. These include: Visual displays, controls, maintainability, conductors, fasteners, test points, test equipment, fuses and circuit breakers, labels and coding, and safety.	met	2.4.2
6	The E300A will display consistent and accurate measurements while exposed to an altitude equivalency of 15,000 feet above sea level.	met	2.5.2
7	The E300A will remain operational and display consistent and accurate measurements while exposed to vibrational stresses.	not met	2.6.2

<u>Item no.</u>	<u>Criteria (source)</u>	<u>Remarks</u>	<u>Applicable subparagraph</u>
8	The E300A will display consistent and accurate measurements during the high temperature operation check.	not met	2.7.2.1
9	The E300A will display consistent and accurate measurements after the high temperature storage.	met	2.7.2.2
10	The E300A will display consistent and accurate measurements during the low temperature operation check.	not met	2.8.2.1
11	The E300A will display consistent and accurate measurements after the low temperature storage.	met	2.8.2.2
12	The E300A will display consistent and accurate measurements while exposed to a high humidity.	met	2.9.2
13	The E300A will not produce emissions in excess of the limits set forth in MIL-STD-461C, Notice 5, paragraph 6.13.	not met	2.10.2.1
14	The E300A will not malfunction when it is subjected to radiated fields as specified in MIL-STD-461C, Notice 5, paragraph 6.20.	partially met	2.10.2.2
15	The E300A will not conduct emissions in excess of the limits set forth in MIL-STD-461C, Notice 5, paragraphs 6.1 and 6.2.	met	2.10.2.3
16	The E300A will not malfunction when it is subjected to conducted emissions as specified in MIL-STD-461C, Notice 5, paragraphs 6.7 and 6.10.	met	2.10.2.4

3.3.2 Significant problems which require corrective action

3.3.2.1 The unit should function properly in response to the test vibration signature in the Z axis.

3.3.2.2 The unit should not suffer LCD illumination failures while operating in a high temperature environment.

3.3.2.3 The unit should operate correctly in a low temperature environment of -25°C.

3.3.2.4 Additional shielding will be needed to reduce the electromagnetic emissions produced by the unit. This shielding also may reduce the susceptibility of the device to electromagnetic emissions.

3.4 REFERENCES

3.4.1 Association for the Advancement of Medical Instruments. Human factors engineering guidelines and preferred practices for the design of medical devices. Arlington, Virginia. AAMI-HE-1988. February.

3.4.2 Department of the Army. 1982. Environmental protection and enhancement. Washington, DC. AR 200-1. June.

3.4.3 Department of the Army. 1987. Maintenance management procedures for medical equipment. Washington, DC. TB 38-750-2. April.

3.4.4 Department of Defense. 1971. EMI characteristics, measurement of. Washington, DC. MIL-STD-462, Notice 3. February.

3.4.5 Department of Defense. 1971. EMI characteristics, requirements for equipment. Washington, DC. MIL-STD-461C, Notice 5. February.

3.4.6 Department of Defense. 1983. Environmental test methods and engineering guidelines. Washington, DC. MIL-STD-810D. July.

3.4.7 Mitchell, G. W., and Adams, J. E. 1988. Technical test and evaluation of aeromedical equipment. Fort Rucker, AL: U.S. Army Aeromedical Research Laboratory. USAARL Letter Report LR-88-16-1-2.

3.4.8 National Fire Protection Association. 1987. Standard for health care facilities. Quincy, Massachusetts. February. NAFF 99.

3.5 ABBREVIATIONS

ac	alternating current
AEST	aeromedical equipment suitability test
AVSCOM	U.S. Army Aviation Systems Command
AWR	airworthiness release
BB	broadband
BP	blood pressure
BPM	beats per minute
CAAF	Cairns Army Airfield
CRT	cathode ray tube
dB	decibel
dc	direct current
ECG	electrocardiograph
EMC	electromagnetic compatibility
EMI	electromagnetic interference
fpm	feet per minute
GFE	government furnished equipment
Gpk	gravity, peak
G(rms)	gravity (root mean square)
Hz	hertz
IAW	in accordance with
ITOP	in-flight test operating procedure
IGE	in-ground effect
MDE® E300A	MDE® Inc. volumetric infusion pump
kg	kilogram
kHz	kilohertz
KIAS	knots indicated airspeed
lb	pound
LCD	liquid crystal display
LED	light emitting diode
LISN	line impedance stabilization network
MAP	mean arterial pressure
MEDEVAC	medical evacuation
MHz	megahertz
MIL-STD	military standard
mL	milliliter
mm	millimeter
mmHg	millimeters of mercury
MSL	mean sea level

NA	not applicable
NAFP	National Association of Fire Prevention
NB	narrowband
NBC	nuclear, biological, and chemical
NiCad	nickel cadmium
NOE	nap-of-the-earth
NVG	night vision goggles
RAM	random access memory
RF	radio frequency
RH	relative humidity
ROM	read only memory
TB	technical bulletin
TFT	technical feasibility testing
T & E	test and evaluation
USAARL	U.S. Army Aeromedical Research Laboratory
V/m	volts per meter

3.6 LIST OF MANUFACTURERS

- 3.6.1 MDE[®] Corporation
127.0 Wentworth Street
Arloca, CA 91331-4309
- 3.6.2 Sikorsky Aircraft
6900 Main Street
Stratford, CT 06601
- 3.6.3 Neurodyne-Dempsey, Inc.
200 Arrowhead Drive
Carson City, NV 89701
- 3.6.4 Tenney Engineering, Inc.
1090 Springfield Road
P.O. Box 3142
Union, NJ 07083
- 3.6.5 Unholtz-Dickey Corporation
6 Brookside Drive
Wallingford, CT 06492
- 3.6.6 Solar Electronics Company
901 North Highland Avenue
Hollywood, CA 90038
- 3.6.7 Tektronix, Inc.
P.O. Box 500
Beaverton, OR 97077
- 3.6.8 Valmedix
32303 Howard Street
Madison Heights, MI 48071

3.7 DISTRIBUTION LIST

Commander, U.S. Army Natick Research,
Development and Engineering Center
ATTN: SATNC-MIL (Documents
Librarian)
Natick, MA 01760-5040

U.S. Army Communications-Electronics
Command
ATTN: AMSEL-RD-ESA-D
Fort Monmouth, NJ 07703

Commander
10th Medical Laboratory
ATTN: Audiologist
APO New York 09180

Naval Air Development Center
Technical Information Division
Technical Support Detachment
Warminster, PA 18974

Commanding Officer, Naval Medical
Research and Development Command
National Naval Medical Center
Bethesda, MD 20814-5044

Deputy Director, Defense Research
and Engineering
ATTN: Military Assistant
for Medical and Life Sciences
Washington, DC 20301-3080

Commander, U.S. Army Research
Institute of Environmental Medicine
Natick, MA 01760

Library
Naval Submarine Medical Research Lab
Box 900, Naval Sub Base
Groton, CT 06349-5900

Director, U.S. Army Human
Engineering Laboratory
ATTN: Technical Library
Aberdeen Proving Ground, MD 21005

Commander
Man-Machine Integration System
Code 602
Naval Air Development Center
Warminster, PA 18974

Commander
Naval Air Development Center
ATTN: Code 602-B (Mr. Brindle)
Warminster, PA 18974

Commanding Officer
Armstrong Laboratory
Wright-Patterson
Air Force Base, OH 45433-6573

Director
Army Audiology and Speech Center
Walter Reed Army Medical Center
Washington, DC 20307-5001

Commander/Director
U.S. Army Combat Surveillance
and Target Acquisition Lab
ATTN: DELCS-D
Fort Monmouth, NJ 07703-5304

Commander, U.S. Army Institute
of Dental Research
ATTN: Jean A. Setterstrom, Ph. D.
Walter Reed Army Medical Center
Washington, DC 20307-5300

Commander, U.S. Army Test
and Evaluation Command
ATTN: AMSTE-AD-H
Aberdeen Proving Ground, MD 21005

Naval Air Systems Command
Technical Air Library 950D
Room 278, Jefferson Plaza II
Department of the Navy
Washington, DC 20361

Director
U.S. Army Ballistic
Research Laboratory
ATTN: DRXBR-OD-ST Tech Reports
Aberdeen Proving Ground, MD 21005

Commander
U.S. Army Medical Research
Institute of Chemical Defense
ATTN: SGRD-UV-AO
Aberdeen Proving Ground,
MD 21010-5425

Commander, U.S. Army Medical
Research and Development Command
ATTN: SGRD-RMS (Ms. Madigan)
Fort Detrick, Frederick, MD 21702-5012

Director
Walter Reed Army Institute of Research
Washington, DC 20307-5100

HQ DA (DASG-PSP-O)
5109 Leesburg Pike
Falls Church, VA 22041-3258

Harry Diamond Laboratories
ATTN: Technical Information Branch
2800 Powder Mill Road
Adelphi, MD 20783-1197

U.S. Army Materiel Systems
Analysis Agency
ATTN: AMXSY-PA (Reports Processing)
Aberdeen Proving Ground
MD 21005-5071

U.S. Army Ordnance Center
and School Library
Simpson Hall, Building 3071
Aberdeen Proving Ground, MD 21005

U.S. Army Environmental
Hygiene Agency
Building E2100
Aberdeen Proving Ground, MD 21010

Technical Library Chemical Research
and Development Center
Aberdeen Proving Ground, MD
21010-5423

Commander
U.S. Army Medical Research
Institute of Infectious Disease
SGRD-UIZ-C
Fort Detrick, Frederick, MD 21702

Director, Biological
Sciences Division
Office of Naval Research
600 North Quincy Street
Arlington, VA 22217

Commander
U.S. Army Materiel Command
ATTN: AMCDE-XS
5001 Eisenhower Avenue
Alexandria, VA 22333

Commandant
U.S. Army Aviation
Logistics School ATTN: ATSQ-TDN
Fort Eustis, VA 23604

Headquarters (ATMD)
U.S. Army Training
and Doctrine Command
ATTN: ATBO-M
Fort Monroe, VA 23651

Structures Laboratory Library
USARTL-AVSCOM
NASA Langley Research Center
Mail Stop 266
Hampton, VA 23665

Naval Aerospace Medical
Institute Library
Building 1953, Code 03L
Pensacola, FL 32508-5600

Command Surgeon
HQ USCENTCOM (CCSG)
U.S. Central Command
MacDill Air Force Base FL 33608

Air University Library
(AUL/LSE)
Maxwell Air Force Base, AL 36112

U.S. Air Force Institute
of Technology (AFIT/LDEE)
Building 640, Area B
Wright-Patterson
Air Force Base, OH 45433

Henry L. Taylor
Director, Institute of Aviation
University of Illinois-Willard Airport
Savoy, IL 61874

Chief, National Guard Bureau
ATTN: NGB-ARS (COL Urbauer)
Room 410, Park Center 4
4501 Ford Avenue
Alexandria, VA 22302-1451

Commander
U.S. Army Aviation Systems Command
ATTN: SGRD-UAX-AL (LTC Gillette)
4300 Goodfellow Blvd., Building 105
St. Louis, MO 63120

U.S. Army Aviation Systems Command
Library and Information Center Branch
ATTN: AMSAV-DIL
4300 Goodfellow Boulevard
St. Louis, MO 63120

Federal Aviation Administration
Civil Aeromedical Institute
Library AAM-400A
P.O. Box 25082
Oklahoma City, OK 73125

Commander
U.S. Army Academy
of Health Sciences
ATTN: Library
Fort Sam Houston, TX 78234

Commander
U.S. Army Institute of Surgical Research
ATTN: SGRD-USM (Jan Duke)
Fort Sam Houston, TX 78234-6200

AAMRL/HEX
Wright-Patterson
Air Force Base, OH 45433

John A. Dellinger,
Southwest Research Institute
P. O. Box 28510
San Antonio, TX 78284

Product Manager
Aviation Life Support Equipment
ATTN: AMCPM-ALSE
4300 Goodfellow Boulevard
St. Louis, MO 63120-1798

Commander
U.S. Army Aviation
Systems Command
ATTN: AMSAV-ED
4300 Goodfellow Boulevard
St. Louis, MO 63120

Commanding Officer
Naval Biodynamics Laboratory
P.O. Box 24907
New Orleans, LA 70189-0407

Assistant Commandant
U.S. Army Field Artillery School
ATTN: Morris Swott Technical Library
Fort Sill, OK 73503-0312

Commander
U.S. Army Aviation Center
Directorate of Combat Developments
Building 507
Fort Rucker, AL 36362

U.S. Army Dugway Proving Ground
Technical Library, Building 5330
Dugway, UT 84022

U.S. Army Yuma Proving Ground
Technical Library
Yuma, AZ 85364

AFFTC Technical Library
6510 TW/TSTL
Edwards Air Force Base,
CA 93523-5000

Commander
Code 3431
Naval Weapons Center
China Lake, CA 93555

Aeromechanics Laboratory
U.S. Army Research and Technical Labs
Ames Research Center, M/S 215-1
Moffett Field, CA 94035

Sixth U.S. Army
ATTN: SMA
Presidio of San Francisco, CA 94129

Commander
U.S. Army Aeromedical Center
Fort Rucker, AL 36362

Strughold Aeromedical Library
Document Service Section
2511 Kennedy Circle
Brooks Air Force Base, TX 78235-5122

Dr. Diane Damos
Department of Human Factors
ISSM, USC
Los Angeles, CA 90089-0021

U.S. Army White Sands
Missile Range
ATTN: STEWS-IM-ST
White Sands Missile Range, NM 88002

U.S. Army Aviation Engineering
Flight Activity
ATTN: SAVTE-M (Tech Lib) Stop 217
Edwards Air Force Base, CA 93523-5000

Ms. Sandra G. Hart
Ames Research Center
MS 262-3
Moffett Field, CA 94035

Commander, Letterman Army Institute
of Research
ATTN: Medical Research Library
Presidio of San Francisco, CA 94129

Commander
U.S. Army Medical Materiel
Development Activity
Fort Detrick, Frederick, MD 21702-5009

Commander
U.S. Army Health Services Command
ATTN: HSOP-SO
Fort Sam Houston, TX 78234-6000

U. S. Army Research Institute
Aviation R&D Activity
ATTN: PERI-IR
Fort Rucker, AL 36362

Commander
U.S. Army Safety Center
Fort Rucker, AL 36362

U.S. Army Aircraft Development
Test Activity
ATTN: STEBG-MP-P
Cairns Army Air Field
Fort Rucker, AL 36362

Commander U.S. Army Medical Research
and Development Command
ATTN: SGRD-PLC (COL Schnakenberg)
Fort Detrick, Frederick, MD 21702

MAJ John Wilson
TRADOC Aviation LO
Embassy of the United States
APO New York 09777

Netherlands Army Liaison Office
Building 602
Fort Rucker, AL 36362

British Army Liaison Office
Building 602
Fort Rucker, AL 36362

Italian Army Liaison Office
Building 602
Fort Rucker, AL 36362

Directorate of Training Development
Building 502
Fort Rucker, AL 36362

Chief
USAHEL/USAAVNC Field Office
P. O. Box 716
Fort Rucker, AL 36362-5349

Commander U.S. Army Aviation Center
and Fort Rucker
ATTN: ATZQ-CG
Fort Rucker, AL 36362

Chief
Test & Evaluation Coordinating Board
Cairns Army Air Field
Fort Rucker, AL 36362

MAJ Terry Newman
Canadian Army Liaison Office
Building 602
Fort Rucker, AL 36362

German Army Liaison Office
Building 602
Fort Rucker, AL 36362

LTC Patrice Cottebrune
French Army Liaison Office
USAAVNC (Building 602)
Fort Rucker, AL 36362-5021

Australian Army Liaison Office
Building 602
Fort Rucker, AL 36362

Dr. Garrison Rapmund
6 Burning Tree Court
Bethesda, MD 20817

Commandant, Royal Air Force
Institute of Aviation Medicine
Farnborough Hampshire GU14 6SZ UK

Commander
U.S. Army Biomedical Research
and Development Laboratory
ATTN: SGRD-UBZ-I
Fort Detrick, Frederick, MD 21702

Defense Technical Information
Cameron Station, Building 5
Alexandra, VA 22304-6145

Commander, U.S. Army Foreign Science
and Technology Center
AIFRTA (Davis)
220 7th Street, NE
Charlottesville, VA 22901-5396

Director,
Applied Technology Laboratory
USARTL-AVSCOM
ATTN: Library, Building 401
Fort Eustis, VA 23604

Commander, U.S. Air Force
Development Test Center
101 West D Avenue, Suite 117
Eglin Air Force Base, FL 32542-5495

Aviation Medicine Clinic
TMC #22, SAAF
Fort Bragg, NC 28305

Dr. H. Dix Christensen
Bio-Medical Science Building, Room 753
Post Office Box 26901
Oklahoma City, OK 73190

Commander, U.S. Army Missile
Command
Redstone Scientific Information Center
ATTN: AMSMI-RD-CS-R
/ILL Documents
Redstone Arsenal, AL 35898

Director
Army Personnel Research Establishment
Farnborough, Hants GU14 6SZ UK

U.S. Army Research and Technology
Laboratories (AVSCOM)
Propulsion Laboratory MS 302-2
NASA Lewis Research Center
Cleveland, OH 44135

COL John F. Glenn
U.S. Army Medical Research
& Development Command
SGRD-ZC
Fort Detrick, Frederick, MD 21702-5012

Dr. Eugene S. Channing
7985 Schooner Court
Frederick, MD 21701-3273

LTC Gaylord Lindsey (5)
USAMRDC Liaison at Academy
of Health Sciences
ATTN: HSHA-ZAC-F
Fort Sam Houston, TX 78234

Dr. A. Kornfield, President
Biosearch Company
3016 Revere Road
Drexel Hill, PA 29026

NVEOD
AMSEL-RD-ASID
(Attn: Trang Bui)
Fort Belvoir, VA 22060

CA Av Med
HQ DAAC
Middle Wallop
Stockbridge Hants S020 8DY UK

Commander and Director
USAE Waterways Experiment Station
ATTN: CEWES-IM-MI-R
Alfrieda S. Clark, CD Dept.
3909 Halls Ferry Road
Vicksburg, MS 39180-6199

Mr. Peter Seib
Human Engineering Crew Station
Box 266
Westland Helicopters Limited Yeovil,
Somerset BA20 2YB UK

Brazilian Army Liaison Office
Building 602
Fort Rucker, AL 36362

Dr. Christine Schlichting
Behavioral Sciences Department
Box 900, NAVUBASE NLON
Groton, CT 06349-5900

COL C. Fred Tyner
U.S. Army Medical Research
& Development Command
SGRD-ZB
Fort Detrick, Frederick, MD 21702-5012

COL MacWillie, Director
Directorate of Combat Developments
ATZQ-CD
Building 515
Fort Rucker, AL 36362

**END
FILMED**

DATE:

8-93

DTIC