

Report Documentation Page

Form Approved
OMB No. 0704-0188

Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.

1. REPORT DATE DEC 1995		2. REPORT TYPE		3. DATES COVERED 00-00-1995 to 00-00-1995	
4. TITLE AND SUBTITLE Versatile sample handling system for scanning tunneling microscopy studies of molecular beam epitaxy				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Naval Research Laboratory, Code 6177, 4555 Overlook Avenue SW, Washington, DC, 20375				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution unlimited					
13. SUPPLEMENTARY NOTES					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT	18. NUMBER OF PAGES	19a. NAME OF RESPONSIBLE PERSON
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified			

Versatile sample handling system for scanning tunneling microscopy studies of molecular beam epitaxy

L. J. Whitman^{a)} and P. M. Thibado
Code 6177, Naval Research Laboratory, Washington, DC 20375

F. Linker and J. Patrin
Park Scientific Instruments, 1171 Borregas Avenue, Sunnyvale, California 94089

(Received 22 December 1995; accepted 1 March 1996)

The ongoing development of short-period semiconductor superlattices for electronic and optoelectronic applications requires atomic-scale control of epitaxial growth, especially at the interfaces.¹ Given this requirement, there is a critical need for *in situ* characterization on the atomic scale as provided by scanning tunneling microscopy (STM).² Here we describe a sample handling system designed to integrate a modified commercial STM³ into a multichamber ultrahigh vacuum (UHV) molecular beam epitaxy (MBE) facility. The system uses a simple, yet versatile, sample holder design that enables quick and easy sample transfers between multiple chambers, including two Riber MBE and two surface characterization chambers interconnected by Riber UHV ModuTracTM.

The chief component of the sample handling system is a molybdenum sample holder, shown in Fig. 1(a). The essential elements of the holder are its square shape and two pairs of clearance holes oriented at right angles to each other: one pair through the front near the top of the holder, and the other running from top to bottom along the edges. As will be described below, the holes are for mounting the holder on the UHV manipulators and STM sample stage. As long as these four holes are present, the basic holder can be modified in a number of ways depending on the sample mounting and heating requirements. The version shown in Fig. 1 is a "window frame" design with four tapped holes through the face of the holder that are used to mount a semiconductor wafer sample. The interior cutout permits direct radiative heating of the wafer sample along with infrared transmission-based thermometry during MBE.⁴ Alternately, the holder can be left solid if heating of the wafer via thermal contact with the holder is desired. Note that mounting the wafer sample on the front face of the holder allows grazing-angle access for reflection high energy electron diffraction (RHEED).

The sample holder mounts onto a commercial STM (with a modified sample stage) by inserting the two vertical pins on the STM stage into the long holes running the length of the holder [Fig. 1(b)]. The pins on the STM are of unequal length, with the longer pin having raised points near the top and bottom, and the shorter pin having a raised point near the top end only [see Fig. 2(c)]. Two identical molybdenum spring clips (0.25 mm thick) push on the front of the sample holder, holding the inside of the long holes against the three raised points on the pins and thereby achieving a stable con-

tact. In addition to the standard 1 cm of horizontal sample motion available, this modified design permits the sample to be stably mounted up to 5 mm above the base, resulting in a viewable area of $\sim 0.5 \text{ cm}^2$. To perform MBE on an STM-compatible sample, the sample holder is mounted on a standard 5-cm-diam "moly-block" [Fig. 1(c)]. The sample holder is held in place by two 0.25-mm-thick molybdenum compression clips screwed to the moly-block (the only required modifications of the moly-block). Screws made of tantalum are used to avoid seizing in the moly-block.




FIG. 1. (a) Schematic drawing of a "window frame" design sample holder. Samples up to $17 \text{ mm} \times 17 \text{ mm}$ (dashed square) can be mounted onto the front by four 0-80 screws. (b) Front view of the modified sample stage on the STM. The sample holder mounts onto two pins, one shorter than the other, and is held in place by molybdenum spring clips. The raised area near the top of the long pin is shown as a white spot. Outlines of the sections of the pins hidden behind the clips are shown as dashed lines. The dash-dot-dash lines show the outline of a sample holder when mounted. See Fig. 2(c) for a side view. (c) Moly-block used for MBE. The sample holder is held onto the face by molybdenum compression clips attached to the moly-block by tantalum screws. All components are drawn approximately to scale.

^{a)}Electronic mail: Lloyd.Whitman@nrl.navy.mil




FIG. 2. Illustrations (side view) of the sample transfers between the manipulator and (a) the load-lock transfer bay or sample storage bay; (b) the ModuTrac transfer arm; (c) the STM. The dashed lines show the outline of the sample holder in the mounted positions.

In addition to the STM and the moly-block, the sample holder also mounts on an UHV sample manipulator, a load-lock transfer bay, a sample storage bay, and a ModuTrac transfer arm (all except the last are on the STM chamber), as shown in Fig. 2. The sample holder mounts onto the sample manipulator via two horizontal pins that are inserted into the holder from the back, and is locked into place by a spring clip that pushes down on the top edge. The load-lock transfer bay and sample storage bay both have a slotted-plate design, holding four and seven sample holders, respectively, and are attached to magnetically coupled linear-rotary transporters. (Sample holders rest upright in slots, held in place by gravity.) The ModuTrac transfer arm, also a magnetically coupled linear-rotary transporter, has a simple pair of pins and spring clips on the end, oriented perpendicular to the axis of rotation. These pins and clips engage the long holes in the sample holder in a manner similar to those on the STM stage.

Transfer of the sample between the various components within the multichamber facility relies on perpendicular motions between the sample manipulators and the sample holder mounts, as indicated in Fig. 2. Transfer from the manipulator to the load-lock transfer bay or sample storage bay occurs by first positioning the empty bay slot under the sample holder, then lowering the manipulator to place the sample holder into the slot, and finally translating the bay




FIG. 3. Top view (inside the ModuTrac) of the motions for transfer of a sample holder between the ModuTrac transfer arm and a moly-block mounted on a cart in the ModuTrac. Note that the track is on the bottom of the ModuTrac chamber.




FIG. 4. Filled-state STM images of a 1- μm -thick GaSb film grown by MBE on a GaSb(001) substrate. (a) 160 nm \times 160 nm area, with each gray level representing a change in height of 0.3 nm (a single GaSb layer). (b) 13 nm \times 13 nm area showing the atomic-scale $2 \times 5/c(2 \times 10)$ reconstruction.

laterally in order to disengage the sample holder from the manipulator. To transfer the holder onto the manipulator from one of the bays, the reverse procedure is used. Transfers between the manipulator and ModuTrac transfer arm are performed in a similar fashion. Transfers between the manipulator and STM occur via motion of the manipulator only, utilizing the 25 mm of in-plane motion permitted by the manipulator micrometers.

The essential components required to transfer the sample holder from the STM chamber to the moly-block (and subsequently to a MBE chamber) are illustrated in Fig. 3. Following transfer from the manipulator to the ModuTrac transfer arm, the sample is rotated 90° and then pulled out of the STM chamber to the far side of the ModuTrac. The moly-block cart, which runs on a track down the center of the ModuTrac, is then positioned in front of the transfer arm. The sample holder (on the end of the transfer arm) is then pushed into the moly-block, engaging the holder between the compression clips. Finally, the cart is moved off to the left, disengaging the sample holder from the transfer arm. Transfer from the moly-block to the STM chamber occurs with the reverse procedure. Note that the simple two-clip mounting system on the moly-block provides considerable latitude in the alignment of the sample holder during these transfers.

The effectiveness of the sample holder design for *in situ* study of MBE-grown films is demonstrated by the images shown in Fig. 4 of a GaSb(001) homoepitaxial surface.⁵ Both the large-scale surface morphology and the atomic-scale structure can be readily observed after cooling the sample to ambient temperature.

ACKNOWLEDGMENTS

The authors thank B. R. Bennett and B. V. Shanabrook for preparing the GaSb(001) surface. This work was supported in part by the Office of Naval Research and an NRL/NRC Postdoctoral Fellowship (P.M.T.).

¹B. R. Bennett, B. V. Shanabrook, R. J. Wagner, J. L. Davis, and J. R. Waterman, *Appl. Phys. Lett.* **63**, 949 (1993); B. V. Shanabrook, B. R. Bennett, and R. J. Wagner, *Phys. Rev. B* **48**, 17,172 (1993).

²P. M. Thibado, B. R. Bennett, M. E. Twigg, B. V. Shanabrook, and L. J. Whitman, *Appl. Phys. Lett.* **67**, 3578 (1995).

³AutoProbe VP, Park Scientific Instruments, 1171 Borregas Ave., Sunnyvale, CA 94089.

⁴B. V. Shanabrook, J. R. Waterman, J. L. Davis, R. J. Wagner, and D. S. Katzer, *J. Vac. Sci. Technol. B* **11**, 994 (1993).

⁵P. M. Thibado, B. R. Bennett, B. V. Shanabrook, and L. J. Whitman (unpublished).