

AD-E400 149

TECHNICAL REPORT ARLCD-TR-77052

ANALYSIS OF IMPURITIES IN COMPOSITION B BY THIN LAYER CHROMATOGRAPATY

T. H. CHEN


H. GILFORD

APRIL 1978



US ARMY ARMAMENT RESEARCH AND DEVELOPMENT COMMAND LARGE CALIBER WEAPON SYSTEMS LABORATORY DOVER. NEW JERSEY

APPROVED FOR PUBLIC RELEASE: DISTRIBUTION UNLIMITED.



The findings in this report are not to be construed as an official Department of the Army position.

DISPOSITION

Destroy this report when no longer needed. Do not return to the originator.

The citation in this report of the names of commercial firms or commercially available products or services does not constitute official endorsement or approval of such commercial firms, products, or services by the US Government.

SECURITY CLASSIFICATION	or rule and forms of a constant	
REPOR	T DOCUMENTATION PAGE	READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER	2. GOVT ACCESSION	NO. 3. RECIPIENT'S CATALOG NUMBER
Trebnical Report	ARLCD-TR-77052	
A. TITLE (and Subtitle)	<i>F</i>	5 TYPE OF REPORT & PERIOD COVERED
		a F
	PURITIES IN COMPOSITION B	Mecha : Al replan
I BY THIN LAYER (CHROMATOGRAPHY.	S. FERFORMING ORG. REPURT NURSEN
7. AUTHOR(a)		S. CONTRACT OR GRANT NUMSER(s)
T. H. Chen	(IRISBEEL C	1 5 1 1 0 0 11 13
M./Gilford	(10) 207 [19	1 AD-E444 149
9. PERFORMING ORGANIZA Feltman Research		10 PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
Picatinny Arsenal	-	
Dover, NJ 07801		
11. CONTROLLING OFFICE		REPORT DATE
ARRADCOM, STIN		Apr 170 78
ATTN: DRDAR-T	rss	TO STORES OF PAGES
Dover NJ 07801	AME & ADDRESSII ditionent from Controlling Office	196-
ARRADCOM, LCW		
ATTN: DRDAR-L	.CE	Unclassified
Dover, NJ 07801		15A DECLASSIFICATION/DOWNGRADING SCHEDULE
	lic release, distribution unlimited	
	ENT (of the eletract entered in State 20, if different	
17. DISTRIBUTION STATEM	ENT (of the eletract entered in State 20, if different	from Report)
17. DISTRIBUTION STATEM 18. SUPPLEMENTARY NOTE 19. KEY WORDS (Centinus on Composition B TNT)	ENT (of the abetract entered in State 29, if different ES reverse side if necessary and identify by block must Analysis Identification	from Report) Net) RAP
17. DISTRIBUTION SYNTEM 18. SUPPLEMENTARY NOTE 19. KEY WORDS (Centinue on Composition B TNT RDX	ENT (of the element entered in Stock 20, if different coveres side if necessary and identify by block must Analysis Identification Thin layer chromatography	From Report) Net) RAP XM650
17. DISTRIBUTION STATEM 18. SUPPLEMENTARY NOTE 19. KEY WORDS (Centinus on Composition B TNT)	ENT (of the abetract entered in State 29, if different ES reverse side if necessary and identify by block must Analysis Identification	hom Report) Net) RAP XM650 Malfunction
17. DISTRIBUTION STATEM 18. SUPPLEMENTARY NOTE 19. KEY WORDS (Centinus on Composition B TNT RDX Impurities	ENT (of the abetract entered in State 29, If different services also if necessary and identify by block must Analysis Identification Thin layer chromatography Ethylenediamine	RAP XM650 Malfunction M116E2-XM201 Cannon m-nitrotoluene
19. KEY WORDS (Continue on Composition B TNT RDX Impurities 21. A thin layer consociated impurities TNT levels in the chromatographic of diamine reagent.	ENT (of the abetract entered in State 29, If different entered in State 29	RAP XM650 Malfunction M110E2-XM201 Cannon m-nitrotoluene eveloped to identify a-TNT d involved reducing RDX and le before two-dimensional tified after applying ethylene- quenching, and R _f values
19. KEY WORDS (Continue on Composition B TNT RDX Impurities A thin layer cassociated impurit a TNT levels in the chromatographic of diamine reagent, were compared to	ENT (of the abstract entered in Stack 20, if different in the stack of	RAP XM650 Malfunction M110E2-XM201 Cannon m-nitrotoluene veloped to identify g-TNT d involved reducing RDX and le before two-dimensional tified after applying ethylene- quenching, and R _f values impurities known to be R St INCLASSIFIED
19. KEY WORDS (Continue on Composition B TNT RDX Impurities A thin layer consociated impurities TNT levels in the chromatographic of diamine reagent, were compared to	ENT (of the abstract entered in Stack 20, if different in the stack of	RAP XM650 Malfunction M110E2-XM201 Cannon m-nitrotoluene veloped to identify a-TNT d involved reducing RDX and le before two-dimensional tified after applying ethylene- quenching, and R _f values impurities known to be

20. ABSTRACT (Continued)

present in Military Grade TNT. Six fragments from a Composition B warhead (which did not premature but was taken from the same lot as those suspected of premature functioning) and a control sample of Composition B obtained from the Holston Army Ammunition Plant were analyzed and compared. The analyses revealed the presence of m-nitrotoluene impurity in only the fragment samples suspected of premature initiation.

TABLE OF CONTENTS

	Page No
Introduction	1
Experimental	1
Apparatus and Reagents Procedure	1 2
Results and Discussion	đ
Conclusions	S
Recommendations	5
References	6
Distribution List	13
Tables	
1 R _f values and colors developed by ethylenediamine reagent of standard substances	7
2 R _f values and colors developed by ethylenediamine reagent of components in Composition B	9
Figures	
1 Two dimensional separation of components in the concentrated mother liquor of a suspected Composition B sample	10
2 Differential visible absorption spectrum of Composition B samples (suspected - controlled) in ethanol	11
in ethanor	·~

A STATE OF THE PARTY OF THE PAR

INTRODUCTION

An investigation was conducted to determine causes of malfunctions which took place during the test firing of a particular group of XM650 RAP rounds (ref. 1). The warheads initiated prematurely, one possible cause of which could be the presence of impurities in the Composition B warhead charges. The time available for identifying such impurities and determining their levels of concentration was limited to approximately five weeks. One approach was to analyze suspect samples of Composition B using thin layer chromatography (TLC) procedures and to confine the analysis to α-TNT associated impurities. Work has appeared (ref. 2, 3, and 4) in which TLC procedures were used to identify impurities in production samples of α -TNT. In the work reported here, RDX is a main constituent of the samples to be analyzed since the nominal formulation of Composition B is 60 parts RDX, 40 parts α -TNT and 1 part wax. The presence of RDX and α -TNT at such high levels interfered with the separation and visualization of the compounds to be identified. A method was developed to remove RDX and α-TNT. An important aspect of this method involved dissolving RDX and α -TNT in acetone and then displacing acetone with ethanol. This resulted in an alcoholic solution of impurities with RDX and α -TNT levels greatly reduced.

Six samples of Composition B from fragments of a suspect warhead charge were received for analysis. A control sample of Composition B, from a lot manufactured at the Holston Army Ammunition Plant, was used for comparison.

EXPERIMENTAL

Apparatus and Reagents

Standard substances for identification purposes were practical grade. Other reagents and solvents were reagent grade except in later experiments where industrial grade absolute ethyl alcohol was used

The spray reagent was a 10 Å ethylenediamine solution in dimethylsuifoxide, 65: 35 by volume (ref. 5). A 100 μ l Hamilton syringe was used to deliver 15 μ l aliquots of preconcentrated solution of impurities to the TLC sheets.

A Perkin-Elmer 202 Ultraviolet-Visible Spectrophotometer was used to record visible spectra of the extracted impurities in absolute ethyl alcohol.

Procedure

Preconcentration of Impurities and RDX Removal

A sample of Composition B (5 to 7 gram piece) was crushed remotely into small granular pieces. A 300 mg portion of the sample was weighed into a 2 ml test tube and approximately 1.3 ml of acetone were added. The mixture was placed in a dibutylphthalate bath and the bath temperature was raised slowly to 60°C to prevent excessive loss of solvent. The sample was mixed frequently to completely dissolve it. More acetone was added as needed. The solution was then evaporated to apparent dryness. During the evaporation, the solid RDX which had accumulated on the inner wall of the test tube was scraped to the bottom of the test tube with a stirring rod to minimize the entrapment of impurities in the solid mass. The mixture was centrifuged for 3 to 4 minutes to separate the mother liquer (about 0.1 ml) from RDX. The mother liquor was decanted to a second 2-ml test tube. About 0.3 mi absolute ethanol was added drop-by-drop to the first test tube to wet the RDX and the combination was mixed thoroughly. This solution was centrifuged and decented into the second test tube. Three alcohol washings were conducted. The combined solution of impurities (about 1 ml) was heated slowly to 78°C to prevent bumping and again evaporated to apparent dryness. The mixture was centrifuged and the concentrated alcohol solution (about 0.1 ml) was used to spot the TLC sheet.

Preparation of Impurity Solution for Recording Visible Spectrum

After sampling for spotting the TLC sheet, the cencent ated supernatant alcohol solution was used to record the visible spectrum. About 1 ml ethanol was added to the remaining supernatant solution and the combination was thoroughly mixed. The mixture was centrifuged and the supernatant solution carefully withdrawn with a dropping pipette and transferred to a 10 ml volumetric flask. The washing was repeated twice. To prevent scattering in the spectrophotometric recording, care was exercised to assure that no solid material was transferred. Additional alcohol was added as needed and the solution was then used to record the absorption spectrum in the visible region in a 1 cm silica cell.

Separation Procedure

A 15 μ I aliquot of the preconcentrated alcohol solution of impurities was placed at a spot on the lower left corner of the TLC sheet 3.3 cm from each edge. Care was taken to avoid sampling of solid residue into the syringe. The syringe was clamped on a stand and the needle tip adjusted to 1.5 mm above the TLC sheet for reproducing spotting without unnecessarily broadening the spot.

After the alcohol evaporated, the sheet was placed in the first chromatographic tank (24.2 cm H x 27.3 cm L x 7.4 cm W) which contained 100 ml benzene-cyclohexane-ethylacetate (50: 45: 5 by volume). To get a solvent vapor-saturated atmosphere, the solvent mixture was used to completely wet two solvent saturation pads, 1 each of which had been positioned against a narrow side of the tank. The pads were cut to the approximate dimensions of the tank sides. The tank atmosphere equilibrated for at least one hour. The sample was chromatographed for 55 ± 5 minutes at 23 ± 4°C by the ascending technique and the solvent front was allowed to travel to approximately 4 mm from the upper edge of the sheet. The sheet was removed from the tank and exposed to the atmosphere for about 15 minutes. The sheet was turned 90° to the direction of the previous development and chromatographed again in a tank (22.6 cm H x 22.0 cm L x 6.3 cm W) containing an equilibrated benzene-cyclohexane (75: 25 by volume) atmosphere. The sheet was removed when the sample had traveled to approximately 4 mm from the upper edge of the sheet (50 ± 5 minutes).

The separated impurities were located by spraying the TLC sheet with ethylenediataine reagent and viewing it under short wavelength ultraviolet illumination. Fluorescence quenching revealed colorless or faintly colored spots. The colored spots faded with time, especially in the presence of light so the developed colors were noted within 15 minutes of spraying. The rates of color development and fading varied considerably among the impurities.

¹Gelman Instrument Company, Catalog Number 51334.

RESULTS AND DISCUSSION

The preliminary experiments, using preconcentrated acetone solution of impurities in Composition B samples, established that RDX in the saturated acetone solution, 7.3 g RDX/100 g acetone at 20°C (ref. 6), caused extensive streaking and interfered with the separation and identification of various components. This problem was solved by displacing acetone with absolute ethanol, which reduced the RDX solubility by seventy fold, 0.105 g RDX/100 g ethanol at 20°C (ref. 6). This procedure removed 99.9 and 99 percent of the RDX and α -TNT, respectively, from the 300 mg sample of Composition B. The 1 percent wax in the Composition B did not interfere with the separation of the various components.

All six samples of suspected Composition B (SCB) and that of controlled Composition B (CCB) were analyzed by two-dimensional TLC to separate and tentatively identify the impurities. The CCB and SCB samples were separated into twelve and thirteen components respectively, the only difference being the presence in the SCB of a relatively non-polar component identified as m-nitrotoluene.

A typical chromatogram of the SCB samples is shown in figure 1. It should be emphasized that the identifications of the various components were based on the comparison of Rf values, colors developed by ethylene-diamine reagent, and fluorescence quenching of standard substances (table 1) and unknowns (table 2). No further identifications were made. The identities of ten of the thirteen components in Composition B are listed in the last column of table 2.

The relative amounts of impurities in SCB and CCB samples could not be determined from the TLC results. However, as shown in figure 2, the SCB samples contained more colored impurities than samples of CCB. The concentrated alcohol solution of impurities absorbs intensely below 370 nm but absorbs weakly above 400 nm. The major difference between SCB and CCB samples is the slightly greater absorption of the SCB in the 400-510 nm region.

The TLC results pertain primarily to α -TNT impurities since the solvent systems used for chromatographic developments were originally intended for α -TNT analysis (ref. 4). Moreover, the limited time available for the present work did not permit further refinements of the procedures used or the development of procedures for analyzing RDX impurities in Composition B.

CONCLUSIONS

An analytical method employing thin layer chromatography was developed to identify impurities associated with $\alpha\textsc{-}TNT$ in Composition B. It was essential that the RDX and $\alpha\textsc{-}TNT$ levels in the solution to be applied to the TLC plate be reduced using ethanol to displace the acetone initially used to dissolve the Composition B sample. 99.9 percent of the RDX and 99 percent of the $\alpha\textsc{-}TNT$ originally present in the Composition B sample were removed, leaving $\alpha\textsc{-}TNT$ associated impurities in the mother liquor. The procedure involved a two-dimensional development of mother liquor spotted on fluorescent silica gel plates. Impurities were identified after the plates were sprayed with 10 M ethylenedlamine in dimethylsulfoxide. Color development, fluorescence quenching, and Rf values were compared to established values for twenty-one compounds known present in $\alpha\textsc{-}TNT$ as it is conventionally manufactured for military use. The analyses of Composition B samples revealed the presence of m-nitrotoluene impurity in the suspected samples but not in the controlled sample.

RECOMMENDATIONS

Additional work in this area would relate to developing procedures by which RDX associated impurities could be identified, making it possible to conduct a complete qualitative and possibly semi-quantitative analysis of Composition B samples. Further work is also indicated to identify α -TNT associated impurities so that greater resolution of separation for the various compounds could be effected. The entire effort would be extensive and would involve choosing appropriate solvent systems and plate materials for chromatographic development and using the appropriate visualization agents to study the color forming reactions of aromatic nitramines.

REFERENCES

- 1. "Preliminary Malfunction Investigation Report. M110E2-XM201 Cannon Firing XM650E4 HE RAP with XM188 Zone 9 Propelling Charge," (Fred Menke Chief Investigator, 1976).
- 2. C. D. Chandler, J. A. Kohlbeck and W. T. Bolleter, "Continuous TNT Process Studies, III. Thin layer Chromatographic Analysis of Oxidation Products from Nitration," J. Chromatog., 64, 123 (1972).
- 3. S. K. Yasuda, "Identification of Impurities in α-Trinitrotoluene by Thin-layer Chromatography," J. Chromatog., 13, 78 (1964).
- 4. J. A. Kohlbeck, C. D. Chandler, and W. T. Bolleter, "Application of Thin-layer Chromatography for Improving the Continuous TNT Process," J. Chromatog., 46, 173 (1970).
- 5. D. J. Glover and E. G. Kayser, "Quantitative Spectrophotometric Analysis of Polynitroaromatic Compounds by Reaction with Ethylene-diamine," Anal. Chem., 40, 2055 (1968).
- AMCP 706-177, "Engineering Design Handbook, Explosives Series, Properties of Explosives of Military Interest," Headquarters, U.S. Army Materiel Command, 22 March 1967, pp 72, 354-355.

Table 1

R_f values and colors developed by ethylenediamine reagent of standard substances

Solvent: benzene-cyclohexane-ethyl acetate (50:45:5 by volume) Concentration: 10 mg/ml or saturated solution in ethanol Sample size for TLC: 15 μ l

		Rel. R		
No.	$\mathbf{\underline{R}_f}$	Ref. a-TNT	Color	Substance 1
1	0.57~0.61	1.02-1.12	Faint Pink-Violet	3-MNT
2	0.58-0.¢0	1.03-1.10	Colorless	2-MNT
3	0.58-0.60	1.05-1.09	Faint Orange-Pink	4-MNT
\$	0.59	1.09	Yellow	2,5-DNT
5	0.53-0.56	1.00	Da: k Reddish Brown	2,4,6-TNT
6	0.54	0.95-0.98	Orange-Yellow	1,4-DNB
7	0.50-0.53	0.93-0.98	Colorless	2,6-DNT
S	0.43-0.52	0.81-0.96	Colorless	3,5-DNT
9	0.51	0 93	Reddish Brown	2,4,6~TNB
10	0.44-0.49	0.82-0.91	Pale Violet	2,4-DNT
11	0 43-0.48	0.31-0.86	Light Orange-Yellow	2,3,5-TNT
12	0.44-0.47	0.83-6.87	Orange	2.3,6-TNT
13	0.40-0.45	0.74-0.83	Orange	2,4,5-TNT
14	0.45-0.47	0.82-0.84	Colorless	1,3-DNB
15	0.34-0.42	0.64-0.77	Yellow	3,4-DNT
16	0.34-0.41	0.64-0.75	Yellow	2,3-DNT
17	0.30-0.37	0.55-0.65	Pink	Hexanitrobibenzyl

¹MNT, DNT, and TNT designate mono-, di-, and trinitrotoluene, respectively. DNB and TNB designate dinitro- and trinitro-benzene, respectively.

Table 1 (Continued)

<u>No</u> .	<u>R</u> f	Ref. R_f Ref. α -TNT	Color	Substance
18	0.33-0.36	0.59-0.65	Orange-Ye!low	1,2-DNB
19	0.30-0.32	0.57	Reddish Orange	2,4,6-Trinitro- benzaldehyde
20	0.23-0.27	0.42-0.47	Yeilow (slightly orange)	3,4,5-TNT
21	0.21-0.24	0.39-0.43	Yellow	2,3,4-TNT
22	0.065- 0.087	0.12-0.16	Brown	2,4,6-Trinitro- berzoic acid
23	0.039- 0.055	0.071-0.10	Colorless	RDX


Table 2

R_f values and colors developed by ethylenediamine reagent of components in Composition B

Solvents: benzene-cyclohexane-ethyl acetate (50:45:5 by volume) Sample size for TLC: 15 μ l

<u>No.</u>	<u>R</u> f	Rel. R _f Ref. α-TNT	Color	Identification
1	0.58-0.64	1.10-1.14	Violet	3-MNT
2	0.50-0.58	1.00	Dark Reddish Brown	2,4,6-TNT
3	0.42-0.53	0.83-0.91	Faint Violet	2,4-DNT
4	0.43-0.50	0.83-0.89	Colorless	3,5-DNT
5	0.38-0.47	0.76-0.82	Yellow (slightly orange)	2,3,5-TNT alone; or 2,3,5-TNT mixed with either 2,4,5-TNT or 2,3,6-TNT or both ¹
6	0.30-0.38	0.57-0.69	Pink	Hexanitrobibenzyl
7	0.28-0.36	0.53-0.61	Faint Tan	2,4,6-Trinitro- benzaldehyde
8	0.23-0.33	0.43-0.58	Purple	Unidentified
9	0.16-0.22	0.31-0.40	Yellow	2,3,4-TNT
10	0.11-0.16	0.19-0.28	Tan	Unidentified
11	0.067- 0.11	0.12-0.19	Tan	2,4,6-Trinitro- benzoic acid
12	0.032- 0.050	0.058- G.095	Colorless	RDX
13	0.00	0.00	Light Brown	Residue unidentified


This identification is uncertain. On the basis of the color developed, it appears to be 2,3,5-TNT or 2,3,5-TNT mixed with 2,3,6-TNT. On the basis of relative R_f value, it appears to be 2,4,5-TNT.



Identification: see table 2

Sample: 300 mg **C** Seg. DP; 15 µl ethanol solution for spotting Location of Components: by the color developed with 10 M ethylenediamine solution in dimethylsulfoxide (all except 4 and 12) and by fluorescence quenching (4 and 12)

Figure 1. Two dimensional separation of components in the concentrated mother liquor of a suspected Composition B sample



officered between the control of the

Samples: Suspected, 300 m.j FC "D" segment Controlled, 300 mg HOL 050-11

Figure 2. Differential visible absorption spectrum of Composition B samples (suspected - controlled) in ethanol

DISTRIBUTION LIST

Director

Chemical Systems Laboratory

ARRADCOM

ATTN: Nr. A. Hillsmeier

Mr. D. Renard (9)

Aberdeen Proving Ground, MD 21010

Commander

US Army Armament Research and Development Command

ATTN: DRDAR-CG, MG B. L. Lewis

DRDAR-TD, Dr. R. Weigle DRDAR-TDS, Mr. V. Lindner DRDAR-LC, COL P. B. Kenyon DRDAR-LCM, Mr. L. Saffien

DRDAR-LCE, Mr. T. C. Castorina (24)

DRDAR-TSS (5)

Dover, NJ 07801

Director

Ballistic Research Laboratory

USA ARRADCOM

ATTN: Technical Library

Aberdeen Proving Ground, MD 21010

Commander

Radford Army Ammunition Plant

ATTN: Dr. W. T. Bolleter

Hercules, Inc.

Radford, VA 24141

Commander

Badger Army Ammunition Plant

Baraboo, WI 53913

Commander

Indiana Army Ammunition Plant

Charlestown, IN 47111

Commar.der

Holston Army Ammunition Plant

Kingsport, TN 37660

Commander
Lone Star Army Ammunition Plant
ATTN: Technical Library
Texarkana, TX 75501

Commander Milan Army Ammunition Plant Milan, TN 38358

Commander
Tooele Army Ammunition Plant
Tooele, UT 84074

Commander
Iowa Army Ammunition Plant
Silas Mason, Mason & Hanger, Inc.
ATTN: Technical Library
Middletown, IA 52638

Commander
Joliet Army Ammunition Plant
Joliet, IL 60436

Commander Longborn Army Ammunition Plant Marshall, TX 75670

Defense Documentation Center (12) Cameron Station Alexandria, VA 22314

Commander Louisiana Army Ammunition Plant Shreveport, LA 71130

Commander Newport Army Ammunition Plant Newport, IN 47966

Commander
Volunteer Army Ammunition Plant
Chattanooga, TN 37401

Commander Kansas Army Ammunition Plant Parsons, KS 67357

Commander
US Army Research Office
ATTN: Dr. G. Wyman
Box CM, Duke Station
Durham, NC 27706

Commander
Naval Surface Weapons Center
White Oak Laboratory
ATTN: Technical Library
Dr. J. Hoffsommer
White Oak, Silver Spring, MD 20910

Director
DARCOM Field Safety Activity
ATTN: DRXOS-ES
Charlestown, IN 47111

Commander
US Army Training and Doctrine Command
ATTN: ATEN-ME
Ft Monroe, VA 23651

Commander
US Naval Sea Systems Command
ATTN: SEA-0332, Dr. A. B. Amstar
Washington, DC 20362

Lawrence Livermore Laboratory ATTN: Technical Library P.O. Box 808 Livermore, CA 94550

Commander
Fort Detrick
US Army Medical Bioengineering
R&D Laboratory
ATTN: Dr. D. H. Rosenblatt
Frederick, MD 21701

Weapon System Concept Team/CSL ATTN: DRDAR-ACW Aberdeen Proving Ground, MD 21010

Technical Library ATTN: DRDAR-TSB-S Aberdeen Proving Ground, MD 21005

Benet Weapons Laboratory Technical Library ATTN: DRDAR-LCB-TL Watervliet, NY 12189

Commander
US Army Armament Materiel Readiness Command
ATTN: DRSAR-LEP-L
Rock Island, IL 61299

Director
US Army TRADOC Systems Analysis Activity
ATTN: ATAA-SL (Tech Lib)
White Sands Missile Range, NM 88002