

AD-A080 906

ARIZONA UNIV TUCSON DEPT OF CHEMISTRY

F/G 20/9

A HIGH POWER INDUCTIVELY COUPLED PLASMA TORCH AND IMPEDANCE MAT--ETC(U)

APR 79 D L WINDSON, D R MEINE, M S DENTON

N00014-78-C-0513

UNCLASSIFIED

TR-22

NL

ALL INFORMATION CONTAINED HEREIN IS UNCLASSIFIED

END
DATE
FILMED

3 - 80

001

⑥ LEVEL II

OFFICE OF NAVAL RESEARCH
Contract N00014-75-C-0513
Task No. NR 051-549
Technical Report No.22

ADA 080906


A HIGH POWER INDUCTIVELY COUPLED PLASMA
TORCH AND IMPEDANCE MATCHING NETWORK

D. L. Windsor, D. R. Heine and M. B. Denton

Department of Chemistry
University of Arizona
Tucson, Arizona 85721

DDC FILE COPY

Prepared for Publication
in
Applied Spectroscopy



Reproduction in whole or in part is permitted for any purpose
of the United States Government.

Approved for Public Release: Distribution Unlimited

80 2 11 081

| REPORT DOCUMENTATION PAGE | | READ INSTRUCTIONS BEFORE COMPLETING FORM |
|--|-------------------------------------|--|
| 1. REPORT NUMBER 22 | 2. GOVT ACCESSION NO. (14) TP-32 | 3. RECIPIENT'S CATALOG NUMBER technical |
| 4. TITLE (and Subtitle) A HIGH POWER INDUCTIVELY COUPLED PLASMA TORCH AND IMPEDANCE MATCHING NETWORK | | 5. TYPE OF REPORT INTERIM neptis |
| 7. AUTHOR(s) D. L. Windsor, D. R. Heine, M. B. Denton | | 6. PERFORMING ORG. REPORT NUMBER |
| 9. PERFORMING ORGANIZATION NAME AND ADDRESS Department of Chemistry University of Arizona Tucson, AZ 85721 | | 8. CONTRACT OR GRANT NUMBER(s) 151 N00014-75-C-0513 |
| 11. CONTROLLING OFFICE NAME AND ADDRESS Office of Naval Research Arlington, VA 22217 11/13 APR 79 | | 10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS NR 051-549 |
| 14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office) 12/15 | | 12. REPORT DATE April 13, 1979 |
| | | 13. NUMBER OF PAGES 3 |
| | | 15. SECURITY CLASS. (of this report) UNCLASSIFIED |
| | | 15a. DECLASSIFICATION/DOWNGRADING SCHEDULE |
| 16. DISTRIBUTION STATEMENT (of this Report) Approved for Public Release; Distribution Unlimited | | |
| 17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report) | | |
| 18. SUPPLEMENTARY NOTES Prepared for publication in Applied Spectroscopy. | | |
| 19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Inductively Coupled Plasma Torch Inductively Coupled Plasma Matching Network | | |
| 20. ABSTRACT (Continue on reverse side if necessary and identify by block number) An improved dismantlable inductively coupled plasma torch and matching network are described. This torch and network provide convenience, improved power handling capabilities, and a high degree of flexibility. | | |

A HIGH POWER ICP TORCH AND IMPEDANCE MATCHING NETWORK

D. L. Windsor*, D. R. Heine and M. B. Denton

Department of Chemistry
University of Arizona
Tucson, Arizona 85721

Index Headings: Inductively Coupled Plasma Torch;
Inductively Coupled Plasma Matching Network

* Present Address: Proctor and Gamble Company
Miami Valley Laboratories
Cincinnati, Ohio 45247

| | | |
|---------------------------------|---------------|-------------------------------------|
| ACCESSION for | | |
| NTIS | White Section | <input checked="" type="checkbox"/> |
| DDC | Buff Section | <input type="checkbox"/> |
| UNANNOUNCED | | <input type="checkbox"/> |
| JUSTIFICATION | | |
| BY | | |
| DISTRIBUTION/AVAILABILITY CODES | | |
| Dist. | AVAIL | and/or SPECIAL |
| A | | |

A HIGH POWER ICP TORCH AND IMPEDANCE MATCHING NETWORK

4 The relative ease of operation, power handling capabilities and overall performance of an inductively coupled plasma-optical emission spectrometer (ICP-OES) depends upon numerous factors, not the least important of which is the design of the torch. A variety of such designs has appeared in the literature (1-6).^c Most workers have employed three concentric fused silica and/or pyrex tubes, referred to as an inner sample gas tube, a "flaired mouth" plasma gas tube, and an outer coolant tube. These tubes are either fused into one or two-piece units through glass blowing techniques or mounted in some type of dismantlable base assembly.

Precise alignment of the glass-blown torch can be achieved during construction through the use of an assembly jig. However, this type of torch is not easily repaired should problems be encountered. The dismantlable base approach can be difficult to align if the three tubes are held only near the bottom. This problem is aggravated by the fact that rarely is commercial quartz and pyrex glass tubing perfectly straight or cylindrical.

A torch design which is both dismantlable and self-aligning has been successfully used in this laboratory for several years. A unique feature is the use of nylon or teflon slip fit spacers between the sample and plasma gas tubes and the plasma and coolant gas tubes. 4

Properly implemented, this approach also provides a significant increase

in tangential flow velocity of the coolant gas, which improves heat transfer from the coolant tube, making possible the use of higher powers and longer coolant tubes. The increase in tangential flow is achieved through the use of threads cut into the coolant-plasma spacer with a 60° "v" tool; two cuts of 1.57 thread/cm, (4 threads/inch) offset by 180° and .2.5 mm deep, are made. This spacer is then placed as shown in Figure 1.

The spacer which fits between the sample and plasma gas tubes contains 12 verticle slots which are 1.4 mm wide and 1.2 mm deep. This permits laminar flow of the plasma or auxiliary gas. The dimensions of the glass components are listed in Table I. The sample tube is constricted at one end to 1.5 mm i.d., and the other end contains a 12/5 ball joint for easy attachment of the sample gas line.

A detailed drawing of the base is shown in Figure 2. In operation this assembly is securely fastened to a moveable mount, allowing precise positioning of the observation region. All glass inserts can be easily changed with the base mounted.

Data obtained using this torch for the analysis of carbon and hydrogen in samples from a gas-liquid chromatograph have shown that absolute precision can be increased more than twofold by increasing the coolant tube from 25 cm to 30 cm. (7). This extended the top of the coolant tube from 5 cm to 10 cm above the load coil.

An impedance coupling network which matches the output impedance of the high power rf power supply (50 ohms) to the input impedance of the load coil is shown in Figure 3.

The capacitors are the vacuum variable type, (Jennings Radio - UCSL 750) able to withstand 5,000 volts, and range from 5-750 pf. The circuit is constructed from 4.76 mm o.d., 3.24 mm i.d. (3/16" i.d.) copper tube which carries cooling water during operation.


This matching network-torch combination readily operates in excess of 5 kW forward power with less than 50 watts of reflected power. At forward power levels near 2 kW, the reflected power is less than 5 watts.

.


TABLE 1. Dimensions of glass used in the dismantlable torch.

| TUBE | LENGTH | O.D. | I.D. | TYPE |
|---------------|---------|-------|-------|--------|
| short coolant | 25 cm | 20 mm | 18 mm | quartz |
| long coolant | 30 cm | 20 mm | 18 mm | quartz |
| Plasma | | | | |
| mouth | 2.5 cm | 15 mm | 13 mm | vycor |
| stem | 22.5 cm | 11 mm | 9 mm | vycor |
| Sample | 30 cm | 7 mm | 5 mm | pyrex |


THIS PAGE IS BEST QUALITY REPRODUCTION
THE COPY FURNISHED TO DOD



THIS PAGE IS BEST QUALITY PHOTOGRAPH
AND ONLY SUBMITTED TO DCS



THIS PAGE IS BEST QUALITY PAPER
PRINTED ON 100% RECYCLED PAPER



Acknowledgement

This research was partially supported by the Office of Naval Research and an Alfred P. Sloan Foundation Research Fellowship awarded to M. Bonner Denton.

Figure Captions

- Figure 1. Dismantleable ICP torch assembly capable of operating at over 5 kW forward power. Spacers are used to provide precise alignment and improved gas flow patterns.
- Figure 2. Detailed drawing of the dismantleable ICP torch base (all dimensions are in inches).
- Figure 3. 27.12 Mhz impedance matching network for use with high powers.

Information Cited

1. K. Ohls and K. Krefta, ICP Information Newsletter, 1, 168 (1976).
2. F.E. Lichte and S.R. Koirtzmann, ICP Information Newsletter, 1, 200 (1976).
3. C.D. Allemand and R.M. Barnes, Appl. Spectrosc., 31, 434 (1977).
4. J.M. Mermet and C. Trassy, Appl. Spectrosc., 31, 237 (1977).
5. Ramon Barnes, ed. ICP Information Newsletter, 1, 49 (1975).
6. J.L. Genna, R.M. Barnes and C.D. Allemand, Anal. Chem., 49, 1450 (1977).
7. David Windsor, "Use of the Inductively-Coupled Plasma for the Elemental Analysis of Organic Compounds," Ph.D. Dissertation, (University of Arizona, 1977) pp. 120-124.

TECHNICAL REPORT DISTRIBUTION LIST, CEN

| | <u>No.</u> <u>Copies</u> | | <u>No.</u> <u>Copies</u> |
|--|-----------------------------|--|-----------------------------|
| Office of Naval Research 800 North Quincy Street Arlington, Virginia 22217 Attn: Code 472 | 2 | Defense Documentation Center Building 5, Cameron Station Alexandria, Virginia 22314 | 12 |
| ONR Branch Office 536 S. Clark Street Chicago, Illinois 60605 Attn: Dr. George Sandoz | 1 | U.S. Army Research Office P.O. Box 1211 Research Triangle Park, N.C. 27709 Attn: CRD-AA-IP | 1 |
| ONR Branch Office 715 Broadway New York, New York 10003 Attn: Scientific Dept. | 1 | Naval Ocean Systems Center San Diego, California 92152 Attn: Mr. Joe McCartney | 1 |
| ONR Branch Office 1030 East Green Street Pasadena, California 91106 Attn: Dr. R. J. Marcus | 1 | Naval Weapons Center China Lake, California 93555 Attn: Dr. A. B. Amster Chemistry Division | 1 |
| ONR Area Office One Wallis Plaza, Suite 601 San Francisco, California 94102 Attn: Dr. P. A. Miller | 1 | Naval Civil Engineering Laboratory Port Hueneme, California 93401 Attn: Dr. R. W. Drisko | 1 |
| ONR Branch Office Building 114, Section D 666 Summer Street Boston, Massachusetts 02210 Attn: Dr. L. H. Peebles | 1 | Professor K. E. Woehler Department of Physics & Chemistry Naval Postgraduate School Monterey, California 93940 | 1 |
| Director, Naval Research Laboratory Washington, D.C. 20390 Attn: Code 6100 | 1 | Dr. A. L. Slafkosky Scientific Advisor Commandant of the Marine Corps (Code RD-1) Washington, D.C. 20380 | 1 |
| The Assistant Secretary of the Navy (R,E&S) Department of the Navy Room 4E736, Pentagon Washington, D.C. 20350 | 1 | Office of Naval Research 800 N. Quincy Street Arlington, Virginia 22217 Attn: Dr. Richard S. Miller | 1 |
| Commander, Naval Air Systems Command Department of the Navy Washington, D.C. 20360 Attn: Code 310C (H. Rosenwasser) | 1 | Naval Ship Research and Development Center Annapolis, Maryland 21401 Attn: Dr. G. Bosmajian Applied Chemistry Division | 1 |
| | | Naval Ocean Systems Center San Diego, California 91232 Attn: Dr. S. Yagamoto, Marine Sciences Division | 1 |

Encl 1

TECHNICAL REPORT DISTRIBUTION LIST, 051C

| | <u>No. Copies</u> | | <u>No. Copies</u> |
|--|-----------------------|--|-----------------------|
| Dr. M. B. Dorton University of Arizona Department of Chemistry Tucson, Arizona 85721 | 1 | Dr. K. Wilson University of California, San Diego Department of Chemistry La Jolla, California | 1 |
| Dr. R. A. Osteryoung Colorado State University Department of Chemistry Fort Collins, Colorado 80521 | 1 | Dr. A. Zirino Naval Undersea Center San Diego, California 92132 | 1 |
| Dr. B. R. Kowalski University of Washington Department of Chemistry Seattle, Washington 98105 | 1 | Dr. John Duffin United States Naval Postgraduate School Monterey, California 93940 | 1 |
| Dr. S. P. Perone Purdue University Department of Chemistry Lafayette, Indiana 47907 | 1 | Dr. G. M. Hieftje Department of Chemistry Indiana University Bloomington, Indiana 47401 | 1 |
| | | Dr. Victor L. Rehn Naval Weapons Center Code 3813 China Lake, California 93555 | 1 |
| Dr. D. L. Venezky Naval Research Laboratory Code 6130 Washington, D.C. 20375 | 1 | Dr. Christie G. Enke Michigan State University Department of Chemistry East Lansing, Michigan 48824 | 1 |
| Dr. H. Freiser University of Arizona Department of Chemistry Tucson, Arizona 85721 | | Dr. Kent Eisentraut, MBT Air Force Materials Laboratory Wright-Patterson AFB, Ohio 45433 | |
| Dr. Fred Saalfeld Naval Research Laboratory Code 6110 Washington, D.C. 20375 | 1 | Walter G. Cox, Code 3632 Naval Underwater Systems Center Building 148 Newport, Rhode Island 02840 | |
| Dr. E. Chernoff Massachusetts Institute of Technology Department of Mathematics Cambridge, Massachusetts 02139 | 1 | | |

USE TO BEST QUALITY PAPER
BY 10/10/68 IN 1000