

AD-A091 666


ARMY MEDICAL BIOENGINEERING RESEARCH AND DEVELOPMENT --ETC F/G 7/4
EVALUATION OF METHODS FOR THE ANALYSIS OF SMALL MOLECULAR WEIGHT--ETC(U)
APR 80 W F COWEN, Y H CHANG, M W SMITH

UNCLASSIFIED

USAMBRDL-TR-8011

NL


For I
AD
201666



END
DATE
FILMED
12-80
DTIC

LEVEL

AD



TECHNICAL REPORT 8011

EVALUATION OF METHODS FOR THE ANALYSIS
OF SMALL MOLECULAR WEIGHT END-PRODUCTS
OF WASTEWATER OZONATION

William F. Cowen, Ph.D.
Yun H. Chang, CPT, M.S.
Michael W. Smith
William J. Baczynskyj

US ARMY MEDICAL BIOENGINEERING RESEARCH and DEVELOPMENT LABORATORY
Fort Detrick
Frederick, MD 21701

APRIL 1980



Approved for public release;
distribution unlimited



US ARMY MEDICAL RESEARCH and DEVELOPMENT COMMAND
Fort Detrick
Frederick, MD 21701

80 11 10 060

AD A091666

DDC FILE COPY

NOTICE

Disclaimer

The findings in this report are not be be construed as an official Department of the Army position unless so designated by other authorized documents.

Disposition

Destroy this report when it is no longer needed. Do not return it to the originator.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER TECHNICAL REPORT 8011	2. GOVT ACCESSION NO. AD-A092666	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) EVALUATION OF METHODS FOR THE ANALYSIS OF SMALL MOLECULAR WEIGHT END-PRODUCTS OF WASTEWATER OZONATION		5. TYPE OF REPORT & PERIOD COVERED Technical Report 1978-1979
7. AUTHOR(s) William F./Cowen Yun H./Chang		6. PERFORMING ORG. REPORT NUMBER
Michael W./Smith William J./Baczynskyj		8. CONTRACT OR GRANT NUMBER(s) 11 Apr 80
9. PERFORMING ORGANIZATION NAME AND ADDRESS U.S. Army Medical Bioengineering Research and Development Laboratory, ATTN: SGRD-UBG Fort Detrick, Frederick, MD 21701		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS 62720A 3E762720A835/00/128
11. CONTROLLING OFFICE NAME AND ADDRESS U.S. Army Medical Research and Development Command ATTN: SGRD-SI Fort Detrick, Frederick, MD 21701		12. REPORT DATE April 1980
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office) 14 SAME FDL-TP-1012		13. NUMBER OF PAGES 69
		15. SECURITY CLASS. (of this report) UNCLASSIFIED
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Acetic Acid Glyoxal Oxalic Acid Dimethylglyoxal Glyoxylic Acid Ozonation Formaldehyde Methylglyoxal Wastewater Analysis Formic Acid		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) The treatment of hospital laboratory wastewaters by ozonation produces a number of organic products: methanol and acetone are converted to formaldehyde, formic acid, and acetic acid and phenols are converted to oxalic acid, glyoxylic acid, and glyoxal. Methods for detecting these residual organic compounds at levels of < 5 mg/L in laboratory wastewater were evaluated and compared: Acetic and Formic Acids - A method of analysis that provides for isolation of formic and acetic acids from water as lipophilic tetra-n-butylammonium hydroxide salts, treatment with benzyl bromide to form benzyl esters, and detection of		

DD FORM 1 JAN 73 1473 EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED
SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)

20. ABSTRACT (Continued)

benzyl esters by FID-GC was tested. In a modified procedure, pentafluorobenzyl bromide was used instead of benzyl bromide; results showed no advantage over the benzyl bromide method, which easily detected acetic and formic acid levels as low as 0.25 mg/L. Recovery was quantitative, with a precision of 1 to 4% of the mean.

Oxalic Acid - Attempts to quantitatively derivatize 0- to 10-mg/L solutions of oxalic acid with pentafluorobenzyl bromide were unsuccessful.

Glyoxals - A method for FID-GC detection of glyoxals in aqueous solution was modified by heating the reaction mixture at 80°C for 30 minutes, by using methylene chloride as a solvent, and by choosing a more polar GC column. Although the method gave reasonable precision for methylglyoxal and dimethylglyoxal at about 2-mg/L levels, the precision for glyoxal was not acceptable.

Formaldehyde - Formaldehyde was readily determined by colorimetry for concentrations up to 8 mg/L and by fluorimetry for levels of 0.05 to 2 mg/L.

Glyoxylic Acid - At least 0.1 mg/L could be detected by a resorcinol procedure, with a precision of 7% relative standard deviation.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)

ACKNOWLEDGMENT

The authors thank Dr. E.S.K. Chian and Dr. P.K. Kuo for furnishing or suggesting several of the procedures that are evaluated in this report and Alan Rosencrance for performing the infrared and NMR analyses in the report.

Final technical editing was performed by Dr. David H. Rosenblatt and Dr. W. Dickinson Burrows of the U.S. Army Medical Bioengineering Research and Development Laboratory. Editorial assistance was provided by Franklin Research Center under Contract No. DAMD17-79-C-9129, Marsha Hall, project supervisor.

Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By _____	
Distribution/	
Availability Codes	
Avail and/or	
Dist	Special
A	

TABLE OF CONTENTS

	Page
I. INTRODUCTION	3
II. OBJECTIVE	3
III. ACETIC AND FORMIC ACIDS	3
A. Literature Review	3
B. Materials and Methods	5
C. Results and Discussion	8
IV. OXALIC ACID	22
A. Literature Review	22
B. Materials and Methods	24
C. Results and Discussion	27
V. GLYOXALS	37
A. Literature Review	37
B. Materials and Methods	38
C. Results and Discussion	39
VI. FORMALDEHYDE	43
A. Literature Review	43
B. Materials and Methods	44
C. Results and Discussion	44
VII. GLYOXYLIC ACID	48
A. Literature Review	48
B. Materials and Methods	49
C. Results and Discussion	49
VIII. CONCLUSIONS	53
A. Acetic and Formic Acids	53
B. Oxalic Acid	53
C. Glyoxals	53
D. Formaldehyde	54
E. Glyoxylic Acid	54
IX. LITERATURE CITED	55
APPENDIX A	58
1. Synthesis of Dipentafluorobenzyl Oxalate	58
2. Mass Spectrometry of PFB Acetate, Formate, and Oxalate	58

LIST OF TABLES

	Page
1. Synthetic Laboratory Reverse Osmosis Permeate Formulation	7
2. Benzyl Method of Recovery of Formic and Acetic Acids from RO-Lab Wastewater	10
3. Kinetic Study of Benzyl Formate and Benzyl Acetate Derivatization Reaction at Room Temperature	12
4. Relative Response of Formate and Acetate Peaks to Various Amounts of PFB-Br Used to Derivatize a Spiked RO-Lab Wastewater	16
5. PFB Method of Recovery of Formic and Acetic Acids from RO-Lab Wastewater	20
6. Relative Peak Areas from Analysis of 20 mL of 9.83 mg/L Oxalic Acid Standard, at 30° to 50°C Rotary Evaporation Bath Temperatures	31
7. Dimethyl Oxalate Method of Recovery of Oxalic Acid from Ozonated RO-Lab Wastewater	32
8. Recoveries of Glyoxals from RO-Lab Wastewater	40
9. Production of Glyoxals During Ozonation of RO-Lab Wastewater	40
10. Precision of Colorimetric Determination of Formaldehyde with Diluted O ₃ /RO-Lab Permeate	47
11. Accuracy of Colorimetric Determination of Formaldehyde with Spiked O ₃ /RO-Lab Permeate	47
12. Formaldehyde Formation During Two Ozonations of RO-Lab Wastewater .	48
13. Precision of Glyoxylic Acid Method (Low-Level Sample)	51
14. Recovery of Glyoxylic Acid Added to Ozonated RO-Lab Permeate	51
15. Production of Glyoxylic Acid During Ozonation of Synthetic RO-Lab Wastewater	52

LIST OF FIGURES

	Page
1. Structures of Suspected Small Molecular Weight End-Products in Ozonated Hospital Laboratory Wastewater	4
2. FID-GC of a 2.5-mg/L Standard of Acetic and Formic Acids, Treated with Benzyl Bromide	9
3. Standard Curves for Benzyl Bromide Derivatization of Formic and Acetic Acids	11
4. Low-level Standard Curves for Benzyl Bromide Derivatization of Formic and Acetic Acids	13
5. Hall Detector-GC of a 2.5-mg/L Standard of Acetic and Formic Acids	14
6. Standard Curves for Pentafluorobenzyl Bromide Derivatization of Formic and Acetic Acids	15
7. Hall Detector-GC of a RO-Lab Wastewater Sample Spiked with 4.35 mg/L of Formic Acid and 4.87 mg/L of Acetic Acid and Derivatized with 45 μ L of 10% (v/v) PFB-Br in Acetone	17
8. Hall Detector-GC of a RO-Lab Wastewater Sample Spiked with 4.35 mg/L of Formic Acid and 4.87 mg/L of Acetic Acid and Derivatized with 10 μ L of 10% PFB-BR in Acetone	18
9. Hall Detector-GC of a RO-Lab Wastewater Sample Treated as in Figure 8 and Using N ₂ Gas Stripping to Remove CO ₂ After Acidification by Ion Exchange	19
10. Low-level Standard Curves for PFB-Br Derivatization of Formic and Acetic Acids	21
11. Ozonation of Synthetic RO-Lab Wastewater, with Analyses of Formic and Acetic Acids by the Benzyl Bromide Method of Bethge and Lindstrom	23
12. Diazomethane Generator	26
13. Gas Chromatogram of TMS Oxalate Reaction Mixture	28
14. Standard Curve for Oxalic Acid by the TMS Method	29
15. Serial Dilution of TMS Oxalate Reaction Mixture with Pyridine-TMCS-BSTFA Reagent	30
16. Gas Chromatogram of Reaction Mixture from a 5-mg/L Oxalic Standard with Diazomethane	33

LIST OF FIGURES (Cont.)

	Page
17. Standard Curves for Dimethyl Oxalate	34
18. Standard Curves for Oxalic Acid with the 30°C Rotary Evaporation Method	35
19. Low-level Standard Curve for Oxalic Acid	36
20. Gas Chromatogram of Quinoxaline Mixture on OV-225	41
21. Standard Curves for Glyoxals	42
22. Standard Curve for Formaldehyde in a Test Using Spectrophotometric Detection	45
23. Standard Curve for Formaldehyde in a Test Using Fluorescence Detection	46
24. Standard Curve for Glyoxylic Acid	50
A-1. NMR Spectrum of Dipentafluorobenzyl Oxalate	59
A-2. Infrared Spectrum of Dipentafluorobenzyl Oxalate	60

I. INTRODUCTION

As part of the mission of the U.S. Army Medical Bioengineering Research and Development Laboratory (USAMBRDL) to develop water quality criteria for reuse of Army wastewaters, an in-house capability to determine concentrations of organic compounds in treated wastewaters is required.

One of the candidate methods for treating hospital laboratory wastewaters is ozonation,¹ which was shown by Chian and Kuo² to convert the major components, methanol and acetone, to formaldehyde, formic acid, and acetic acid. Similarly, ozonation of phenols by Gould and Weber³ was shown to produce oxalic acid, glyoxylic acid, and glyoxal. Since phenols may also be present in hospital laboratory wastewaters, their ozonation products were considered in the analytical requirements list. Figure 1 shows the structures of the compounds suspected as the major products from ozonation of hospital laboratory wastewaters. In general, ozonated organic compounds are converted to more oxidized forms such as aldehydes, ketones, and acids before loss of organic carbon as carbon dioxide. One form of evidence for this is the ratio of chemical oxygen demand to total organic carbon (COD/TOC), which has been shown by McCarthy and co-workers¹ to decrease during the ozonation of synthetic laboratory wastewater.

Because the level of residual TOC in highly ozonated wastewaters would be expected to lie at or below 5 mg/L, relatively sensitive methods are required for the detection of residual organic compounds. For example, a 1-mg/L residual TOC concentration due to oxalic acid would be equivalent to 3.75 mg/L of oxalic acid.

II. OBJECTIVE

The objective of these studies was to evaluate analytical methods capable of detecting levels of < 5 mg/L. The sources of methods were: (a) Army contractors (Dr. P.K. Kuo and Dr. E.S.K. Chian, University of Illinois, Champaign/Urbana, Illinois), (b) chemical literature, and (c) methods modified during the present studies from those of sources (a) and (b). By documenting the successful methods found from these studies, it is hoped that this report will be useful to researchers as a manual of analytical methods in future ozonation studies.

III. ACETIC AND FORMIC ACIDS

A. Literature Review

Acetic and formic acids have been determined by direct aqueous injection onto a gas-liquid chromatography column packed with 15% SP-1220/1% H₃PO₄ on 100/120 mesh Chromosorb WAW, with thermal conductivity detection.⁴ Acetic acid has been determined by the same method, using a graphitized carbon packing (Carbopack B/3% Carbowax 20 M/0.5% H₃PO₄)⁵ or Chromosorb 101





Figure 1. Structures of suspected small molecular weight end-products in ozonated hospital laboratory wastewater.

coated with 3% FFAP liquid phase.⁵ The sensitivity of thermal conductivity (TC) detectors is generally less than that of flame ionization detectors (FID), so that detectability by TC would not be expected to match the levels reached by FID detectors. In the case of aqueous formic acid, however, the FID response would be expected to be minimal, because there is only one C-H bond. By judicious selection of derivatives of these acids, the FID response per mole of acid can be increased, and better separation of the two compounds can be achieved.

Benzyl esters of formic, acetic, and higher homologous fatty acids have been used by Bethge and Lindstrom⁶ for quantification of these acids. They isolated the acids from water by first forming lipophilic tetra-n-butyl-ammonium hydroxide (TBAOH) salts of the acids and then removing the water by rotary evaporation. The remaining salts were dissolved in acetone and treated with benzyl bromide to form the benzyl esters. The detection limit was related to the volume of the sample aliquot taken for analysis; starting with a 10-mL aliquot, they were able to detect acids present at the 5×10^{-5} M level in the original sample. This would be equivalent to about 3 mg/L of acetic acid. Davis⁷ used 18-crown-6 (crown ether) as a catalyst with the potassium salts of acetic and higher homologous acids to form pentafluorobenzyl esters with pentafluorobenzyl bromide (PFB-Br) in an organic solvent. However, no procedure was given for extracting the acetic acid from a water solution. The pentafluorobenzyl esters were detected with a highly sensitive electron capture detector. The lower limit of sensitivity was 10 ng because the unreacted pentafluorobenzyl bromide peak and artifact peaks made work at lower attenuations impossible. In addition, Davis⁷ did not attempt to form the pentafluorobenzyl ester of formic acid.

Thus, the literature indicated that a good method existed for the isolation of formic and acetic acids from water in a form that would react in an organic solvent to displace bromide from benzyl bromide. By replacing benzyl bromide with pentafluorobenzyl bromide, the resulting derivative would be susceptible to detection by gas chromatography (GC) detectors more sensitive and specific than FID, namely electron capture or electrolytic conductivity detectors. In addition, there appeared to be a chance that oxalic acid could be determined simultaneously with the acetic and formic acids.

B. Materials and Methods

1. Standard Acetic and Formic Acids

Reagent grade formic acid and glacial acetic acid were diluted to approximately 1-g/L concentrations in glass-distilled water. These stock solutions were titrated to pH 7 (formic) or pH 7.5 (acetic) with standard sodium hydroxide on a Metrohm Herisau Dosimat-pH meter combination. The sodium hydroxide was standardized against reagent grade potassium acid phthalate. Appropriate volumes of the two stock solutions were diluted to 100 mL with glass-distilled water to give a combined 100-mg/L standard of the two acids. This standard was then diluted as required for standard curves.

2. Benzyl Formate and Acetate

The following procedure is essentially that presented by Bethge and Lindstrom.⁶ To 10 mL of sample or standard were added 100 μ L of approximately 500-mg/L propionic acid in water as an internal standard. The mixture of acids was passed through a 12-mm diameter x 100-mm length column of 20-50 mesh AG-50W x 8 ion exchange resin (Bio-Rad Laboratories) in the hydrogen form, followed by two 5-mL water washes. The solution passing through the column was collected in a small beaker and titrated to pH 8 with 0.03 N tetra-n-butylammonium hydroxide. The resulting solution was evaporated under aspirator vacuum on a rotary evaporator with a 50°C water bath. A special flask designed by Junk et al.⁸ was used for the evaporations. This flask was made by joining a tapered centrifuge tube to the bottom of a 100-mL round-bottom flask. To form the tetra-n-butylammonium salt derivatives of the acids in the residue, 0.5 mL of redistilled acetone and an appropriate excess of 10% (v/v) benzyl bromide in acetone were added to the residue in the rotary evaporation flask. The flask was stoppered and allowed to stand at least 15 minutes at room temperature.

The benzyl esters formed in the reaction were determined by injecting 1 μ L of the solution onto the gas chromatograph. The GC conditions were:

GC:	Hewlett-Packard 5830A, FID detector
Column:	6 ft x 1/8 in O.D. stainless steel, packed with 3% butane-1,4-diol succinate on 80/100 mesh Chromosorb WHP (Supelco, Inc.)
Carrier:	Helium, 25 cc/minute
Column temperature:	110°C isothermal
FID temperature:	190°C
Injection port temperature:	160°C

3. Pentafluorobenzyl (PFB) Formate and Acetate

To 10 mL of sample or standard were added 100 μ L of approximately 500-mg/L n-butyric acid as an internal standard. This solution was passed through a 12-mm diameter x 100-mm length column of 20-50 mesh AG-50W x 8 ion exchange resin (Bio-Rad Laboratories) in the hydrogen form, followed by two 5-mL portions of water washes. The solution passing through the column was collected in a small beaker and titrated to pH 8 with 0.03 N tetra-n-butylammonium hydroxide. The resulting solution was evaporated under aspirator vacuum at 50°C in the special flasks described previously. The lipophilic tetra-n-butylammonium salts in the residue were dissolved in 0.5 mL of redistilled acetone added directly to the evaporation flask. An appropriate excess of 10% PFB-Br in acetone was added to the flask, which was then stoppered and allowed to stand at least 20 minutes before a 1- μ L sample was taken for GC analysis.

The GC conditions were:

GC:	Hewlett-Packard 5830A, equipped with a Tracor Model 700 Hall electrolytic conductivity detector.
Transfer line temperature:	250°C
Detector resin column:	65% IRN-77 pump side 35% IRN-150 cell side
Detector solvent:	50:50 1-propanol:water (v/v) at 0.8 cc/minute flow
Detector mode:	Quartz tube, in halogen mode with hydrogen reactor gas at 50 cc/minute, 850°C oven temperature
Column:	6 ft x 1/8 in O.D. stainless steel, packed with 3% OV-225 on 100/120 mesh Gas Chrom Q (Applied Science Laboratories, Inc.)
Carrier:	Helium, ultra high purity, at 43 cc/minute
Column temperature:	80°C (4 minutes), 8°C/minute to 175°C (1 minute)
Injection port temperature:	150°C

4. Ozonated Laboratory Wastewater

Table 1 shows the composition of the synthetic reverse osmosis permeate of laboratory wastewater (RO-Lab) from an Army hospital. This wastewater formulation was also used by McCarthy et al.¹ The wastewater was ozonated

Table 1. Synthetic Laboratory Reverse Osmosis Permeate Formulation

<u>Chemical</u>	<u>Concentration^a</u> <u>(mg/L except where noted)</u>
Diethyl ether	1.8
Methanol	289.1
Urea	0.9
Glycerol	2.1
Ethanol	3.3
Formaldehyde (37% sol.)	3.5 µL/L
o-Phenylphenol	0.4
o-Benzyl-p-chlorophenol	0.3
2,4-Xylenol	0.2
2-Propanol	0.1
Acetone	32.8

a. Tap water was used for all batches of this synthetic wastewater.

for 8 hours in the Life Systems Modified Torricelli Ozone Contactor (LMTOC), designed by Life Systems, Inc., Cleveland, Ohio, under Army contract No. DAMD 17-76-C-6041.⁹

The following conditions were used in the ozonation:

Gas flow:	Air at 20 standard cubic ft/hr
O ₃ dose:	1.5 weight %
UV light:	On
pH control:	None
Sample period:	Hourly
Temperature:	Ambient

This formulation and the ozonation protocol were used whenever a source of ozonated wastewater was required for testing analytical methods, at various periods in fall 1977 and spring 1978. No preservative was added to the samples, which were stored at 4°C until analyzed. Hourly samples were collected in glass bottles with Teflon-lined caps, whereas samples for spiking tests were collected in plastic containers or glass jars.

C. Results and Discussion

1. Benzyl Formate and Acetate

Figure 2 shows the separation of these esters from the benzyl bromide and internal standard (benzyl propionate) peaks on the 3% butane-1,4-diol succinate column. The only problem seen was the presence of an artifact peak following the internal standard (I.S.) peak, ascribed to benzyl alcohol by Bethge and Lindstrom.⁶ However, the GC integrator was able to reproducibly split the peak areas, as shown by the least-squares regression of relative peak area (area standard peak/area I.S. peak) versus concentration (Figure 3). For these standard curves, 36 µL of 10% (v/v) benzyl bromide in acetone was added to the standards and blanks. This value was chosen from a titration of a spiked (see below) sample of RO-Lab wastewater after passage through the H⁺ form ion exchange resin, with equation (1) used to calculate the required amount of reagent.

$$V = 36 \mu\text{L} = (0.673 \text{ mL TBAOH}) (0.03 \text{ N TBAOH}) \times \frac{171.04 \text{ mg/mmol}}{1.438 \text{ mg}/\mu\text{L}} \quad (1)$$

x 10 x 1.5

where

factor 1.438 mg/µL = the density of benzyl bromide

171.04 mg/mmol = the molecular weight of benzyl bromide

factor 10 = the correction for use of 10% reagent benzyl bromide

factor 1.5 = the 1.5-fold molar excess chosen for the ratio of
benzyl bromide/titrated acids

0.673 mL TBAOH = volume of TBAOH reagent required for titration of acids
in the sample.




Figure 2. FID-GC of a 2.5-mg/L standard of acetic and formic acids, treated with benzyl bromide.

Because the RO-Lab waste was made up in tap water, the titrated acids included carbonic acid from bicarbonate in the tap water as well as formic and acetic acids added as spikes.

The standard curves in Figure 3 showed sample correlation coefficients (r^2) of 0.99981 (formic) and 0.99945 (acetic) for the range of 1 to 10 mg/L; hence, there was excellent linearity of FID response relative to the internal standard.

Addition of 4.35 mg/L of formic acid and 4.87 mg/L of acetic acid to the RO-Lab wastewater provided a sample with a mid-range level of acids for recovery and precision tests. Analysis of one spiked sample with only 8 μ L of 10% benzyl bromide (calculated to be 1.5 molar excess for a 10-mg/L combined acid standard plus internal standard in distilled water) showed approximate concentrations of 0.7 mg/L (formic) and 4.8 mg/L (acetic). Hence, the carbonic acid in these samples was probably competing effectively with the formic acid for the benzyl bromide, but the acetic acid was not affected. When 36 μ L were added, as required by the titration calculation of equation (1), the recoveries of both acids were comparable. Table 2 shows the data collected from seven replicate aliquots from a single spiked RO-Lab waste. The coefficients of variation were 4.34 and 1.35% for formic and acetic acids, respectively, at the 4- to 5-mg/L level of these acids.

Table 2. Benzyl Method of Recovery of Formic and Acetic Acids from RO-Lab Wastewater

Spiked Sample No.	Concentration Found (mg/L)	
	Formic Acid	Acetic Acid
1	4.04	5.22
2	4.56	5.39
3	4.53	5.42
4	4.32	5.38
5	4.32	5.29
6	4.43	5.33
7	4.18	5.27
Mean value (\bar{x})	4.34	5.33
Standard deviation (s)	0.187	0.072
Coefficient of variation ($s/\bar{x} \times 100$)	4.3%	1.4%
mg/L acid added (z) ^a	4.35	4.87
Average % recovery ($\bar{x}/z \times 100$)	100%	109%

a. No formic or acetic acids were present in the unspiked wastewater.




Figure 3. Standard curves for benzyl bromide derivatization of formic and acetic acids.

As a test to see whether the derivatization reaction conditions (at least 15 minutes, room temperature) were adequate, the data in Table 3 were collected from a spiked RO-Lab sample. It appeared that the peak area ratio was constant by 15 minutes and did not change significantly for periods up to an hour.

Table 3. Kinetic Study of Benzyl Formate and Benzyl Acetate Derivatization Reaction at Room Temperature

Reaction Time (min)	Peak Areas			Relative Peak Areas ^a	
	Formate	Acetate	Propionate	Formate	Acetate
0	0	0	0	0	0
15	124,900	123,000	100,600	1.242	1.223
30	127,600	124,500	102,300	1.247	1.217
45	130,500	126,200	102,900	1.268	1.226
60	137,900	132,200	107,600	1.282	1.233

a. Area formate or acetate/area propionate.

In order to find a detection limit for this procedure, standards of 0, 0.1, 0.25, and 0.5 mg/L formic and acetic acids were run. The data were analyzed by peak height at an attenuation of 2^4 on the Hewlett-Packard recorder. Figure 4 indicates that 0.25 mg/L of these acids could be easily distinguished from the blank responses.

2. Pentafluorobenzyl (PFB) Formate and Acetate

Because the volatilities of the PFB esters were close to the volatility of PFB-Br and the PFB esters were more volatile than the benzyl esters, a lower column temperature and different liquid phase was chosen for the PFB esters. Figure 5 shows the separation of esters from the PFB-Br and side-reaction peaks as detected by the halogen-specific Hall electrolytic conductivity detector. GC/MS studies indicated that the large peak just preceding the PFB-Br peak was pentafluorobenzyl alcohol. The peak between the PFB-acetate and PFB-butyrate was not identified, but may have been an ether of structure PFB-O-PFB, as suggested by Davis.⁷

Because the electronic integration system was not able to reliably skim the acetate and formate peaks from the PFB-Br peak, these peaks were manually skimmed and their peak heights were measured relative to the internal standard peak height. Figure 6 shows the resulting calibration curve for distilled water standards permitted to react with 10 μ L of 10% PFB-Br (1.5 molar excess for a 10-mg/L standard with internal standard). Fairly good correlation coefficients (r^2) of 0.99807 and 0.99176 were obtained for the least-squares regressions of formic and acetic acids, respectively.




Figure 4. Low-level standard curves for benzyl bromide derivatization of formic and acetic acids.




Figure 5. Hall detector-GC of a 2.5-mg/L standard of acetic and formic acids.




Figure 6. Standard curves for pentafluorobenzyl bromide derivatization of formic and acetic acids.

Titration of a spiked RO-Lab wastewater demonstrated a requirement for 45 μL of 10% PFB-Br for a 1.5 molar ratio to titrated acids. Figure 7 shows a chromatogram from such a sample. The PFB alcohol and other artifact peaks were very large, probably because of the PFB-Br reaction with carbonic acid (derived by ion exchange) from the bicarbonate in the tap water used to make up the wastewater. Carbonic acid would react with TBAOH in the titration of the acids to form $(\text{TBA})_2\text{CO}_3$, which would then react with PFB-Br to form $(\text{PFB})_2\text{CO}_3$ and this would possibly hydrolyze to form the alcohol (PFB-OH) and HCO_3^- in the presence of moisture. The relative (to internal standard) peak heights were 2.33 and 3.10 for formate and acetate, respectively, in Figure 7. By comparison, a spiked sample with only 10 μL PFB-Br used to make the acid derivatives is shown in Figure 8. It can be seen that the PFB-Br peak was almost entirely converted to PFB alcohol, and the relative peak heights for formate and acetate were 1.35 and 3.41, respectively. The formate peak, but not the acetate peak, was reduced relative to butyrate; thus, the reaction to form the alcohol competed with the reaction to form PFB-formate but not with that to form PFB-acetate. This is also apparent from the data in Table 4, which show relative peak heights for three different volumes of 10% PFB-Br added to prepare the acid derivatives in a separate test from that shown in Figure 8.

Table 4. Relative Response of Formate and Acetate Peaks to Various Amounts of PFB-Br Used to Derivatize a Spiked RO-Lab Wastewater^a

Volume of 10% PFB-Br Used (μL)	Relative Peak Height	
	PFB-Formate	PFB-Acetate
10	0.83	2.40
20	2.50	2.43
30	1.97	2.47

a. Spike added: 4.35 mg/L formic acid, 4.87 mg/L acetic acid added to raw RO-Lab wastewater containing no acetic or formic acids.

An attempt was made to see whether the interference from bicarbonate in the sample could be eliminated by stripping a spiked RO-Lab wastewater sample with nitrogen gas for 5 minutes after acidification with the ion exchange column. This treatment reduced the TBAOH titer by about 18% and resulted in a chromatogram similar to that shown in Figure 9, with 10 μL of PFB-Br. Comparison with Figure 8 shows that the PFB alcohol was still present, but the formate peak had increased relative to the internal standard peak. Relative peak height ratios of 2.10 (formate) and 2.67 (acetate) were recorded for this sample, in approximate agreement with the values of Table 4 for 20 and 30 μL PFB-Br used. Consequently, the competitive interference from bicarbonate could be reduced by stripping out the carbon dioxide in equilibrium with the carbonic acid of the sample before TBAOH titration, or more PFB-Br could be




Figure 7. Hall detector-GC of a RO-Lab wastewater sample spiked with 4.35 mg/L of formic acid and 4.87 mg/L of acetic acid and derivatized with 45 μ L of 10% (v/v) PFB-Br in acetone. Peak identification is the same as in Figure 5.




Figure 8. Hall detector-GC of a RO-Lab wastewater sample spiked with 4.35 mg/L of formic acid and 4.87 mg/L acetic acid and derivatized with 10 μ L of 10% (v/v) PFB-Br in acetone. Peak identification is the same as in Figure 5.




Figure 9. Hall detector-GC of a RO-Lab wastewater sample treated as in Figure 8 and using N_2 gas stripping to remove CO_2 after acidification by ion exchange. Peak identification is the same as in Figure 5.

used, in accord with equation (1) to cover the extra acids present from tap water. In the latter case, the chromatograms (see Figure 7) were difficult to quantify accurately because of the large by-product peaks. Recovery data (Table 5), for example, showed relatively poor precision for acetic acid and poor recoveries (accuracy) for both acids, compared to the benzyl bromide procedure (Table 2). The data in Table 5 were collected from derivatization using 45 μ L of 10% (v/v) PFB-Br, in accord with equation (1).

Table 5. PFB Method of Recovery of Formic and Acetic Acids from RO-Lab Wastewater

Spiked Sample No.	Concentration Found (mg/L)	
	Formic Acid	Acetic Acid
1	2.86	4.91
2	3.06	5.29
3	3.06	4.40
4	3.02	5.50
5	3.60	7.97
6	3.28	5.88
Mean value (\bar{x})	3.15	5.66
Standard deviation (s)	0.26	1.24
Coefficient of variation (s/ \bar{x} x 100)	8.2%	21.9%
mg/L acid added (z) ^a	4.35	4.87
Average % recovery (\bar{x}/z x 100)	72%	116%

a. No formic or acetic acids were present in the unspiked wastewater.

Figure 10 shows the results of a low-level standard curve run with distilled water standards and 2 μ L of 10% (v/v) PFB-Br added to derivatize the acids. A Hall detector range of 30 with attenuation of 1 to 4 was used to detect these levels. From the figure, it was apparent that 0.25-mg/L levels could be detected.

The PFB-Br procedure was tested unsuccessfully several times with oxalic acid in the combined standard. This lack of success may have been due to the poor nucleophilicity expected of the oxalate anion, because of mutual electron withdrawing by one COO⁻ group from the other. A dipentafluorobenzyl oxalate diester was prepared with PFB-Br and tetra-n-butylammonium oxalate as shown by infrared, NMR, and mass spectra (Appendix A). However, the retention time of this diester was very close to that of the high-boiling reaction by-product eluting just ahead of PFB-butyrate, making quantification difficult. The oxalate diester peak never showed a proportional area increase as a function of oxalic acid concentration, as in the case of acetic and formic acids in the procedure outlined previously. Appendix A shows the mass spectrum of the PFB-oxalate, as well as GC/MS data for PFB-acetate and PFB-formate.




Figure 10. Low-level standard curves for PFB-Br derivatization of formic and acetic acids.

3. Test of Benzyl Acetate and Formate Method with Ozonated RO-Lab Wastewater

Because there appeared to be no difference in sensitivity of the benzyl and pentafluorobenzyl methods, and because the benzyl method was much easier to quantify, the benzyl procedure was used to measure the concentrations of acetic and formic acids in an ozonated RO-Lab wastewater. Figure 11 shows the changes in formic and acetic acids and TOC as a function of time. As expected from the high initial methanol concentration, the formic acid concentration increased to a maximum early in the ozonation. Acetic acid, primarily from acetone breakdown, did not reach a peak until later in the run. The only point which did not appear to fall on a smooth curve was the 180-minute point for both acids. No explanation is currently available for this decrease in concentration from the values on either side of 180 minutes. The samples from 60 through 300 minutes were diluted 1:10 for this analysis, and this may not have been the optimal dilution for all of these samples. The final concentrations of 1.48 mg/L of formic acid and 0.5 mg/L of acetic acid accounted for 0.4 and 0.2 mg/L of a final TOC, which could only be measured as < 5 mg/L on a Dohrman DC-50 carbon analyzer.

IV. OXALIC ACID

A. Literature Review

Oxalic acid has been determined fluorimetrically by a rather complex procedure involving extraction with tri-n-butyl phosphate, co-precipitation with calcium sulfate, reduction to glyoxylic acid, then reaction with resorcinol to form a fluorescent complex.¹⁰ A linear fluorescence curve was found between 0 and 5 mg/L oxalic acid.

Oxalic acid has also been determined by GC after derivatization to its bis-trimethylsilyl (TMS) ester, dimethyl ester, or diisopropyl ester. Rumsey and Noller¹¹ used hydrogen chloride and methanol to form the methyl ester in the range of about 45 to 630 mg/L. Mee and Stanley¹² used the same reaction in the range of 2 to 25 µg of oxalic acid (derivatized) injected onto the GC column. This would be a concentration of 2,000 to 25,000 mg/L in the injected sample, for a 1-µL injection. More recently, Chian and Kuo² used a small-scale diazomethane procedure similar to that of Schlenk and Gellerman¹³ to esterify oxalic acid, which was concentrated from a water sample by rotary evaporation, in the range of about 2 to 30 mg/L. Gould and Weber³ reacted dried residues of wastewater samples with 2-propanol and sulfuric acid, then extracted with chloroform, and analyzed the extract on the GC for diisopropyl oxalate.

For trimethylsilylation, Butts and Rainey¹⁴ made the ammonium salt of oxalic acid and then reacted it with bis-(trimethylsilyl) trifluoroacetamide (BSTFA) in dimethylformamide. No quantification was attempted. Von Nicolai and Zilliken¹⁵ concentrated oxalic acid from urine by rotary evaporation and then permitted the acid to react with a mixture of pyridine-BSTFA-trimethylchlorosilane (TMCS) to form the trimethylsilyl (TMS) ester. Quantification by




Figure 11. Ozonation of synthetic RO-Lab wastewater, with analyses of formic and acetic acids by the benzyl bromide method of Bethge and Lindstrom.⁶

GC was attempted in the range of 2 to 16 μ g of oxalic acid injected as its derivative in a 1- μ L injection or 2,000- to 16,000-mg/L solution concentration.

B. Materials and Methods

1. Trimethylsilylation of Oxalic Acid

The derivatizing reagent of Von Nicolai and Zilliken¹⁵ was used with the ammonium salt of oxalic acid as proposed by Butts and Rainey.¹⁶ Because of the volatility of oxalic acid (melting point 101°C, sublimation beginning at 100°C), it was expected that losses of the acid might be experienced during rotary evaporation of the water; this would not be the case with the ammonium salt.

A 10-mL sample or standard was passed through a 12-mm diameter x 100-mm length column of AG-50W x 8 resin (Bio-Rad Laboratories), 20-50 mesh, in the NH_4^+ form. The column was then washed with two 5-mL portions of distilled water. All solutions were collected from the column in a 50-mL round-bottom flask, and this solution was evaporated under aspirator vacuum on a rotary evaporator with a 56°C water bath. The dry flask containing the oxalate was purged with nitrogen gas to remove moisture and then 200 μ L of 5:3:0.1 (v/v/v) pyridine:BSTFA:TMCS reagent was added. The flask was stoppered and heated for 2 minutes at 80°C and then cooled. A 20- μ L volume of 10 μ L/10 mL dodecane in pyridine was added as an internal standard, and 2 μ L was injected onto the GC.

The GC conditions were:

GC:	Hewlett-Packard 5830A, FID detector
Column:	6 ft x 2 mm inner diameter glass, packed with 10% OV-1 on 80/100 mesh Chromosorb WHP (Alltech Associates)
Carrier:	Nitrogen at 24 cc/minute
Column temperature:	150°C (6 minutes), 4°C/minute to 125°C (2 minutes)
FID temperature:	220°C
Injection port temperature:	200°C

All silylation reagents were silylation grade (Pierce Chemical Co.). Oxalic acid was dried by placing ACS grade oxalic acid dihydrate in an oven at 98° to 99°C for 2 hours and then storing over P_2O_5 in a desiccator. Purity by pH titration with standard sodium hydroxide was found to be 100.2%.

2. Methylation of Oxalic Acid

The procedure evaluated in this study was recommended by Chian and Kuo,² and is derived from that proposed by Schlenk and Gellerman.¹³ However, Chian and Kuo used a butanediol succinate (polyester) GC column packing, whereas the procedure reported below used an OV-225 (cyanopropyl phenyl silicone) column packing. In addition, ether was used as a solvent instead of methylene chloride for the derivatization.

A 20-mL aliquot of sample or standard was passed through a 12-mm diameter x 100-mm length column of 20-50 mesh AG-50W x 8 ion exchange resin (Bio-Rad Laboratories) in the hydrogen form, followed by 5 mL of distilled water. The effluent from the resin was evaporated at 30°C on a rotary evaporator under water aspirator vacuum in a 50-mL round-bottom flask. Five mL of laboratory grade ether, which had been redistilled in glass before use, were added, and the exit tube of the diazomethane generator (Figure 12) was dipped into the ether to bubble in CH_2N_2 gas (performed in a fume hood) until the solution turned yellow, indicating an excess of CH_2N_2 . Diazomethane was generated from a mixture of 2 mL of potassium hydroxide (60 g/100 mL H_2O) plus approximately 0.3 g of Diazald (Aldrich Chemical Co., N-methyl-N-nitroso-p-toluenesulfonamide) in the second tube. The first tube contained ether to scrub the helium sweep gas. Rubber stoppers and plastic and glass tubing were used for flexible connections between the helium source and between the tubes.

The ether solution was poured into a 10-mL Kuderna-Danish (Kontes No. K-570050) concentrator tube, the round-bottom flask was washed with 2 to 3 mL of ether, and the washings were added to the concentrator tube (final volume 6 to 8 mL).

An evaporative concentrator microreflux column (Kontes No. K-569251) was placed on the tube and heated in a Kuderna-Danish tube heater block maintained at 70°C, with bumping protection provided by capillary tubes that bubbled nitrogen gas into the tubes. The solutions were concentrated to a minimum of 0.5 mL (caution was taken not to allow further evaporation). Finally, 30 μL of 10 $\mu\text{L}/10\text{ mL}$ hexadecane in methylene chloride was added as an internal standard, and 2 μL were injected onto the gas chromatograph with conditions as follows:

GC:	Hewlett-Packard 5830A, FID detector
Column:	6 ft x 1/8 in O.D. stainless steel, packed with 3% OV-225 on 100/120 mesh Gas Chrom Q (Applied Science Laboratories, Inc.)
Carrier:	Nitrogen at 30 cc/minute
Column temperature:	90°C (4 minutes), 15°C/minute to 150°C (3 minutes)
FID temperature:	200°C
Injection port temperature:	150°C




Figure 12. Diazomethane generator.

C. Results and Discussion

1. (TMS)₂ Oxalate

Figure 13 shows a typical gas chromatogram of (TMS)₂ oxalate, from analysis of a 2.0-mg/L standard of oxalic acid. The peak was well separated from the reagent peaks. Figure 14 shows a calibration curve prepared from distilled water standards. Since the line of best fit for the data did not appear to intersect the x-axis at zero, a constant loss of material in the procedure was indicated. To check the FID response, a 20-mg/L oxalic acid sample was derivatized and the resulting (TMS)₂ oxalate was diluted serially with derivatizing reagent and injected onto the GC. Figure 15 shows that the response curve still did not intersect at zero; thus, an irreversible column adsorption or breakdown of the (TMS)₂ oxalate was suspected. The test was repeated after the column had been conditioned with Silyl-8 (Pierce Chemical Co.) silylating reagent at 175°C (300 µL), but the results were the same. Silylation of a 6-mg/L standard with 1/4-strength reagent mixture (diluted with pyridine) gave an area ratio of 1.321 relative to dodecane. After the mixture stood overnight, the ratio dropped to 0.205. However, injection of 3 µL of Silyl-8 at 180°C onto the GC column, followed by another injection of the standard, showed an increase in the ratio to 0.871. A repeated (3 µL) GC column conditioning with Silyl-8 before sample injection gave a ratio of 1.089. Finally, preconditioning with 12 µL of Silyl-8 showed an area ratio of 1.528 and 1.476 for duplicate injections of the standard. Further tests indicated that injection of TMS reagent after an injection of 10-mg/L standard oxalate (peak area = 282,700 units) would result in a "ghost" peak of TMS oxalate (area = 42,910 units), so that column adsorption of (TMS)₂ oxalate was reversible or the column breakdown of the ester was reversible by injection of excess reagent. It was apparent that at low levels this method could lead to very confusing results, and further evaluation was not attempted.

2. Dimethyl Oxalate

Because of suspected breakdown¹⁶ of dimethyl oxalate on the 1,4-butanediol succinate polyester column used by Chian and Kuo,² a 3% OV-225, polar silicone column was chosen to determine this ester for these studies.

Although methylene chloride was suggested by Chian¹⁷ for use as a solvent for the diazomethane reaction, the solubility of oxalic acid in this solvent was found to be low. Consequently, ether was tried, because this solvent readily dissolved the oxalic acid, and Schlenk and Gellerman¹⁸ had used ether (with 10% methanol) as a solvent for diazomethane's reaction with fatty acids.




Figure 13. Gas chromatogram of TMS oxalate reaction mixture (2 mg/L oxalic acid standard).




Figure 14. Standard curve for oxalic acid by the TMS method.




Figure 15. Serial dilution of TMS oxalate reaction mixture (20 mg/L oxalic acid standard) with pyridine-TMCS-BSTFA reagent.

Figure 16 shows the separation of dimethyl oxalate from the solvent peak and the hexadecane internal standard under the GC conditions described previously. Early work at low levels, with rotary evaporation temperatures of 50°C or more, showed chromatograms without the oxalate peak. Figure 17 presents a typical curve for 50°C rotary evaporated distilled water standards, plotted without the blank value, to show that the data (least-squares regression) extrapolated to 66 µg of oxalic acid on the x-axis for a relative area of zero. This indicated a loss of 66 µg from all the standards. When oxalic acid standards in ether were converted to derivatives and analyzed on the GC without the rotary evaporation step, the least-squares regression curve extrapolated to a value of 5 µg on the x-axis, which was much closer to the expected value of zero (Figure 17). Hence, most of the loss must have occurred during rotary evaporation at 50°C. Table 6 indicates that the results were low or erratic at 40° or 50°C, whereas the evaporation gave higher and more consistent results at 30°C.

Table 6. Relative Peak Areas from Analysis of 20 mL
of 9.83 mg/L Oxalic Acid Standard, at 30° to 50°C
Rotary Evaporation Bath Temperatures

<u>Temperature (°C)</u>	<u>Relative Peak Area of Dimethyl Oxalate</u>
30	1.118
30	1.115
40	0.766
40	1.111
50	0.726
50	0.798

Even at 30°C, the rotary evaporation still caused small losses of oxalic acid (Figure 18). The two extrapolated curves in Figure 18 show 1.3 mg/L on the x-axis, so that a lower limit of detection would have to be somewhere above 1.3 mg/L. A low-level standard curve (Figure 19) indicated that about 2 mg/L was the lower limit of detection.

Recovery of oxalic acid from synthetic RO-Lab wastewater that had been ozonated for 4 hours is given in Table 7. The recovery at the level tested (9.9 mg/L) was 97%, and the coefficient of variation was 7%.

Table 7. Dimethyl Oxalate Method of Recovery of Oxalic Acid
from Ozonated RO-Lab Wastewater

<u>Spiked Sample No.</u>	<u>Concentration Found (mg/L)</u>
1	10.00
2	9.68
3	9.68
4	10.59
5	8.32
6	9.55
7	<u>9.62</u>
Mean value (\bar{x})	9.63
Standard deviation (s)	0.68
Coefficient of variation (s/ \bar{x} x 100)	7%
mg/L acid added (z) ^a	9.9
Average % recovery (\bar{x}/z x 100)	97%

a. Oxalic acid in unspiked wastewater = 0 mg/L by analysis.

Analysis of samples from an 8-hour ozonation of RO-Lab wastewater, taken at 60, 180, 300, and 420 minutes, showed no oxalic acid at these times during the ozonation; hence, the concentration was less than 2 mg/L in this wastewater at the intervals sampled. Two large peaks at retention times greater than that of dimethyl oxalate were seen. Qualitative standards of malonic, succinic, maleic, and fumaric acids were run with diazomethane, but these did not match the unknown peak retention times.




Figure 16. Gas chromatogram of reaction mixture from a 5-mg/L oxalic standard with diazomethane.




Figure 17. Standard curves for dimethyl oxalate: (a) from aqueous standards using rotary evaporation at 50°C, then methylation in ether for GC; (b) from standards in ether, methylated for GC.




Figure 18. Standard curves for oxalic acid with the 30°C rotary evaporation method, 1 March and 24 March 1978.





Figure 19. Low-level standard curve for oxalic acid.

V. GLYOXALS

A. Literature Review

Glyoxal, OHC-CHO, has been determined colorimetrically by the formation of bis(2,4-dinitrophenylhydrazone) (DNP). Gould and Weber³ measured the loss in color of a standard solution of 2,4-DNP as a measure of the amount of glyoxal reacted (molar extinction coefficient = 13,200 L/mole·cm at 365 nm for DNP). Johnson et al.¹⁸ used a mildly basic solution of diethanolamine in pyridine to develop the color of the phenylhydrazone formed. This method was used by Gilbert¹⁹ to follow the production of glyoxal during ozonation of organic compounds at 1-mM/L initial concentrations.

Gas chromatographic methods have been developed to determine glyoxal, methylglyoxal (pyruvaldehyde), and dimethylglyoxal (2,3-butanedione) in a single chromatogram, using the reaction of glyoxals with ortho-phenylenediamine to form quinoxalines (equation 2):



R = H (in glyoxal and methylglyoxal) or CH₃ (in methyl- and dimethylglyoxal)

Nojima et al.²⁰ used 4-chloro-o-phenylenediamine to detect glyoxals in air pollution studies, with electron capture (⁶³Ni) gas chromatography for quantification. Chian and Kuo²¹ developed a method for glyoxals in aqueous solution by FID-GC detection. They reacted a 10-mL water sample with 1 mL of 0.185 M methanolic-o-phenylenediamine solution and 1 mL of 10% aqueous sodium bisulfite for 30 minutes at room temperature. Benzene was added to extract the quinoxalines and the extract was evaporated to 0.5 mL in a micro Kuderna-Danish evaporator and analyzed on a 3% OV-17 GC column. Yields of quinoxaline (Q), methylquinoxaline (MQ), and dimethylquinoxaline (DMQ) (from glyoxal, methylglyoxal, and dimethylglyoxal, respectively) were 81%, 32%, and 33%.

Extraction efficiencies were 27%, 85%, and 83%, whereas the solvent evaporation efficiencies were 95%, 98%, and 94%, respectively. Detection limits were reported as 0.08, 0.05, and 0.2 mg/L. Precision estimates were 2%, 4.3%, and 2.5% (relative standard deviation) at the 10-mg/L level.

B. Materials and Methods

The method of Chian and Kuo¹¹ was modified by heating the reaction mixture at 80°C for 30 minutes in sealed (Teflon-lined caps) test tubes instead of shaking 30 minutes at room temperature. In addition, methylene chloride was used as a solvent instead of benzene, and a more polar (OV-225 instead of OV-17) GC column was chosen.

Ten mL of combined glyoxal standard or sample were pipetted into a 30-mL test tube. One mL each of 0.185 M o-phenylenediamine (Aldrich Chemical Co., recrystallized from ether after treatment with Norit carbon and filtration) in methanol and 10% (wt/vol) aqueous sodium bisulfite were added, and the tubes were capped and shaken to mix the solutions. The mixture was heated in an 80°C water bath for 30 minutes, cooled, and shaken with 5 mL of methylene chloride for 1 minute. Four mL of the lower phase was withdrawn with a pipette and evaporated to 0.5 mL in a 10-mL micro Kuderna-Danish concentrator tube with reflux column and nitrogen gas ebullator. The heating block temperature was set at 80°C. A 1/1,000 (v/v) solution of hexadecane in hexane was added to each concentrate as an internal standard, and 2 µL of the concentrate was injected onto the GC with the following conditions:

GC:	Hewlett-Packard 5830A, FID detector
Column:	6 ft x 1/8 in O.D. stainless steel, packed with 3% OV-225 on 100/120 mesh Gas Chrom Q (Applied Science Labora- tories, Inc.)
Carrier:	Helium at 30 cc/minute
Column temperature:	130°C (2 minutes), 10°C/minute to 200°C (2 minutes)
FID temperature:	220°C
Injection port temperature:	200°C

Glyoxal was purchased as a 40% aqueous solution (Pfaltz and Bauer Chemical Co.), which was found to have a density of 1.26 g/mL. Dimethylglyoxal (2,3-butanedione) was obtained as a 99% pure liquid (density = 0.9809 g/mL) from Aldrich Chemical Co. Methylglyoxal (density = 1.045 g/mL) was obtained as a 99% pure liquid from Pfaltz and Bauer. A stock standard was made up with 250 µL of 40% glyoxal solution, 100 µL of methylglyoxal, and 100 µL of dimethylglyoxal diluted to 100 mL with water. The concentrations of this stock solution were calculated to be:

1,260 mg/L glyoxal
1,000 µL/L methylglyoxal = 1,045 mg/L
1,000 µL/L dimethylglyoxal = 981 mg/L

The concentrations of methyl- and dimethylglyoxal in this standard were considered as 1,000 mg/L to within 5% error, based on the densities of the pure liquids, and subsequent concentrations are reported as mg/L in this report.

For recovery studies, 500 mL of RO-Lab wastewater were spiked with 10 mL of 1/10 diluted stock glyoxal combined standard to give concentrations of 2.47 mg/L of glyoxal and 1.96 mg/L of methylglyoxal and dimethylglyoxal.

C. Results and Discussion

Substitution of OV-225 for OV-17 as the liquid phase for gas chromatography did not change the order of elution of quinoxalines, but did shift the peak due to unreacted o-phenylenediamine from its position between quinoxaline (Q) and methylquinoxaline (MQ), as reported by Chian and Kuo,²¹ to a position after the dimethylquinoxaline (DMQ) peak. This peak was generally very large and tailed significantly, but caused no problem because it was at the end of the chromatogram (Figure 20).

A standard curve for low levels of glyoxals (0.25 to 6.3 mg/L range) is shown in Figure 21. The curves were linear and showed that the sensitivity followed the order glyoxal < methylglyoxal < dimethylglyoxal. The lower limit of detection was no greater than 0.25 mg/L in these tests; no attempts were made to try lower concentrations. The data in Figure 21 were obtained with methylene chloride as an extractant for the quinoxalines. When pentane was used, the relative areas for Q, MQ, and DMQ were reduced by about 86%, 45%, and 22%, respectively. The amount of o-phenylenediamine extracted also decreased when pentane replaced methylene chloride. Chian and Kuo²¹ reported an extraction efficiency of only 27% for Q with benzene, which resembles the reduced recovery seen with pentane. Hence, methylene chloride was a much better solvent for quinoxaline extraction, and the extra amount of o-phenylenediamine extracted did not interfere with determination of the quinoxalines on the OV-225 column.

Table 8 gives the recoveries of glyoxals added to synthetic RO-Lab permeate wastewater. At approximately 2-mg/L concentrations, the recoveries averaged 79 to 88%. The variation in glyoxal recovery was much higher than that of methyl- or dimethylglyoxal. A test with distilled water standards showed better precision for glyoxal (11% coefficient of variation, for four replicates at 3.15 mg/L) and approximately the same precision values as those in Table 8 for the other glyoxals. The low recoveries for all the glyoxals and poor precision for glyoxal were probably due to the low level of the spike, which was only about 10 times the estimated limit of detection.

The effect of temperature on the derivatization reaction was not studied here. Chian and Kuo²¹ indicated that heating was not necessary in tests with 10 mL of 1,000 ppm of glyoxals. An 80°C temperature was used here to ensure a complete reaction even for low levels of glyoxals.

The results of an 8-hour ozonation of RO-Lab wastewater are presented in Table 9. Glyoxal concentration was less than 0.25 mg/L in all samples. Methylglyoxal increased to about 3 mg/L after 360 minutes and disappeared by 420 minutes. Dimethylglyoxal showed 0.4 mg/L in the 60-minute sample and undetectable concentrations thereafter. Hence, of the three glyoxals, only methylglyoxal was a significant reaction product in the ozonation of this wastewater.

Table 8. Recoveries of Glyoxals from RO-Lab Wastewater

Spiked Sample No.	Concentration Found (mg/L)		
	Glyoxal	Methylglyoxal	Dimethylglyoxal
1	1.19	1.46	1.70
2	1.02	1.55	1.83
3	3.17	1.73	1.72
4	1.53	1.50	1.81
5	1.82	1.74	1.81
6	2.37	1.62	1.54
7	2.64	1.68	1.72
Mean value (\bar{x})	1.96	1.61	1.73
Standard deviation (s)	0.79	0.11	0.10
Coefficient of variation ($s/\bar{x} \times 100$)	40.0%	6.9%	5.8%
mg/L added (z) ^a	2.47	1.96	1.96
Average % recoveries ($\bar{x}/z \times 100$)	79%	82%	88%

a. No glyoxals were detected in the unspiked wastewater.

Table 9. Production of Glyoxals During Ozonation of RO-Lab Wastewater

O ₃ Time (min)	mg/L		
	Glyoxal	Methylglyoxal	Dimethylglyoxal
0	<0.25	<0.25	<0.25
60	<0.25	0.59	0.41
120	<0.25	1.05	<0.25
180	<0.25	1.79	<0.25
240	<0.25	2.45	<0.25
300	<0.25	3.18	<0.25
360	<0.25	3.26	<0.25
420	<0.25	<0.25	<0.25
480	<0.25	<0.25	<0.25




Figure 20. Gas chromatogram of quinoxaline mixture on OV-225.





Figure 21. Standard curves for glyoxals.

VI. FORMALDEHYDE

A. Literature Review

Formaldehyde has been determined in low ppm amounts by gas chromatography and colorimetric/fluorimetric methods. Because of the low FID detector response to formaldehyde, TC detectors have had to be used, with limited sensitivity. Hammarstrand²² let formaldehyde react with dimedone (5,5-dimethyl-1,3-cyclohexanedione) to form a derivative from one molecule of formaldehyde and two of dimedone, thus increasing the FID response. Kallio et al.²³ demonstrated that the 2,4-dinitrophenylhydrazone of formaldehyde could be detected in less than 100- to 200-pg amounts by electron capture detection on an SE-30 column. However, the derivatives were prepared in milligram amounts, not at the low levels expected in wastewaters.

Colorimetric analyses of formaldehyde have often been performed with chromotropic acid in hot, strong sulfuric acid. A milder reaction was found to be the Hantzsch reaction, in which ammonia reacts with formaldehyde and a β -diketone or a β -keto ester to form a dihydropyridine derivative. An example is the method used by Nash²⁴ with acetylacetone (equation 3):



The product, diacetyldihydrolutidine (DDL), was found to have a molar extinction of 8,000 L/mole \cdot cm at 412 nm. Belman²⁵ applied this reaction to the assay of formaldehyde at sub-ppm levels, using excitation at 410 nm and fluorescence at 510 nm. In an extensive literature review by Bartos and Pesez,²⁶ the uses of ethyl acetoacetate, 2,4-pentanedione, and dimedone for fluorimetric analysis of formaldehyde were described. In addition, J-acid (6-amino-1-naphthol-3-sulfonic acid) has been used to form a fluorescent derivative of formaldehyde. These reagents have been utilized in the determination of alpha glycolic compounds via analysis of the formaldehyde generated on reaction of the compound with permanganate or periodate.²⁷ Because of their specificity and ease, the methods of Belman and Nash were chosen for use with ozonated RO-Lab wastewaters.

B. Materials and Methods

Color-developing reagent was prepared by combining 150 g of ammonium acetate and 3 mL of glacial acetic acid with 2 mL of redistilled acetylacetone and diluting the mixture to 1 liter with water.

Ten mL of color-developing reagent was mixed with 10 mL of formaldehyde standard or 10 mL of a sample containing less than 8 mg/L of formaldehyde. The mixture was heated in glass test tubes with Teflon-lined screw caps at 58°C for 10 minutes and then allowed to cool.

Fluorescence was read on a Perkin-Elmer Model 204 fluorescence spectrophotometer with 1-cm quartz cuvettes, using $\lambda(\text{excitation}) = 440 \text{ nm}$, $\lambda(\text{emission}) = 510 \text{ nm}$.

Absorbance was read on a Bausch & Lomb Spectronic 700 UV-visible spectrophotometer at 410 nm with 1-cm or 5-cm quartz cells.

Standard formaldehyde solution was prepared by weighing approximately 6.0 g of 40% (wt/vol) formaldehyde (ACS grade) solution and transferring to a 100-mL volumetric flask. Ten mL of 1 N sulfuric acid (H_2SO_4) were added and the solution allowed to stand 10 minutes. Then 10 mL of 1 N sodium hydroxide were added and the mixture was diluted to 100 mL. This stock solution was standardized by measuring 25 mL of the solution into an Erlenmeyer flask, adding 2 drops of thymolphthalein indicator solution (1.0 g dye/L ethanol) and neutralizing with dilute acid or base. In another flask, 25 mL of sodium sulfite solution (125 g/L in water) was similarly neutralized with thymolphthalein indicator. The formaldehyde and sodium sulfite solutions were combined and titrated with standard 1 N sulfuric acid (1 mL 1 N H_2SO_4 is equivalent to 30.03 mg formaldehyde).

C. Results and Discussion

Figure 22 shows the linear standard curve obtained from 0 to 8 mg/L of formaldehyde, analyzed after the Hantzsch reaction with a spectrophotometer. A standard curve from 0 to 0.4 mg/L formaldehyde, analyzed by fluorimetry, is shown in Figure 23. Standards of 0.005 and 0.010 mg/L did not show a detectable response at the settings used (10 x 1), but the response to 0.050 mg/L was readily seen. The interference from a 10-mg/L glyoxylic acid standard was equivalent to only 0.004 mg/L of formaldehyde with this method.

It should be noted that Belman²⁵ reported an excitation maximum at 410 nm, whereas the maximum with the Perkin-Elmer unit was found to be at 440 nm, perhaps because of mechanical problems with the instrument wavelength settings. The lower limit of 0.01 mg/L formaldehyde by fluorimetric analysis reported by Belman²⁵ was not achieved here. A practical limit of detection should probably be about 0.02 mg/L. With a 5-cm cell, spectrophotometric analysis showed an absorbance of 0.095 for a 0.1-mg/L standard of formaldehyde.




Figure 22. Standard curve for formaldehyde in a test using spectrophotometric detection.




Figure 23. Standard curve for formaldehyde in a test using fluorescence detection.

The precision of the spectrophotometric method was tested on replicates of an ozonated synthetic laboratory wastewater RO permeate (RO-Lab) that was diluted by 1/10 with water for analysis. Table 10 shows that the coefficient of variation was only 4%.

Table 10. Precision of Colorimetric Determination of Formaldehyde with Diluted O₃/RO-Lab Permeate

<u>Sample No.</u>	<u>Formaldehyde (mg/L)</u>
1	4.49
2	5.09
3	4.54
4	4.60
5	4.57
6	4.63
7	4.55
Mean \bar{x}	4.64
Variance	4.16×10^{-2}
Standard deviation	2.04×10^{-1}
Coefficient of variation	4%

Addition of a formaldehyde spike to the ozonated RO-Lab solution (Table 11) showed that an average of 83% of the expected value was found. As in the precision test, all solutions were diluted by 1/10 before an analysis to get the concentrations into the range of 0 to 8 mg/L.

Table 11. Accuracy of Colorimetric Determination of Formaldehyde with Spiked O₃/RO-Lab Permeate^a

<u>Sample No.</u>	<u>Formaldehyde (mg/L)</u>
1	6.08
2	6.19
3	6.08
4	6.02
5	6.08
6	6.08
7	6.13
Mean \bar{x}	6.09
Variance	2.795×10^{-3}
Standard deviation	5.29×10^{-2}
Spike recovered	1.45
Percent recovery	83%
Coefficient of variation	1%

a. Initial concentration = 4.64 mg/L; added spike = 1.74 mg/L.

The progress of ozonation of synthetic RO-Lab wastewater was followed with fluorimetric analyses for formaldehyde, although the levels were much higher than expected and dilutions of 1/20 or 1/40 were necessary to permit work within the linear standard curve from 0 to 2.0 mg/L. These analyses could have been performed much more easily with the spectrophotometric methods outlined previously. Table 12 shows the progress of two ozonations on separate days. In both runs, the formaldehyde peaked at around 180 minutes. Hence, the chemical analyses appeared able to track the progress of the ozonation with respect to formaldehyde.

Table 12. Formaldehyde Formation During Two Ozonations of RO-Lab Wastewater

Time (min)	Formaldehyde (mg/L)	
	Test 1	Test 2
0	0.8	0.8
60	53.6	16.8
120	56.0	49.2
180	59.2	56.4
240	44.0	34.8
300	24.2	23.6
360	0.2	10.6
420	<0.1	<0.1
450	0.8	No data
480	<0.1	1.0

VII. GLYOXYLIC ACID

A. Literature Review

Gas chromatography of glyoxylic acid requires the formation of volatile derivatives. Among those that were synthesized for possible GC applications have been: (a) O-benzyl oxime trimethylsilyl esters, (b) O-methyl oxime trimethylsilyl esters, and (c) 1,3-diphenylimidazolidine-2-carboxylic acid trimethylsilyl ester.²⁸ Kallio and Linko²⁹ synthesized the 2,4-dinitro-phenylhydrazone of glyoxylic acid and then made the methyl ester of that product. However, in these cases, there was no attempt to analyze for glyoxylic acid in dilute aqueous (1 to 10 ppm) solution.

A colorimetric method using oxidation of glyoxylic acid phenylhydrazone to form 1,5-diphenylformazancarboxylic acid was proposed by Kramer et al.,³⁰ but formaldehyde was an interfering compound.

A relatively specific and sensitive procedure was reported by Zarembski and Hodgkinson,³¹ who permitted glyoxylic acid to react with resorcinol and measured the product by its fluorescence. Their method appeared to be the best candidate for evaluation in the studies reported here.

B. Methods and Materials

Potassium Carbonate-Bicarbonate Buffer (pH 9.6). Two hundred mL of 0.1 M K_2CO_3 and 800 mL of 0.1 M $KHCO_3$ were mixed.

Glyoxylic Acid Stock. A 1,000-mg/L solution was prepared with glyoxylic acid in distilled water. This solution was diluted to a 100-mg/L working solution.

Reagents. HCl, resorcinol, K_2CO_3 , and ascorbic acid were prepared as aqueous solutions.

Procedure. Each sample (2 mL) or standard was placed in a screw-cap test tube capable of holding 20 mL. The cap was lined with a Teflon disk. To each sample or standard was added 1 mL of 3 N HCl, 0.1 mL of 5% (wt/vol) resorcinol, and 1.5 mL concentrated HCl. The tubes were capped and heated in boiling water for 5 minutes; then they were cooled in an ice bath. Five mL of isobutanol were added and the tube was shaken for 2 minutes and allowed to stand for phase separation. The isobutyl alcohol phase was transferred to a 25-mL stoppered tube, 10 mL of 25% (wt/vol) K_2CO_3 were added, and the tube was shaken again for 2 minutes. After settling, the aqueous phase was transferred to a 25-mL volumetric flask, 1 mL of 10% (wt/vol) ascorbic acid was added, and the volume was made up to 25 mL with potassium carbonate-bicarbonate buffer. After mixing, the solution was allowed to stand for 30 minutes. The fluorescence was measured on a Perkin-Elmer Model 204 spectrofluorimeter with 1-cm quartz cuvettes at an excitation wavelength of 490 nm and an emission wavelength of 530 nm.

C. Results and Discussion

A linear calibration curve was found for the range of 0 to 5 mg/L glyoxylic acid (Figure 24). The detection limit was no higher than 0.1 mg/L for the method.

The precision of the glyoxylic acid procedure was checked with replicate analyses of an ozonated RO-Lab permeate, as shown in Table 13. The coefficient of variation was about 15% at the low level tested. When a spike was added to increase the concentration to 0.88 mg/L, the precision as expressed by coefficient of variation was about 8% (Table 14). The percent recovery of the spike was 96%.




Figure 24. Standard curve for glyoxylic acid.

Table 13. Precision of Glyoxylic Acid Method
(Low-Level Sample)

<u>Sample No.</u>	<u>Glyoxylic Acid (mg/L)</u>
1	0.16
2	0.21
3	0.14
4	0.14
5	0.16
6	0.16
7	<u>0.16</u>
Mean \bar{x}	0.16
Variance	0.0005
Standard deviation	0.0234
Coefficient of variation	15%

Table 14. Recovery of Glyoxylic Acid Added to
Ozonated RO-Lab Permeate^a

<u>Sample No.</u>	<u>Glyoxylic Acid (mg/L)</u>
1	0.89
2	0.84
3	0.74
4	0.84
5	0.92
6	0.92
7	<u>0.79</u>
Mean \bar{x}	0.85
Variance	0.0045
Standard deviation	0.0674
Spike recovered	0.69
Percent recovery	96%
Coefficient of variation	8%

a. Unspiked concentration = 0.16 mg/L; added spike = 0.72 mg/L.

Analysis of samples from ozonation of RO-Lab permeate showed very low levels of glyoxylic acid (Table 15). A maximum of only 0.4 mg/L was seen after 180 to 240 minutes of ozonation. These values may be low due to decreases in glyoxylic acid during 12-day storage at 4°C. The effect of storage time was not tested in these studies.

Table 15. Production of Glyoxylic Acid During Ozonation
of Synthetic RO-Lab Wastewater

<u>Time (min)</u>	<u>Glyoxylic Acid (mg/L)</u>
0	<0.1
60	0.2
120	0.2
180	0.4
240	0.4
300	0.2
360	<0.1
420	<0.1
480	<0.1
485	0.1

VIII. CONCLUSIONS

A. Acetic and Formic Acids

The method of Bethge and Lindstrom⁶ was found to be effective for analysis of acetic and formic acids. There was no advantage to the use of pentafluorobenzyl bromide over benzyl bromide, even though the Hall detector showed very good sensitivity to the PFB esters. The key disadvantage of the PFB-Br system was the more difficult GC resolution of the PFB-acetate and formate esters from the unreacted PFB-Br. In both the PFB-Br and benzyl bromide procedures, the carbonate alkalinity of the samples interfered more with the esterification of formic acid than acetic acid. The interference could be reduced by purging the samples with inert gas after passage of the samples through acidic ion exchange resin or by adding an excess of benzyl bromide reagent to cover the carbonic acid formed with the ion exchange resin. Although Bethge and Lindstrom⁶ suggested computation of the amount of benzyl bromide reagent needed for each sample from equation (1), a more practical procedure would be to spike a control sample (containing no ozonated products) with 10 mg/L of each of the acids to be determined, plus internal standard, then to titrate this sample to calculate the amount of reagent to be added to all samples and standards, with equation (1). Since 0- to 10-mg/L calibration curves prepared with one level (1.5 molar excess over the 10-mg/L level) of reagent were linear, the presence of excess benzyl bromide reagent did not interfere at the lower concentrations. Acetic and formic acid levels as low as 0.25 mg/L were easily reached when the benzyl bromide method was used. Recoveries were essentially quantitative and the precision was 1 to 4% of the mean (relative standard deviation) at 4 to 5 mg/L of added formic or acetic acids.

B. Oxalic Acid

Attempts to quantitatively derivatize 0- to 10-mg/L solutions of oxalic acid with pentafluorobenzyl bromide were not successful. Although a TMS diester was easily formed, it apparently was not very stable during gas chromatography and TMS esterification was discarded as a method for low levels of oxalic acid. The dimethyl ester, however, was readily prepared and chromatographed, although care still had to be taken to avoid losses of oxalic acid during rotary evaporation of water. A temperature of 30°C was found to be a good compromise temperature for rotary evaporation. A lower level of detection of about 2 mg/L was seen with the procedure described. The spike recovery was 97% with a 7% coefficient of variation at 9.9 mg/L added oxalic acid.

C. Glyoxals

The modification of Chian and Kuo's method²¹ showed linear standard curves for glyoxal, methylglyoxal, and dimethylglyoxal. The sensitivity increased with the number of methyl groups in the series. Although the method

gave reasonable precision for methylglyoxal and dimethylglyoxal at about 2 mg/L levels, the precision for glyoxal was not acceptable at this level. The spike recoveries for all glyoxals were in the range of 79 to 88%. A multiple extraction procedure should be tried to obtain more consistent recoveries of the quinoxalines formed in the reaction step.

D. Formaldehyde

Formaldehyde was readily determined by colorimetry²⁴ for concentrations up to 8 mg/L and by fluorimetry²⁵ for levels up to 2 mg/L. As little as 0.05 mg/L was detected by the fluorimetric method. Glyoxylic acid did not interfere with the formaldehyde determination, and recovery of an added spike of less than 2 mg/L was about 83%, with a precision of about 1% relative standard deviation using the colorimetric method. For most studies, the extra sensitivity of the fluorimetric method may not be necessary; colorimetry with 1- to 5-cm cells would be recommended in such cases.

E. Glyoxylic Acid

The resorcinol procedure of Zarembski and Hodgkinson³¹ was found to work well for 0- to 5-mg/L levels of glyoxylic acid. At least 0.1 mg/L could be detected. The precision was shown to be about 7% relative standard deviation, and a spike of 0.8 mg/L was recovered at 96% of theoretical. Because of the multiple extractions used, this method was more time-consuming than most colorimetric procedures and it did require the use of a spectrofluorimeter.

IX. LITERATURE CITED

1. McCarthy, J.J., W.F. Cowen, E.S.K. Chian, and B.W. Peterman, "Evaluation of an Air Stripping-Ozone Contactor System," Tech. Report No. 7707, U.S. Army Medical Bioengineering Research and Development Laboratory, Ft. Detrick, Frederick, MD (June 1977).
2. Chian, E.S.K. and P.K. Kuo, "Fundamental Study on the Post Treatment of RO Permeates from Army Wastewaters," Second Annual Summary Report to U.S. Army Research and Development Command, Contract No. DAMD 17-75-C-5006 (October 1976).
3. Gould, J.P. and J.W. Weber, Jr., "Oxidation of Phenols by Ozone," J. Water Pollut. Control Fed. 48: 47-60 (1976).
4. Supelco, Inc., Newsletter, "Analysis of VFA's from Anaerobic Fermentation," Chromatogr. Lipids IX: (May 1975).
5. Van Huyssteen, J.J., "The Determination of Short-Chain Fatty Acids in Aqueous Solution by Gas-Liquid Chromatography," Water Res. 4: 645-659 (1970).
6. Bethge, P.O. and K. Lindstrom, "Determination of Organic Acids of Low Relative Molecular Mass (C_1 to C_4) in Dilute Aqueous Solution," Analyst 99: 137-142 (1974).
7. Davis, B., "Crown Ether Catalyzed Derivatization of Carboxylic Acids and Phenols with Pentafluorobenzyl Bromide for Electron Capture Gas Chromatography," Anal. Chem. 49: 832-834 (1977).
8. Junk, G.A., J.J. Richard, M.D. Grieser, D. Witiak, J.L. Witiak, M.D. Arguello, R. Vick, H.J. Svec, J.S. Fritz, and G.V. Calder, "Use of Macroreticular Resins in the Analysis of Water for Trace Organic Contaminants," J. Chromatogr. 99: 745-762 (1974).
9. See, G.G., P.Y. Yang, and K.K. Kacholia, "Research and Development of an Ozone Contactor System," Final Report to U.S. Army Medical Research and Development Command, Life Systems, Inc., Cleveland, OH, Contract No. DAMD 17-76-C-6041 (1976).
10. Zarembski, P.M. and A. Hodgkinson, "The Fluorimetric Determination of Oxalic Acid in Blood and Other Biological Materials," Biochem. J. 96: 717-721 (1965).
11. Rumsey, T.S. and C.H. Noller, "A Study of the Quantitative Measurement of Certain Metabolic Acids by Gas-Liquid Chromatography," J. Chromatogr. 24: 325-334 (1966).
12. Mee, J.M.L. and R.W. Stanley, "A Rapid Gas-Liquid Chromatographic Method for Determining Oxalic Acid in Biological Materials," J. Chromatogr. 76: 242-243 (1973).

13. Schlenk, H. and J.L. Gellerman, "Esterification of Fatty Acids with Diazomethane on a Small Scale," Anal. Chem. 32: 1412-1414 (1960).
14. Butts, W.C. and W.T. Rainey, Jr., "Gas Chromatography and Mass Spectrometry of the Trimethylsilyl Derivatives of Inorganic Anions," Anal. Chem. 43: 538-542 (1971).
15. Von Nicolai, N. and F. Zilliken, "Gas Chromatographische Bestimmung von Oxalsäure, Malonsäure und Bernsteinsäure aus Biologischem Material," J. Chromatogr. 92: 431-434 (1974).
16. Ackman, R.G., M.A. Bannerman, and F.A. Vanderhevel, "Decomposition of Short-Chain Dicarboxylic Acid Esters during Separation on Polyester Gas Chromatography Media," Anal. Chem. 32: 1209 (1960).
17. Chian, E.S.K., Personal Communication (20 December 1977), Georgia Institute of Technology, Atlanta, Georgia 30332.
18. Johnson, D.P., F.E. Critchfield, and J.E. Ruch, "Spectrophotometric Determination of Trace Quantities of Alpha-dicarbonyl Compounds," Anal. Chem. 34: 1389-1391 (1962).
19. Gilbert, E., "Reactions of Ozone with Organic Compounds in Dilute Aqueous Solution: Identification of Their Oxidation Products," in Ozone/Chlorine Dioxide Oxidation Products of Organic Materials, Rice, R.C. and J.A. Cotruvo (eds.), Ozone Press International, Cleveland, OH, pp. 227-242 (1978).
20. Nojima, K., K. Fukaya, S. Fukui, and S. Kanno, "The Formation of Glyoxals by the Photochemical Reaction of Aromatic Hydrocarbons in the Presence of Nitrogen Monoxide," Chemosphere 5: 247-252 (1974).
21. Chian, E.S.K. and P.K. Kuo, "Fundamental Study on the Post Treatment of RO Permeates from Army Wastewaters," Third Quarter Report to U.S. Army Research and Development Command, Contract No. DAMD 17-75-C-5006 (1977).
22. Hammarstrand, K., "A Sensitive GC Method for Formaldehyde Determinations," presented at the 25th Pittsburgh Conference (1974), available from Varian Instrument Division, Palo Alto, CA, March 1974.
23. Kallio, H., R. Linko, and J. Kaitaranta, "Gas-Liquid Chromatographic Analysis of 2,4-Dinitrophenyl Hydrazones of Carbonyl Compounds," J. Chromatogr. 65: 355-360 (1972).
24. Nash, T., "The Colorimetric Estimation of Formaldehyde by Means of the Hantzsch Reaction," Biochem. J. 55: 416-421 (1953).
25. Belman, S., "The Fluorimetric Determination of Formaldehyde," Anal. Chim. Acta 29: 120-126 (1963).
26. Bartos, J. and M. Pesez, "Organic Analysis by Luminescence Methods," Talanta 19: 93-124 (1972).

27. Pesez, M. and J. Bartos, "Elements De Fluorimetrie Organique Fonctionnelle," Talanta 14: 1097-1108 (1967).
28. Chalmers, R.A. and R.W.E. Watts, "Derivatives for the Identification and Quantitative Determination of Some Keto and Aldo-carboxylic Acids by Gas-Liquid Chromatography," Analyst 97: 951-957 (1972).
29. Kallio, H. and R.R. Linko, "Gas-Liquid Chromatographic Analysis of 2,4-Dinitrophenylhydrazones of Keto Acid Methyl Esters," J. Chromatogr. 76: 229-232 (1973).
30. Kramer, D.N., N. Klein, and R.A. Basalice, "Quantitative Determination of Glyoxylic Acid," Anal. Chem. 31: 250-252 (1959).
31. Zarembski, P.M. and A. Hodgkinson, "The Fluorimetric Microdetermination of Glyoxylic Acid in Blood, Urine and Bacterial Extracts," Biochem. J. 96: 218-223 (1965).

APPENDIX A

1. Synthesis of Dipentafluorobenzyl Oxalate

Oxalic acid dihydrate (0.5 g) was titrated to pH 8 with 0.386 M tetra-n-butylammonium hydroxide (TBAOH) in a small volume of water. The water was removed by rotary evaporation at 43°C. The residue was a yellow oil.

To the residue was added 270 μ L of pure pentafluorobenzyl bromide, and the mixture was allowed to stand about 15 minutes until white crystals formed. About 10 mL of benzene was added to dissolve the crystals, and the TBA bromide salts were extracted with two portions of water. The benzene was removed from the organic layer by rotary evaporation, and the crystals were recrystallized from ethyl ether and air dried (melting point = 133°-134°C).

Figure A-1 shows the single nuclear magnetic resonance (NMR) proton peak expected for the four equivalent $-\text{CH}_2-$ protons. An infrared spectrum (Figure A-2) showed a carbonyl band at $1,770\text{ cm}^{-1}$, but no OH absorption in the $2,500\text{--}3,000\text{ cm}^{-1}$ range, indicating that a free-COOH group was not present, as expected for an ester.

2. Mass Spectrometry of PFB Acetate, Formate, and Oxalate

The mass spectra were recorded at 70 electron volts ionizing potential on a DuPont 21-490B gas chromatograph/mass spectrometer and 21-094 data collection system.

a. Dipentafluorobenzyl Oxalate

The crystals of dipentafluorobenzyl oxalate were introduced into the mass spectrometry source through the probe.

The major fragments were:

$(\text{C}_6\text{F}_5\text{CH}_2)^+$, m/e 181, as a base peak;

$(\text{C}_6\text{F}_5\text{CH}_2\text{CH}_2)^+$, m/e 195;

$(\text{C}_6\text{F}_5\text{CH}_2(\text{OCO})_2)^+$, m/e 269;

and $(\text{C}_6\text{F}_5\text{CH}_2\text{OCO})_2^+$, m/e 450, as a parent ion.




Figure A-1. NMR spectrum of dipentafluorobenzyl oxalate.




Figure A-2. Infrared spectrum of dipentafluorobenzyl oxalate.

b. Pentafluorobenzyl Formate

The GC conditions were:

GC column: 10 ft x 2 mm I.D. glass, packed with 10%
OV-1 on 100/120 mesh Gas Chrom Q
(Applied Science Laboratories, Inc.)

Carrier: Ultra high purity grade helium, 20 psi

Column temperature: 100°C to 250°C (2 minutes),
12°C/minute

Injection port temperature: 150°C

Detector temperature: 225°C

The major fragments were:

$(C_6F_5CH_2)^+$, m/e 181, as a base peak, and
 $(C_6F_5CH_2-OCHO)^+$, m/e 226, as a parent ion.

c. Pentafluorobenzyl Acetate

The GC conditions were:

GC column: 10 ft x 2 mm I.D. glass, 10% OV-1 on
100/120 mesh Gas Chrom Q (Applied Science
Laboratories, Inc.)

Carrier: Ultra high purity grade helium, 20 psi

Column temperature: 100°C to 250°C (2 minutes),
8°C/minute

Injection port temperature: 150°C

Detector temperature: 220°C

The major fragments were:

$(C_6F_5CH_2)^+$, m/e 181;
 $(C_6F_5CH_2CH_2)^+$, m/e 195;
and $(C_6F_5CH_2-O-CO-CH_3)^+$, m/e 240, as a parent ion.

DISTRIBUTION LIST

No. of
Copies

5	US Army Medical Research and Development Command ATTN: SGRD-SI Fort Detrick, Frederick, MD 21701
12	Defense Technical Information Center (DTIC) ATTN: DTIC-DDA Cameron Station Alexandria, VA 22314
1	Commandant Academy of Health Sciences, US Army ATTN: AHS-COM Fort Sam Houston, TX 78234
2	USAMBRDL Technical Library

DATA
FILM
2-