

AD-A116 804

ARIZONA UNIV TUCSON DEPT OF CHEMISTRY
PAIRED ION CHROMATOGRAPHY OF CERTAIN MONOSACCHARIDES. (U)
SEP 81 T GNANASAMBANDAN, H FREISER
TR-5

F/6 7/3

N00014-81-K-0576

UNCLASSIFIED

NL

[RE]
AL-1
01104



END
DATE
FILMED
6-82
DTIC

AD A116804

DTIC FILE COPY

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		1. READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER 5	2. GOVT ACCESSION NO. AD-A116804	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) "Paired Ion Chromatography of Certain Mono-saccharides"		5. TYPE OF REPORT & PERIOD COVERED Technical -- Interim
7. AUTHOR(s) T. Gnanasambandan and H. Freiser		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS University of Arizona Department of Chemistry Tucson, Arizona 85721		8. CONTRACT OR GRANT NUMBER(s) N00014-81-K-0576
11. CONTROLLING OFFICE NAME AND ADDRESS Materials Sciences Division Office of Naval Research Arlington, Virginia 22217		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE
		13. NUMBER OF PAGES 5
		15. SECURITY CLASS. (of this report) unclassified
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) This document has been approved for public release and sale; its distribution is unlimited		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES to be published in Analytical Chemistry		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) monosaccharide chromatographic analyses, fructose, ribose, and xylose		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) A paired ion dye-impregnated reversed phase chromatographic column is shown to give well defined chromatographic zones for monosaccharides as their borate complexes, permitting determination of sub-microgram levels of fructose, xylose, and ribose.		

DD FORM 1473

EDITION OF 1 NOV 65 IS OBSOLETE
S/N 0102-LF-014-6601

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

OFFICE OF NAVAL RESEARCH
Contract N00014-81-K-0576

TECHNICAL REPORT NO. 5

"Paired Ion Chromatography
of Certain Monosaccharides"

by

T. Gnanasambandan and H. Freiser

submitted to

Analytical Chemistry

University of Arizona
Department of Chemistry
Tucson, Arizona 85721

Reproduction in whole or in part is permitted for
any purpose of the United States Government

This document has been approved for public release
and sale; its distribution is unlimited

Paired ion Chromatography of Certain Monosaccharides

Sir:

Several chromatographic methods (1,2,3) have been reported for separation of sugars but this still remains a formidable task for determination of trace levels. The packings which have been most successful in separating sugars are those with amine groups bound to inert silica or synthetic supports such as U-Bondapak carbohydrate (Waters Assoc.), Lichrosorb-NH₂ (Merck) and Aminex resins (Bio-Rad Lab) (4). These methods involve elution with water/acetonitrile system or water at 80-85°C.

In the course of our study in reverse phase ion-pair HPLC (5) we found that sugars such as fructose, can be ionized on column and paired with a chromophore to effect improved separation and detection.

One of the interesting aspects of this separation is that the non-ionic species react on the column to form suitable species and separated. Generally the analyte is subject to chemical pre-treatment before separation. This is the first report of an on-column reaction in reverse phase HPLC.



Accession For	
NTIS	CRS
DTIC	TR
Unannounced	
Justification	
By	
Distribution	
Availability Codes	
Dist	Special

A

Experimental

A modular chromatographic unit consisting of an Altex Pump (Model 100A), a variable-wave length uv-vis detector (Schoeffel Instruments No. 770) and a Fisher Recordall Series 500 strip chart recorder was used. Absorbance measurements were made at 651nm. A 25cm Zorbax ODS column, described by the manufacturer as having 15% carbon loading, was employed.

The mobile phase used was methanol-(0.1M aqueous borate buffer PH 7.5) in 2:98 V/V% containing 1×10^{-4} M methylene blue (as chloride).

Before the chromatographic separations were undertaken, the column was conditioned by passing sufficient mobile phase to "saturate" it with methylene blue (6). The flow rate used for both column saturation and subsequent chromatography was 1.0ml per min. (1% W/V) solutions of the fructose, xylose and ribose (separately) in the mobile phase were added through 10 μ l loop sampling valve (rotary valve injector SP-419-0410). As may be seen from the data, reasonable separation of the three sugars tested was obtained at sub-microgram detection limits ($\sim 10^{-7}$ g).

Table 1. Capacity Factors of Monosaccharides

Sugars	k'
Fructose	3.5
Ribose	2.9
Xylose	1.2
Mobile phase:	0.1M Borate Buffer pH7.5 Methylene blue $10^{-4}M$ 2% MeOH/H ₂ O V/V
Flow rate:	1.0ml per min.
Sample:	10 μ l 1.0% wt/V Solutions

The use of borate complex formation to produce anionic sugar species results in a dramatic improvement in selectivity. The formation constant for the glucose-borate complex is far smaller ($\sim 10^2$) than the corresponding fructose complex ($\sim 10^4$) (7). As a result, glucose does not give any detectable response under the experimental conditions employed here, making the determination of fructose entirely free from glucose interference. Changing these conditions to accomodate the lower stability of the glucose complex, e.g. by increasing borate ion concentration, would very likely permit determination of glucose also.

Although further work is required to clearly establish the optimum conditons necessary to separate all sugars, it would appear that reactions occur on the column to give suitable products for effective separation and detection. Further work along these lines is under way in this laboratory.

Literature cited

1. Ladish, M.R., Tsao, G.T., J. Chromatogr., 1978, 166, 18-100.
2. Palla, G., Anal. Chem., 1981, 53, 1966-1967.
3. Linden, J.L., Lawhead, C.L., J. Chromatogr. 1975, 105, 125-233.
4. Conrad, E.C., Palmer, J.K., Food Technol (Chicago), 1976, 84-92.

5. Gnanasambandan T., Freiser, H., Anal. Chem., 1981, 53, 909-911.
6. Gnanasambandan, T., Freiser, H., Anal. Chem., in press.
7. Roy, G.L., Laferriere and Edwards, J.O., J. Inorganic Nuclear Chem., 1957, 4, 106.

Acknowledgment. This work was supported by the Office of Naval Research.

T. Gnanasambandan

H. Freiser

Department of Chemistry
University of Arizona
Tucson, Arizona
April 19, 1982

MAILING LIST FOR OFFICE OF NAVAL RESEARCH REPORTS

Mr. John Boyle
Materials Branch
Nabal Ship Engineering Center
Philadelphia, PA 19112

Dr. H. Chernciff
Department of Mathematics
Massachusetts Institute of Technology
Cambridge, MA 02139

Walter G. Cox, Cod 3632
Naval Underwater Systems Center
Building 148
New Port, RI 02840

Defense Technical Information Center (12)
Building 5, Cameron Station
Alexandria, VA 22314

Dr. M.B. Denton
Department of Chemistry
University of Arizona
Tucson, AZ 85721

The Assistant Secretary of the Navy
(RExS)
Department of The Navy
Room 4E736, Pentagon
Washington, DC 20350

Department of Physics & Chemistry
Naval Postgraduate School
Monterey, CA 93940

Dr. John Duffin
United States Naval Postgraduate School
Monterey, CA 93940

Dr. Kent Eisentraut, MBT
Air Force Materials Laboratory
Wright-Patterson AFB, OH 45433

Dr. Christie G. Enke
Michigan State University
Department of Chemistry
East Lansing, MI 48824

Dr. G.M. Hieftje
Department of Chemistry
Indiana University
Bloomington, IN 47401

Mr. James Kelley
DTNSRDC Code 2803
Annapolis, MD 21402

Dr. Bruce Kowalski
Chemistry Department
University of Washington-Seattle
Seattle, WA 98105

Dr. Rudolph J. Marcus
Office of Nabal Research
Scientific Liaison Group
American Embassy
APO San Francisco, CA 96503

Dr. George H. Morrison
Department of Chemistry
Cornell University
Ithaca, NY 14853

Naval Civil Engineering Laboratory
Attn: Dr. R.W. Drisko
Port Hueneme, CA 93401

Naval Weapons Center
Attn: Dr. A.B. Amster
Chemistry Division
China Lake, CA 93555

Naval Ocean Systems Center
Attn: Dr. S. Yamamoto
Marine Sciences Division
San Diego, CA 91232

Naval Ocean Systems Center
Attn: Mr. Joe McCartney
San Diego, CA 92152

Commander, Naval Air Systems Command
Attn: Code 310C (H. Rosenwasser)
Department of the Navy
Washington, DC 20390

Director, Naval Research Laboratory
Attn: Code 6100
Washington, DC 20390

Naval Ship Research & Development Center
Attn: Dr. G. Bosmajian
Applied Chemistry Division
Annapolis, MD 21401

Office of Naval Research (2)
Attn: Dr. Richard S. Miller
800 N. Quincy Street
Arlington, VA 22217

Office of Naval Research
Attn: Code 472
800 N. Quincy Street
Arlington, VA 22217

ONR Branch Office
Attn: Dr. George Sandoz
536 S. Clark Street
Chicago, IL 60605

ONR Eastern/Central Regional Office
Attn: Dr. L. H. Peebles
Building 114, Section D
666 Summer Street
Boston, MA 02210

ONR Western Regional Office
1030 East Green Street
Pasadena, CA 91106

Dr. R. A. Osteryoung
Department of Chemistry
SUNY, Buffalo
Buffalo, NY 14214

Dr. S. P. Perone
Department of Chemistry
Purdue University
West Lafayette, IN 47907

Dr. Victor L. Rehn
Naval Weapons Center
Code 3813
China Lake, CA 93555

Dr. Freed Saalfeld
Chemistry Division, Code 6100
Naval Research Laboratory
Washington, DC 20375

U.S. Army Research Office
Attn: CRD-AA-IP
P. O. Box 1211
Research Triangle Park, NC 27709

Dr. A. L. Slafkosky
Scientific Advisor
Commandant of the Marine Corps
(Code RD-1)
Washington, DC 20380

Dr. D. L. Venezky
Naval Research Laboratory
Code 6130
Washington, DC 20375

Dr. Isiah H. Warner
Texas A&M University
Department of Chemistry
College Station, TX 77840

Dr. K. Wilson
Department of Chemistry
University of California, San Diego
La Jolla, CA 92037

Dr. A. Zirino
Naval Underseas Center
San Diego, CA 92132

James R. Delaney
Contract Specialist
Office of Naval Research
422 University District Bldg.
1107 NE 45th Street JD-27
Seattle, WA

9/23/81