

AD-A130 159

MUTAGENIC POTENTIAL OF THE HOLSTON COMPOUNDS: VIRGIN  
DMSO DMSO RECYCLE SO..(U) LETTERMAN ARMY INST OF  
RESEARCH PRESIDIO OF SAN FRANCISCO CA

1//


UNCLASSIFIED

L J SAUERS ET AL. JUN 83 LAIR-149

F/G 6/20

NL

						END							
						DATE							
						FINED							
						7 83							
						DTIC							



MICROCOPY RESOLUTION TEST CHART  
NATIONAL BUREAU OF STANDARDS-1963-A

(12)

INSTITUTE REPORT NO. 149

**MUTAGENIC POTENTIAL OF THE HOLSTON COMPOUNDS:**

Virgin DMSO  
DMSO Recycle Solvent  
DMSO Evaporator Sludge

LEONARD J. SAUERS, MS, SP5  
THOMAS P. KELLNER, BA, SP4  
and  
JOHN T. FRUIN, DVM, PhD, COL VC

TOXICOLOGY GROUP,  
DIVISION OF RESEARCH SUPPORT

DTIC  
JUL 6 1983  
A

JUNE 1983

Toxicology Series 57

LETTERMAN ARMY INSTITUTE OF RESEARCH  
PRESIDIO OF SAN FRANCISCO, CALIFORNIA 94129

83 07 06 127

DTIC FILE COPY

AD A130159


**Mutagenic Potential of the Holston Compounds: Virgin DMSO, DMSO Recycle Solvent, DMSO Evaporator Sludge (Toxicology Series 57)--Sauers et al**

Reproduction of this document in whole or in part is prohibited except with the permission of the Commander, Letterman Army Institute of Research, Presidio of San Francisco, California 94129. However, the Defense Technical Information Center is authorized to reproduce the document for United States Government purposes.

Destroy this report when it is no longer needed. Do not return it to the originator.

Citation of trade names in this report does not constitute an official endorsement or approval of the use of such items.

This material has been reviewed by Letterman Army Institute of Research and there is no objection to its presentation and/or publication. The opinions or assertions contained herein are the private views of the author(s) and are not to be construed as official or as reflecting the views of the Department of the Army or the Department of Defense. (AR 360-5)

  
(Signature and date)

This document has been approved for public release and sale; its distribution is unlimited.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER Institute Report No.149	2. GOVT ACCESSION NO. A730159	RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) Mutagenic Potential of the Holston Compounds: Virgin DMSO, DMSO Recycle Solvent, and DMSO Evaporator Sludge	5. TYPE OF REPORT & PERIOD COVERED Final 3 Jan - 20 Mar 83	
7. AUTHOR(s) Leonard J. Sauers, MS, SP5 Thomas P. Kellner, BA, SP4 John T. Fruin, DVM, PhD, COL VC	6. PERFORMING ORG. REPORT NUMBER	
9. PERFORMING ORGANIZATION NAME AND ADDRESS US Army Medical Research and Development Command Letterman Army Institute of Research Presidio of San Francisco, CA 94129	8. CONTRACT OR GRANT NUMBER(s)	
11. CONTROLLING OFFICE NAME AND ADDRESS US Army Medical Research and Development Command Fort Detrick Frederick, MD 21701	10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS APC: TL01	
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)	12. REPORT DATE June 1983	
	13. NUMBER OF PAGES 47	
	15. SECURITY CLASS. (of this report) UNCLASSIFIED	
16. DISTRIBUTION STATEMENT (of this Report)  THIS DOCUMENT HAS BEEN APPROVED FOR PUBLIC RELEASE AND SALE: ITS DISTRIBUTION IS UNLIMITED.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Mutagenicity, Toxicology, Ames Assay, Holston Compounds, Virgin DMSO, DMSO Recycle Solvent, DMSO Evaporator Sludge, Frameshift Mutagen, Dimethyl Sulfoxide		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) The mutagenic potential of the Holston Compounds (Virgin DMSO*, DMSO Recycle Solvent, and DMSO Evaporator Sludge) was assessed by using the Ames Salmonella/ Mammalian Microsome Mutagenicity Assay. Tester strains TA 98, TA 100, TA 1535, TA 1537, and TA 1538 were exposed to doses ranging from 0.1 ml of a 100% to 0.1 ml of a 0.032% solution. Negative mutagenic responses were observed for the Virgin DMSO and the DMSO Recycle Solvent. Mutagenic potential was observed for the DMSO Evaporator Sludge. *DMSO = Dimethyl Sulfoxide		

1/

ABSTRACT

The mutagenic potential of the Holston Compounds (Virgin DMSO, DMSO Recycle Solvent, and DMSO Evaporator Sludge) was assessed by using the Ames Salmonella/Mammalian Microsome Mutagenicity Assay. Tester strains TA 98, TA 100, TA 1535, TA 1537, and TA 1538 were exposed to doses ranging from 0.1 ml of a 100% to 0.1 ml of a 0.032% solution. Negative mutagenic responses were observed for the Virgin DMSO and the DMSO Recycle Solvent. Mutagenic potential was observed for the DMSO Evaporator Sludge.

\* DMSO = Dimethyl Sulfoxide

KEY WORDS: Mutagenicity, Toxicology, Ames Assay, Holston Compounds, Virgin DMSO, DMSO Recycle Solvent, DMSO Evaporator Sludge, Frameshift Mutagen, Dimethyl Sulfoxide



## PREFACE

TYPE REPORT: Ames Assay GLP Study Report

TESTING FACILITY: US Army Medical Research and Development Command  
Letterman Army Institute of Research  
Presidio of San Francisco, CA 94129

SPONSOR: US Army Medical Research and Development Command  
US Army Medical Bioengineering Research  
and Development Laboratory  
Fort Detrick, Frederick, MD 21701

PROJECT: DMSO Recrystallization Solution  
TLO1

GLP STUDY NUMBER: 83001

STUDY DIRECTOR: COL John T. Fruin, DVM, PhD, VC  
Diplomate, American College of  
Veterinary Preventive Medicine

PRINCIPAL INVESTIGATOR: SP5 Leonard J. Sauer, MS

REPORT AND DATA MANAGEMENT: A copy of the final report, study protocols, raw data, retired SOPs, chemical data, and an aliquot of each test compound will be retained in the LAIR Archives.

TEST SUBSTANCE: The Holston Compounds (Virgin DMSO, DMSO Recycle Solvent, and DMSO Evaporator Sludge).

INCLUSIVE STUDY DATES: 3 January - 20 March 1983

OBJECTIVE: To determine the mutagenic potential of the Holston Compounds using the Ames Assay. Tester strains TA 98, TA 100, TA 1535, TA 1537, and TA 1538 were used. The plate incorporation method was followed. The test substances was diluted in reagent grade dimethyl sulfoxide (DMSO) and this diluent was checked for sterility.

#### ACKNOWLEDGMENTS

The authors wish to thank SP4 Larry Mullen, BS and John Dacey for their assistance in performing the research.



Signatures of Principal Scientists involved  
in the Study

We, the undersigned, believe the study number 83001 described  
in this report to be scientifically sound and the results and  
interpretation to be valid. The study was conducted to comply, to  
the best of our ability, with the Good Laboratory Practice  
Regulations outlined by the Food and Drug Administration.

*Leonard J. Sauers 25 Apr 83*  
LEONARD J. SAUERS, MS / DATE  
SP5, USA  
Principal Investigator

*John T. Fruin 25 Apr 83*  
JOHN T. FRUIN, DVM, PhD / DATE  
COL, VC  
Study Director

*Thomas P. Kellner / 25 Apr 83*  
THOMAS P. KELLNER, BA / DATE  
SP4, USA  
Toxicology Group



DEPARTMENT OF THE ARMY  
LETTERMAN ARMY INSTITUTE OF RESEARCH  
PRESIDIO OF SAN FRANCISCO, CALIFORNIA 94129

REPLY TO  
ATTENTION OF:

SGRD-ULZ-QA

6 May 1983

MEMORANDUM FOR RECORD

SUBJECT: Report of GLP Compliance

I hereby certify that in relation to LAIR GLP study 83001 the following inspections were made:

12 Jan 83  
26 Jan 83  
27 Jan 83  
16 Mar 83  
19 Mar 83

The report and raw data for this study were audited on 5 May 1983.

Routine inspections with no adverse findings are reported quarterly, thus these inspections are also included in the April 83 report to management and the Study Director.

NELSON R. POWERS, Ph.D.  
CPT, MSC  
Quality Assurance Officer

## TABLE OF CONTENTS

Abstract.....	i
Preface.....	iii
Acknowledgments.....	iv
Signatures of Principal Scientists.....	v
Report of Quality Assurance Unit.....	vi
Table of Contents.....	vii
BODY OF REPORT	
INTRODUCTION	
Rationale for using the Ames Assay.....	1
Description of Test, Rationale for strain selection.....	1
Description of Strains, History, Methods, and Data.....	2
Objective of Study.....	3
METHODS	
Rationale for Dosage Levels and Response Tabulations.....	3
Test Format.....	3
Statistical Analysis.....	4
Chemical Analysis.....	4
RESULTS and DISCUSSION.....	5
CONCLUSIONS.....	6
RECOMMENDATION.....	6
REFERENCES.....	7
APPENDICES	
Appendix A (Chemical Analysis for Holston Compounds).....	10
Appendix B (Tables 1 through 12).....	16
Appendix C (Figures 1 through 6).....	34
DISTRIBUTION LIST.....	40

MUTAGENIC POTENTIAL OF THE HOLSONT COMPOUNDS: Virgin DMSO, DMSO Recycle Solvent,  
DMSO Evaporator Sludge-Sauers et al

Rationale for using the Ames Assay

The Ames Salmonella/Mammalian Microsome Mutagenicity Test is one of a standard bank of tests used by our laboratory for the assessment of the mutagenic potential of a test substance. It is a short-term screening assay, which we use for the prediction of potential mutagenic agents in mammals. It is inexpensive when compared to in vivo tests, yet is highly predictive and reliable in its ability to detect mutagenic activity and therefore carcinogenic probability (1). It relies on basic genetic principles and allows for the incorporation of a mammalian microsomal enzyme system to increase sensitivity through enzymatically altering the test substance into an active metabolite. It has proven highly effective in assessing human risk (1).

Description of Test (Rationale for the selection of strains)

The test was developed by Bruce Ames, Ph.D. from the University of California-Berkeley. The test involves the use of several different genetically altered strains of Salmonella typhimurium, each with a specific mutation in the histidine operon (2). The test substance demonstrates mutagenic potential if it is able to revert the mutation in the bacterial histidine operon to the wild type and reestablish prototrophic growth within the test strain. This reversion also can occur spontaneously due to a random mutational event. If, after adding a test substance, the number of revertants is significantly greater than the spontaneous reversion rate, then the test substance physically altered the locus involved in the operon's mutation and is able to induce point mutations (2).

In order to increase the sensitivity of the test system, other mutations in the Salmonella are used (2). To insure a higher probability of uptake of test substance, the genome for the lipopolysaccharide layer (LP) is mutated and, therefore, larger molecules are allowed to enter the bacteria. Each strain has another induced mutation which causes loss of excision repair mechanisms. A

mammalian microsomal enzyme system is incorporated since many chemicals are not by themselves mutagenic but have to be activated by an enzymatic process. These microsomal enzymes are obtained from livers of rats induced with Aroclor 1254; the enzymes allow for the expression of the metabolites which would occur in the mammalian system. This activated rat liver microsomal enzyme homogenate is termed S-9.

Description of Strains (History of the strains used method to monitor the integrity of the organisms, and data pertaining to current and historical control and spontaneous reversion rates)

The test consists of using five different strains of *Salmonella typhimurium* that are unable to grow in absence of histidine because of a specific mutation in the histidine operon. This histidine requirement is verified by attempting to grow the tester strains on minimal glucose agar (MGA) plates, both with and without histidine. The dependence on this amino acid is shown when growth occurs only in its presence. The plasmids in strains TA 98 and TA 100 contain an ampicillin resistant R factor. Strains deficient in this plasmid demonstrate a zone of inhibition around an ampicillin impregnated disc. The alteration of the LP layer allows uptake by the *Salmonella* of larger molecules. If a crystal violet impregnated disc is placed onto a plate containing any one of the bacterial strains, a zone of growth inhibition will occur because the LP layer is altered. The absence of excision repair mechanisms can be determined by using ultraviolet (UV) light. These mechanisms function primarily by repairing photodimers between pyrimidine bases. Exposure of bacteria to UV light will activate the formation of these dimers and cause cell lethality, since excision of these photodimers can not be made. The genetic mutation resulting in UV sensitivity also induces a dependence by the *Salmonella* to biotin. Therefore, this vitamin must be added. In order to prove that the bacteria are responsive to the mutation process, positive controls are run with known mutagens. If after exposure to the positive control substance, a revertant count is obtained which is greater than twice the spontaneous reversion rate, then the bacteria are adequately responsive. Sterility controls are performed to determine the presence of contamination. Sterility of the test compound is also confirmed in each first dilution. Verification of the tester strains occurs simultaneously with the running of each assay. The value of the spontaneous reversion rate is obtained by using the same inoculum of bacteria that is used in the assay (3).

Strains were obtained directly from Dr. Ames, University of California-Berkeley, propagated and then maintained at  $-80^{\circ}\text{C}$  in our laboratory. Before any substance was tested, quality controls were performed on the bacterial strains to establish the presence of their special features and also to determine the spontaneous reversion rate (2). Records are maintained of all the data to determine if deviations from the set trends have occurred. These records are kept in the archives of the Quality Assurance Unit.

In this series of tests for the detection of mutagenic potential of different agents, we compared the spontaneous reversion values with our own historical values and those cited by Ames et al (2). Our conclusions are based on the spontaneous reversion rate compared to the experimentally induced rate of mutation. When operating effectively, these strains detect substances that cause base pair mutations (TA 1535, TA 100) and frameshift mutations (TA 1537, TA 1538, and TA 98).

### Objective of Study

The objective of the study is to determine the mutagenic potential of the Holston Compounds using the Ames Assay. Tester strains TA 98, TA 100, TA 1535, TA 1537, and TA 1538 were used. The plate incorporation method was followed. The test substances were diluted in reagent grade dimethyl sulfoxide (DMSO) and this diluent was checked for sterility.

### METHODS (3)

#### Rationale for Dosage Levels and Dose Response Tabulations

To insure readable and reliable results, a sublethal concentration of the test substance had to be determined. This toxicity level was found by using MGA plates, various concentrations of the substance, and approximately  $10^8$  cells of TA 100 per plate, unless otherwise specified. Top agar containing trace amounts of histidine and biotin were placed on MGA plates. TA 100 was used because it is the most sensitive strain. Strain verification was confirmed on the bacteria, along with a determination of the spontaneous reversion rate. After incubation, the growth was observed on the plates. (The auxotrophic Salmonella will replicate a few times and potentially express a mutation. When the histidine and biotin supplies are exhausted, only those bacteria that reverted to the prototrophic phenotype will continue to reproduce and form macrocolonies; the remainder of the bacteria comprises the background lawn. The minimum toxic level is defined as the lowest serial dilution at which decreased macrocolony formation, below that of the spontaneous revertant rate, and an observable reduction in the density of the background lawn occurs.) A maximum dose of 1 mg/plate is used when no toxicity is observed. The densities were recorded as normal, slight, and no growth.

### Test Format

After we validated our bacterial strains and determined the optimal dosage of the test substance, we began the Ames Assay. In the actual experiment, 0.1 ml of the particular strain of Salmonella ( $10^8$  cells) and the specific dilutions of the test substance were added to 2 ml of molten top agar, which contained trace amounts of histidine and biotin. Since survival is better from cultures which have just passed the log phase, the Salmonella strains are used 16 hours

(maximum) after initial inoculation into nutrient broth. The dose of the test substance spanned a 1000-fold, decreasing from the minimum toxic level by a dilution factor of 5. All the substances were tested with and without S-9 microsome fraction. The optimal titer of the S-9 was determined and 0.5 ml was added to the molten top agar. After all the ingredients were added, the top agar was mixed, then overlaid on minimum glucose agar plates. These plates contained 2% glucose and Vogel Bonner "E" Concentrate (4). The water used in this medium and all reagents came from a polymetric system. Plates were incubated, upside down in the dark at 37° C for 48 hours. Plates were prepared in triplicate and the average revertant counts were recorded. The corresponding number of revertants obtained was compared to the number of spontaneous revertants; the conclusions were recorded statistically. A correlated dose response is considered necessary to declare a substance as a mutagen. Commoner (5), in his report, "Reliability of Bacterial Mutagenesis Techniques to Distinguish Carcinogenic and Non-Carcinogenic Chemical," and McCann et al (1) in their paper, "Detection of Carcinogens as Mutagen in the Salmonella/Mammalian Microsome Mutagenicity Test: Assay of over 300 Chemicals," have concurred on the test's ability to detect mutagenic potential.

#### Statistical Analysis

Quantitative evaluation was ascertained by the method of Ames et al (2). They assumed that a compound which causes twice the spontaneous reversion rate and a correlated dose response is mutagenic.

#### Chemical Analysis

Information on the chemical analysis on the Holston Compounds appears in Appendix A. The stability of the Holston Compounds under these test conditions has not been determined but assumed to be stable at room temperature.

## RESULTS AND DISCUSSION

On 19 January 1983, the toxicity level determination was performed for the Virgin DMSO, DMSO Evaporator Sludge, and DMSO Recycle Solvent. For this experiment, all sterility, strain verification, and negative controls were normal (Appendix B, Table 1). Other results appear in Tables 2-12, Appendix B and Figures 1-6, Appendix C. No toxicity was observed after exposure of the tester strains to the compounds at the highest dose used (0.1 ml of a 100% solution) (Tables 2-4). It was observed that the DMSO Recycle Solvent precipitated when added to the top agar. The DMSO Recycle Solvent has a low solubility in water, the major component of the top agar solution. It is speculated that only about 0.5% of the solvent went into solution (personal communication, Thomas Kawakami, PhD, Laboratory for Energy-Related Health Research, 28 January 1983).

On 26 January 1983, the Ames Assay was performed on the test substances. In this assay, normal results were observed for all sterility and strain verification controls (Table 5). Normal results were also observed for all positive and negative controls (Table 6). Following exposure of the bacteria to the Virgin DMSO and DMSO Recycle Solvent, no incidences of mutagenicity were observed (Tables 7,8). Following exposure of the bacteria to the DMSO Evaporator Sludge, a doubling of the spontaneous reversion rate was induced at the 100% solution level in TA 98, TA 1537, and TA 1538. Increased reversion counts were seen for TA 100 and TA 1535 (Table 9). No evidence of mutagenicity was observed at the 20% solution level. For a positive mutagenic response in the Ames Assay, a test compound must induce a doubling of the spontaneous reversion rate and a correlated dose response. It was speculated that components in the DMSO Evaporator Sludge were mutagenic, but were in such small concentrations that a dose response could not be seen because of the wide difference between dilutions. A second Ames Assay was performed on 17 March 1983, using 0.1 ml per plate volumes of 100%, 80%, 60%, 40%, 20%, and 1% solutions of the DMSO Evaporator Sludge. For this experiment, all strain verifications and sterility controls were normal (Table 10). Normal results were observed for all positive and negative controls (Table 11). Following exposure of the bacteria to the DMSO Evaporator Sludge, mutagenic responses were seen for TA 98 at the 100%, 80%, and 60% solutions without S-9, and at the 100% and 80% solutions with S-9. A positive response was seen for TA 1537 at the 100%, 80%, and 60% solution with and without S-9, and at the 40% dose level with S-9. Mutagenicity was evident for TA 1538 at the 100% through 40% solution with and without S-9, and at the 20% solution without S-9 (Table 12). For each of the strains exhibiting a mutagenic response, a graph has been constructed to illustrate the correlated dose response (Figures 1-6). A deviation from a definitive dose response is observed for TA 1537. This can be attributed to the insensitivity of this particular strain and the closely spaced dilutions.



#### CONCLUSION

On the basis of the Ames Assay, the Virgin DMSO and DMSO Recycle Solvent are not mutagenic at the levels tested. The DMSO Evaporator Sludge shows characteristics of a frameshift mutagen at the levels tested but does not require the presence of metabolic activation for mutagenic induction.

Since the DMSO Recycle Solvent precipitated in the top agar, it is plausible that components within the solution are mutagenic but do not exhibit a response, since they were in such small concentrations when exposed to the Salmonella.

#### RECOMMENDATION

Components of the DMSO Evaporator Sludge should be identified and tested to determine the cause of the mutagenic response.

REFERENCES

1. McCann JE, Choi E, Yamasaki E, Ames BN. Detection of carcinogens as mutagens in the Salmonella/microsome test: Assay of 300 chemicals. Proc Nat Acad Sci, USA 1975;72:5135-5139.
2. Ames BN, McCann J, Yamasaki E. Methods for detection of carcinogens and mutagens with Salmonella/mammalian microsome mutagenicity test. Mutation Res 1975;31:347-364.
3. LAIR SOP OP-STX-1, Ames Salmonella/mammalian microsome mutagenicity test.
4. Vogel HJ, Bonner DM. Acetylornithinase of E. coli: Partial purification and some properties. J Biol Chem 1956;218:97-106.
5. Commoner B. Reliability of the bacterial mutagenesis techniques to distinguish carcinogenic and non-carcinogenic chemicals. EPA 600/1 76-002, 1976.

	page
Appendix A, Chemical Analysis of Holston Compounds.....	10
Appendix B, Tables 1-12.....	16
Appendix C, Figures 1 6.....	33

APPENDICES

Toxicity Test Sample Composition<sup>a</sup>

Concentration by HPLC, g/1

Sample	<sup>b</sup> RDX	<sup>c</sup> HMX	<sup>d</sup> TAX	<sup>e</sup> SEX	<sup>h</sup> %H <sub>2</sub> O	%DMSO
Virgin DMSO <sup>f</sup>	0	0	0	0	0.63	99.37
DMSO Recycle Solvent <sup>g</sup>	24.188	39.542	0.263	0	35.48 <sup>i</sup>	58.64 <sup>j</sup>
DMSO Evaporator Sludge <sup>f</sup>	0.548	0.942	3.521	0	5.35 <sup>i</sup>	94.19 <sup>j</sup>

Calculated Data In Weight Percent<sup>a</sup>

Sample	RDX	HMX	TAX	SEX	H <sub>2</sub> O	DMSO
Virgin DMSO	0	0	0	0	0.63	99.37
DMSO Recycle Solvent	2.22	3.64	0.02	0	35.48	58.64
DMSO Evaporator Sludge	0.05	0.09	0.32	0	5.35	94.19

<sup>a</sup> Data supplied by sponsor

<sup>b</sup> RDX: Cyclotrimethylenetrinitramine

<sup>c</sup> HMX: Cyclotetramethylenetetranitramine

<sup>d</sup> TAX: 1-Acetylhexahydro-3,5-Dinitro-1,3,5-Triazine

<sup>e</sup> SEX: Octahydro-1-(N)-Acetyl-3,5,7-Trinitro-1,3,5,7-Tetrazine

<sup>f</sup> At ambient temperature.


<sup>g</sup> Analysis of equilibrium liquid at 40°C.

<sup>h</sup> By Karl Fisher

<sup>i</sup> Water content calculated by difference.


<sup>j</sup> DMSO content by gas chromatography using virgin DMSO sample as the standard.

CHART SPEED 8.0 CM/MIN  
ATTEN: 64 ZERO: 5% 1 MIN/TICK  
STAT: 0 INJECT  
PS: 1.0



TITLE: 10:28 31 MAR 83  
CHANNEL NO: 1 SAMPLE: "SLUDGE" METHOD: THIOLS  
PEAK NO PEAK NAME RESULT AREA % TIME (MIN) TIME OFFSET AREA COUNTS SEP CODE W1/2 (SEC)  
1 65.3886 5.220 8169 BB ? 15.60  
2 34.6114 11.720 4324 VV ? 27.00  
TOTALS: 100.0000 12493  
DETECTED PKS: 13 REJECTED PKS: 11  
MULTIPLIER: 1.00000  
NOISE: 0.0 OFFSET: -0  
SAVED FILE: RDX006  
NOTES:  
COL: 2 M GLASS - 5% OV-17. 80-100 MESH  
CARRIER: NITROGEN - 20 ML/MIN  
INJ: 150° C DET: 370° C  
TEMP PROG: 100° TO 260° C @ 20°/MIN  
SOLVENT: DMSO SAMP SIZE: 0.5 UL  
DETECTOR: FID SENSITIVITY: 10-10  
RUN LENGTH: 15 MINUTES

# Gas Chromatograph Analysis of DMSO Evaporator Sludge



```


TITLE:                                     9:26      31 MAR 82
CHANNEL NO: 1      SAMPLE: DMSO                      METHOD: THIOLS
PEAK  PEAK  RESULT  TIME  TIME  AREA  SEP  W1/2
NO    NAME  AREA %  (MIN)  OFFSET  COUNTS  CODE  (SEC)
TOTALS:                0.0000                        0

DETECTED PKS:      4      REJECTED PKS:      4
MULTIPLIER: 1.00000
NOISE:      2.4  OFFSET:      -4
SAVED FILE: RDX003
ERRORS:
NO PEAKS
NOTES:
COL: 2 M GLASS - 5% OV-17, 80-100 MESH
CARRIER: NITROGEN - 20 ML/MIN
INJ: 150° C      DET: 270° C
TEMP PROG: 100° TO 260° C @ 20° C/MIN
SOLVENT: DMSO      SAMP SIZE 1UL
DETECTOR: FID -      SENSITIVITY: 10-10
RUN LENGTH: 15 MINUTES

```

## APPENDIX A (cont.)

CHART SPEED 0.8 CM/MIN  
ATTEN: 16 ZERO: 5% 1 MIN/TICK  
LS: 1.0



TI: ON/11: OFF

141

RECALC  
TITLE:

12:13 29 MAR 83

CHANNEL NO: 1

SAMPLE: TP013

METHOD: THIOLS

PEAK NO	PEAK NAME	RESULT AREA %	TIME (MIN)	TIME OFFSET	AREA COUNTS	SEP CODE	HI/2 (SEC)
1		12.4597	5.346		35291	BV	12.50
2		19.7905	6.224		56055	VV	10.70
3		17.5359	7.481		49669	VV	6.75
4		4.2547	8.685		12051	VB	3.10
5		44.8299	11.789		126977	VV	53.15
6		1.1294	12.636		3199	VB	6.40

TOTALS:

100.0000

283242

DETECTED PKS: 15 REJECTED PKS: 9

MULTIPLIER: 1.00000

NOISE: 1.2 OFFSET: -5

NOTES:

COL: 2 M GLASS - 5% DV-17, 80-100 MESH  
CARRIER: NITROGEN - 20 CC/MIN  
INJ: 150° C DET: 270° C  
TEMP PROG: 100° TO 260° @ 20°/MIN  
SOLVENT: DMSO ? SAMP SIZE: 1 UL  
DETECTOR: FID - SENSITIVITY 10-10  
RUN LENGTH: 15 MINUTES

Gas Chromatograph of DMSO Recycle Solvent

APPENDIX A (concluded)

## LIST OF TABLES

		Date	Page
Table 1	Strain Verification for Toxicity Level Determination	19 Jan 83	16
Table 2	Toxicity Level Determination (Virgin DMSO)	19 Jan 83	17
Table 3	Toxicity Level Determination (Recycle DMSO)	19 Jan 83	18
Table 4	Toxicity Level Determination (Evaporator Sludge)	19 Jan 83	19
Table 5	Strain Verification Control	26 Jan 83	20
Table 6	Number of Revertants/Plate	26 Jan 83	21
Table 7	Number of Revertants/Plate (Virgin DMSO)	26 Jan 83	22
Table 8	Number of Revertants/Plate (Recycled DMSO)	26 Jan 83	24
Table 9	Number of Revertants/Plate (Evaporator Sludge)	26 Jan 83	26
Table 10	Strain Verification Control	17 Mar 83	28
Table 11	Number of Revertants/Plate	17 Mar 83	29
Table 12	Number of Revertants/Plate (Evaporator Sludge)	17 Mar 83	30



TABLE 1  
STRAIN VERIFICATION FOR TOXICITY LEVEL DETERMINATION

Strains	Histidine Requirement	Ampicillin Resistance	UV	Sensitivity to Crystal Violet	Sterility Control	Response (1)
100	NG	G	NG	NG 18mm	NG	+
1537	NG	NG - 25mm	NG	NG 15mm	NG	+
WT	G	NA	G	NA	NA	+

STERILITY CONTROL

His-Bio Mix Initial: NG End: NG MGA Plate: NG

Top Agar Initial: NG End: NG

Diluent: NG Nutrient Broth: NG

Test Compound (a) Virgin-NG (b) Recycle-NG (c) Sludge-NG (d) NA (e) NA

G = Growth NG = No Growth NT = Not Tested NA = Not Applicable WT = Wild Type

Spontaneous Revertants: TA 100, No S-9 86, 99, 104 average = 96

(1) + = expected response - = unexpected response

Study Number: 83001 Date: 19 Jan 83 By: Sauers, Kellner

TABLE 2

TOXICITY LEVEL DETERMINATION

Substance assayed: Virgin DMSO      Substance dissolved in: DMSO (reagent grade)  
 Study Number: 83001      Date: 19 Jan 83      Performed by: Sauers, Kellner, Mullen, Dacey

TA 100 REVERTANT PLATE COUNT

Test Compound Concentration	Plate #1	Plate #2	Plate #3	Average	Background Lawn (1)
100% solution	77	74	84	78	NL
10% solution	79	93	108	93	NL
1% solution	98	93	110	100	NL
0.1% solution	111	99	110	107	NL
0.01% solution	95	108	118	107	NL
0.001% solution	77	81	79	79	NL
0.0001% solution	70	95	89	85	NL
0.00001% solution	81	70	76	76	NL

(1) NG = No Growth      ST = Slight Growth      NL = Normal Lawn

TABLE 3

TOXICITY LEVEL DETERMINATION

Substance assayed: DMSO Recycle Solvent Substance dissolved in: Reagent grade DMSO  
 Study Number: 83001 Date: 19 Jan 83 Performed by: Sauers, Kellner, Mullen, Dacey

TA 100 REVERTANT PLATE COUNT

Test Compound Concentration	Plate #1	Plate #2	Plate #3	Average	Background Lawn (1)
100% solution	132	116	123	124	NL
10% solution	106	84	91	94	NL
1% solution	94	96	93	94	NL
0.1% solution	92	107	77	92	NL
0.01% solution	81	86	95	87	NL
0.001% solution	82	89	99	90	NL
0.0001% solution	83	97	111	97	NL
0.00001% solution	89	90	96	92	NL

(1) NC = No Growth      ST = Slight Growth      NL = Normal Lawn

TABLE 4

TOXICITY LEVEL DETERMINATION

Substance assayed: Evaporator Sludge Substance dissolved in: Reagent grade DMSO  
 Study Number: 83001 Date: 19 Jan 83 Performed by: Sauers, Kellner, Mullen, Dacey

TA 100 REVERTANT PLATE COUNT

Test Compound Concentration	Plate #1	Plate #2	Plate #3	Average	Background Lawn (1)
100% solution	153	164	159	159	NL
10% solution	117	119	89	108	NL
1% solution	105	96	106	102	NL
0.1% solution	107	87	106	100	NL
0.01% solution	88	79	88	85	NL
0.001% solution	84	99	82	88	NL
0.0001% solution	83	87	89	86	NL
0.00001% solution	83	90	94	89	NL

(1) NG = No Growth      ST = Slight Growth      NL = Normal Lawn

TABLE 5  
STRAIN VERIFICATION CONTROL

Strains	Histidine Requirement	Ampicillin Resistance	UV	Sensitivity to Crystal Violet	Sterility Control	Response (1)
98	NG	G	NG	NG 15mm	NG	+
100	NG	G	NG	NG 20mm	NG	+
1535	NG	NA	NG	NG 15mm	NG	+
1537	NG	NG - 29mm	NG	NG 13mm	NG	+
1538	NG	NA	NG	NG 15mm	NG	+
WT	G	NA	G	NA	NA	+

STERILITY CONTROL

His-Bio Mix Initial: NG End: NG Diluent: NG

Top Agar Initial: NG End: NG MCA Plate: NG

S-9 Mix Initial: NG End: NG Nutrient Broth: NG

Test Compound (a) Virgin-NG (b) Recycle-NG (c) Sludge-NG (d) NA (e) NA (f) NA

G = Growth NG = No Growth NT = Not Tested NA = Not Applicable WT = Wild Type

Study Number: 83001 By: Sauers, Kellner (1) + = expected response

Date: 26 Jan 83 - = unexpected response

TABLE 6  
NUMBER OF REVERTANTS/PLATE

Compd.	Amount of Compd. Added	S-9 Added	98	100	Strain No. 1535	1537	1538
AF	2 ug/plate	yes	(256, 497, 339) 364	(278, 335, con) 307			(408, 645, 368) 474
BP	2 ug/plate	yes	( 92, 97, 67) 85	(586, 412, 627) 542		( 81, 86, 70) 79	(115, 142, 112) 123
AA	2 ug/plate	yes	(709, 758, 934) 800	(999, 999, 999) 999		(116, 68, 82) 89	(999, 867, 999) 955
MING	2 ug/plate	no		(999, 999, 999) 999			
	20 ug/plate	no			(999, 999, 999) 999		

Spontaneous Reversion Rate

before	yes	( 26, 23, 31) ( 36, 37, 39) 32	( 98, 95, 96) (128, 134, 134) 114	( 20, 11, con) ( 14, 20, 13) 16	( 5, 4, 7) ( 3, 10, 10) 7	( 18, 16, 18) ( 17, 28, 28) 21
after	no	( 11, 21, 26) ( 36, 25, 23) 24	( 95, 102, 92) (129, 113, 119) 108	( 22, 23, 19) ( 22, 22, 24) 22	( 7, 6, 10) ( 7, 8, 6) 7	( 14, 24, 12) ( 20, 14, 19) 17

con - plate value was disregarded due to contamination  
999 - signifies a revertant count greater than 1000

Study Number: 83001

Date: 26 Jan 83 By: Sauers, Kellner, Mullen, Dacey

TABLE 7

NUMBER OF REVERTANTS/PLATE

Compd.	0.1 ml Solution Added	S-9 Added	Strain Number				
			98	100	1535	1537	1538
Virgin DMSO	100% solution	no	( 26, 19, 15) 20	( 88, 99, 95) 94	( 22, 16, 15) 18	( 6, 6, 3) 5	( 14, 15, 18) 16
		yes	( 28, 34, 26) 29	(122, 92, 112) 109	( 18, 10, 11) 13	( 6, 3, 6) 5	( 15, 18, 17) 17
Virgin DMSO	20% solution	no	( 16, 16, 17) 16	( 96, 106, 95) 99	( 21, 14, 20) 18	( 2, 5, 10) 6	( 15, 13, 14) 14
		yes	( 32, 32, 24) 29	( 99, 88, 101) 96	( 17, 16, 10) 14	( 4, 6, 7) 6	( 16, 13, 19) 16
Virgin DMSO	4% solution	no	( 28, 23, 14) 22	(105, 89, con) 97	( 19, 18, 19) 19	( 11, 5, 6) 7	( 11, 12, 14) 12
		yes	( 23, 17, 24) 21	( 89, 99, 101) 96	( 17, 15, 23) 18	( 6, 2, 2) 3	( 22, 17, con) 20
Virgin DMSO	0.8% solution	no	( 21, 17, 19) 19	(102, 98, 95) 98	( 20, 11, 26) 19	( 5, 6, 3) 5	( 11, 14, 14) 13
		yes	( 29, 29, 29) 29	(104, 99, 120) 108	( 18, 19, 14) 17	( 9, 3, 7) 6	( 29, 19, 17) 22

con - plate value was disregarded due to contamination

Study Number: 83001

Date: 26 Jan 83

By: Sauers, Kellner, Mullen, Dacey

TABLE 7 (cont.)  
NUMBER OF REVERTANTS/PLATE

Compd.	Q.1 ml Solution Added	S-9 Added	Strain Number			
			98	100	1535	1537 1538
Virgin DMSO	0.16% solution	no	( 17, 18, 14) 16	(120,102,101) 108	( 25, 24, 31) 27	( 3, 7, 6) ( 15, 18, 14) 5 16
		yes	( 25, 32, 27) 28	(101,119,104) 108	( 10, 38, 14) 21	( 3, 6, 5) ( 23, 20, 24) 5 22
Virgin DMSO	0.032% solution	no	( 20, 19, 13) 17	(121,106, 89) 105	( 20, 12, 16) 16	( 4, 3, 5) ( 9, 12, 15) 4 12
		yes	( 31, 22, 25) 26	( 86,103,127) 105	( 18, 12, 14) 15	( 6, 3, 4) ( 20, 22, 24) 4 22

Study Number: 83001 Date: 26 Jan 83 By: Sauers, Kellner, Mullen, Dacey



TABLE 8  
NUMBER OF REVERTANTS/PLATE

Compd.	0.1 ml Solution Added	S-9 Added	98	100	Strain Number		
					1535	1537	1538
Recycled DMSO	100% solution	no	( 20, 15, 16) 17	(116, 98, 88) 101	( 12, 16, 12) 13	( 6, 7, 7) 7	( 13, 17, 16) 15
		yes	( 18, 16, 21) 18	(101, 118, 104) 108	( 12, 9, 11) 11	( 7, 8, 4) 6	( 20, 16, 12) 16
Recycled DMSO	20% solution	no	( 25, 21, 25) 24	(107, 107, 118) 111	( 14, 25, 18) 19	( 11, 8, 5) 8	( 22, 23, 21) 22
		yes	( 34, 30, 25) 30	(111, 117, 99) 109	( 27, 11, 15) 18	( 5, 9, 11) 8	( 26, 25, 21) 24
Recycled DMSO	4% solution	no	( 18, 19, 19) 19	(112, 105, 140) 119	( 18, 18, 35) 24	( 7, 5, 7) 6	( 20, 16, 13) 16
		yes	( 33, 34, con) 34	(114, 119, 86) 106	( 26, 19, 20) 22	( 4, 8, 9) 7	( 24, 20, 21) 22
Recycled DMSO	0.8% solution	no	( 24, 26, 21) 24	(101, 114, 114) 110	( 20, 21, 23) 21	( 7, 7, son) 7	( 22, 20, 21) 21
		yes	( 20, 21, 29) 23	( 91, 88, 109) 96	( 23, 21, 20) 21	( 1, 4, 6) 4	( 29, 20, 21) 23

con - plate value was disregarded due to contamination

Study Number: 83001

Date: 26 Jan 83

By: Sauers, Kellner, Mullen, Dacey

TABLE 8 (cont.)  
NUMBER OF REVERTANTS/PLATE

Compd.	Solution Added	S-9 Added	Strain Number			
			98	100	1535	1537
Recycled DMSO	0.16% solution	no	(25, 21, 17) 21	(111, 120, 95) 109	(27, 27, 22) 25	(4, 14, 6) 8
		yes	(42, 34, 33) 36	(104, 90, 99) 98	(10, 17, 27) 18	(5, 2, 6) 4
Recycled DMSO	0.032% solution	no	(16, 19, 23) 19	(110, 99, 110) 106	(25, 22, 14) 20	(7, 4, 5) 5
		yes	(38, 18, con) 28	(123, 107, 118) 116	(13, 11, 13) 12	(4, 4, 4) 4

con - plate value was disregarded due to contamination

Study Number: 83001      Date: 26 Jan 83      By: Sauers, Kellner, Mullen, Dacey

TABLE 9  
NUMBER OF REVERTANTS/PLATE

Compd.	Solution Added	S-9 Added	98	100	Strain Number		
					1535	1537	1538
Evaporator sludge	100% solution	no	(115, 84, 121) 107	(140, 137, 154) 144	(42, 27, 38) 36	(21, 27, 14) 21	(123, 103, 96) 107
		yes	(53, 58, 78) 63	(126, 159, 142) 142	(29, 26, 25) 27	(23, 21, 20) 21	(66, 55, 74) 65
Evaporator sludge	20% solution	no	(33, 29, 29) 30	(110, 109, 156) 125	(22, 22, 21) 22	(14, 5, 6) 8	(30, 21, 29) 27
		yes	(38, 33, 27) 33	(113, 86, 105) 101	(18, 15, 20) 18	(10, 5, 8) 8	(33, 29, 26) 29
Evaporator sludge	4% solution	no	(24, 20, 26) 23	(95, 87, 97) 93	(22, 19, 22) 21	(4, 3, 4) 4	(19, 12, 14) 15
		yes	(36, 22, 27) 28	(126, 120, 108) 118	(7, 15, 13) 12	(4, 5, 3) 4	(19, 23, 23) 22
Evaporator sludge	0.8% solution	no	(24, 18, 15) 19	(94, 120, 103) 106	(15, 18, 18) 17	(5, 6, 4) 5	(15, 26, 22) 21
		yes	(35, 24, 00n) 30	(125, 104, 111) 113	(17, 16, 12) 15	(5, 6, 9) 7	(27, 27, 15) 23

con - plate value disregarded due to contamination

Study Number: 83001 Date: 26 Jan 83 By: Sauers, Kellner, Mullen, Dacey

TABLE 9 (cont.)

NUMBER OF REVERTANTS/PLATE

Compd.	0.1 ml Solution Added	S-9 Added	Strain Number			
			98	100	1535	1537
Evaporator sludge	0.16% solution	no	( 14, 23, 26) 21	(111,101,114) 109	( 22, 22, 23) 22	( 2, 5, 6) 4
		yes	( 29, 21, 30) 27	(123,102,108) 111	( 19, 15, 18) 17	( 5, 4, 3) 4
Evaporator sludge	0.032% solution	no	( 16, 18, 23) 19	(107, 98,109) 105	( 27, 19, 13) 20	( 7, 5, 3) 5
		yes	( 39, 26, 23) 29	(125,107,125) 119	( 12, 18, 23) 18	( 7, 5, 8) 7

Study Number: 83001

Date: 26 Jan 83

By:

Sauers, Kellner, Mullen, Dacey

TABLE 10

STRAIN VERIFICATION CONTROL

Strains	Histidine Requirement	Ampicillin Resistance	UV	Sensitivity to Crystal Violet	Sterility Control	Response (1)
98	NG	G	NG	NG (15 mm)	NG	+
100	NG	G	NG	NG (18 mm)	NG	+
1535	NG	NA	NG	NG (16 mm)	NG	+
1537	NG	NG (28 mm)	NG	NG (13 mm)	NG	+
1538	NG	NA	NG	NG (19 mm)	NG	+
WT	NA	NA	G	G	NA	+

STERILITY CONTROL

His-Bio Mix      Initial: NG      End: NG      Diluent: NG  
 Top Agar          Initial: NG      End: NG      MGA Plate: NG  
 S-9 Mix            Initial: NG      End: NG      Nutrient Broth: NG  
                          Evaporator  
 Test Compound    (a) sludge - NG(b)    NA    (c) NA    (d) NA    (e) NA    (f) NA

G = Growth    NG = No Growth    NT = Not Tested    NA = Not Applicable    WT = Wild Type

Study Number: 83001      By: Sauers, Kellner      (1) + = expected response

Date: 17 Mar 83      - = unexpected response

TABLE 11

NUMBER OF REVERTANTS/PLATE

Compd.	Amount of Compd. Added	S-9 Added	98	100	Strain No. 1535	1537	1538
AF	2 ug/plate	yes	(439,344,359) 381	(289,385,148) 274	(553,516,452) 507		
BP	2 ug/plate	yes	( 88,102, 76) 89	(440,422,487) 450	( 74, 66, 70) 70	( 97,102, 97) 99	
AA	2 ug/plate	yes	(999,999,999) 999	(999,999,999) 999	(243,176,186) 202	(965,999,999) 988	
MUNG	2 ug/plate	no		(999,999,999) 999			
	20 ug/plate	no			(999,999,999) 999		
<u>Spontaneous Reversion Rate</u>							
	before	yes	( 30, 34, 39) ( 27, 26, 34) 32	( 97,104,111) ( 99,103,112) 104	( 8, 14, 12) ( 9, 12, 13) 11	( 5, 4, 3) ( 5, 2, 7) 4	( 21, 19, 17) ( 22, 20, 25) 21
	after	no	( 17, 30, 15) ( 14, 20, 29) 21	( 95, 82,101) (104,109, 97) 98	( 22, 13, 15) ( 10, 14, 18) 15	( 4, 4, 4) ( 4, 9, 6) 6	( 9, 17, 12) ( 19, 17, 10) 14

999 - signifies a revertant count greater than 1000

Study Number: 83001

Date: 17 Mar 83 By: Sauers, Kellner, Mullen, Dacey

TABLE 12  
NUMBER OF REVERTANTS/PLATE

	Compd.	0.1 ml Solution Added	S-9 Added	Strain Number				
				98	100	1535	1537	1538
1	Evaporator sludge	100% solution	no	(122, 72, 183) 126	(153, 149, 144) 149	(17, 18, 26) 20	(27, 14, 23) 21	(119, 140, 131) 130
2			yes	(72, 71, 76) 73	(139, 141, 120) 133	(27, 21, 29) 26	(21, 14, 20) 18	(82, 71, 54) 69
3	Evaporator sludge	80% solution	no	(77, 41, 129) 82	(158, 119, 129) 135	(20, 19, 12) 17	(33, 18, 20) 24	(85, 107, 108) 100
4			yes	(71, 59, 65) 65	(127, 126, 130) 128	(24, 23, 22) 23	(16, 17, 15) 16	(62, 84, 65) 70
5	Evaporator sludge	60% solution	no	(59, 76, 55) 63	(112, 117, 117) 115	(14, 10, 17) 14	(14, 17, 17) 16	(54, 60, 80) 65
6			yes	(52, 50, 46) 49	(136, 116, 118) 123	(23, 18, 32) 24	(6, 12, 10) 9	(48, 54, 62) 55
7	Evaporator sludge	40% solution	no	(46, 47, 34) 42	(116, 142, 136) 131	(14, 17, 13) 15	(11, 7, 12) 10	(37, 53, 46) 45
8			yes	(36, 35, 44) 38	(106, 103, 94) 101	(22, 20, 15) 19	(14, 7, 11) 11	(51, 42, 40) 44

Study Number: 83001 Date: 17 Mar 83 By: Sauers, Kellner, Mullen, Dacey

12  
cont.

Sauers--31

TABLE 12 (cont.)

NUMBER OF REVERTANTS/PLATE

	Compd.	0.1 ml Solution Added	S-9 Added	98	100	Strain Number		
						1535	1537	1538
3	Evaporator sludge	20% solution	no	( 40, 36, 39) 38	(122, 117, 124) 121	( 14, 19, 21) 18	( 6, 4, 8) 6	( 34, 34, 30) 33
0			yes	( 31, 32, 27) 30	(127, 101, 116) 115	( 19, 13, 14) 15	( 5, 6, 5) 5	( 43, 34, 45) 41
1	Evaporator sludge	1% solution	no	( 18, 12, 14) 15	( 84, 74, 118) 92	( 7, 17, 23) 16	( 6, 4, 8) 6	( 9, 13, 7) 10
2			yes	( 26, 33, 27) 29	(103, 123, 96) 107	( 7, 12, 10) 10	( 7, 8, 4) 6	( 13, 25, 18) 19

Study Number: 83001      Date: 17 Mar 83      By: Sauers, Kellner, Mullen, Dacey



LIST OF FIGURES

	Page
FIGURE 1 Dose Response for TA 98 with S-9.....	34
FIGURE 2 Dose Response for TA 98 without S-9.....	35
FIGURE 3 Dose Response for TA 1537 with S-9.....	36
FIGURE 4 Dose Response for TA 1537 without S-9.....	37
FIGURE 5 Dose Response for TA 1538 with S-9.....	38
FIGURE 6 Dose Response for TA 1538 without S-9.....	39

FIGURE 1




FIGURE 2




FIGURE 3




FIGURE 4




FIGURE 5





FIGURE 6



OFFICIAL DISTRIBUTION LIST

Commander  
US Army Medical Research  
and Development Command  
ATTN: SGRD-RMS/Mrs. Madigan  
Fort Detrick, Frederick MD 21701

Defense Technical Information Center  
ATTN: DTIC-DDA (12 copies)  
Cameron Station  
Alexandria VA 22314

Director of Defense Research and Engineering  
ATTN: Assistant Director, Environmental  
and Life Sciences  
Washington DC 20301

The Surgeon General  
ATTN: DASG-TLO  
Washington DC 20314

HQ DA (DASG-ZXA)  
WASH DC 20310

Commandant  
Academy of Health Sciences  
ATTN: HSHA-CDM  
Fort Sam Houston TX 78234

Assistant Dean  
Institute and Research Support  
Uniformed Services University  
of Health Sciences  
6917 Arlington Road  
Bethesda MD 20014

Commander  
US Army Environmental Hygiene Agency  
Aberdeen Proving Ground MD 21070

US Army Research Office  
ATTN: Chemical and Biological Sciences  
Division  
P.O. Box 1221  
Research Triangle Park NC 27709

Biological Sciences Division  
Office of Naval Research  
Arlington VA 22217

Director of Life Sciences  
USAF Office of Scientific Research (AFSC)  
Bolling AFB  
Washington DC 20332

Director  
Walter Reed Army Institute of Research  
Washington DC 20307

Commander  
US Army Medical Research Institute  
of Infectious Diseases  
Fort Detrick, Frederick MD 21701

Commander  
US Army Research Institute  
of Environmental Medicine  
Natick MA 01760

Commander  
US Army Institute of Surgical Research  
Brooke Army Medical Center  
Fort Sam Houston TX 78234

Commander  
US Army Medical Bioengineering  
Research and Development Laboratory  
Fort Detrick, Frederick MD 21701

Commander  
US Army Aeromedical Research Laboratory  
Fort Rucker AL 36362

Commander  
US Army Research Institute  
of Chemical Defense  
Aberdeen Proving Ground  
Edgewood Arsenal MD 21010

Commander  
Naval Medical Research Institute  
National Naval Medical Center  
Bethesda MD 20014

Commander  
USAF School of Aerospace Medicine  
Aerospace Medical Division  
Brooks Air Force Base TX 78235