


AD-A166 597 THE USE OF A PYROCHEMILUMINESCENT NITROGEN ANALYZER AS 1/1  
A CAPILLARY GC (GA (U) AIR FORCE WRIGHT AERONAUTICAL  
LABS WRIGHT-PATTERSON AFB OH E W PITZER FEB 86  
UNCLASSIFIED AFWAL-TR-85-2077 F/G 7/4 NL






MICROCOPY RESOLUTION TEST CHART  
NATIONAL BUREAU OF STANDARDS-1963-A

AD-A166 597

AFWAL-TR-85-2077

2



# THE USE OF A PYROCHEMILUMINESCENT NITROGEN ANALYZER AS A CAPILLARY G.C. DETECTOR

Edward W. Pitzer

Experimental Support Branch  
Aero Propulsion Laboratory

February 1986

Final Report for Period March 1985 - May 1985

Approved for Public Release, Distribution Unlimited.

DTIC FILE COPY

DTIC  
ELECTE  
APR 15 1986  
S D  
E

AERO PROPULSION LABORATORY  
AIR FORCE WRIGHT AERONAUTICAL LABORATORIES  
AIR FORCE SYSTEMS COMMAND  
WRIGHT-PATTERSON AIR FORCE BASE, OHIO 45433-6563


86 4 14 047

NOTICE

When Government drawings, specifications, or other data are used for any purpose other than in connection with a definitely related Government procurement operation, the United States Government thereby incurs no responsibility nor any obligation whatsoever; and the fact that the government may have formulated, furnished, or in any way supplied the said drawings, specifications, or other data, is not to be regarded by implication or otherwise as in any manner licensing the holder or any other person or corporation, or conveying any rights or permission to manufacture use, or sell any patented invention that may in any way be related thereto.


This report has been reviewed by the Office of Public Affairs (ASD/PA) and is releasable to the National Technical Information Service (NTIS). At NTIS, it will be available to the general public, including foreign nations.

This technical report has been reviewed and is approved for publication.



EDWARD W. PITZER  
PROJECT SCIENTIST, EXPERIMENTAL  
SUPPORT BRANCH, FUELS AND  
LUBRICATION DIVISION

FOR THE COMMANDER



ARTHUR V. CHURCHILL  
CHIEF, FUELS BRANCH  
FUELS AND LUBRICATION DIVISION  
AERO PROPULSION LABORATORY



ROBERT D. SHERRILL  
CHIEF, FUELS AND LUBRICATION DIVISION  
AERO PROPULSION LABORATORY

"If your address has changed, if you wish to be removed from our mailing list, or if the addressee is no longer employed by your organization please notify AFWAL/PO5X, W-PAFB, OH 45433 to help us maintain a current mailing list".

Copies of this report should not be returned unless return is required by security considerations, contractual obligations, or notice on a specific document.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE

ADA 166597

## REPORT DOCUMENTATION PAGE

1a. REPORT SECURITY CLASSIFICATION Unclassified			1b. RESTRICTIVE MARKINGS N/A	
2a. SECURITY CLASSIFICATION AUTHORITY N/A			3. DISTRIBUTION/AVAILABILITY OF REPORT Approved for public release; distribution unlimited.	
2b. DECLASSIFICATION/DOWNGRADING SCHEDULE N/A				
4. PERFORMING ORGANIZATION REPORT NUMBER(S)			5. MONITORING ORGANIZATION REPORT NUMBER(S) AFWAL-TR-85-2077	
6a. NAME OF PERFORMING ORGANIZATION Aero Propulsion Laboratory		6b. OFFICE SYMBOL (If applicable) AFWAL/POSF		7a. NAME OF MONITORING ORGANIZATION Experimental Support Branch Aero Propulsion Laboratory
6c. ADDRESS (City, State and ZIP Code) AFWAL/POSX Air Force Wright Aeronautical Laboratories Wright-Patterson AFB OH 45433-6563			7b. ADDRESS (City, State and ZIP Code) AFWAL/POSX Air Force Wright Aeronautical Laboratories Wright-Patterson AFB OH 45433-6563	
8a. NAME OF FUNDING/SPONSORING ORGANIZATION Aero Propulsion Laboratory		8b. OFFICE SYMBOL (If applicable) AFWAL/POSF		9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER N/A
8c. ADDRESS (City, State and ZIP Code) AFWAL/POSX Air Force Wright Aeronautical Laboratories Wright-Patterson AFB OH 45433-6563			10. SOURCE OF FUNDING NOS.	
			PROGRAM ELEMENT NO. 62203F	PROJECT NO. 3048
			TASK NO. 05	WORK UNIT NO. 91
11. TITLE (Include Security Classification) See Block 16				
12. PERSONAL AUTHOR(S) E.W. Pitzer				
13a. TYPE OF REPORT Final		13b. TIME COVERED FROM Mar 85 TO May 85		14. DATE OF REPORT (Yr., Mo., Day) February 1986
15. PAGE COUNT 17				
16. SUPPLEMENTARY NOTATION The Use of a Pyrochemiluminescent Nitrogen Analyzer as a Capillary GC. Dector				
17. COSATI CODES			18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number)	
FIELD	GROUP	SUB. GR.	Pyrochemiluminescence, Gas Chromatograph (GC), Chemiluminescent Detector (CLD), Flame Ionization Detector, (FID), Parts per million nitrogen (ppm-N), Relative Response Factor (RRF)	
07	01	03		
21	04	05		
19. ABSTRACT (Continue on reverse if necessary and identify by block number)				
<p>A study was conducted to determine the feasibility of interfacing a pyrochemiluminescent nitrogen analyzer to a high resolution glass capillary gas chromatograph. The instruments used were a Varian Model 3700 gas chromatograph and an Antek Model 703C pyrochemiluminescent nitrogen analyzer.</p> <p>Known standards of nitrogen containing compounds were prepared in shale derived jet propulsion fuel that had been previously determined to have no nitrogen containing species. These solutions were then chromatographed on the GC/CLD system. The resulting data were reduced using relative response factors with an internal standard. The nitrogen content of the standard solutions were validated by analyses on a separate Antek nitrogen analyzer used to determine total nitrogen content.</p> <p>The method of GC/CLD was deemed feasible. However, it is recommended that appropriate design changes be made to the Antek Model 703C nitrogen analyzer to accommodate the type of flows common with capillary gas chromatographs. <i>Keywords:</i></p>				
20. DISTRIBUTION/AVAILABILITY OF ABSTRACT UNCLASSIFIED/UNLIMITED <input checked="" type="checkbox"/> SAME AS RPT. <input type="checkbox"/> DTIC USERS <input type="checkbox"/>			21. ABSTRACT SECURITY CLASSIFICATION UNCLASSIFIED	
22a. NAME OF RESPONSIBLE INDIVIDUAL EDWARD W. PITZER			22b. TELEPHONE NUMBER (Include Area Code) 513/255-3887	22c. OFFICE SYMBOL AFWAL/POSX

## FOREWORD

This research was conducted in the Experimental Support Branch, Fuels and Lubrication Division, Aero Propulsion Laboratory, Air Force Wright Aeronautical Laboratories, Wright-Patterson Air Force Base, Ohio under Project No. 3048, Task No. 304805, Project No. 30480591 by Edward W. Pitzer as Project Scientist.

This research was conducted during the period March 1985 - May 1985.

Accession For	
NTIS GRA&I	<input checked="checked" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By _____	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A-1	



## TABLE OF CONTENTS

SECTION		PAGE
I	INTRODUCTION	1
II	EXPERIMENTAL	3
	1. Overview	3
	2. Sample Preparation	3
	3. Methodology	3
III	RESULTS AND RECOMMENDATIONS	11
	REFERENCES	14

## LIST OF ILLUSTRATIONS

FIGURE		PAGE
1	Quartz Combustion Tube	2
2	Block Diagram of the Nitrogen Analyzer	2
3	GC/CLD Chromatogram	7
4	GC/FID Chromatogram	9
5	Third Order Calibration Curve	10
6	Measurement of Signal-to-Noise Ratio	12

## LIST OF TABLES

TABLE		PAGE
1	Chromatographic Operating Conditions	4
2	C.L.D. Chromatographic Data	5
3	Antek Calibration Curves	8
4	Effect of Basicities on Response Factors	9

# LIST OF ABBREVIATIONS

GC	Gas Chromatograph
CLD	Chemiluminescent Detector
FID	Flame Ionization Detector
ppm-N	Parts per million Nitrogen
RRF	Relative Response Factor
C	Degrees Celsius
RT	Retention Time
Kb	Base Dissociation Constant
ng-N	nanograms Nitrogen
nm	nanometers

## SECTION I



## INTRODUCTION

The Antek Model 703C nitrogen analyzer is a very reliable, and highly quantitative method for determining the total bound nitrogen content in a wide variety of sample matrices (References 1,2,3). The excellent quantitative capability is attributable to the high temperature (1000+ Deg. C) and abundant oxygen supply in the quartz combustion tube which allows for total oxidation of the sample. Figure 1 shows a typical combustion tube used with Antek nitrogen analyzers.

Nitrogen containing organic species introduced to the combustion tube are converted, via oxidative pyrolysis, to carbon dioxide, water, and nitric oxide as shown below (Reference 2).



The resultant nitric oxide is then routed to a reaction chamber which is supplied with ozone from an internal ozone generator. In this reaction chamber the nitric oxide reacts with the ozone to form a metastable nitrogen dioxide, which upon relaxation to its stable state, emits a photon of light in the wavelength range of 700-900 nm (Reference 2).



The Antek Model 703C measures these photons of light with a photomultiplier tube and bandpass filter. Figure 2 shows a block diagram of the nitrogen analyzer system used in this study.

The objective of this study was to determine if the eluents of a capillary gas chromatograph could be detected by a nitrogen analyzer to take advantage of its nitrogen specificity and excellent quantitative capability of the analyzer.




Figure 1. Quartz Combustion Tube




Figure 2. Block Diagram of the Nitrogen Analyzer

## SECTION II

## EXPERIMENTAL

## 1. OVERVIEW

A Varian 3700 gas chromatograph, equipped with a split/splitless injector was interfaced with an Antek Model 703C nitrogen analyzer via a heated transfer line. The fused silica column used for the chromatography was routed through the transfer line, through a high temperature septum, and inserted into the quartz combustion tube of the nitrogen analyzer. Initially, the column end extended four inches into the combustion tube (Figure 1). Inspection of the column end after initial operations revealed that the polyimide coating of the column had been totally oxidized along a section of the column end from 1.25 to 3.5 inches. This made it apparent that the column eluents were being released downstream of the optimum oxidation zone of the combustion tube. Accordingly, the column end was inserted only 1.25 inches into the combustion tube which resulted in increased response. Table 1 lists the operating conditions and parts of the GC/CLD system.

## 2. SAMPLE PREPARATION

Two solutions of nitrogen containing compounds were prepared, gravimetrically, in shale derived jet propulsion fuel that had been previously determined to have no nitrogen containing species. These solutions contained a total of 19 different compounds. Each solution contained approximately the same amount of 2,5-dimethylpyrrole used as an internal standard. The solutions were then gravimetrically diluted in the same shale fuel to concentrations in the 5 to 10 ppm-N range. Table 2 lists the compounds used, their concentrations, and the resultant chromatographic data from this study.

## 3. METHODOLOGY

A previous study (in this laboratory) on the column chosen for this research, revealed that the minimum operating pressure for the column was 46 PSIG (Reference 4). The column was operated at this minimum pressure to avoid introducing the sample to the combustion tube too quickly.

TABLE 1  
CHROMATOGRAPHIC OPERATING CONDITIONS

Detector	Pyrochemiluminescent (Antek Model 703C)
Detector Temp.	1050 Deg. C
Gas Chromatograph	Varian Model 3700 with Split/Splitless Injector
Injector	Split Mode, Fritted Glass Liner
Injector Temp.	270 Deg. C
Liquid Sample	5 uL (Manual Injection), Split 1:2
Detector Make-up	100 ml/min Argon
Det. Pyrolysis 02	260 ml/min
Det. Ozonolysis 02	120 ml/min
Column	DB-1 (Fused Silica Bonded Phase MeSi) WCOT, 60 M, 0.32 mm ID, 1.0 uM Film Thickness, 80,000 Effective Theoretical Plates
Column Flow Rate	8.8 ml/min Argon at 46 PSIG and 40 Deg. C
Column Temp. Profile	40 to 240 Deg. C at 5 Deg. C/min
Data Processing	Antek Attenuation 5 to a Hewlett Packard 3357 Computer via an HP18652A A/D Converter

TABLE 2  
C.L.D. CHROMATOGRAPHIC DATA

Compound -----	R.T. (min) -----	[ppm-N] -----	R.R.F. (mean) -----	Std.Dev. -----	95%CL -----
Pyrrolidine	5.46	13.07	1.85333	0.20	0.50
1-Methyl- pyrrolidine	5.50	9.88	1.04798	0.03	0.06
Piperidine	7.40	8.36	1.29789	0.00	0.01
1-Methyl- piperidine	7.60	6.59	1.03673	0.02	0.04
2-Methyl- piperidine	8.55	4.63	1.09290	0.02	0.06
2-Methyl- pyridine	8.69	5.35	0.98190	0.06	0.15
3-Methyl- piperidine	9.37	6.64	1.79694	0.14	0.36
2,6-Dimethyl- pyridine	9.65	6.03	1.01336	0.01	0.01
4-Methyl- pyridine	10.13	11.50	1.02576	0.04	0.10
2,6-Dimethyl- pyridine	11.16	4.82	1.01336	0.01	0.01
2,5-Dimethyl- pyrrole	12.15	8.96/9.17	1.00000	(Internal Standard)	
2,2,6,6-Tetra- methylpiperidine	14.14	6.76	1.48471	0.30	0.75
3,5-Dimethyl- pyridine	14.32	9.54	1.03570	0.03	0.07
3,4-Dimethyl- pyridine	15.19	9.72	1.05710	0.03	0.09
2,3,6-Tri- methylpyridine	15.66	6.27	1.06642	0.05	0.13
2,4-Dimethyl- 3-Ethylpyrrole	18.83	4.81	1.30819	0.01	0.01
1-Methylindole	24.36	4.22	1.06968	0.04	0.10
2-Methylindole	27.37	5.91	1.13143	0.10	0.25
5-Methylindole	27.56	4.57	1.00503	0.02	0.06

Figure 3 is a GC/CLD chromatogram of one of the dilute standard solutions analyzed. When compared to Figure 4, a GC/FID chromatogram of the same solution, the selectivity of the GC/CLD detector for nitrogen compounds is readily apparent.

A separate Antek nitrogen analyzer was used to validate the nitrogen content of the dilute standard solutions. The analyzer had been previously calibrated with series of dilutions of 2-methylpiperidine in toluene. Table 3 lists the results of four different calibration curves investigated for the concentration range of the dilute solutions. The data in this table are listed as the back-calculated amount of the standard solution over the relative percent error. For example, the 20.33 ng-N standard was predicted by the first order curve to be 19.95 ng-N which is a relative error of 2.13%. The results were compared by each curve's relative ability to predict its own calibration points. The obviously superior third order curve (Figure 5) was used to analyze the dilute standard solutions with the following results.

	<u>ppm-N(known)</u>	<u>ppm-N(pred.)</u>	<u>Rel % Error</u>
Solution 1	74.57	74.28	0.38
Solution 2	72.23	71.61	0.86

The GC/CLD response of each compound analyzed, relative to the internal standard 2,5-dimethylpyrrole, is listed in Table 2.




Figure 3. GC/CLD Chromatogram

TABLE 3  
ANTEK CALIBRATION CURVES

PREDICTED ng-N / RELATIVE % ERROR					
ng-N	Area Counts	Exponential	First Order	Second Order	Third Order
8.75	1091	7.86 / 10.20	9.95 / 13.73	9.85 / 12.57	8.81 / 0.69
16.01	2073	14.93 / 6.74	16.13 / 0.73	16.07 / 0.69	15.99 / 0.12
20.33	2680	19.30 / 5.29	19.95 / 2.13	19.91 / 2.31	20.31 / 0.34
38.20	5342	38.48 / 0.72	36.69 / 3.96	36.75 / 3.80	38.25 / 0.13
78.35	12065	86.90 / 10.91	78.97 / 0.79	79.18 / 1.06	78.34 / 0.01
212.48	33293	239.72 / 12.85	212.49 / 0.00	212.41 / 0.03	212.48 / 0.00

Exponential:  $ng-N = \text{Area Counts} \exp(a)$ ;  $a = \ln(ng-n/\text{Area Counts})$   
 $a = -4.93333$

First Order: Coefficient of Correlation = 0.999928

Second Order: Coefficient of Correlation = 0.999930

Third Order: Coefficient of Correlation = 1.000000




Figure 4. GC/FID Chromatogram




Figure 5. Third Order Calibration Curve

## SECTION III

## RESULTS AND RECOMMENDATIONS

The nitrogen concentrations listed in Table 2 by no means indicate that nitrogen compounds do not exist at levels much lower than these in shale derived jet fuels. The Antek 703C nitrogen analyzer was not designed to receive eluents from capillary GC columns. Only packed column GC/CLD work has been published to date (Reference 2). Accordingly, a wide bore column, large injection volume, and relatively concentrated samples were used in order to judge the merit of the particular GC/CLD system reported herein.

The minimum detectable limit was investigated by preparing a solution containing a 1.43 ppm-N concentration of one of the higher boiling compounds studied (2-methylindole). Figure 6 shows the measurements of the signal and noise levels which resulted in a signal-to-noise ratio of approximately 6. This indicates that the minimum detectable limit for 2-methylindole, on the system described in this research, is in the range of 0.5 ppm-N.

It would seem obvious that a much smaller combustion tube, and slower flow rates of both pyrolysis and ozonolysis oxygen would result in a much lower minimum detectable limit.

If the Antek nitrogen analyzer totally oxidizes the sample molecules, any relative response factors generated, other than unity response factors, must be due to injection or column effects. The response factors listed in Table 2 show two definite trends: relative response factors greater than one for non-aromatic nitrogen compounds, and relative response factors of near one for aromatic nitrogen compounds.

Although the activity of the column was not determined as a part of this study, it is safe to assume that basic nitrogen compounds would preferentially adhere to any acidic sites on the column, injector, or syringe. Furthermore, it is known that aromatic nitrogen compounds are weaker bases than non-aromatic nitrogen compounds due to a resonance effect (Reference 5). Table 4 lists some of the compounds studied, their relative response factors, and their base dissociation constants in




Figure 6. Measurement of Signal-to-Noise Ratio

water (Reference 6). Even though the compounds were not analyzed in water, the comparison of basicities to relative response factors show strong evidence that the non-aromatic compounds were being preferentially adsorbed.

The GC/CLD system investigated in this study appears to have considerable potential as a totally specific, highly quantitative GC nitrogen detector system.

It is recommended that the combustion tube be drastically reduced in size, and the flow rates be reduced accordingly. The combination of these design changes, the use of on-column injection, and proper deactivation techniques for the injection system and the column should yield a very sensitive, highly quantitative nitrogen specific detector for use with capillary gas chromatographs.

TABLE 4  
EFFECT OF BASICITIES ON RESPONSE FACTORS

<u>Compound Name</u>	<u>Response Factor</u>	<u>Kb</u>
Pyrrolidine	1.85333	$1.86 \times 10^{-3}$
Piperidine	1.29789	$1.33 \times 10^{-3}$
2,2,6,6-Tetra- methylpiperidine	1.48471	$1.18 \times 10^{-3}$
3,5 Dimethyl- pyridine	1.03570	$1.41 \times 10^{-8}$
4-Methyl- pyridine	1.02576	$1.05 \times 10^{-8}$
2-Methyl- pyridine	0.98190	$4.78 \times 10^{-9}$

# REFERENCES

1. H. V. Drushel, "Analytical Chemistry", Vol. 49, No. 7, 932-938, June 1977.
2. H. A. Hernandez, "American Laboratory", 72-76, June 1981.
3. R. C. Striebig, Nitrogen Analysis by Chemiluminescence, Aero Propulsion Laboratory, WPAFB, Ohio 45433, Report No. AFWAL-TR-84-2089, 1984.
4. Dues/Steward, RD Report Book, AFSC Form 79, No. 60016, AFWAL/POSF, WPAFB, OH 1984.
5. T. W. G. Solomons, "Organic Chemistry", 2nd Ed., 1980, Wiley and Sons.
6. R. C. Weast (Editor), C.R.C. Handbook of Chemistry and Physics, 60th Ed., 1979, CRC Press Inc.

END

Dtic

5-86