

MICROCOPY RESOLUTION TEST CHART NATIONAL BUREAU OF STANDARDS-1963-A



OFFICE OF NAVAL RESEARCH

K 1

Contract N00014-86-K-0659

Technical Report No. 3

SOLID STATE C-13 NMR CHARACTERIZATION OF TEXTILES: QUALITATIVE AND QUANTITATIVE ANALYSIS

by

Ronald F. Colletti and Lon J. Mathias

Prepared for Publication in

J. Appl. Polym. Sci., in press



Department of Polymer Science University of Southern Mississippi Southern Station Box 10076 Hattiesburg, MS 39406-0076

Reproduction in whole or in part is permitted for any purpose of the United States Government.

This document has been approved for public release and sale; its distribution is unlimited.

JECONITY CEASSIFICATION OF THIS PAGE		_ ــــــــــــــــــــــــــــــــــــ			
	REPORT DOCU	MENTATION	PAGE		
1a. REPORT SECURITY CLASSIFICATION		16. FRESTRICTIVE MARKINGS			
NONE		NONE			
2a. SECURITY CLASSIFICATION AUTHORITY NONE		3. DISTRIBUTION/AVAILABILITY OF REPORT			
2b. DECLASSIFICATION / DOWNGRADING SCHEDULE NONE		UNLIMITED			
4. PERFORMING ORGANIZATION REPORT NUMBER(S)		5. MONITORING ORGANIZATION REPORT NUMBER(S)			
Technical Report #3		ONR N00014-86-K-0659			
6a. NAME OF PERFORMING ORGANIZATION University of Southern	6b. OFFICE SYMBOL (If applicable)				
Mississippi		Office of Naval Research			
6c. ADDRESS (City, State, and ZIP Code)	7b. ADDRESS (City, State, and ZIP Code)				
University of Southern Missis	<b>[</b>				
6c ADDRESS (Gry. State. and ZIP Code) University of Southern Mississippi Polymer Science Department Southern Station Box 10076 Hattiesburg, MS 39406-0076		800 North Quincy Avenue			
Hattiesburg, MS 39406-0076		Arlington, VA 22217			
8a. NAME OF FUNDING/SPONSORING	9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER				
ORGANIZATION	(If applicable)				
Office of Naval Research	·				
8c. ADDRESS (City, State, and ZIP Code)		10. SOURCE OF FUNDING NUMBERS			
800 N. Quincy Avenue		PROGRAM	PROJECT	TASK	WORK L
Arlington , VA 22217		ELEMENT NO.	NO.	NO.	ACCESSI
		<u> </u>	<u></u>	<u> </u>	
11. TITLE (Include Security Classification)					
Solid State C-13 NMR Characte	erization of Tex	tiles: Qual	itative and	Quantitativ	e Analy
12. PERSONAL AUTHOR(S) Ronald F. Colletti and Lon J.	Mathian				
المتنبية والمستور والمتناق المتناسب والبراء والمستورة والمتنا والمتناق والم				- 1 10 2 2	
13a. TYPE OF REPORT 13b. TIME CO Technical FROM 8-		14. DATE OF REPO 87-8-27	RI (Year, Month,	Day) 15. PAGE	COUNT
16. SUPPLEMENTARY NOTATION		0, 0 1,			
10. JOFF ELMERTARY HOTATION					
J. Appl. Polym. Sci.					
17. COSATI CODES	18. SUBJECT TERMS (	Continue on reverse	if necessary and	d identify by bloc	k number)
FIELD GROUP SUB-GROUP					
19. ADS MACT (Continue on reverse if necessary	and identify by block r	number)			
cross-polarization mag	Ic engle enin	ning (CP/M	AS) carbon	13 NMR w	a g
Cross-polarization magic angle spinning (CP/MAS) carbon 13 NMR was used to qualitatively characterize textiles of various					
compositions. Rapid and accurate determination of					
cotton/polyester blends was accomplished with minimal sample					
	ative identi				
various fiber blends in				s techniqu	
demonstrates ability to characterize as-obtained fibers and					
fabrics and should ext					
prepregs.	•		•		
· · · · · · · · · · · · · · · · · · ·					
1 '11'					
1					
20. DISTRIBUTION/AVAILABILITY OF ABSTRACT		21. ABSTRACT SE	CURITY CLASSIFIC	ATION	
☑UNCLASSIFIED/UNLIMITED ☑ SAME AS RPT. ☐ DTIC USERS					
22a. NAME OF RESPONSIBLE INDIVIDUAL		226: TELEPHONE (	Include Area Code	) 22c. OFFICE SY	MBOL
Lon J. Mathias		(601)266-48	000		
DD FORM 1473, 84 MAR 83 AP	Redition may be used ur	itil exhausted.	SECURITY	CLASSIFICATION (	OF THIS PA

All other editions are obsolete.

# SOLID STATE C-13 NMR CHARACTERIZATION OF TEXTILES: QUALITATIVE AND QUANTITATIVE ANALYSIS

RONALD F. COLLETTI and LON J. MATHIAS\*

Department of Polymer Science

University of Southern Mississippi

Hattiesburg, MS 39406-0076

Synopsis

Cross-polarization magic angle spinning (CP/MAS) carbon 13

NMR was used to qualitatively and quantitatively characterize textiles of various compositions. Rapid and accurate determination of cotton/polyester blends was accomplished with minimal sample preparation. Qualitative identification was demonstrated for various fiber blends.

# QUALITY Dist Special A 1

## INTRODUCTION

Only within the past decade has the textile industry replaced wet chemical analysis for composition of fabrics and fibers with modern instrumental techniques. The advantages of

these new techniques are that they are quick, generally non-destructive, less operator intensive, do not require difficult sample handling procedures, and are almost always more precise and more accurate. Previous instrumental analysis in the textile industry has included X-ray diffraction, infrared spectroscopy (IR), and electron spectroscopy for chemical analysis (ESCA)<sup>1</sup>. Nuclear magnetic resonance spectroscopy (NMR) has also been utilized to study textiles<sup>2</sup>, and the advent of powerful solid state NMR spectrometer systems allows characterization with little or no preparation.

Specific NMR applications have included both qualitative and quantitative analyses. Proton broad line NMR was used to examine segment mobility in synthetic fibers<sup>3</sup> and the adsorption of water in oriented fibers<sup>4</sup>. Cellulose and regenerated cellulose were also studied with wide line proton NMR<sup>5</sup> as well as with solid state carbon-13 NMR<sup>6,7</sup>. Cotton/polyester blends have been examined quantitatively by IR<sup>1</sup>, X-ray diffraction<sup>8</sup>, density measurements<sup>9</sup> and ESCA<sup>10</sup>.

The technique presented here allows rapid, accurate, and non-destructive compositional determination of cotton/polyester blends and general qualitative identification of textile fibers.

#### EXPERIMENTAL

Textile samples were obtained from clothing and cut into approximately 5cm by 2cm pieces. The actual composition of the fabric was taken to be that on the garment label. Each sample was rolled and placed in a rotor and tamped for best packing and spinning. If a sample did not spin well, a small amount of fine silica gel was put into the spinner to help balance it.

The carbon-13 NMR spectra were recorded at 50.32 MHz on a Bruker MSL-200 Spectrometer equipped with a magic angle spinning probe and cross polarization circuitry to satisfy the Hartmann-Hahn condition. Spectra were acquired with a contact time of 5 ms, an acquisition time of 50 ms and a recycle time of 3 s. The 90 degree pulse was set at 5 us and the magic angle was set using potassium bromide. Between 1,000 and 3,000 FID's were collected for each sample. Adamantane was used to set the Hartmann-Hahn relationship and as the chemical shift reference. Fused alumina rotors with Kel-F caps were spun at 4.9 KHz. Total scanning time was between 30 - 60 m, with about 15 m necessary to pack the sample and start spinning. All measurements were taken at a temperature of 300 K.

#### RESULTS AND DISCUSSION

The quantitative aspects of the  $^{13}\mathrm{C}$  magic angle spinning NMR experiment have already been shown 11,12. The determination of the cotton/polyester content of commercial blends was accomplished by a peak height ratio method. The peaks at 104 and 162 ppm, for cotton (C-1 of cellulose) 13,14 and for the polyester carbonyl, respectively, were chosen since they are free of overlap with other peaks. Cotton/polyester blends examined included nominal ratios of 85:15, 75:25, 63:37, 60:40, 55:45, 53:47, 50:50 and 35:65. The ratio of the nominal cotton/polyester content was then plotted against the ratio of the cotton and polyester peak heights. Analysis by linear regression produced a straight line plot with a slope of 0.71, a y-intercept of -0.13 and a correlation coefficient of 0.993 (Figure 1). Eliminating the highest point and recalculating the straight line gave a slope of 0.67, y-intercept of -0.08 and a correlation coefficient of 0.963.

A plot of the ratio of the respective peak areas was also generated. Linear regression produced a slope of 0.72, y-intercept of 0.35 and a correlation coefficient of 0.969. Peak areas apparently do not correlate as well, perhaps due to differences in noise and baseline linearity. More careful setting of the integration regions and baseline, and a greater number of acquisitions to reduce noise may give more precise integration values. In this particular system, however, use of

4

peak heights gives excellent results, is much easier, and is more rapid.

The spectra of 100% cotton, 100% polyester, a 35:65 and a 50:50 cotton/polyester blend are all shown in Figure 2. The peak height ratio does not exactly equal the actual blend ratio possibly due to variations in cross-polarization efficiency between cotton and polyester, and certainly due to the differences in the ratios of crystalline to amorphous content for each polymer component.

In addition to the clothing samples used to establish the relationship in Figure 1, additional samples were obtained from a fabric store. Cotton/polyester blends with nominal ratios of 86:14, 84:16, 73:27, 72:28, and 20:80 were evaluated but did not correlate at all with the earlier data. In fact, one sample labeled as containing only 20% cotton was clearly almost all cotton; the cotton content calculated from measured peak heights corresponded to a ratio of 85:15. We are suggesting that these samples are mislabeled, either deliberately or through careless handling. In any event, fabric purchasers are not receiving the fabrics they think they are getting. It is also clear, especially from the non-zero intercept obtained from the peak height data, that careful determination of known sample blends would greatly improve the accuracy of this technique.

The utility of <sup>13</sup>C CP/MAS for qualitative analysis of textile fabrics was demonstrated for a number of clothing samples. These included fabrics nominally containing 100% rayon,

100% Qiana<sup>T.M.</sup>nylon, 100% shetland wool, and blends of 95:5 cotton/polyurethane, 50:50 cotton/Creslan<sup>T.M.</sup>acrylic, 82:12:6 cotton/acrylic/viscose rayon, 44:44:12 cotton/polyester/rayon, and 80:20 Arnel<sup>T.M.</sup>nylon. In all the spectra, peak assignments were made for each component of the blend by comparison with the spectra of individual fiber component and representative model compounds. Interestingly, it was found that the spectra of the garment labelled 80:20 Arnel<sup>T.M.</sup>/nylon did not contain peaks for nylon. Instead it showed peaks for polyester; compare Figure 3 and the polyester spectrum in Figure 2. The spectra of Qiana<sup>T.M.</sup> nylon and the 50:50 cotton/Creslan<sup>T.M.</sup>acrylic blend are shown in Figures 4 and 5, respectively.

These solid state carbon NMR examples demonstrate the value of this technique for rapid quantitative analysis of textile blends, and for quick, positive identification of fabrics and blends. This technique has great promise for routine screening of fabric samples to confirm the identity and composition of domestic and especially imported textile products.

### **ACKNOWLEDGEMENT**

The authors wish to gratefully acknowledge a grant from the Department of Defense to purchase our Bruker MSL-200 NMR spectrometer. This work was supported in part by the Office of Naval Research.

# REFERENCES

- 1. R. T. O'Conner, Instrumental Analysis of Cotton Cellulose and Modified Cotton Cellulose, Marcel Dekker, New York (1972).
- 2. B. Focher and A. Jus, Tinctoria, 72, 11 (1975).
- 3. W. O. Statton, Am. Dyestuff. Reptr., 54, 314 (1965).
- 4. R. E. Dehl, J. Chem. Phys., 48, 831 (1968).
- 5. R. A. Pittman and V. W. Tripp, J. Polym. Sci., Part A-2, 8, 969 (1970).
- 6. W. I. Earl and D. L. VanderHart, Macromolecules, 14, 570 (1981).
- 7. A. Hirai, F. Horii and R. Kitamaru, Bull. Inst. Chem. Res.,

  Kyoto Univ., 63, 340 (1985).
- 8. Te-Hsiung Wang, Hsin Hsien Wei, 20, 31 (1978).
- 9. Te-Hsiung Wang, Hsin Hsien Wei, 18, 18 (1976).
- 10. D. M. Soignet, Charact. Met. Polym. Surf., 2, 73 (1977).
- L. B. Alemany, D. M. Grant, R. J. Pugmire, T. D. Alger and K.
   W. Zilm, J. Am. Chem. Soc., 105, 2133 (1983).
- L. B. Alemany, D. M. Grant, R. J. Pugmire, T. D. Alger and K.
   W. Zilm, J. Am. Chem. Soc., 105, 2142 (1983).


- R. H. Atalla, J. C. Gast, D. J. Sindorf, V. T. Bartuska, G.
   E. Maciel, J. Am. Chem. Soc., 102, 3251 (1980).
- 14. W. I. Earl, D. L. VanderHart, J. Am. Chem. Soc., 102, 3251 (1980).


# LIST OF FIGURES


- Figure 1. Plot of peak height vs. cotton/polyester blend ratio.
- Figure 2. Carbon CP/MAS spectrum of A: 100% cotton; B: 50:50 and C: 35:65 cotton/polyester blend; and D: 100% polyester.
- Figure 3. Carbon CP/MAS spectrum of a nominal 80:20 Arnel // nylon blend.


y<u>e</u>


- Figure 4. Carbon CP/MAS spectrum of Qiana nylon.
- Figure 5. Carbon CP/MAS spectrum of 50:50 cotton/Creslan<sup>T.M.</sup> acrylic.











0-87 / /