



DTIC FILE COPY

2

AD-A223 183

Laboratory Note No. 90-79

## Combined Primate Transport and Restraint Device


S.B. Reynolds,  
H. Zwick  
and  
D.O. Robbins

DIVISION OF OCULAR HAZARDS RESEARCH

### DISTRIBUTION STATEMENT A

Approved for public release  
Distribution Unlimited

April 1990



LETTERMAN ARMY INSTITUTE OF RESEARCH PRESIDIO OF SAN FRANCISCO CALIFORNIA 94129


This document has been approved for public release and sale; its distribution is unlimited.

Destroy this report when it is no longer needed. Do not return to the originator.

Citation of trade names in this report does not constitute an official endorsement or approval of the use of such items.

The experimental studies of the author described in this report were reviewed and approved by the Institutional Review Committee/Animal Care and Use Committee at Letterman Army Institute of Research. The Manuscript was peer reviewed for compliance prior to submission for publication. In conducting the research described here, the author adhered to the "Guide for the Care and Use of Laboratory Animals," DHEW Publication (NIH) 85-23.

This material has been reviewed by Letterman Army Institute of Research and there is no objection to its presentation and/or publication. The opinions or assertions contained herein are the private views of the author(s) and are not to be construed as official or as reflecting the views of the Department of the Army or the Department of Defense. (AR 360-5)

 20 Apr '90  
DONALD G. CORBY (date)  
COL, MC  
Commander

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE

## REPORT DOCUMENTATION PAGE

Form Approved  
OMB No. 0704-0188

1a. REPORT SECURITY CLASSIFICATION UNCLASSIFIED			1b. RESTRICTIVE MARKINGS		
2a. SECURITY CLASSIFICATION AUTHORITY			3. DISTRIBUTION / AVAILABILITY OF REPORT  UNLIMITED DISTRIBUTION		
2b. DECLASSIFICATION / DOWNGRADING SCHEDULE					
4. PERFORMING ORGANIZATION REPORT NUMBER(S)  Laboratory Note - 90-97			5. MONITORING ORGANIZATION REPORT NUMBER(S)		
6a. NAME OF PERFORMING ORGANIZATION Letterman Army Institute of Research		6b. OFFICE SYMBOL (If applicable) SGRD-ULY-OH	7a. NAME OF MONITORING ORGANIZATION US Army Medical Research and Development Command		
6c. ADDRESS (City, State, and ZIP Code) Letterman Army Institute of Research Division of Ocular Hazards Presidio of San Francisco, CA 94129-6800			7b. ADDRESS (City, State, and ZIP Code) Fort Detrick Frederick, MD 21701-5012		
8a. NAME OF FUNDING / SPONSORING ORGANIZATION		8b. OFFICE SYMBOL (If applicable)	9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER		
8c. ADDRESS (City, State, and ZIP Code)			10. SOURCE OF FUNDING NUMBERS		
			PROGRAM ELEMENT NO.	PROJECT NO. 3E162 777A878✓	TASK NO.
					WORK UNIT ACCESSION NO. 161
11. TITLE (Include Security Classification)  Combined Primate Transport and Restraint Device					
12. PERSONAL AUTHOR(S) Scottie B. Reynolds, Harry Zwick and Dave O. Robbins					
13a. TYPE OF REPORT Laboratory Note		13b. TIME COVERED FROM Mar 89 TO Jun 89		14. DATE OF REPORT (Year, Month, Day) 1990, March	
15. PAGE COUNT 9					
16. SUPPLEMENTARY NOTATION					
17. COSATI CODES			18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number)		
FIELD	GROUP	SUB-GROUP			
			Primate, Device, Transport, Restraint, Animal, Laboratory Animals		
19. ABSTRACT (Continue on reverse if necessary and identify by block number)  This paper describes a practical, inexpensive alternative to the standard poling and chairing method of transferring awake primates. Employing this device, a technician working alone can retrieve, restrain and position an awake primate without endangering either the animal or the technician. The device, basically a Plexiglas box, is placed against an animal's home cage. The animal is trained to avoid the squeeze mechanism of its home cage by entering the box. We describe modifications made to the box to adapt the standard one-ring animal collar and to accommodate animals of various sizes. We also present a top panel that can be angularly adjusted to allow for the animal's normal head carriage, thus avoiding the strain of holding the chin up for long periods of time.					
20. DISTRIBUTION / AVAILABILITY OF ABSTRACT <input checked="" type="checkbox"/> UNCLASSIFIED/UNLIMITED <input type="checkbox"/> SAME AS RPT. <input type="checkbox"/> DTIC USERS			21. ABSTRACT SECURITY CLASSIFICATION UNCLASSIFIED		
22a. NAME OF RESPONSIBLE INDIVIDUAL Donald G. Corby, COL, MC Commanding			22b. TELEPHONE (Include Area Code) (415) 561-6300		22c. OFFICE SYMBOL SGRD -ULY-Z

# ABSTRACT

This paper describes a practical, inexpensive alternative to the standard poling and chairing method of transferring awake primates. Employing this device, a technician working alone can retrieve, restrain and position an awake primate without endangering either the animal or the technician. The device, basically a Plexiglas box, is placed against an animal's home cage. The animal is trained to avoid the squeeze mechanism of its home cage by entering the box. We describe modifications made to the box to adapt the standard one-ring animal collar and to accommodate animals of various sizes. We also present a top panel that can be angularly adjusted to allow for the animal's normal head carriage, thus avoiding the strain of holding the chin up for long periods of time.



Accession For	
NTIS CRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution /	
Availability Codes	
Dist	Avail and/or Special
A-1	

## **COMBINED PRIMATE TRANSPORT AND RESTRAINT DEVICE**

### **INTRODUCTION**

The visual performance and retinal anatomy of humans and rhesus monkeys are remarkably similar; therefore, the rhesus monkey is an excellent human prototype for ocular studies that require long-term, daily periods of visual, behavioral, and physiological response evaluation. The major drawback in the regular use of nonhuman primates is that they are unpredictable and can be very difficult and sometimes dangerous to handle.

In studies that require daily periods of physiological experimental manipulation by awake, task-oriented rhesus monkeys (1), a technique is needed to transfer the animals cautiously and expeditiously from their home cage to a transport and restraint device, thereby making the daily administration of anesthesia unnecessary. In the past a poling method was used to transfer awake primates. The poling method required that a permanent, lightweight, plastic collar be fitted around an anesthetized animal's neck. The collar had a metal fastening ring that could be hooked with a spring-loaded pole while the animal was in its home cage. By maneuvering the pole and employing positive reinforcement, an animal could be conditioned to go from its home cage to a standard temporary restraining chair with little difficulty.

The poling technique works well with average-sized, young monkeys. Older, larger animals that have been pole-trained can become difficult to transfer and may require two handlers. For older, larger animals that have not been previously pole-trained, initiation of pole training is nearly impossible. In this instance, the likelihood of harm to the laboratory technician as well as the animal is high. Therefore, we have explored several alternative techniques for transporting unanesthetized primates to experimental settings.

### **APPARATUS**

A new, simple, and inexpensive alternative to pole transferring and the standard primate chair was recent -

ly developed and described (2). This new transfer apparatus is designed to allow easy removal of fully conscious monkeys from their home cages to the capture box. The transfer procedure requires a minimum of cooperation from the animal, uses no form of anesthesia, and significantly decreases the possibility of escape or injury to the animal or handler. The top panel of the box provides head restraints for experimental tasks which require that the animal's visual line of fixation remain constant. The capture box also serves as the transport device to and from the laboratory.

The apparatus, basically a portable Plexiglas box with a vertically sliding door on the front side, is mounted onto a portable hydraulic lift platform and aligned with the vertically sliding door of a standard, nonhuman primate cage. With both vertical doors open, the animal is then trained to avoid the squeeze mechanism of the home cage by entering the Plexiglas box. Once the animal is inside the Plexiglas box, two aluminum rods are used to hook the animal's two-ring collar and draw it up against the top panel of the box. The animal's head protrudes through a 14 cm diameter hole cut in the center of the panel. The rods are then secured with base plates and the collar is locked into position, thus positioning the animal by its collar and maintaining the correct line of fixation and visual pathway to the viewing screen. Employing this practical restraint and transport device, a technician working alone can retrieve, restrain, and position an awake rhesus.

#### **MODIFICATIONS**

This paper describes modifications made to the Plexiglas box (1) to capture and hold the standard one-ring animal collar, (2) to accommodate a variety of animal sizes, and (3) to allow the top panel to angularly adjust to the animal's normal angle of view. The modified Plexiglas box is shown in Figure 1 with overall dimensions that can easily accommodate the larger-sized rhesus monkeys.




Figure 1. Diagram of the Plexiglas restraint device. The box is constructed with 3/8" Plexiglas all around and 3/8" diameter aluminum rods form the floor.

The front panel of the box has a 26 cm wide center opening that runs the full height of the box. This allows the animals maximum height clearance while being transferred to and from their home cage, and provides more room for the task manipulation keys and levers in the behavioral task procedure. Before transport, the front opening is closed by inserting the removable, vertical panel into slots located on both sides of the opening (Figure 2).




Figure 2. Removable Front Panel

The top panel of the Plexiglas box is adjustable for the animal's height and angle of view (Figure 3), and is constructed with a center opening 12 cm wide. Thin strips of Plexiglas (1.5 cm X .5 cm) are used to construct slots for the removable plates and animal collar. All Plexiglas strips used to construct the slots have smooth, rounded edges. The slot to hold the animal's restraint collar is formed by overlapping (1.5 cm) both sides of the opening with strips that are attached to the bottom of the top panel with screws. These strips are 23 cm long and leave a 10.5 cm opening in front for pulling up and capturing the animal's collar. The top of this slot is formed by two strips

that run the full length of the opening and overlap it 1.5cm on each side. Both strips have a 1 cm indentation at the collar capture opening which provides an 11 cm clearance to allow passage of only the animal's head. By constructing the opening in this fashion, the commercially available single-ring, Plexiglas animal collar can be utilized. The front end of this slot is also used to hold a plate mounted with a drinking spout. The spout delivers 1 ml squirts of fruit juice as positive reinforcement during the behavior task procedure. The topmost slot is formed with two strips identical to the strips forming the bottom of the animal's collar slot. This slot holds the removable, sliding top cover plate that seals the center opening during capture and transport (Figure 4).




Figure 3. The Plexiglas adjustable, restraining top is diagramed from three different views. It is constructed with a center neck opening 12 cm wide x 28 cm long, and slotted for the removable plates (juicer plate, animal collar, top cover plate).

When the animal is first moved from its home cage to the restraint device, the top opening of the box is completely covered by the top sliding cover and the juicer plate (Figure 4a). Initially, the tightening bolt on the juicer plate is loosened and the plate is pulled out of the center opening only far enough to secure the latching pole to the animal's collar ring. The juicer plate is then removed as the ring and collar are pulled up and forward into the juicer slot (Figure 4b). To eliminate any possibility of escape, only the animal's head is allowed to pass through the indented front opening (for much larger animals the back cover plate can be opened slightly). After the collar is in the juicer slot, the tightening bolt on the cover plate is loosened. By pressing down on the collar ring with the pole and pushing back, the collar and animal are moved back into the collar restraint slot. The juicer plate is replaced in the slot in front of the collar ring and tightened; the collar is secured (Figure 4c) and the pole is removed from the ring.




Figure 4. Center opening of the top panel (enlargement in Figure 1). (A) Covered with the top sliding cover 1 and juicer plate 2. (B) Animal collar 3 replaces juicer plate. (C) Collar pushed back into collar restraint slot and juicer plate replaced.

The angularly adjustable top has screw-in tightening knobs fitted through the vertical slots in both side panels. This allows the height to be varied by 14 cm and the animal's angle of view to be adjusted from 0 to 15 degrees, thereby relieving the animal of possible head and neck discomfort (Figure 5). The vertical adjustment can be checked from the see-through sides to assure that a comfortably crouched animal's shoulders barely touch the restraining top panel. There is 8.5 cm of space between the back of the restrained collar and the back of the box. When the animal is crouching, this space allows for the natural curvature of the spine. Each animal should be checked for comfortable positioning and for the absence of neck stretching or body cramming.




FIGURE 5. Side Panel with Slots for Tightening Knobs

We have tested this apparatus on two mature rhesus monkeys. Both animals wear one-ring collars and were previously pole transferred, but have now become too large to pole transfer. The procedure and apparatus described above were successful in allowing both animals to be transferred from their home cage and utilized in experimental behavioral situations.

**REFERENCES**

1. Robbins DO, Zwick H, and Holt GC. A method for producing foveal retinal exposures in an awake, task-oriented rhesus monkey. Behav Res Method and Instr 1973;5(6):457-461.
2. Robbins DO, Zwick H, Leedy M, Stearns G. Acute restraining device for rhesus monkeys. Lab Anim Sci 1986;36(1):68-70.

OFFICIAL DISTRIBUTION LIST

Commander  
USAMRDC  
ATTN:SGRD-PLC/MG Russell  
ATTN:SGRD-PLC/COL Lam  
SGRD-RMS/Ms. Madigan  
Fort Detrick, MD 21701-5000

Director  
Defense Technical Information  
ATTN:DTIC-DDA (2 copies)  
Cameron Station  
Alexandria, VA 22314

Commander  
US Army SMO  
ATTN:AMXCM-EO/MAJ Dedmond  
2800 Powder Mill Road  
Adelphi, MD 20783

Commander  
USAMSAA  
ATTN: DRXSY-CSD/P. Baers  
ATTN: DRXST-GWD/F. Campbell  
Aberdeen Proving Ground  
Maryland 21010

Commander  
ATTN: AFWAL/MLPJ/G. Kepple  
Wright Patterson AFB  
Ohio 45433

Commander  
US AEHA  
ATTN: HSHB-RL/D. Sliney  
Aberdeen Proving Ground  
Maryland 21010

Dr. John Ewen  
PO Box 1925  
Washington, DC 20013

Commander  
HQ, USAMMDA  
ATTN: SGRD-UMA/Channing  
Fort Detrick  
Frederick, MD 21701-5012

Commander  
ATTN:AFAMRL/HEF/R. Susnik  
Wright Patterson AFB  
Ohio 45433

Headquarters  
Department of the Army  
ATTN: DASG-TLO  
Washington, DC 20310

Commander  
CACDA/ATZL-OPS-SE  
ATTN: MAJ J.C. Johnson  
Fort Leavenworth  
Kansas 66027

Director  
NADC  
ATTN: Code 601B/Dr. Chisum  
Warminster, PA 18974-5000

Commander  
NMRDC  
ATTN: Code 43  
National Naval Med Center  
Bethesda, MD 20814

Commander  
USAF SAM  
ATTN: RZW/Dr. Farrer  
ATTN: RZW/LTC Cartledge  
Brooks AFB, Texas 78235

Official Distribution List, cont.

Director  
AMSAA  
ATTN: AMXSY-CR/Mr. Brand  
Aberdeen Proving Ground  
Maryland 21005

Commander  
USA Aviation Systems Command  
ATTN: AMCPM-ALSE/H. Lee  
4300 Goodfellow Blvd  
St. Louis, MO 63120

Director  
Defense Intelligence Agency  
ATTN: DT-5A/Hal Hock  
Pentagon  
Washington, DC 20301

Commander  
USA Aeromedical Research Lab  
ATTN: COL La Mothe  
Ft. Rucker, AL 36330-5000

Director  
EWL/RSTA Center  
ATTN: AMSEL-EW-C/J. Allen  
Ft Monmouth, NJ 07703-5303

Director  
USA HEL  
ATTN: AMXHE-IS/D. Giordano  
Aberdeen Proving Ground  
Maryland 21005-5001

Director  
US Army AMMRC  
ATTN: AMXMR-O/Fitzpatrick  
Watertown, MA 02172-0001

Director  
Armed Forces Medical  
Intelligence Center  
ATTN: AFMIC-SA/MAJ Downs  
Fort Detrick  
Maryland 21701-5004

Director  
DTD Directorate  
ATTN: EOGWCM-CCM/Kasperek  
White Sands Missile Range  
New Mexico 88002-5519

Commander  
HQ TRADOC  
ATTN: ATCD-ML/J. Gray  
Ft. Monroe, VA 23651-5000

Commander  
LAIR  
ATTN: SGRD-ULZ (1 copy)  
SGRD-IR (10 copies)  
PSF, CA 94129-6800