

Technica! Report 917

# Application of Adaptive Nulling to Electromagnetic Hyperthermia for Improved Thermal Dose Distribution in Cancer Therapy



A.J. Fenn

3 July 1991

# **Lincoln Laboratory**

MASSACHUSETTS INSTITUTE OF TECHNOLOGY

LEXINGTON, MASSACHUSETTS



Prepared for the Department of the Air Force under Contract F19628-90-C-0002.

Approved for public release: distribution is unlimited.

91-11490

This report is based on studies performed at Lincoln Laboratory, a center for research operated by Massachusetts Institute of Technology. The work was sponsored by the Department of the Air Force under Contract F19628-90-C-0002.

This report may be reproduced to satisfy needs of U.S. Government agencies.

The ESD Public Affairs Office has reviewed this report, and it is releasable to the National Technical Information Service, where it will be available to the general public, including foreign nationals.

This technical report has been reviewed and is approved for publication.

FOR THE COMMANDER

がある。 のでは、 のでは

というと というという というとう かんとうしょうしょ こうしょうしょ

Hugh L. Southall, Lt. Col., USAF

Chief. ESD Lincoln Laboratory Project Office

Hugh L. Southall

Non-Lincoln Recipients

PLEASE DO NOT RETURN

Permission is given to destroy this document when it is no longer needed.

# 5 September 1991

# **ERRATA**

Document: Technical Report 917

Application of Adaptive Nulling to Electromagnetic Hyperthermia

for Improved Thermal Dose Distribution in Cancer Therapy

3 July 1991

Errors in the production process require the following corrections:

The figure on page 59 is Figure A-6, and belongs with the caption on page 61.

The figure on page 60 is Figure A-7, and belongs with the caption on page 62.

The figure on page 61 is Figure A-4, and belongs with the caption on page 59.

The figure on page 62 is Figure A-5, and belongs with the caption on page 60.

Publications Office
MIT Lincoln Laboratory
PO Box 73
Lexington, MA 02173-9108

# MASSACHUSETTS INSTITUTE OF TECHNOLOGY LINCOLN LABORATORY

# APPLICATION OF ADAPTIVE NULLING TO ELECTROMAGNETIC HYPERTHERMIA FOR IMPROVED THERMAL DOSE DISTRIBUTION IN CANCER THERAPY

A.J. FENN Group 61

**TECHNICAL REPORT 917** 

3 JULY 1991

Approved for public release; distribution is unlimited.


LEXINGTON MASSACHUSETTS

#### ABSTRACT

Adaptive nulling is applied to the problem of generating a therapeutic thermal dose distribution in electromagnetic hyperthermia treatment of cancer. A system design concept for implementing adaptive hyperthermia is introduced. With the proposed design concept, it may be possible to maximize the applied electric field at a tumor position in the target body and simultaneously minimize or reduce the electric field at target positions where undesired high temperature regions (hot spots) occur. In a clinical situation, either a gradient search algorithm or sample matrix inversion algorithm would be used to rapidly form the adaptive null (or nulls) prior to any significant tissue heating. Auxiliary short-dipole field probes are used in effecting the desired electric field nulls. An adaptive null formed at the surface of the target has a finite width and extends into the target. This finite null width can allow for noninvasively positioned auxiliary probes in an adaptive hyperthermia system. The allowed spacing or resolution between a deep null and focus is fundamentally equal to the hyperthermia antenna half-power beamwidth in the tissue. A closer spacing between the desired null position and focus is achieved by reducing the null depth via the signal-to-noise ratio at the auxiliary probe.

Analysis of an annular phased-array antenna embedded in an infinite homogeneous medium shows the potential merit of combining adaptive nulling with conventional near-field focusing used in hyperthermia. The analysis is based on a well-known moment-method theory for conducting thin-wire antennas in a homogeneous conducting medium. The theory and software used to compute the moment-method received voltage at a short-dipole probe due to a transmitting dipole array are documented. Computer simulation results are presented for a fully adaptive eight-element annular array operating at 120 MHz and embedded in homogeneous muscle tissue. Multiple simultaneous nulls are used in adaptively generating a desired radiation pattern in the muscle tissue. The power received at a movable short-dipole E-field sensor is then used as the power source in a second computer simulation that performs a transient thermal analysis of an elliptical target body surrounded by a constant-temperature water bolus. The computer simulations show that adaptive nulling can prevent undesired high-temperature regions from occurring while

simultaneously heating a deep-seated tumor site.



# TABLE OF CONTENTS

	Abstract	iii
	List of Illustrations	vii
	Acknowledgments	viii
1.	INTRODUCTION	1
2.	THEORY	7
	2.1 Noninvasive Adaptive Hyperthermia System Concept	7
	2.2 Adaptive Transmit Array Formulation	12
	2.3 Moment-Method Formulation	16
3.	SIMULATION RESULTS	27
	3.1 Electric Field for Array in Homogeneous Tissue	28
	3.2 Temperature Distribution in Elliptical Phantom	40
4.	CONCLUSION	53
AF	PPENDIX A - SYSTEM DEGRADATION DUE TO INSUFFICIENT NUMBER OF	
ΑŪ	UXILIARY PROBES	55
ΑF	PPENDIX B - SOFTWARE DOCUMENTATION	63
RE	EFERENCES	99

# LIST OF ILLUSTRATIONS

Figure No.		Page
1	A clinical hyperthermia phased-array antenna system.	2
2	General concept for focused near-field adaptive nulling.	4
3	Two applications of focused near-field adaptive nulling: (a) Radar system testing and (b) Hyperthermia treatment of cancer.	5
4	Noninvasive adaptive hyperthermia system concept.	8
5	(a) Transverse cross section at the prostate level. An eight-element hyperthermia ring array of dipole elements is used to irradiate a target body with sufficient power to elevate the temperature of the target body to a therapeutic level. (b) Elliptical phantom target used to model the prostate (target)-level cross section. Shown are four noninvasive auxiliary E-field probes that are used in forming null zones with widths controlled by the individual null depths. The null zones extend into the body to reduce the internal electric field at specified locations.	9
6	Simulated two-dimensional thermal pattern at time $t=20$ minutes in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\text{C}$ constant-temperature water bolus. The incident RF power distribution is at 120 MHz and the initial temperature of the phantom is $25^{\circ}\text{C}$ . Temperature contour levels are given in $2^{\circ}\text{C}$ steps. (a) Before nulling; hot spots on the left and right sides of the target are present. (b) After nulling; hot spots are eliminated.	10
7	Adaptive transmit phased-array antenna near-field focusing concept.	11
8	Implementation of the SMI algorithm for an adaptive hyperthermia system. Four performance measures are used in evaluating the adaptive response of the system: electric-field distribution, covariance matrix eigenvalues, adaptive transmit weights, and interference cancellation.	13
9	Transmit array and receive antenna probe. Open-circuit mutual impedance between array elements is denoted $Z_{m,n}$ . The open-circuit voltage at the probe is computed from the array terminal currents and from $Z_n^j$ , the open-circuit mutual impedance between the <i>n</i> th array element and the <i>j</i> th probe antenna.	17
10	Equivalent circuit model for receive antenna probe.	17
<b>4</b> U	Equitation officer model for receive amenia proce.	- '

# LIST OF ILLUSTRATIONS (Continued)

Figure No.		Page
11	Equivalent circuit model used in thermal analysis. The inhomogeneous body of interest is modeled with a grid of nodes interconnected by resistors. Each node has an associated power source and capacitor.	24
12	The moment-method simulation's use with the transient thermal analysis. RF power distribution over the target region is converted into a time-dependent temperature distribution.	26
13	Geometry for eight-element ring array and four E-field auxiliary sensors.	30
14	Simulated two-dimensional quiescent radiation pattern at 120 MHz for uniformly illuminated eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon_r = 73.5, \sigma = 0.5$ ). A fictitious elliptical target region is indicated by the dashed curve. Radiation contour levels are in 10-dB steps.	31
15	Simulated two-dimensional quiescent radiation pattern at 120 MHz for uniformly illuminated eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon_{\tau} = 73.5$ , $\sigma = 0.5$ ). A fictitious elliptical target region is indicated by the dashed curve. Radiation contour levels are given in 1-dB steps.	32
16	Simulated one-dimensional quiescent radiation pattern cut ( $z=0$ ) at 120 MHz for uniformly illuminated eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon_{\tau}=73.5, \sigma=0.5$ ).	33
17	Simulated one-dimensional quiescent radiation pattern cut $(x=0)$ at 120 MHz for uniformly illuminated eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon_{\tau}=73.5, \sigma=0.5$ ).	34
18	Simulated two-dimensional adaptive radiation pattern at 120 MHz for eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon_{\tau}=73.5, \sigma=0.5$ ). Radiation contour levels are given in 10-dB steps. Four auxiliary sensors are used in forming the adaptive pattern. The quiescent focus is at $(0.0)$ .	36
19	Simulated one-dimensional quiescent and adaptive radiation patterns in the $z=0$ cut at 120 MHz for eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon_{\tau}=73.5, \sigma=0.5$ ). Four auxiliary sensors are used in forming the adaptive pattern.	37

# LIST OF ILLUSTRATIONS (Continued)

Figure No.		Page
20	Simulated one-dimensional quiescent and adaptive radiation patterns in the $x=0$ cut at 120 MHz for eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon_r=73.5, \sigma=0.5$ ). Four auxiliary sensors are used in forming the adaptive pattern.	<b>3</b> 8
21	Transmit array weights before and after adaptive nulling. A dynamic range of about 5 dB is evident for the adaptive weights. Four auxiliary sensors are used in the adaptive process. (a) Amplitude; (b) Phase.	39
22	Channel covariance matrix eigenvalues (degrees of freedom) used in the adaptive process with four auxiliary sensors.	40
23	Simulated two-dimensional thermal pattern at time $t=20$ minutes before nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\text{C}$ constant-temperature water bolus. The incident RF power distribution, from Figure 14. is at 120 MHz. Temperature contour levels are given in $2^{\circ}\text{C}$ steps. Hot spots on the left and right sides of the target are observed.	42
24	Simulated one-dimensional ( $z=0$ ) thermal pattern at time $t=20$ minutes before nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\mathrm{C}$ constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. Hot spots on the left and right sides of the target are observed.	43
25	Simulated one-dimensional $(x = 0)$ thermal pattern at time $t = 20$ minutes before nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\text{C}$ constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. No hot spots are present.	44
26	Simulated two-dimensional thermal pattern at time $t=20$ minutes (with adaptive nulling at four auxiliary sensors in effect) in elliptical phantom muscle-tissue target surrounded with $10^{\circ}$ C constant-temperature water bolus. The adapted incident RF power distribution, from Figure 18, is at 120 MHz. Temperature contour levels are given in $2^{\circ}$ C steps. Hot spots on the left and right sides of the target are eliminated.	45

# LIST OF ILLUSTRATIONS (Continued)

Figure No.		Page
27	Simulated one-dimensional $(z=0)$ thermal patterns at time $t=20$ minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}$ C constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. Hot spots on the left and right sides of the target are clearly eliminated by the adaptive nulling process. Four auxiliary sensors are used in the adaptive process.	46
28	Simulated one-dimensional $(x=0)$ thermal patterns at time $t=20$ minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\text{C}$ constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. No undesired hot spots are present. Four auxiliary sensors are used in the adaptive process.	47
29	E-field probe-sample-spacing convergence check for simulated two-dimensional thermal pattern at time $t=20$ minutes (before adaptive nulling) in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\text{C}$ constant-temperature water bolus. The quiescent incident RF power distribution is at $120 \text{ MHz}$ . Temperature contour levels are given in $2^{\circ}\text{C}$ steps. The grid spacing is one-half the spacing of that in Figure 23. Hot spots on the left and right sides of the target are present, as previously observed for the coarser probe-sample spacing.	48
30	E field probe-sample-spacing convergence check for simulated two-dimensional thermal pattern at time $t=20$ minutes (with adaptive nulling in effect) in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\text{C}$ constant-temperature water bolus. The adapted incident RF power distribution is at 120 MHz. Temperature contour levels are given in $2^{\circ}\text{C}$ steps. The grid spacing is one-half the spacing of that in Figure 26. Hot spots on the left and right sides of the target are eliminated, as previously observed for the	
	coarser probe-sample spacing.	49

# LIST OF ILLUSTRATIONS

(Continued)

	Page
E-field probe-sample-spacing convergence check for simulated one-dimensional $(z=0)$ thermal patterns at time $t=20$ minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}C$ constant-temperature water bolus. The quiescent incident RF power distribution is at 120 MHz with 0.9525-cm sample spacing. Hot spots on the left and right sides of the target are eliminated by the adaptive nulling process, as previously observed in Figure 27. Four auxiliary sensors are used in the adaptive process.	50
E-field probe-sample-spacing convergence check for simulated one-dimensional $(x=0)$ thermal patterns at time $t=20$ minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\text{C}$ constant-temperature water bolus. The quiescent incident RF power distribution is at 120 MHz with 0.9525-cm sample spacing. No undesired hot spots are present. Four auxiliary sensors are used in the adaptive process.	51
Geometry for eight-element ring array and two E-field auxiliary sensors.	55
Simulated two-dimensional adaptive radiation pattern at 120 MHz for eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon_{\tau}=73.5, \sigma=0.5$ ). Radiation contour levels are given in 10-dB steps. Two auxiliary sensors are used in forming the adaptive pattern. The quiescent focus is at $(0.9)$ .	57
Simulated two-dimensional thermal pattern at time $t=20$ minutes (with adaptive nulling at only two auxiliary sensors in effect) in elliptical phantom muscle-tissue target surrounded with $10^{\circ}$ C constant-temperature water bolus. The adapted incident RF power distribution, from Figure 18, is at 120 MHz. Temperature contour levels are given in $2^{\circ}$ C steps. Initial hot spots (before nulling) on the left and right sides of the target (see Figure 23) are redistributed to the top and bottom (anterior and posterior positions).	58
Simulated one-dimensional $(z=0)$ thermal patterns at time $t=20$ minutes before and after nulling in elliptical phantom muscle-tissuc target surrounded with $10^{\circ}$ C constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. Hot spots on the left and right sides of the target are eliminated by the adaptive nulling process. Two auxiliary sensors are used in the adaptive process.	59
	( $z=0$ ) thermal patterns at time $t=20$ minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\mathrm{C}$ constant-temperature water bolus. The quiescent incident RF power distribution is at 120 MHz with 0.9525-cm sample spacing. Hot spots on the left and right sides of the target are eliminated by the adaptive nulling process, as previously observed in Figure 27. Four auxiliary sensors are used in the adaptive process.  E-field probe-sample-spacing convergence check for simulated one-dimensional ( $x=0$ ) thermal patterns at time $t=20$ minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\mathrm{C}$ constant-temperature water bolus. The quiescent incident RF power distribution is at 120 MHz with 0.9525-cm sample spacing. No undesired hot spots are present. Four auxiliary sensors are used in the adaptive process.  Geometry for eight-element ring array and two E-field auxiliary sensors.  Simulated two-dimensional adaptive radiation pattern at 120 MHz for eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue: $\epsilon=73.5$ , $\sigma=0.5$ ). Radiation contour levels are given in 10-dB steps. Two auxiliary sensors are used in forming the adaptive pattern. The quiescent focus is at (0.0).  Simulated two-dimensional thermal pattern at time $t=20$ minutes (with adaptive nulling at only two auxiliary sensors in effect) in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\mathrm{C}$ constant-temperature water bolus. The adapted incident RF power distribution, from Figure 18, is at 120 MHz. Temperature contour levels are given in $2^{\circ}\mathrm{C}$ steps. Initial hot spots (before nulling) on the left and right sides of the target (see Figure 23) are redistributed to the top and bottom (anterior and posterior positions). Simulated one-dimensional ( $z=0$ ) thermal patterns at time $t=20$ minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}\mathrm{C}$ constant-temperature water bol

# LIST OF ILLUSTRATIONS

(Continued)

Figure No.		Page
A-5	Simulated one-dimensional $(x = 0)$ thermal patterns at time $t = 20$ minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with $10^{\circ}$ C constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. Two undesired hot spots are present in the adaptive pattern. Two auxiliary sensors are used in the adaptive process.	60
A-6	Transmit array weights before and after adaptive nulling. A dynamic range of about 7 dB is evident for the adaptive weights. Two auxiliary sensors are used in the adaptive process. (a) Amplitude: (b) Phase.	61
A-7	Channel covariance matrix eigenvalues (degrees of freedom) used in the adaptive process with two auxiliary sensors.	62

# **ACKNOWLEDGMENTS**


The author wishes to express his gratitude to K.P. Lawton for software support and to D.H. Temme, R.P. Rafuse, D.M. Nathanson, and J.R. Johnson for technical discussions. The encouragement of G.N. Tsandoulas is sincerely appreciated. Technical discussions with numerous individuals outside Lincoln Laboratory are gratefully acknowledged.

#### 1. INTRODUCTION

The successful treatment of deep-seated malignant tumors within a patient is often a difficult task. The objective of the treatment is to reduce in size or completely remove the tumor mass by one or more modalities available at the treatment facility. Common modalities are surgery, chemotherapy, and x-ray therapy [1]. A modality used alone or in conjunction with one of the above modalities is "tissue heating," or hyperthermia [1,2,3,4]. Hyperthermia can be considered as a form of high fever within the body; a controlled thermal dose distribution is required for hyperthermia to have a therapeutic value. Typical localized-hyperthermia temperatures required for therapeutic treatment of cancer are in the 43-45°C range. Normal tissue should be kept at temperatures below 43°C during the treatment. The most difficult aspect of implementing hyperthermia, with either radio-frequency (RF) waves or acoustic (ultrasound) waves, is producing sufficient heating at depth. Multiple-applicator RF hyperthermia arrays are commonly used to provide a focused nearfield main beam at the tumor position. A focal region should be concentrated at the tumor with minimal energy delivered to surrounding normal tissue. As the hyperthermia antenna beamwidth is proportional to the wavelength, a small focal region suggests that the RF wavelength be as small as possible. However, due to propagation losses in tissue, the RF depth of penetration decreases with increasing transmit frequency. One of the major problems in heating a deep-seated tumor with a hyperthermia antenna is the formation of undesired "hot spots" in surrounding tissue. This additional undesired heating often produces pain, burns, and blistering in the patient, which requires terminating the treatment immediately. The patient does not receive anesthetics during the hyperthermia treatment in order to provide direct verbal feedback of any pain. Thus, techniques for reducing hot spots are necessary in hyperthermia treatment.

Several journals have published special issues on acoustic and electromagnetic hyperthermia treatment of cancer [5,6]. Electromagnetic transmitting arrays in the frequency band 60–2000 MHz are used to localize heating of malignant tumors within a target body. Many studies have been conducted to produce improved therapeutic field distributions with hyperthermia phased arrays [7–20]. Phase control can be used to synthesize therapeutic RF radiation patterns without adaptive control of the transmit array weights [7–12]. Array transmit weights can be adaptively controlled to maximize the tumor temperature (or RF power delivered to the tumor) while minimizing the surrounding tissue temperature (or RF power delivered to the surrounding tissue) [13–20]. All other studies require invasive techniques to optimize the radiation pattern [13–20]. This report addresses the potential benefit of using adaptive nulling with noninvasive auxiliary sensors to reduce the field intensity at selected positions in the target body while maintaining a desired focus at a tumor.

It is useful to describe certain features of a clinical hyperthermia system used to treat deep-seated tumors. A hyperthermia annular phased-array antenna system (Model BSD-2000, SIGMA-60 applicator [21], BSD Medical Corporation, Salt Lake City, Utah) is shown in Figure 1. By fully surrounding the patient with an annular phased array, it is possible to obtain constructive interference (or signal enhancement) deep within the target volume. This hyperthermia system uses a 60-cm array diameter with eight uniformly spaced dipole elements operating over the frequency



# DIPOLE PAIR

BSD MEDICAL CORP. MODEL BSD-2000		
OPERATING FREQ.	60–120 MHz	
NO. CHANNELS	4	
POWER	500 W/CHANNEL	
ANTENNA	8-DIPOLE RING	

Figure 1. A clinical hyperthermia phased-array antenna system.

band 60-120 MHz [22]. The eight dipoles are fed as four active pairs of elements. There are four high-power amplifiers which drive the dipole pairs with up to 500 W average power per channel. Each of the four active channels has an electronically controlled variable-phase shifter for focusing the array. Temperature and electric-field probe sensors (both invasive and noninvasive) are used to monitor the treatment. A cool-water (5-40°C) bolus between the patient and the phased array is used to prevent excess heating of the skin surface. The bolus is filled with circulating distilled water, which has a very low propagation loss. Computer simulation models for this hyperthermia antenna system with inhomogeneous targets have been developed [7,9]. To the author's knowledge, the hyperthermia system in Figure 1 has not been used to form adaptive nulls, although the system may have sufficient temperature and electric-field sensing capability to implement adaptive nulling with an appropriate transmit-weight control algorithm. A candidate adaptive nulling algorithm would be a gradient search based on minimizing the signal received by electric-field sensors at the desired null positions while maximizing the signal at the focal point. The reason for considering this algorithm is that the BSD-2000 system is currently capable of measuring the applied electricfield amplitude but not the phase. Note: the E-field detectors used in the BSD-2000 system are semiconductor diode detectors. A potentially better approach to performing adaptive hyperthermia is to use a channel-covariance-matrix-based algorithm to control the weights. If the electric-field phase and amplitude are measured (by modification of the E-field probe hardware to form in-phase and quadrature signals) the sample matrix inversion (SMI) algorithm [23] could be used to form the nulls. The concept of adaptive nulling as it applies to hyperthermia will now be described.

Recently, the author has been investigating a technique called focused near-field adaptive nulling [24-30], for which the general concept is shown in Figure 2. A calibration probe antenna is used to focus the main beam at approximately one diameter D of the phased array. Additional (auxiliary) probe antennas are located in the sidelobe region of the phased-array quiescent radiation pattern. Adaptive nulls are formed in the direction of these auxiliary probes. Computer simulations show that a radar system can be tested by using the focused near-field nulling technique [29]. Experimental measurements of focused near-field adaptive nulling are available [30]. During the radar system testing study, similarities between the radar application and electromagnetic hyperthermia were noticed by the author. These similarities are that 1) phased-array near-field focusing is used to form the main beam and 2) uncontrolled or high sidelobes are potentially deleterious to system performance. For the radar system, interference signals entering through uncontrolled sidelobes can reduce the signal-to-noise ratio (SNR); adaptive nulling [31] is commonly used to counteract this possible degradation. In the hyperthermia system, high-transmit antenna sidelobes can give rise to hot spots in the target tissue. These high-temperature regions could possibly be alleviated with adaptive nulling of the electric-field sidelobes in the target. Much of the software developed for the radar system testing analysis may be readily modified to handle adaptive hyperthermia analysis.

There is continued interest at Lincoln Laboratory in applying techniques used in Department of Defense (DoD)-sponsored research to other disciplines, such as medical applications. This study applies radar system testing technology (currently in the form of software simulation) to the medical application of hyperthermia. Figure 3 depicts the similarities between radar system near-field




Figure 2. General concept for focused near-field adaptive nulling.

testing and electromagnetic hyperthermia treatment of cancer. In Figure 3(a), a radar phased-array antenna system is positioned in an anechoic test chamber with radiating sources that simulate typical radar signals. Notice that the test antenna is focused in the near field at one aperture diameter. Computer simulations have shown that the near-field signals (jammer, radar clutter, and target) received and processed by an adaptive radar system are equivalent to fielded radar conditions [29]. In Figure 3(b), a hyperthermia (transmit) phased array is used to therapeutically illuminate a target body by focusing its electromagnetic energy on a malignant tumor deep within the body. Nulls are formed adaptively to reduce the electromagnetic energy delivered to potential hot spots. The nulls are shown invasive to the target; however, noninvasive null positions on the surface of the target are also possible as shown in this report.

An example of a deep-seated tumor is cancer of the prostate [1]. The tumor volume often has a decreased blood flow which aids in heating the tumor, compared to normal tissue for which heat is carried away by normal blood flow. Note: a blood-flow thermal analysis [32,33] is beyond the scope of this report; we would like to show initially that with a computer simulation of RF-illuminated phantom muscle tissue [34] adaptive nulling can improve the thermal dose distribution. The thermal analysis used in this report is similar to the method presented in Zhang, et al. [35], where blood flow effects are ignored. In practice, undesired high-temperature regions away from the focus can occur inside the target volume. For example, scar tissue, which has a decreased blood flow rate, will tend to heat up more rapidly than normal tissue. In the proposed adaptive hyperthermia concept, electric-field nulls are used to reduce the power delivered to potential hot

spots. Computer simulations show that noninvasive field probes can be used to eliminate hot spots interior to the phantom target tissue.




Figure 3. Two applications of focused near-field adaptive nulling: (a) Radar system testing and (b) Hyperthermia treatment of cancer.

An antenna analysis code (called WIRES) capable of analyzing a wide variety of antenna or radar cross section problems has been developed by J.H. Richmond [36–38]. The software can analyze complex geometries involving multiple-connected and/or isolated wires and is well-suited to analyzing low-frequency phased arrays radiating into an infinite homogeneous conducting medium. WIRES is a moment-method code that uses the electric field integral equation (EFIE) to enforce the boundary condition of the tangential electric field being zero at the surface of the antenna of interest. The moment-method basis and testing functions used in this code are piecewise sinusoidal. For a complete description of the theory and capabilities of the WIRES code, the reader is referred to the user's manual [37]. The WIRES software has been modified by the author to include near-field focusing and adaptive nulling capabilities. This modified code has been used to analyze the

near-field and far-field adaptive nulling performance of thin-wire phased arrays in free space [27,29]. A new version of the thin-wire code that can analyze adaptive hyperthermia arrays in an infinite homogeneous conducting medium has recently been written and is documented in this report.

This report is organized in the following manner. Section 2 discusses the concept for a noninvasive adaptive hyperthermia system and describes in detail the theory used in performing the focused near-field adaptive nulling simulation. Pertinent theory for wave propagation in a conducting medium is reviewed. A brief description of the transient thermal analysis approach used to compute the temperature distribution in a target is given. Section 3 presents simulation results for a fully adaptive eight-element hyperthermia ring array operating at 120 MHz. The calculated received RF-power distributions, before and after adaptive nulling, at a short-dipole field probe with both the field probe and the hyperthermia array embedded in homogeneous muscle tissue, are presented. The RF power distribution is then fed into transient thermal analysis software which computes the two-dimensional temperature distribution in an elliptical muscle-tissue target surrounded with a constant-temperature water bolus. The temperature distribution obtained by using multiple independent adaptive nulls clearly indicates the elimination of undesired hot spots. Conclusions are made in Section 4. Performance degradation due to insufficient number of auxiliary null positions is described in Appendix A. Appendix B contains data files and source code for the focused near-field adaptive nulling computer simulations.

#### 2. THEORY

## 2.1 NONINVASIVE ADAPTIVE HYPERTHERMIA SYSTEM CONCEPT

The concept of a noninvasive adaptive-nulling hyperthermia system [39,40] is shown in Figure 4. Theoretically, to generate the desired field distribution in a clinical adaptive hyperthermia system, receiving sensors are positioned as close as possible to the focus (tumor site) and to where high temperatures are to be avoided (such as near the spinal cord and scar tissue). As simulations in this report will show, for an annular array configuration the receiving sensors can be located noninvasively on the surface (skin) of the target. Initially, the hyperthermia array is focused to produce the required field intensity at the tumor. An invasive probe is required to achieve the optimum focus at depth. To avoid undesired hot spots, it is necessary to minimize the power received at the desired null positions and to constrain the array weights to deliver a required amount of transmitted or focal region power. The adaptive array weights (with gain g and phase g) are controlled by either the SMI algorithm or a gradient search algorithm to rapidly (within seconds) form the nulls before a significant amount of target heating takes place. With this adaptive technique, it should be possible to avoid unintentional hot spots and maintain a therapeutic thermal dose distribution at the tumor.

Figure 5(a) shows a specific example of an eight-element hyperthermia ring phased array with a target cross section at the prostate level [9]. Figure 5(b) shows an elliptical phantom target which is used to model the prostate-level cross section. A noninvasive adaptive nulling system is achieved by placing auxiliary sensors  $1.2...., N_{aux}$  on the target skin as shown. The null zones centered at each auxiliary probe naturally extend into the elliptical target region to eliminate undesired hot spots. (An example of a two-dimensional near-field null zone is shown in Figure 10(b) of Fenn [28]). The width of each null zone is directly related to the depth of each null (sometimes referred to as the amount of cancellation) is directly related to the SNR at the sensor position. A low SNR produces a small amount of nulling, and a high SNR produces a large amount of nulling. The resolution or minimum spacing between the focus and null position is normally equal to the half-power beamwidth of the antenna. The resolution is enhanced somewhat by using weak nulls whenever the separation between the null and focus is closer than the half-power beamwidth.

Figure 6 summarizes the thermal simulation for a hyperthermia ring array transmitting into a homogeneous elliptical target region. The tumor site is assumed to be at the center of the ellipse. The thermal distribution in Figure 6(a) contains two undesired hot spots to the left and right of the focus. In Figure 6(b), adaptive nulling is applied and the hot spots are eliminated. The details of the computer simulations are contained in Section 3.

In the near-field nulling technique described here, it is assumed that the hyperthermia phasedarray antenna is focused (as it normally is) in the near field and that a main beam and possibly sidelobes are formed in the target. In this report, it is assumed that phase focusing is used to produce the desired quiescent main beam. Figure 7 shows a primary (calibration) probe antenna located at a desired focal point of an array. The array can maximize the signal received by the




Figure 4. Noninvasive adaptive hyperthermia system concept.




Figure 5. (a) Transverse cross section at the prostate level. An eight-element hyperthermia ring array of dipole elements is used to irradiate a target body with sufficient power to elevate the temperature of the target body to a therapeutic level. (b) Elliptical phantom target used to model the prostate (target)-level cross section. Shown are four noninvasive auxiliary E-field probes that are used in forming null zones with widths controlled by the individual null depths. The null zones extend into the body to reduce the internal electric field at specified locations.




Figure 6. Simulated two-dimensional thermal pattern at time t=20 minutes in elliptical phantom muscle-tissue target surrounded with  $10^{\circ}$  C constant-temperature water bolus. The incident RF power distribution is at 120 MHz and the initial temperature of the phantom is 25° C. Temperature contour levels are given in 2° C steps. (a) Before nulling; hot spots on the left and right sides of the target are present. (b) After nulling; hot spots are eliminated.

calibration probe by adjusting its phase shifters so that the observed element-to-element phase variation is removed. One way to do this in a numerical simulation is to choose a reference path length as the distance from the focal point to the phase center of the array. This distance is denoted  $r_F$  and the distance from the focal point to the *n*th array element is denoted  $r_n^F$ . The voltage received at a probe (located at the focal point) due to the *n*th array element is computed here using the method of moments as described in Section 2.3. To maximize the received voltage at the focal-point probe output, it is necessary to apply the phase conjugate of the observed signals, due to the array elements, at the transmit array. The resulting near-field radiation pattern will have a main beam and sidelobes. The main beam will be pointed at the array focal point, and sidelobes will exist at angles away from the main beam. Auxiliary probe antennas can then be placed at the desired null positions in the quiescent sidelobe region. These sidelobes are where tissue hot spots are likely to occur; they are nulled by the adaptive nulling algorithms described below.




Figure 7. Adaptive transmit phased-array antenna near-field focusing concept.

#### 2.2 ADAPTIVE TRANSMIT ARRAY FORMULATION

Consider the hyperthermia array and probe geometry shown in Figure 7. Typically, the hyperthermia array contains N identical antenna elements. The input signal to each of the N array elements is obtained from the weighted signal distributed by the power divider. The number of adaptive channels is denoted M, and for a fully adaptive array M=N. Note: the report will concentrate on fully adaptive arrays (arrays in which all transmit array elements have adaptive weight control); however, the computer simulation program is developed for both fully adaptive and sidelobe canceller arrays. In this report, ideal transmit weights (a complex voltage gain vector) are assumed in the computer simulation, with  $\mathbf{w}=(w_1,w_2,\cdots,w_N)^T$  denoting the adaptive channel weight vector as shown in Figure 4. (Superscript T means transpose). To generate adaptive nulls, the transmit weights (phase and gain) are controlled by either the SMI algorithm or a gradient search algorithm. The SMI algorithm has the flexibility to operate in either open- or closed-loop feedback modes 41; the gradient search algorithm operates only in a feedback mode.

## 2.2.1 Sample Matrix Inversion Algorithm

With the SMI algorithm [23], the fundamental quantities required to fully characterize the incident field for adaptive nulling purposes are the adaptive channel cross correlations. To implement this algorithm it is necessary to know the complex received voltage at the auxiliary probes. For example, the moment-method formulation (described in Section 2.3) allows computation of the complex-received voltage at the auxiliary probes.

The block diagram in Figure 8 shows how the SMI algorithm is used in the hyperthermia analysis. The diagram indicates four performance measures used to quantify the computer simulations: electric-field distribution, covariance matrix eigenvalues, adaptive transmit weights, and interference cancellation. The calculation of these performance measures is described in detail below.

Referencing Figure 7, let a spheric. wavefront be incident at the *i*th probe antenna, due to each array element (radiating one at a time with a unity-amplitude reference signal), which results in a set of probe-received complex voltages denoted  $v_1^i, v_2^i, \dots, v_N^i$ . The cross correlation  $R_{mn}^i$  of the received voltages due to the *m*th and *n*th adaptive transmit channels at the *i*th probe is given by

$$R_{mn}^{t} = \boldsymbol{E}(v_{m}v_{n}^{\bullet}) \tag{1}$$

<sup>&</sup>lt;sup>1</sup>The effects of transmit-weight quantization and transmit-weight random errors are addressed in the software.

 $V_1$   $V_2$  • • •  $V_N$  TRANSMIT ARRAY




Figure 8. Implementation of the SMI algorithm for an adaptive hyperthermia system. Four performance measures are used in evaluating the adaptive response of the system: electric-field distribution, covariance matrix eigenvalues, adaptive transmit weights, and interference cancellation.

where \* means complex conjugate and  $E(\cdot)$  means mathematical expectation. Because  $v_m$  and  $v_n$  represent voltages of the same waveform but at different times,  $R_{mn}^i$  is also referred to as an autocorrelation function. Note: for convenience, in Equation (1) the superscript i in  $v_m$  and in  $v_n$  has been omitted.

In the frequency domain, assuming the transmit waveform has a band-limited white noise power spectral density (as commonly assumed in jamming of a radar system), Equation (1) can be expressed as the frequency average

$$R_{mn}^{i} = \frac{1}{B} \int_{f_{1}}^{f_{2}} v_{m}(f) v_{n}^{*}(f) df \qquad , \tag{2}$$

where  $B = f_2 - f_1$  is the nulling bandwidth and f is the frequency. It should be noted that  $v_m(f)$  takes into account the transmit wavefront shape, which is spherical for the hyperthermia application. For the special case of a continuous wave (CW) transmit waveform, as normally used in hyperthermia, the cross correlation reduces to

$$R_{mn}^i = v_m(f_o)v_n^*(f_o) \qquad , \tag{3}$$

where  $f_o$  is the transmit frequency of the hyperthermia array.

Let the channel or interference covariance matrix be denoted R. (Note: in hyperthermia, interference is used to refer to the signals received at the auxiliary probes. The undesired hot spots can be thought of as interfering with the therapy.) If there are  $J_{aux}$  independent desired null positions or auxiliary probes, the  $J_{aux}$ -probe covariance matrix is the sum of the covariance matrices observed at the individual probes; that is,

$$R = \sum_{i=1}^{J_{aux}} R_i + I \qquad . \tag{4}$$

where  $R_{\tau}$  is the sample covariance matrix observed at the *i*th probe and I is the identity matrix used to represent the thermal noise level of the receiver. (It should be noted that both  $N_{aux}$  and  $J_{zux}$  can be used interchangeably to describe the number of auxiliary nulling sensors or received-interference signal positions. The symbol  $J_{aux}$  is perhaps more descriptive and can equivalently be thought of as the number of positions where jamming signals are nulled).

Prior to generating an adaptive null, the adaptive channel weight vector,  $\boldsymbol{w}$ , is chosen to synthesize a desired quiescent radiation pattern. When nulling is desired, the optimum set of transmit weights to form an adaptive null (or nulls), denoted  $\boldsymbol{w}_a$ , is computed [23] by

$$\boldsymbol{w}_a = \boldsymbol{R}^{-1} \boldsymbol{w}_q \qquad . \tag{5}$$

where  $^{-1}$  means inverse and  $\boldsymbol{w}_q$  is the quiescent weight vector. During array calibration, the normalized quiescent transmit weight vector, with element 1 radiating, is chosen to be  $\boldsymbol{w}_q = (1,0,0,\cdots,0)^T$ ; that is, the transmit channel weight of element 1 is unity and the remaining transmit channel weights are zero. Similar weight settings are used to calibrate the remaining transmit elements. For a fully adaptive annular array focused at the origin in homogeneous tissue, the normalized quiescent weight vector is simply  $\boldsymbol{w}_q = (1,1,1,\cdots,1)^T$ . Commonly, the weight vector is constrained to deliver a required amount of power to the hyperthermia array or to the tumor. For simplicity in the software used to analyze the hyperthermia array, the weights are constrained such that

$$\sum_{m=1}^{M} |w_m|^2 = 1 \qquad , {(6)}$$

where  $w_n$  is the transmit weight for the *n*th element. Note: in the computer simulations, the electric field due to the normalized weight vector is scaled appropriately to deliver the required amount of power to the tissue so that a desired focal-region temperature level is achieved after t minutes.

The summation of power received at the probes is given by

$$p = \boldsymbol{w}^{\dagger} \boldsymbol{R} \boldsymbol{w} \tag{7}$$

where  $^{\dagger}$  means complex conjugate transpose. The interference-plus-noise-to-noise ratio, denoted INR, is computed as the ratio of the auxiliary probe array output power [defined in Equation (7)] with the transmit signal present to the probe array output power with only receiver noise present; that is,

$$INR = \frac{\boldsymbol{w}^{\dagger} \boldsymbol{R} \boldsymbol{w}^{\dagger}}{\boldsymbol{w}^{\dagger} \boldsymbol{w}} \qquad . \tag{8}$$

The adaptive array cancellation ratio, denoted C, is defined here as the ratio of the summation of probe-received power after adaptation to the summation of probe-received power before adaptation; that is.

$$C=\frac{p_a}{p_a} \qquad .$$

Substituting Equation (7) in Equation (9) yields

$$C = \frac{\boldsymbol{w_a}^{\dagger} \boldsymbol{R} \boldsymbol{w_a}}{\boldsymbol{w_q}^{\dagger} \boldsymbol{R} \boldsymbol{w_q}} \qquad . \tag{10}$$

Next, the covariance matrix defined by the elements in Equations (2) or (3) is Hermitian (that is,  $\mathbf{R} = \mathbf{R}^{\dagger}$ ), which, by the spectral theorem, can be decomposed in eigenspace [42] as

$$\mathbf{R} = \sum_{k=1}^{M} \lambda_k \mathbf{e}_k \mathbf{e}_k^{\dagger} \qquad , \tag{11}$$

where  $\lambda_k, k = 1, 2, \dots, M$  are the eigenvalues of R, and  $e_k, k = 1, 2, \dots, M$  are the associated eigenvectors of R. The interference covariance matrix eigenvalues  $(\lambda_1, \lambda_2, \dots, \lambda_M)$  are a convenient quantitative measure of the use of the adaptive array degrees of freedom [43]. The amplitude spread between the largest and smallest eigenvalues is a quantitative measure of the dynamic range of the interference (hot spot) signals.

## 2.2.2 Gradient Search Algorithm

Under conditions where only the probe-received voltage amplitude is measured (as in the BSD-2000 system previously mentioned), it is appropriate to consider a gradient search algorithm to minimize the interference power at selected positions. The gradient search is used to iteratively control the transmit weights so that the RF signal received by the probe array is minimized. The transmit array weights (gain and phase) are adaptively changed in small increments and the probe array output power is monitored to determine weight settings that reduce the output power most rapidly to a null. The mathematical formulation for the gradient search is developed in a straightforward manner [44] and will not be considered further here. Only the sample covariance matrix inversion algorithm has been implemented in the focused near-field nulling simulation presented in this report.

#### 2.3 MOMENT-METHOD FORMULATION

This section describes a method of moments formulation to compute the probe-received voltages in Equation (2) due to the transmitting hyperthermia phased-array antenna in an infinite homogeneous conducting medium. The medium is described by the three parameters  $\mu$ ,  $\epsilon$ , and  $\sigma$ , which are discussed in the next section. The formulation given here is analogous to that developed under array-receiving conditions for an adaptive radar [29]. The software used to analyze a hyperthermia array is based on the receive-array analogy but the theory presented below is given in the context of a transmit array. Referring to Figure 9, assume that each transmit element is fed with a generator having a known impedance  $Z_L$ . Let  $v_{n,j}^{o.c.}$  represent the open-circuit voltage at the jth probe due to the nth transmit-array element. Here, the jth probe can denote either the focal point calibration probe or one of the auxiliary probes used to null a sidelobe. The number of auxiliary probes is denoted by  $J_{aux}$ . The receive probe is assumed to be terminated in an impedance  $Z_r$ .




Figure 9. Transmit array and receive antenna probe. Open-circuit mutual impedance between array elements is denoted  $Z_{m,n}$ . The open-circuit voltage at the probe is computed from the array terminal currents and from  $Z_n^j$ , the open-circuit mutual impedance between the nth array element and the jth probe antenna.




Figure 10. Equivalent circuit model for receive antenna probe.

Next, let Z denote the open-circuit mutual impedance matrix (with dimensions  $N \times N$  for the N-element array). The open-circuit mutual impedance between array elements m and n is denoted  $Z_{m,n}$ . It is assumed that multiple interaction between the hyperthermia array and the probe sensor can be neglected. Thus, the hyperthermia array terminal current vector i can be computed in terms of the transmit weights w as

$$\boldsymbol{i} = [\boldsymbol{Z} + \boldsymbol{Z}_L \boldsymbol{I}]^{-1} \boldsymbol{w} \qquad . \tag{12}$$

Next, let  $Z_n^j$  be the open-circuit mutual impedance between the *j*th probe and the *n*th array element. The induced open-circuit voltage  $v_{n,j}^{o.c.}$  at the *j*th receive probe, due to the *n*th array element transmit current  $i_n$ , can then be expressed as

$$v_{n,i}^{o.c.} = Z_n^j \cdot i_n \qquad . \tag{13}$$

In matrix form, the induced open-circuit probe-voltage matrix  $v_{probe}^{o.c.}$  is

$$\boldsymbol{v_{probe}^{o.c.}} = \boldsymbol{Z_{probe,array}} \boldsymbol{i} \tag{14}$$

or

$$\boldsymbol{v}_{probe}^{o.c.} = \boldsymbol{Z}_{probe,array} [\boldsymbol{Z} + Z_L \boldsymbol{I}]^{-1} \boldsymbol{w} , \qquad (15)$$

where  $Z_{probe,array}$  is a rectangular matrix of order  $J_{aux} \times N$  for the open-circuit mutual impedance between the probe array and the hyperthermia array. Note: the jth row of the matrix  $Z_{probe,array}$  is written as  $(Z_1^j, Z_2^j, \dots, Z_N^j)$ , where  $j = 1, 2, \dots, J_{aux}$ . The receive voltage matrix is then computed by the receiving circuit equivalence theorem for an antenna [50]. The receive-antenna equivalent circuit is depicted in Figure 10, where it is readily determined that

$$v_{probe}^{rec} = v_{probe}^{o.c.} \frac{Z_r}{Z_{in} + Z_r} \qquad , \tag{16}$$

where  $Z_{in}$  is the input impedance of the probe. It should be noted that the  $v_{probe}^{rec}$  matrix is a column vector of length  $J_{aux}$  and  $v_j^{rec}$  is the jth element of the matrix. The probe-receive current matrix is given by

$$i_{prob\epsilon}^{rec} = v_{prob\epsilon}^{o.c.} \frac{1}{Z_{in} + Z_r} \tag{17}$$

The jth element of the column vector  $i_{probe}^{rec}$  is denoted  $i_j^{rec}$ ,  $j=1,2,\cdots,J_{aux}$ . Finally, the power received by the jth probe is

$$p_j^{rec} = \frac{1}{2} Re(v_j^{rec} \cdot i_j^{rec^*}) \qquad , \tag{18}$$

where Re means real part. Substituting Equations (16) and (17) into Equation (18) yields

$$p_j^{rec} = \frac{1}{2} |v_j^{o.c.}|^2 \frac{Re(Z_r)}{|Z_{in} + Z_r|^2} \qquad . \tag{19}$$

The total interference power received by the auxiliary probe array is given by

$$p^{rec} = \sum_{j=1}^{J_{aux}} p_j^{rec} \qquad .$$

The incident electric field E is related to the open-circuit voltage  $v^{o.c.}$  by the effective height h of the probe antenna [48] as

$$v^{o.c.} = hE (21)$$

If the length  $L_p$  of the probe antenna is approximately  $0.1\lambda$  or less, the current distribution is triangular and the effective height is  $h = 0.5L_p$  [48]. Thus, for a short-dipole probe the open-circuit voltage can be expressed as

$$v^{o.c.} = \frac{L_p}{2}E \qquad . ag{22}$$

It then follows from Equation (22) that the E field for a short-dipole probe at position (x, y, z) is given by

$$E(x, y, z) = \frac{2v^{o.c.}(x, y, z)}{L_p} . (23)$$

Finally, the quiescent and adapted E-field radiation patterns are computed using the quiescent and adapted weight vectors  $w_q$  and  $w_a$ , respectively, in Equations (15) and (23).

The moment-method expansion and testing functions are assumed to be sinusoidal. The opencircuit mutual impedances in Equation (15) between thin-wire dipoles in a homogeneous conducting medium are computed based on subroutines from a well-known moment-method computer code [38]. In evaluating  $Z_n^j$  for the jth auxiliary probe, double precision computations are used.

As mentioned previously, the array is calibrated (phased focused) initially using a short dipole at the focal point. To accomplish this numerically, having computed  $v_{focus}^{rec}$ , the receive array weight vector  $\boldsymbol{w}$  will have its phase commands set equal to the conjugate of the corresponding phases in  $v_{focus}^{rec}$ . Transmit antenna radiation patterns are obtained by scanning (moving) a dipole probe with half-length l in the near-field and computing the receive probe-voltage response.

The received voltage matrix for the jth probe (denoted  $\boldsymbol{v}_{j}^{rec}$ ) is computed at K frequencies across the nulling bandwidth. Thus,  $\boldsymbol{v}_{j}^{rec}(f_{1}), \boldsymbol{v}_{j}^{rec}(f_{2}), \cdots, \boldsymbol{v}_{j}^{rec}(f_{K})$  are needed. In this report, the impedance matrix is computed at K frequencies and is inverted K times. The probe covariance matrix elements are computed by evaluating Equation (2) numerically, using Simpson's rule numerical integration. For multiple auxiliary probes, the covariance matrix is evaluated using Equation (3). Adaptive array radiation patterns are computed by superimposing the quiescent radiation pattern with the weighted sum of auxiliary-channel-received voltages.

## 2.3.1 WAVE PROPAGATION IN CONDUCTING MEDIA

To gain insight into the effect of a lossy medium on the propagation of an electromagnetic wave, it is useful to review certain fundamental equations which govern the field characteristics. In a conducting medium, Maxwell's curl equations in time-harmonic form are

$$\nabla \times \boldsymbol{H} = \boldsymbol{J} + j\omega \epsilon \boldsymbol{E} \tag{24}$$

$$\nabla \times \boldsymbol{E} = -j\omega \mu \boldsymbol{H} \qquad , \tag{25}$$

where E and H are the electric and magnetic fields, respectively. J is the conduction current density,  $\omega = 2\pi f$  is the radian frequency,  $\epsilon$  is the permittivity of the medium, and  $\mu$  is the permeability of the medium. The permittivity is expressed as  $\epsilon = \epsilon_{\tau} \epsilon_{o}$ , where  $\epsilon_{\tau}$  is the dielectric constant (relative permittivity) and  $\epsilon_{o}$  is the permittivity of free space. Similarly,  $\mu = \mu_{\tau} \mu_{o}$ , where  $\mu_{\tau}$  is the relative permeability and  $\mu_{o}$  is the permeability of free space. For a medium with electrical conductivity  $\sigma$ , J and E are related as

$$\boldsymbol{J} = \sigma \boldsymbol{E} \qquad . \tag{26}$$

Substituting Equation (26) into Equation (24) yields

$$\nabla \times \boldsymbol{H} = (\sigma + j\omega \epsilon)\boldsymbol{E} \qquad . \tag{27}$$

From Equations (24) and (25), the vector wave equation in terms of E is derived as

$$\nabla^2 E - \gamma^2 E = 0 (28)$$

It is readily shown [49] that

$$\gamma = \pm \sqrt{j\omega\mu(\sigma + j\omega\epsilon)} = \pm j\omega\sqrt{\mu\epsilon}\sqrt{1 - j\frac{\sigma}{\omega\epsilon}} \qquad . \tag{29}$$

The quantity  $\sigma/\omega\epsilon$  is referred to as the loss tangent. It is common to express the complex propagation constant as

$$\gamma = \alpha + j\beta \qquad , \tag{30}$$

where  $\alpha$  is the attenuation constant and  $\beta$  is the phase constant. The constants  $\alpha$  and  $\beta$  are found by setting Equation (29) equal to Equation (30) and then squaring both sides, equating the real and imaginary parts, and solving the pair of simultaneous equations, with the result

$$\alpha = \frac{\omega\sqrt{\mu\epsilon}}{\sqrt{2}} \left\{ \sqrt{1 + (\frac{\sigma}{\omega\epsilon})^2} - 1 \right\}^{1/2} \tag{31}$$

$$\beta = \frac{\omega\sqrt{\mu\epsilon}}{\sqrt{2}} \left\{ \sqrt{1 + (\frac{\sigma}{\omega\epsilon})^2} + 1 \right\}^{1/2} \qquad . \tag{32}$$

The wavelength  $\lambda$  in the lossy dielectric is then computed from

$$\lambda = \frac{2\pi}{3} \qquad . {33}$$

The intrinsic wave impedance  $\eta$  is given by

$$\eta = \sqrt{\frac{j\omega\mu}{\sigma + j\omega\epsilon}} = \sqrt{\frac{\mu}{\epsilon}} \frac{1}{\sqrt{1 - j\frac{\sigma}{\omega\epsilon}}} \qquad (34)$$

The instantaneous power density of the electromagnetic field is given by Poynting's vector, denoted P.

$$P = \frac{1}{2}E \times H^* \qquad , \tag{35}$$

which has units of  $(W/m^2)$ . The time-average power flow density is equal to the real part of the complex Poynting's vector. The time-average power dissipation per unit volume  $P_d$   $(W/m^3)$  is derived from Maxwell's equations, with the result

$$P_d = \frac{1}{2} \boldsymbol{E} \cdot \boldsymbol{J}^* = \frac{1}{2} \sigma |\boldsymbol{E}|^2 \qquad . \tag{36}$$

The specific absorption rate (SAR) [2] is the power dissipated or absorbed per unit mass (W/kg) of the medium (tissue), or

$$SAR = \frac{P_d}{\rho} = \frac{\sigma}{2\rho} |\mathbf{E}|^2 \qquad , \tag{37}$$

where  $\rho$  is the density of the medium in kg/m<sup>3</sup>.

It is convenient to have a simple equation for computing the propagation loss between any two points in the near field of an isolated transmitting antenna. Thus, mutual coupling effects are ignored for the time being. Consider a time-harmonic source radiating a spherical wave into an infinite homogeneous conducting medium. For an isotropic radiator, and suppressing the  $e^{j\omega t}$  time dependence, the electric field as a function of range r can be expressed as

$$E(r) = E_o \frac{e^{-\gamma r}}{r} \qquad , {38}$$

where  $E_o$  is a constant.

For a source at the origin, the amplitude of the electric field at range  $r_1$  is given by

$$|E(r_1)| = E_o \frac{e^{-\alpha r_1}}{r_1} \qquad , {39}$$

and at range  $r_2$  by

$$|E(r_2)| = E_o \frac{e^{-\alpha r_2}}{r_2} . (40)$$

The total propagation loss between ranges  $r_1$  and  $r_2$  is found by taking the ratio of Equations (40) and (39), or

$$\frac{|E(r_2)|}{|E(r_1)|} = \frac{r_1}{r_2} e^{-\alpha(r_2 - r_1)} \qquad . \tag{41}$$

The field attenuation  $A_{\alpha}$  in dB from range  $r_1$  to range  $r_2$  due to the lossy dielectric is simply

$$A_{\alpha} = 20 \log_{10}(e^{-\alpha(r_2 - r_1)}) \tag{42}$$

Similarly, the 1/r attenuation loss  $A_r$  in dB is

$$A_r = 20\log_{10}\frac{r_1}{r_2} . (43)$$

#### 2.3.2 THERMAL MODELING OF AN INHOMOGENEOUS TARGET

A thermal analysis computer program called the transient thermal analyzer (TTA), developed by Arthur D. Little, Inc., has been used at Lincoln Laboratory for space payload analysis since 1969. Recent documentation for TTA exists [51,52]. The important details of the analysis are described in this section. The TTA software is used in the next section to accomplish the thermal modeling of homogeneous muscle tissue surrounded by a constant-temperature water bolus. Thermal modeling techniques have been described in detail [53].

The TTA program uses the finite-difference technique to solve a set of nonlinear energy balance equations. Consider a system of interconnected nodes that model an inhomogeneous volume for which the temperature  $T_i$  of the *i*th node is to be determined. The heat-balance equation, which is solved by TTA, is expressed as


$$\sum_{i=1}^{N} Q_{i,j} - P_i(t) + M_i \frac{dT_i}{dt} = 0 , (44)$$

where  $Q_{i,j}$  is the net outward heat flow from node i in the direction of node j,  $P_i(t)$  is the power into node i at time t, and  $M_i$  is the thermal mass (mass times specific heat) of node i. Figure 11 shows an electric circuit analog which is used to model the two-dimensional thermal characteristics of the material volume. Power  $P_i$  in watts is delivered to the ith node. Capacitor  $C_i$  (with units Joules i0) is used to model the thermal capacitance at the ith node. Resistor  $R_{i,j}$  (with units i0) is used to model the heat resistance between nodes i and j.

With a spacing of  $\Delta l$  between nodes (assuming cubic cells), the values of  $R_{i,j}$ ,  $C_i$  and  $P_i$  are computed as

$$R_{i,j} = \frac{1}{k_{i,j}\Delta l} \qquad , \tag{45}$$

where  $k_{i,j}$  is the thermal conductivity (with units W/m°C) between nodes i and j;



$$\begin{split} \mathbf{R} &= \frac{1}{\mathsf{k}\Delta^f} \left[ {^{\circ}\mathbf{C}} / \mathbf{W} \right] , \, \mathbf{k} = \mathsf{THERMAL} \; \mathsf{CONDUCTIVITY} \\ \mathbf{C} &= \rho C_\mathsf{P} \; \left( \Delta \ell \right)^3 \; \left[ \mathsf{J}/{^{\circ}\mathbf{C}} \right] , \, \rho = \mathsf{DENSITY}, \, \mathsf{C}_\mathsf{P} = \mathsf{SPECIFIC} \; \mathsf{HEAT} \\ \mathbf{P} &= \mathsf{SAR} \; \rho \; \left( \Delta \ell \right)^3 \; \left[ \mathbf{W} \right] , \mathsf{SAR} = \frac{\sigma}{2 \; \rho} \; \left| \mathbf{E} \right|^2 \; \mathsf{(Specific Absorption Rate)} \\ \sigma &= \mathsf{ELECTRICAL} \; \mathsf{CONDUCTIVITY} \\ \left| \mathbf{E} \right| \; = \mathsf{MAGNITUDE} \; \mathsf{OF} \; \mathsf{ELECTRICF FIELD} \end{split}$$

Figure 11. Equivalent circuit model used in thermal analysis. The inhomogeneous body of interest is modeled with a grid of nodes interconnected by resistors. Each node has an associated power source and capacitor.

$$C_i = \rho_i C_{pi} (\Delta l)^3 \qquad , \tag{46}$$

where  $C_{pi}$  is the specific heat at the *i*th node and  $\rho_i$  is the density  $(kg/m^3)$  at the *i*th node; and

$$P_i = (SAR)_i \rho_i (\Delta l)^3 \qquad , \tag{47}$$

where  $(SAR)_i$  is the SAR for the *i*th node, which is given by

$$(SAR)_i = \frac{\sigma_i}{2\rho_i} |E_i|^2 \qquad , \tag{48}$$

where  $\sigma_i$  is the electrical conductivity of the *i*th node and  $|E_i|$  is the magnitude of the electric field delivered by the hyperthermia array to the *i*th node. Note: in substituting Equation (48) into Equation (47), the density  $\rho_i$  cancels. Thus, an equivalent approach to computing the power delivered to the *i*th node is written in terms of the time-average power dissipated per unit volume of the *i*th node (denoted  $P_{di}$ ) as

$$P_{t} = P_{di}(\Delta l)^{3} \qquad . \tag{49}$$

A block diagram showing how TTA is used in the hyperthermia analysis is given in Figure 12. First, the method of moments, controlled by the SMI nulling algorithm, is used to compute the electric field radiation pattern throughout a homogeneous region (muscle tissue). This power distribution is then read into the TTA program, which computes the temperature distribution inside an elliptical muscle-tissue target surrounded with a constant-temperature water bolus.




Figure 12. The moment-method simulation's use with the transient thermal analysis. RF power distribution over the target region is converted into a time-dependent temperature distribution.

#### 3. SIMULATION RESULTS

To demonstrate the effectiveness of focused near-field adaptive nulling in reducing undesired hot spots, several computer simulations are presented. The E-field simulations are for the signal received by a short-dipole probe due to a transmitting phased array embedded in an infinite homogeneous lossy dielectric (muscle tissue). However, the thermal simulations are for an elliptical target (muscle tissue) surrounded by a constant-temperature water bolus. Because the RF wavelengths in the target and water bolus are similar, the E-field simulations are believed to give a reasonable approximation to the field distribution inside the elliptical target.

The E-field calculation in the assumed infinite homogeneous medium introduces additional field attenuation not present in a clinical hyperthermia system with an annular array transmitting through a water bolus into a patient. (As mentioned earlier, the water bolus has very little RF propagation loss.) In addition, the transmit array weights are normalized according to Equation (6). Thus, in this report no attempt is made to compute the absolute E-field strength in volts/meter in the elliptical target. Instead, the peak power in the elliptical target is adjusted (by a scale factor) to produce a desired maximum focal-region temperature  $(T_{max})$  after t minutes. Note: an approximate absolute scale factor could be computed by making an initial computer simulation with an infinite homogeneous water bolus and then matching the target boundary field to the infinite homogeneous muscle tissue simulation.

The computer simulation model is related, in part, to the hyperthermia annular phased-array antenna system shown in Figure 1. This hyperthermia system uses a 60-cm array diameter with eight uniformly spaced dipole elements (fed as four pairs) which operate over the frequency band 60-120 MHz [21]. The simulated array is assumed to be fully adaptive, and with eight adaptive elements seven independent nulls can be formed while simultaneously focusing on a tumor. For comparison, the BSD-2000 hyperthermia system could be used to form up to three independent nulls.

In this section, we assume that the adaptive radiation pattern null-width characteristics in a homogeneous target will be similar to the characteristics observed in an inhomogeneous target. The null width is directly related to the wavelength and there is only a 5 percent change in wavelength between the assumed muscle tissue and water bolus. With this assumption, the transmit array is embedded in homogeneous tissue, which allows direct use of the thin-wire moment-method formulation in Section 2.3. After computing the two-dimensional electric-field distribution in the homogeneous medium, we then consider only an elliptical portion of the homogeneous region and use the ellipse as the homogeneous target. In the thermal analysis, the elliptical target is surrounded with a constant  $10^{\circ}$ C water bolus. The electric field amplitude is scaled to produce a  $46^{\circ}$ C peak temperature, at time t = 20 minutes, at the center of the elliptical phantom. The initial temperature of the phantom is assumed to be  $25^{\circ}$ C (room temperature).

All computer simulations presented in this section are at a 120-MHz operating frequency with four auxiliary nulling sensors: that is,  $N_{aux} = 4$ . Degraded system performance due to an

insufficient number of auxiliary probes is considered in Appendix A. The parameters used in the electrical and thermal analyses are summarized in Table 1. These parameters, obtained in part

TABLE 1
Parameters Used in Electrical/Thermal Analysis

171588-5

PARAMETER	PHANTOM MUSCLE TISSUE	DISTILLED WATER
DIELECTRIC CONSTANT @ 100 MHz	73.5	80.0
ELECTRICAL CONDUCTIV!TY @ 100 MHz	0.5 S/m	0.0001 S/m
DENSITY	970.0 kg/m <sup>3</sup>	1000.0 kg/m <sup>3</sup>
SPECIFIC HEAT	3516.0 J/kg °C	4200.0 J/kg °C
THERMAL CONDUCTIVITY	0.544 W/m °C	0.6019 W/m °C

from Sullivan [9], are for a frequency of 100 MHz; we assume that similar values of the parameters will exist at 120 MHz. Notice that the relative dielectric constants of phantom muscle tissue and distilled water are very similar; however, the electrical conductivities are vastly different. The relevant thermal characteristics—density, specific heat, and thermal conductivity [35]—are very similar for phantom muscle tissue and distilled water.

#### 3.1 ELECTRIC FIELD FOR ARRAY IN HOMOGENEOUS TISSUE

Substituting the values f=120 MHz,  $\sigma=0.5$  S/m, and  $\epsilon_{\tau}=73.5$  into Equation (29) yields  $\gamma_m=10.0+j23.8$  for the muscle tissue. With  $\beta_m=23.8$  radians/m, the wavelength in

the phantom muscle tissue is  $\lambda_m=26.5$  cm. The attenuation constant for the muscle tissue is  $\alpha_m=10.0$  radians/m. Similarly, for distilled water  $\gamma_w=0.0021+j22.5$ , so the wavelength is  $\lambda_w=27.9$  cm. The attenuation constant for the distilled water medium is  $\alpha_w=0.0021$  radians/m. The propagation loss in the phantom muscle tissue is  $20\log_{10}e^{-10.0}$ , or -0.87 dB/cm. Similarly, the propagation loss in the distilled water is found to be -0.0002 dB/cm. Thus, the total loss due to propagation through 15 cm of distilled water is 0.003 dB. For 15 cm of muscle tissue the corresponding loss is 13.1 dB. The wave impedance in the muscle tissue is computed from Equation (34) as  $\eta_m=33.9+j$  14.2  $\Omega$ , and similarly in the distilled water  $\eta_w=42.1+j$  0.004  $\Omega$ .

The geometry used in the simulations is shown in Figure 13. A 60-cm-diameter ring array of eight dipoles uniformly surrounds a fictitious elliptical target zone with major axis 30 cm and minor axis 20 cm. The length of each perfectly conducting center-fed dipole array element at 120 MHz in the infinite homogeneous muscle tissue is  $\lambda/2$ , or 13.25 cm. The array focus is assumed at the origin, and four auxiliary short-dipole sensors with length 1.27 cm  $(0.05\lambda)$  are positioned in (x,y,z) coordinates at (15 cm, 0, 0), (-15 cm, 0, 0), (0, 0, 10 cm), and (0, 0, -10 cm); that is, the auxiliary E-field sensors are located every 90° in azimuth on the perimeter of the target. In rectangular coordinates, each dipole is oriented in the  $\hat{y}$  direction and the feed terminals of each dipole are located at y=0. The moment-method computer simulations were run on a Sun 3/260 workstation. The total CPU time for a complete moment-method run is 19.2 minutes. This CPU time includes computing the quiescent and adaptive radiation patterns on a 41 by 41 grid of points. The CPU time without radiation pattern calculations is 33 seconds.

The two-dimensional radiation pattern in the plane y=0, before nulling, at 120 MHz with uniform amplitude and phase illumination, is shown in Figure 14. The calculated data are collected on a 41 by 41 grid of points over a square region, with side length 76.2 cm, centered at the focus. The spacing between data points is 1.905 cm, or  $0.072\lambda$ , and the contour levels are displayed in 10-dB steps. The E-field data are computed for the case of a 1.27-cm short-dipole observation probe. The positions of the eight dipole radiators are clearly evident by the -20-dB contours surrounding each element. The radiation pattern is symmetric because of the symmetry of the array and the assumed homogeneous medium.

Finer contour levels (1-dB steps) for the quiescent radiation pattern are displayed in Figure 15. Here, it is evident that the focused main beam of the ring array is increasing in amplitude as the observation point moves closer to the focus. Away from the main beam region, the pattern amplitude is seen to increase as the observation position moves away from the focus.

A quiescent radiation pattern cut at z=0 is shown in Figure 16. The large amplitude that occurs at  $\pm 30$  cm is due to the E-field probe's close proximity to the transmitting elements. The large attenuation that occurs from the ar.ay diameter to the focus is due to the 1/r attenuation loss and the loss in the uniform homogeneous muscle tissue. The radiation pattern cut for x=0 in Figure 17 is identical to the pattern in Figure 16 due to the symmetry of the array. In both Figures 16 and 17 the boundary of the fictitious elliptical target zone is indicated. The increasing radiation pattern amplitude near the left and right sides of the elliptical target (Figure 16) will be




Figure 13. Geometry for eight-element ring array and four E-field auxiliary sensors.




Figure 14. Simulated two-dimensional quiescent radiation pattern at 120 MHz for uniformly illuminated eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue:  $\epsilon_{\tau}=73.5, \sigma=0.5$ ). A fictitious elliptical target region is indicated by the dashed curve. Radiation contour levels are in 10-dB steps.




Figure 15. Simulated two-dimensional quiescent radiation pattern at 120 MHz for uniformly illuminated eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue:  $\epsilon_r = 73.5, \sigma = 0.5$ ). A fictitious elliptical target region is indicated by the dashed curve. Radiation contour levels are given in 1-dB steps.

shown to produce hot spots in the thermal distribution. As shown in the next section, because the top (anterior) and bottom (posterior) of the elliptical target are not as strongly illuminated as the left and right sides, no quiescent hot spots occur at the top or bottom. From Figure 16, the




Figure 16. Simulated one-dimensional quiescent radiation pattern cut (z=0) at 120 MHz for uniformly illuminated eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue:  $\epsilon_r = 73.5, \sigma = 0.5$ ).

ring-array half-power beamwidth in the target region is approximately 13 cm, or approximately one-half the wavelength (26.5 cm) in the phantom muscle tissue. <sup>2</sup> The adaptive nulling resolution

<sup>&</sup>lt;sup>2</sup>Sullivan [54] has recently published simulations and measurements of SAR patterns at 70-110 MHz. using a homogeneous  $35 \times 25$ -cm elliptical target surrounded with a layer of fat and a water




Figure 17. Simulated one-dimensional quiescent radiation pattern cut (x=0) at 120 MHz for uniformly illuminated eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue:  $\epsilon_r = 73.5, \sigma = 0.5$ ).

or closest allowed spacing between a deep adaptive null and the main beam is equal to the half-power beamwidth of the antenna [55]. Thus, the closest allowed null position is 13 cm from the focus. Since the target width is 30 cm, two nulls can be formed at  $(x = \pm 15 \text{ cm}, z = 0)$  without disturbing the focus. However, if two deep nulls are formed at  $(x = 0, z = \pm 10 \text{ cm})$  the focus will be compromised. In practice, the water bolus would restrict the placement of short-dipole sensors to the surface of the target. Thus, only weak nulls can be formed at  $(x = 0, z = \pm 10 \text{ cm})$  so that the focus will not be affected by the adaptive nulling process.

Next, adaptive radiation patterns are computed with four auxiliary dipole sensors; their positions are shown in Figure 13. The value of the receiving gain for auxiliary dipoles 1 and 2 is adjusted to produce a greater-than-35-dB SNR. This amount of SNR results in greater than 35 dB of nulling in the direction of auxiliary dipoles 1 and 2. In contrast, the gain values for auxiliary dipoles 3 and 4 are turned down to produce about a 3-dB SNR. Thus, only about 3 dB of nulling will occur at sensor positions 3 and 4 as the adaptive algorithm reduces the interference to the noise level of the receiver. The reason for choosing these null depths will become apparent with the data that follow. Figure 18 shows the two-dimensional radiation pattern after nulling with four auxiliary sensors. Two deep adaptive nulls at  $x = \pm 15$  cm occur as expected, and weak nulling occurs at  $z=\pm 10$  cm, also as expected. The two deep nulls in the z=0 cut are quantified in Figure 19. where greater than 35 dB of interference nulling or pattern reduction occurs at  $x = \pm 15$  cm. The peak level at the focus is adjusted to 0 dB for both the quiescent and adaptive patterns. weak adaptive nulls are in effect in the x=0 radiation pattern cut shown in Figure 20; however, weak nulls are desired in this cut due to temperature considerations, as shown in the next section. The weak nulls in effect in the adaptive patterns reduce variation from the quiescent radiation pattern. Finally, the transmit array weights before and after nulling and the covariance matrix eigenvalues are shown in Figures 21 and 22, respectively. The adaptive transmit weights exhibit a 5-dB dynamic range in Figure 21(a). There are two large eigenvalues and two weak (nonzero) eigenvalues shown in Figure 22. These eigenvalues are directly associated with the two high-SNR auxiliary sensors and the two weak-SNR auxiliary sensors. Note that the 0-dB level in Figure 22 is equal to the noise level. The probe-array output power before and after adaptive nulling is 31.4 dB and 0.9 dB, respectively. This difference in power before and after nulling means that the adaptive cancellation is -30.5 dB. The output file for the moment-method simulation of this section is given in Appendix B.

bolus. A comparison of Sullivan's 110-MHz SAR simulation (Figure 5(a) of his paper) and the results (at 120 MHz) presented in Figures 16 and 17 of this report indicate reasonable agreement over the elliptical target region.




Figure 18. Simulated two-dimensional adaptive radiation pattern at 120 MHz for eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue:  $\epsilon_{\tau}=73.5, \sigma=0.5$ ). Radiation contour levels are given in 10-dB steps. Four auxiliary sensors are used in forming the adaptive pattern. The quiescent focus is at (0,0).




Figure 19. Simulated one-dimensional quiescent and adaptive radiation patterns in the z=0 cut at 120 MHz for eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue:  $\epsilon_r=73.5, \sigma=0.5$ ). Four auxiliary sensors are used in forming the adaptive pattern.




Figure 20. Simulated one-dimensional quiescent and adaptive radiation patterns in the x=0 cut at 120 MHz for eight-element ring array in infinite homogeneous conducting medium (phantom muscle tissue:  $\epsilon_{\tau}=73.5, \sigma=0.5$ ). Four auxiliary sensors are used in forming the adaptive pattern.




Figure 21. Transmit array weights before and after adaptive nulling. A dynamic range of about 5 dB is evident for the adaptive weights. Four auxiliary sensors are used in the adaptive process. (a) Amplitude: (b) Phase.




Figure 22. Channel covariance matrix eigenvalues (degrees of freedom) used in the adaptive process with four auxiliary sensors.

## 3.2 TEMPERATURE DISTRIBUTION IN ELLIPTICAL PHANTOM

In this section, the transient thermal analysis software is used to compute the temperature distribution in an elliptical phantom surrounded with a constant-temperature water bolus. The 41 by 41 two-dimensional E-field radiation pattern data from the previous section are used as the power source for the thermal node network. Two spacings of nodes are considered here: initially, the node spacing  $\Delta x = \Delta z = \Delta l = 1.905$  cm (coarse grid) is used; to check convergence, the node spacing is then decreased by a factor of two to  $\Delta l = 0.9525$  cm (fine grid). The coarser spacing is shown to be adequate.

#### 3.2.1 1.905-cm Thermal Node Spacing

Note: the scale factors used to convert the normalized E-field distributions to a power level that induces a  $46^{\circ}$ C peak temperature at t=20 minutes are 94.1 dB and 96.0 dB for the quiescent and adaptive patterns, respectively. These scale factors are determined by trial and error. From equations given in Section 2.3.2 and the parameter values in Table 1, all resistors in the phantom muscle tissue had a value of  $96.5^{\circ}$ C/W and all resistors in the water bolus had a value of  $87.2^{\circ}$ C/W. The value of the capacitors in the phantom muscle tissue is 23.6 J/ $^{\circ}$ C; values for capacitors in the

water-bolus region are not used in the input to the transient thermal analysis software. Instead, a constant temperature of 10°C is enforced at each water-bolus node. With a 41 by 41 grid, a total of 3280 resistors and 1681 capacitors are used in the thermal simulation. The CPU time required to compute the temperature distribution on a 41 by 41 grid is under four minutes.

Figure 23 shows the two-dimensional temperature distribution produced at time t=20 minutes in the elliptical phantom muscle tissue target without adaptive nulling. To generate Figure 23, the power source used in the transient thermal analysis is the quiescent radiation pattern given in Figure 14. The initial temperature (at time t=0) is 25°C. Notice the occurrence of two hot spots on the left and right sides of the elliptical phantom. The peak temperature on-focus is 46°C, which is achieved by scaling the normalized quiescent E-field as described earlier. The two hot spots are quantified in the z=0 temperature pattern cut shown in Figure 24. While the peak temperature at each hot spot in Figure 24 is only 41°C, the temperature profile for x=0 in Figure 25 shows no hot spots. As any undesired hot spot is a potential source for compromising the therapy session, adaptive nulling is used to reduce the sidelobes corresponding to the hot spots.

Figure 26 shows the simulated two-dimensional thermal distribution at time t=20 minutes, with adaptive nulling at four auxiliary sensors in effect. The focal-spot diameter with adaptive nulling is equivalent to the focal-spot diameter before adaptive nulling. Hot spots on the left and right sides of the target are eliminated. A comparison of the temperature distribution before and after nulling along the major axis (z=0) of the target ellipse is made in Figure 27. Similarly, the temperature distribution before and after nulling along the minor axis of the target ellipse is shown in Figure 28.

#### 3.2.2 Convergence Check: 0.9525-cm Thermal Node Spacing

The convergence of the previous thermal simulations was verified by increasing the density of E-field observation probe positions by a factor of two, with a new spacing between points of 0.9525 cm, still with a 41 by 41 grid. The ring array operates as before at 120 MHz, and there are four auxiliary sensors laid out as shown in Figure 13. As the auxiliary positions are the same, the adaptive weights and covariance matrix eigenvalues in Figures 21 and 22, respectively, remain the same. From the parameter values in Table 1, all resistors in the finer-grid muscle-tissue phantom had a value of 193.0°C/W and all resistors in the water bolus had a value of 174.4°C/W. The value of the capacitors in the phantom muscle-tissue is 2.95 J/°C. Again, a constant temperature of 10°C is enforced at each water-bolus node. The E-field scaling factors to raise the focal-point temperature to 46°C before and after nulling are 76.5 dB and 78.4 dB, respectively. The finer-grid two-dimensional thermal distributions before and after nulling are shown in Figures 29 and 30, respectively. Although the temperature contours are smoother, the general agreement between these patterns and the coarser-grid patterns in Figures 23 and 26 are evident. Similarly, onedimensional thermal pattern cuts with the finer grid are shown in Figures 31 (x axis) and 32 (z axis); good agreement with the coarse-grid patterns (Figures 27 and 28) is observed. In particular, the finer detail in Figure 31 shows that the hot spots are at a 42°C level compared to the 41°C level observed for the coarse grid (Figure 24). Thus, convergence of the coarse-grid thermal patterns is demonstrated.




Figure 23. Simulated two-dimensional thermal pattern at time t=20 minutes before nulling in elliptical phantom muscle-tissue target surrounded with 10°C constant-temperature water bolus. The incident RF power distribution, from Figure 14, is at 120 MHz. Temperature contour levels are given in 2°C steps. Hot spots on the left and right sides of the target are observed.




Figure 24. Simulated one-dimensional (z=0) thermal pattern at time t=20 minutes before nulling in elliptical phantom muscle-tissue target surrounded with  $10^{\circ}$ C constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. Hot spots on the left and right sides of the target are observed.




Figure 25. Simulated one-dimensional (x=0) thermal pattern at time t=20 minutes before nulling in elliptical phantom muscle-tissue target surrounded with  $10^{\circ}C$  constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. No undesired hot spots are present.




Figure 26. Simulated two-dimensional thermal pattern at time t=20 minutes (with adaptive nulling at four auxiliary sensors in effect) in elliptical phantom muscle-tissue target surrounded with 10°C constant-temperature water bolus. The adapted incident RF power distribution, from Figure 18, is at 120 MHz. Temperature contour levels are given in 2°C steps. Hot spots on the left and right sides of the target are eliminated.




Figure 27. Simulated one-dimensional (z=0) thermal patterns at time t=20 minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with 10°C constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. Hot spots on the left and right sides of the target are clearly eliminated by the adaptive nulling process. Four auxiliary sensors are used in the adaptive process.




Figure 28. Simulated one-dimensional (x=0) thermal patterns at time t=20 minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with  $10^{\circ}$  C constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. No undesired hot spots are present. Four auxiliary sensors are used in the adaptive process.




Figure 29. E-field probe-sample-spacing convergence check for simulated two-dimensional thermal pattern at time t=20 minutes (before adaptive nulling) in elliptical phantom muscle-tissue target surrounded with  $10^{\circ}$  C constant-temperature water bolus. The quiescent incident RF power distribution is at 120 MHz. Temperature contour levels are given in 2° C steps. The grid spacing is one-half the spacing of that in Figure 23. Hot spots on the left and right sides of the target are present, as previously observed for the coarser probe-sample spacing.




Figure 30. E-field probe-sample-spacing convergence check for simulated two-dimensional thermal pattern at time t=20 minutes (with adaptive nulling in effect) in elliptical phantom muscle-tissue target surrounded with  $10^{\circ}$  C constant-temperature water bolus. The adapted incident RF power distribution is at 120 MHz. Temperature contour levels are given in  $2^{\circ}$  C steps. The grid spacing is one-half the spacing of that in Figure 26. Hot spots on the left and right sides of the target are eliminated, as previously observed for the coarser probe-sample spacing.




Figure 31. E-field probe-sample-spacing convergence check for simulated one-dimensional (z=0) thermal patterns at time t=20 minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with  $10^{\circ}$  C constant-temperature water bolus. The quiescent incident RF power distribution is at 120 MHz with 0.9525-cm sample spacing. Hot spots on the left and right sides of the target are eliminated by the adaptive nulling process, as previously observed in Figure 27. Four auxiliary sensors are used in the adaptive process.




Figure 32. E-field probe-sample-spacing convergence check for simulated one-dimensional (x=0) thermal patterns at time t=20 minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with 10°C constant-temperature water bolus. The quiescent incident RF power distribution is at 120 MHz with 0.9525-cm sample spacing. No undesired hot spots are present. Four auxiliary sensors are used in the adaptive process.

### 4. CONCLUSION

Focused near-field adaptive nulling has been applied to the problem of generating therapeutic temperature distributions, free of undesired hot spots, for deep-regional cancer treatment in a target body. A noninvasive adaptive hyperthermia system concept has been described and is analyzed in detail. Auxiliary short-dipole nulling sensors are positioned on the surface of the target. Due to the finite width of an electric-field null, a null formed at the target surface extends into the interior region of the target. When the width of the null is properly chosen, undesired hot spots inside the target are eliminated while simultaneously focusing on a deep-seated tumor. The resolution between a deep null and focus is limited by the half-power beamwidth of the hyperthermia array. The resolution can be somewhat enhanced by using weak nulls whenever the spacing between the null and focus is less than one half-power beamwidth. The depth of null is controlled by the SNR at the auxiliary probe position.

A theory for analyzing adaptive nulling with a phased array in an infinite homogeneous conducting medium has been developed. The method of moments is used to compute the electric field at a short-dipole sensor due to a thin-wire dipole ring array. The SMI algorithm is used to adaptively control the transmit array weights and to form radiation pattern nulls. The method used for analyzing the transient thermal behavior of an RF-illuminated target by using an equivalent electric circuit network has been described.

Computer simulations of the electric field and induced temperature distribution for an eightelement dipole ring array at 120 MHz have been presented. Four auxiliary E-field short-dipole
sensors are used in eliminating two widely separated hot spots surrounding the array focus at the
center of an elliptical target with dimensions 30 cm by 20 cm. Calculated data indicate that simultaneous nulling and focusing can yield desirable electric-field distributions for hyperthermia.
Undesired hot spots controlled with an insufficient number of auxiliary nulling probes have been
shown to simply redistribute within the target. Thus, an adequate number of auxiliary probes
must be provided to ensure that no undesired hot spots are formed in deep-regional hyperthermia therapy. Although the computer simulations are implemented for a phantom target, adaptive
nulling may prove beneficial in clinical trials of hyperthermia with living targets. Further adaptivenulling simulations with inhomogeneous targets are desirable, as are experimental measurements of
adaptive hyperthermia in phantom targets. A study of proper number and positioning of the noninvasive sensors for treating arbitrary tumor positions needs to be performed. The adaptive-nulling
technique described here should be applicable to both invasive and noninvasive RF hyperthermia
systems, as well as ultrasound hyperthermia systems.

# APPENDIX A SYSTEM DEGRADATION DUE TO INSUFFICIENT NUMBER OF AUXILIARY PROBES

This section shows simulation results which indicate that an insufficient number of auxiliary nulling probes will degrade adaptive hyperthermia system performance. The simulation parameters are the same as presented in Section 3, except only two auxiliary probes  $(N_{aux} = 2)$  are used to form nulls at the surface of the elliptical muscle-tissue target. The geometry used in the simulations is shown in Figure A-1. A 60-cm diameter ring array of eight dipoles uniformly surround a fictitious elliptical target zone with major axis 30 cm and minor axis 20 cm. The length of each perfectly conducting center-fed dipole array element, at 120 MHz in the infinite homogeneous muscle tissue, is  $\lambda/2$  or 13.25 cm. The array focus is assumed to be at the origin, and two auxiliary short-dipole sensors with length 1.27 cm are positioned at (15 cm, 0) and (-15 cm, 0). In rectangular coordinates, each dipole is oriented in the  $\hat{y}$  direction and the feed terminals of each dipole are located at y = 0. The quiescent radiation patterns and thermal distributions are the same as




Figure A-1. Geometry for eight-element ring array and two E-field auxiliary sensors.

presented in Figures 14, 15, 16, 17, 23, 24, and 25.

With only two auxiliary sensors, the two-dimensional adapted electric-field radiation pattern shown in Figure A-2 indicates that both auxiliary positions are nulled. However, the two-dimensional adaptive thermal distribution in Figure A-3 shows that the quiescent hot spots (see Figure 23) have only been redistributed. The one-dimensional thermal distribution in Figure A-4 shows that the hot spots are reduced on the left and right sides of the elliptical target, and Figure A-5 shows that the hot spots are redistributed to lie along the z axis. Thus, two additional auxiliary probes (probes 3 and 4 in Figure 13) are necessary to avoid the possibility of redistributed hot spots. Finally, the adaptive weights and covariance matrix eigenvalues are shown in Figures A-6 and A-7, respectively. There are only small differences between the weights and eigenvalues compared to the four-auxiliary probe case (see Figures 21 and 22). The transmit weights cover a 7-dB dynamic range and there are only two nonzero eigenvalues. The probe-array power outputs before and after adaptive nulling are 31.4 dB and 0.0 dB, respectively. Thus, the interference cancellation is to the noise level, or -31.4 dB. An output file for the moment-method simulation for this section is given in Appendix B.




Figure A-2. Simulated two-dimensional adaptive radiation pattern at 120 MHz for eightelemen, ring array in infinite homogeneous conducting medium (phantom muscle tissue:  $\epsilon_{\tau}=73.5, \sigma=0.5$ ). Radiation contour levels are given in 10-dB steps. Two auxiliary sensors are used in forming the adaptive pattern. The quiescent focus is at (0,0).




Figure A-3. Simulated two-dimensional thermal pattern at time t=20 minutes (with adaptive nulling at only two auxiliary sensors in effect) in elliptical phantom muscle-tissue target surrounded with 10°C constant-temperature water bolus. The adapted incident RF power distribution, from Figure 18, is at 120 MHz. Temperature contour levels are given in 2°C steps. Initial hot spots (before nulling) on the left and right sides of the target (see Figure 23) are redistributed to the top and bottom (anterior and posterior positions).




Figure A-4. Simulated one-dimensional (z=0) thermal patterns at time t=20 minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with 10°C constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14, is at 120 MHz. Hot spots on the left and right sides of the target are eliminated by the adaptive nulling process. Two auxiliary sensors are used in the adaptive process.




Figure A-5. Simulated one-dimensional (x=0) thermal patterns at time t=20 minutes before and after nulling in elliptical phantom muscle-tissue target surrounded with  $10^{\circ}C$  constant-temperature water bolus. The quiescent incident RF power distribution, from Figure 14. is at 120 MHz. Two undesired hot spots are present in the adaptive pattern. Two auxiliary sensors are used in the adaptive process.




Figure A-6. Transmit array weights before and after adaptive nulling. A dynamic range of about 7 dB is evident for the adaptive weights. Two auxiliary sensors are used in the adaptive process. (a) Amplitude; (b) Phase.




Figure A-7. Channel covariance matrix eigenvalues (degrees of freedom) used in the adaptive process with two auxiliary sensors.

### APPENDIX B SOFTWARE DOCUMENTATION

This appendix lists sample input and output files for the adaptive hyperthermia simulation. Listings of the key FORTRAN programs that implement adaptive nulling in a homogeneous conducting medium are also given. Additional moment-method software not listed here is contained in Richmond [37].

The first data file (four auxiliary sensors) is used to generate the E-field results presented in Section 3.1. The second data file (two auxiliary sensors) was used in calculating the electric-field distribution presented in Appendix A. The corresponding output files give the values for the array mutual coupling, quiescent and adaptive transmit weights, covariance matrix, eigenvalues, and cancellation.

```
***Input data file for adaptive nulling with four auxiliary probes.
***sdipjamhyper data file, filename sdipjamhyper.datacirconicrr4
&DIPOLE NCOLX=8,NROWY=1,HZIN=0.001,HLIN=2.6,
 DXIN=6.888, DYIN=6.888, ARADIN=0.0039,
 ICIRC=1.CRADIN=11.81.HLSIN=0.25.
 IWL=0,FCHZ=120.0E6,BWFHZ=1.0E0,NFREQ=5,IWR=0,
 NPDX=0,NXDUM=0,NYDUM=0,ER3=73.5,SIG3=0.5,TD3=-1.0,
 ZLOAD=50.0.ZCHAR=0.0.NGEN=8.IGEN=0.
 NSCANS=1.
 THSINC=5.0.
 IMUT=1. IBLTSL=0.
 IENORM=1, ICHEB=0, SLLDB=20., EDGTDB=0.,
 IPATRN=1, IPRCOM=1, IANGLP=0,
 NPHCT=0,NTHPT=499,
 THDR=180., THDMIN=-90.
 NCOLXN=121, NROWYN=1, NCOLZN=1,
 RLSYIN=0.0.RLSZIN=0.0,
 NCOLXN=117,NROWYN=1,NCOLZN=1,
 INEAR=1.
 IPOL=2, IGRNDP=0
 ITLTPR=0.ITLTDP=0.
 NFCOLX=1,NFROWY=1,
 IANTX=1, IANTY=0, NPOWER=0,
 IPCONN=1, IPCONF=0, IPCUTF=0, IPCFX=0, IPCFY=0, IPCFZ=0,
 ITEK=0.
 IQUAN=0, IRNERR=0, ELERDB=0.02, ELERDG=0.2, NBMOD=12,
 NRAN=1, NBADWT=32, AWERDB=0.0, AWERDG=0.0,
 NAUX=7, IAUXA(1)=1,2,3,4,5,6,7,
 IATTEN=1.
 NJAMS=7.ISLC=0,AUXADB(1)=8+0.0.
 PWRJDB(1)=40.,40.,15.,15.,3*-99.0,
 YJAMIN(1)=7=0.0,
 XFOCIN=0.0, ZFOCIN=0.0,
 IJAMIN(1)=5.9,-5.9,2*0.0,3*-0.0,
 ZJAMIN(1)=0.0,0.0,4.0,-4.0,3*0.0
 YNIN=0.0,
 NCDLIN#101 NTHPT#101.
 ININ=-15., ZNIN=-15., RLSXIN=30., RLSZIN=30., NCOLXN=41, NCOLZN=41,
MEND
```

```
***Output data file for adaptive nulling with four auxiliary probes.
***sdipjamhyper data file, filename sdipjamhyper.datacirconicrr4
 XJAMIN(1),YJAMIN(1),ZJAMIN(1)=
                                   5.900
                                               0.000
 FGHZ=
           0.1200000
 DX .DY .HL .ARAD=
                                     0.17496
                                                  0.06604
                       0.17496
                                                                0.00010
 IV,VTA=
               0.5916E-02 0.2211E-02
          1
               0.5916E-02 0.2211E-02
 IV.VTA=
 IV, VTA=
               0.5916E-02 0.2211E-02
 IV, VTA=
          4
               0.5916E-02 0.2211E-02
 IV, VTA=
               0.5916E-02 0.2211E-02
          5
 IV.VTA=
          6
               0.5916E-02 0.2211E-02
 IV, VTA=
               0.5916E-02 0.2211E-02
 IV, VTA= 8
               0.5916E-02 0.2211E-02
******* array mutual impedance matrix (first row)*******
 Z(1, 1) = 0.91705E+02 0.12403E+02
 Z(1. 2)=-0.33552E-01 0.19136E+00
 Z(1.
       3)=-0.14488E-01-0.54125E-02
Z(1,
       4)= 0.29420E-02 0.13780E-02
 Z(1,
       5)= 0.15283E-02-0.11371E-02
       6)= 0.29420E-02 0.13780E-02
 Z(1.
       7) =-0.14488E-01-0.54125E-02
 Z(1,
       8)=-0.33552E-01 0.19136E+00
 Z(1.
 NEL, NCOLX, NROWY=
                  8
                       8
 CURRENTS
    1 1.000000
                        0.683E-04
                                       13.
    2 1.000000
                        0.683E-04
                                        13.
    3 1.000000
                        0.683E-04
                                       13.
     4 1.000000
                        0.683E-04
                                        13.
    5 1.000000
                        0.683E-04
                                       13
    6 1.000000
                        0.683E-04
                                       13.
    7 1.000000
                        0.683E-04
                                       13.
     8 1.000000
                        0.683E-04
                                       13.
IC.RWTA= 1 0.35355E+00
IC.RWTA= 2 0.35355E+00
IC,RWTA= 3 0.35355E+00
IC.RWTA=
               0.35355E+00
           4
IC,RWTA=
           5
               0.35355E+00
IC, RWTA=
               0.35355E+00
           6
IC,RWTA= 7
               0.35355E+00
IC.RWTA= 8 0.35355E+00
I,CWTA(I,1)=
               1
                   0.34510E+00-0.76830E-01
I,CWTA(I,1)=
               2 0.34510E+00-0.76830E-01
I.CWTA(I.1)= 3
                   0.34510E+00-0.76830E-01
I,CWTA(I,1) \approx 4 0.34510E+00-0.76830E-01
I,CWTA(I,1)=
I,CWTA(I,1)=
               5
                   0.34510E+00-0.76830E-01
               6
                   0.34510E+00-0.76830E-01
I,CWTA(I,1)=
               7
                  0.34510E+00-0.76830E-01
I,CWTA(I,1)= 8 0.34510E+00-0.76830E-01
******ring array weights before nulling (amp,phase)**********
    1 CWTADB, CWTADG=
I =
                             0.00000
                                         -12.55098
    2 CWTAP" "TADG=
T=
                             0.00000
                                         -12.55099
    3 CWTADB, CWTADG=
I=
                             0.00000
                                         -12.55098
Ī=
    4 CWTADB, CWTADG=
                             0.00000
                                         -12.55097
1=
    5 CWTADB.CWTADG=
                             0.00000
                                         -12.55098
I=
     6 CWTADB, CWTADG=
                             0.00000
                                         -12.55097
    7 CWTADB, CWTADG=
                             0.00000
                                         -12.55098
] =
I= 8 CWTADB,CWTADG=
                             0.00000
                                         -12.55097
Z(1.
      1)= 0.91705E+02 0.12403E+02
      2)=-0.33552E-01 0.19136E+00
```

```
3)=-0.14488E-01-0.54125E-02
Z(1,
Z(1, 4) = 0.29420E-02 0.13780E-02
Z(1, 5) = 0.15283E-02-0.11371E-02
       6)= 0.29420E-02 0.13780E-02
Z(1.
Z(1,
       7)=-0.14488E-01-0.54125E-02
       8)=-0.33552E-01 0.19136E+00
Z(1,
NEL, NCOLI, NROWY=
                     8
                        8 1
II,PWRJ IN POWER=
                        0.10000E+05
                    1
II, PURJ IN POWER=
                       0.10000E+05
                     3 0.31623E+02
II,P&RJ IN POWER=
II, PWRJ IN POWER=
                       0.31623E+02
                         0.12589E-09
II,PWRJ IN POWER=
                     5
II, PWRJ IN POWER=
                     6
                         0.12589E-09
                         0.12589E-09
II,PWRJ IN POWER=
                    7
******COV&riance matrix*****************
                                  36.35
I, J, CNDB, PHASEN= 1 1
I.J.CNDB.PHASEN=
                   1
                                  31.68
                                              91.33
                   1
                                             -117.95
                                  24.99
I, J, CNDB, PHASEN=
                       3
I, J, CNDB, PHASEN= 1
                                  20.88
                                               6.26
                                               -0.14
                                  20.72
I, J, CNDB, PHASEN= 1 5
                                  20.88
                                               6.26
I, J, CNDB, PHASEN=
                        6
                   1
                                             -117.95
I, J, CNDB, PHASEN=
                        7
                                  24.99
                                  31.68
                                              91.33
I, J, CNDB, PHASEN=
                        8
                                              -91.33
                                  31.68
I, J, CNDB, PHASEN=
                    2 1
                                  27.05
                                               0.00
I, J, CNDB, PHASEN=
                                              155.55
                        3
                                  20.23
I, J, CNDB, PHASEN=
                                             177.91
                                  8.12
I, J, CNDB, PHASEN=
                                              -6.14
                                  20.90
I, J, CNDB, PHASEN=
                                   8.28
                                              177.98
                        6
 I, J, CNDB, PHASEN=
                        7
                                  20.24
                                              155.46
 I.J.CNDB.PHASEN=
                                  27.04
                                               0.00
                        8
 I, J, CNDB, PHASEN=
 I.J.CNDB.PHASEN=
                                  24.99
                                              117.95
                                  20.23
                                             -155.55
 I, J, CNDB, PHASEN=
                       2
                                  16.67
                                               0.00
I, J, CNDB, PHASEN=
                                             -155.49
                                  20.26
 I, J, CNDB, PHASEN=
                        4
                                  25.01
                                              117.95
 1, J, CNDB, PHASEN=
                                             -155.41
 I, J, CNDB, PHASEN=
                        6
                                  20.27
                                  16.52
                                               0.00
                        7
                    3
 I, J, CNDB, PHASEN=
                                             -155.46
                        8
                                  20.24
 I, J, CNDB, PHASEN=
                                               -6.26
                                  20.88
 I, J, CNDB, PHASEN=
                        1
                                             -177.91
                                   8.12
 I, J, CNDB, PHASEN=
                                              155.49
                                  20.26
                        3
 I, J, CNDB, PHASEN=
                                   27.08
                                               0.00
 I, J, CNDB, PHASEN=
                                   31.71
                                              -91.33
 I, J, CNDB, PHASEN=
                        5
                                   27.07
                                               0.00
 I.J.CNDB.PHASEN=
                                   20.27
                                              155.41
                        7
 I, J, CNDB, PHASEN=
 I.J.CNDB.PHASEN=
                        3
                                    8.28
                                             -177.98
                                                0.14
                                   20.72
 I, J, CNDB, PHASEN=
                        1
                                   20.90
                                                6.14
 I, J, CNDB, PHASEN=
                                             -117.95
                                   25.01
 I, J, CNDB, PHASEN=
                    5 3
 I, J, CNDB, PHASEN=
                    5
                                   31.71
                                               91.33
                                               0.00
 I, J, CNDB, PHASEN=
                     5
                         5
                                   36.38
                                   31.71
                                               91.33
 I, J, CNDB, PHASEN=
                    5
                         6
                                             -117.95
                         7
                                   25.01
 I, J, CNDB, PHASEN=
                         8
                                   20.90
                                                6.14
 I, J, CNDB, PHASEN=
                     5
                                               -6.26
                         1
                                   20.88
 I, J, CNDB, PHASEN=
                                   8.28
                                             -177.98
                         2
 I, J, CNDB, PHASEN=
                     6
                                   20.27
                                              155.41
 I.J.CNDB.PHASEN=
                                   27 07
                                                0.00
 I.J.CNDR PHACEN-
```

```
31.71
                                          -91.33
I, J, CNDB, PHASEN=
                                27.08
                                           0.00
I, J, CNDB, PHASEN= 6 6
                                20.26
                                          155.49
                     7
I, J, CNDB, PHASEN= 6
I, J, CNDB, PHASEN=
                                 8.12
                                          -177.91
                 6
                      8
I, J, CNDB, PHASEN=
                                          117.95
                                24.99
                  7
                      1
                               20.24
                                         -155.46
I, J, CNDB, PHASEN= 7
                                           0.00
                               16.52
I, J, CNDB, PHASEN= 7
                      3
                                20.27
                                          -155.41
I, J, CNDB, PHASEN=
                  7
                      4
                                          117.95
                                25.01
I, J, CNDB, PHASEN=
                  7
                      5
                               20.26
                                          -155.49
I, J, CNDB, PHASEN=
                      6
                                           0.00
                      7
                               16.67
I, J, CNDB, PHASEN= 7
                                          -155.55
I, J, CNDB, PHASEN= 7
                      8
                                20.23
                               31.68
                                          -91.33
I, J, CNDB, PHASEN= 8
                      1
I,J,CNDB,PHASEN= 8 2
                                            0.00
                               27.04
                                           155.46
I,J,CNDB,PHASEN= 8 3
                              20.24
                                          177.98
                                8.28
I, J, CNDB, PHASEN= 8 4
I,J,CNDB,PHASEN= 8 5
I,J,CNDB,PHASEN= 8 6
                                20.90
                                           -6.14
                                8.12
                                           177.91
I, J, CNDB, PHASEN= 8 7
                                20.23
                                           155.55
                                27.05
                                            0.00
I, J, CNDB, PHASEN= 8 8
*********eigenvalues**************
I,EVLDBN= 1
                   37.437
I,EVLDBN= 2
                   37.195
I,EVLDBN= 3
                   2.771
                     3.103
I,EVLDBN= 4
                     0.000
I,EVLDBN=
                     0.000
I,EVLDBN=
           6
                      0.000
I.EVLDBN= 7
                     0.000
I,EVLDBN= 8
 CPROD1, CPROD2= 0.13839E+04-0.85265E-13 0.10000E+01 0.00000E+00
 INR NORMALIZATION PARAMETER, INRNOR= 1
I= 1 WANDB=
                   -12.318
 I= 2 WANDB=
                    -8.734
    3 WANDB=
                    -7.602
I=
I=
     4 WANDB=
                    -8.734
    5 WANDB=
                   -12.318
T =
    6 WANDB=
                    -8.734
 I=
     7 WANDB=
                     -7.602
 I=
     8 WANDB=
                     -8.734
 I,WAN(I,1)= 1 0.11307E+00-0.21414E+00
 I,WAN(I,1)= 2 0.34937E+00 0.10851E+00
I,WAN(I,1)= 3 0.32832E+00-0.25671E+00
I,WAN(I,1)= 4 0.34937E+00 0.10851E+00
 I, WAN(I,1)= 5 0.11307E+00-0.21414E+00
 I.WAN(I.1) = 6 0.34937E+00 0.10851E+00
 I,WAN(I,1)= 7 0.32832E+00-0.25671E+00
 I.WAN(I,1)= 8 0.34937E+00 0.10851E+00
********adaptive array weights (amp.,phase)*******
 I= 1 WANDB, WANDG=
                         -4.71576
                                       -62.16515
                          -1.13211
                                        17.25469
    2 WANDB, WANDG=
 I=
    3 WANDB, WANDG=
                           0.00000
                                       -38.02192
 T=
                           -1.13210
                                        17.25474
      4 WANDB, WANDG=
 I=
                                       -62.16529
     5 WANDB, WANDG=
                           -4.71576
 I=
                                        17.25474
     6 WANDB, WANDG=
                           -1.13210
 I=
                                        -38.02193
   7 WANDB, WANDG=
                           0.00000
 T=
    8 WANDB, WANDG=
                           -1.13210
                                        17.25469
********cancellation************
 CPROD1, CPROD2= 0.12271E+01 0.12185E-13 0.10000E+01 0.00000E+00
 INR NORMALIZATION PARAMETER, INRNOR= 1
```

```
INR= QUI, ADAP .CANCEL=
                              31.411 0.889 -30.522 DB
                                      1
CANCDB, NRAN, AVECAN= -30.52207
                                             -30.52207
***Input data file for adaptive nulling with two auxiliary probes.
***sdipjamhyper data file, filename sdipjamhyper.datacirconicrr3
ADIPOLE NCOLX=8,NROWY=1,HZIN=0.001,HLIN=2.6,
  DXIN=6.888,DYIN=6.888,ARADIN=0.0039,
  ICIRC=1, CRADIN=11.81, HLSIN=0.25,
  IWL=0,FCHZ=120.0E6,BWFHZ=1.0E0,NFREQ=5,IWR=0,
  MPDX=0, NXDUM=0, NYDUM=0, ER3=73.5, SIG3=0.5, TD3=-1.0,
  ZLOAD=50.0, ZCHAR=0.0, NGEN=8, IGEN=0,
  NSCANS=1,
  THSINC=5.0
  IMUT=1, IBLTSL=0,
  IENORM=1,ICHEB=0,SLLDB=20.,EDGTDB=0.,
  IPATRN=1, IPRCOM=1, IANGLP=0,
  NPHCT=0,NTHPT=499,
  THDR=180., THDMIN=-90.,
  RLSYIN=0.0, RLSZIN=0.0,
  NCOLXN=117, NROWYN=1, NCOLZN=1,
  INEAR=1,
  IPOL=2.IGRNDP=0.
  ITLTPR=0, ITLTDP=0,
  NFCOLX=1,NFROWY=1,
  IANTX=1, IANTY=0, NPOWER=0,
  IPCONN=1, IPCONF=0, IPCUTF=0, IPCFX=0, IPCFY=0, IPCFZ=0,
  ITEK=0,
  IQUAN=0, IRNERR=0, ELERDB=0.02, ELERDG=0.2, NBMOD=12,
  NRAN=1.NBADWT=32.AWERDB=0.0.AWERDG=0.0.
  NAUX=7,IAUXA(1)=1,2,3,4,5,6,7,
  IATTEN=1.
  INRNOR=1.
  NJAMS=7, ISLC=0, AUXADB(1) +8+0.0,
  PWRJDB(1)=40.,40.,5*-99.0,
 YJAMIN(1)=7+0.0.
  XFDCIN=0.0.ZFOCIN=0.0.
  XJAMIN(1)=5.9,-5.9,2*0.0,2*-4.0,1*0.0,
  ZJAMIN(1)=0.0,0.0,4.0,-4.0,3.0,-3.0,1*0.0
  YNIN=0.0,
  NCOLXN=101,NTHPT=101,
 XNIN=-15., ZNIN=-15., RLSXIN=30., RLSZIN=30., NCOLXN=41, NCOLZN=41,
AEND
***Output data file for adadptive nulling with two auxiliary probes.
***sdipjamhyper data file, filename sdipjamhyper.datacirconicrr3
XJAMIN(1), YJAMIN(1), ZJAMIN(1)=
                                                   0.000
FGHZ=
            0.1200000
DI, DY, HL, ARAD=
                         0.1749€
                                        0.17496
                                                      0.06604
                                                                      0.00010
IV,VTA= 1 0.5916E-02 0.2211E-02
IV,VTA= 2 0.5916E-02 0.2211E-02
IV.VTA= 3
              0.5916E-02 0.2211E-02
IV,VTA= 4
              0.5916E-02 0.2211E-02
IV,VTA= 5
               0.5916E-02 0.2211E-02
IV,VTA= 6
IV,VTA= 7
                0.5916E-02 0.22 1E-02
                0.5916E-02 0.2211E-02
IV,VTA= 8
               0.5916E-02 0.2211E-02
******** array mutual impedance matrix (first row)*******
Z(1, 1)= 0.91705E+02 0.12403E+02
Z(1, 2)=-0.33552E-01 0.19136E+00
Z(1, 3)=-0.14488E-01-0.54125E-02
```

```
Z(1, 4)= 0.29420E-02 0.13780E-02
Z(1,
       5) = 0.15283E-02-0.11371E-02
Z(1, 6)= 0.29420E-02 0.13780E-02
Z(1, 7)=-0.14488E-01-0.54125E-02
Z(1, 8)=-0.33552E-01 0.19136E+00
NEL, NCOLI, NROWY = 8
                          8
CURRENTS
    1 1.000000
                          0.683E-04
                                           13.
                                           13.
                          0.683E-04
     2 1.000000
     3 1.000000
                          0.683E-04
                                           13.
                                           13.
                          0.683E-04
     4 1.000000
                          0.683E-04
                                           13.
     5 1.000000
                          0.683E-04
                                           13.
     6 1.000000
                          0.683E-04
                                           13.
     7 1.000000
                                           13.
     8 1.000000
                          0.683E-04
IC,RWTA= 1 0.35355E+00
IC,RWTA= 2 0.35355E+00
IC,RWTA= 3 0.35355E+00
IC,RWTA= 4 0.35355E+00
IC,RWTA= 5 0.35355E+00
IC,RWTA= 6 0.35355E+00
IC,RWTA= 7 0.35355E+00
IC,RWTA= 8 0.35355E+00
I,CWTA(I,1)= 1 0.34510E+00-0.76830E-01
I,CWTA(I,1)= 2 0.34510E+00-0.76830E-01
I,CWTA(I,1)= 3 0.34510E+00-0.76830E-01
I,CWTA(I,1)= 4 0.34510E+00-0.76830E-01
I,CWTA(I,1)= 5 0.34510E+00-0.76830E-01
I,CWTA(I,1)= 6 0.34510E+00-0.76830E-01
I.CWTA(I,1) = 7 0.34510E+00-0.76830E-01
I,CWTA(I,1)= 8 0.34510E+00-0.76830E-01
*****ring array weights before nulling (amp,phase)**********
I= 1 CWTADB, CWTADG= 0.00000
                                           -12.55098
I= 2 CWTADB,CWTADG=
                               0.00000
                                              -12.55099
                               0.00000
                                              -12.55098
I=
    3 CWTADB, CWTADG=
    4 CWTADB, CWTADG=
                                0.00000
                                              -12.55097
1=
                                              -12.55098
                                0.00000
I=
     5 CWTADB, CWTADG=
                                              -12.55097
I= 6 CWTADB,CWTADG=
                                0.00000
                                 0.00000
                                              -12.55098
    7 CWTADB, CWTADG=
I= 8 CWTADB, CWTADG=
                                 0.00000
                                              -12.55097
Z(1, 1) = 0.91705E+02 0.12403E+02
Z(1,
        2)=-0.33552E-01 0.19136E+00
Z(1, 3)=-0.14488E-01-0.54125E-02
Z(1, 4)= 0.29420E-02 0.13780E-02
Z(1, 5)= 0.15283E-02-0.11371E-02
        6)= 0.29420E-02 0.13780E-02
Z(1.
Z(1, 7)=-0.14488E-01-0.54125E-02
        8)=-0.33552E-01 0.19136E+00
Z(1,
NEL , NCOLX , NROWY = 8
                          8 1
                     1 0.10000E+05
II,PWRJ IN POWER=
II,PWRJ IN POWER= 2 0.10000E+05
II.PWRJ IN POWER= 3 0.12589E-09
II,PWRJ IN POWER* 4 0.12589E-09
II,PWRJ IN POWER* 5 0.12589E-09
II,PWRJ IN POWER# 6 0.12589E-09
II,PWRJ IN POWER* 7 0.12589E-09
******COVariance matrix******************
                             36.35
31.68
I,J,CNDB,PHASEN= 1 1
I,J,CNDB,PHASEN= 1 2
                                                  0.00
                                                  91.33
```

I, J, CNDB, PHASEN=	1	3	24.98	-117.94
I, J, CNDB, PHASEN=	1	4	20.88	6.28
I, J, CNDB, PHASEN=	1	5	20.72	-0.14
I, J, CNDB, PHASEN=	1	6	20.88	6.28
I, J, CNDB, PHASEN=	1	7	24.98	-117. <del>94</del>
I, J, CNDB, PHASEN=	1	8	31.68	91.33
I, J, CNDB, PHASEN=	2	1	31.68	-91.33
I, J, CNDB, PHASEN=	2	2	27.05	0.00
I, J, CNDB, PHASEN=	2	3	20.24	155.47
I, J, CNDB, PHASEN	2	4	8.24	177 .97
I, J, CNDB, PHASEN=	2	5	20.90	-6.17
I, J, CNDB, PHASEN=	2	6	8.24	177.97
I, J, CNDB, PHASEN=	2	7	20.24	155.47
I, J, CNDB, PHASEN=	2	8	27.04	0.00
I, J, CNDB, PHASEN=	3	1	24.98	117.94
I, J, CNDB, PHASEN=	3	2	20.24	-155.47
I, J, CNDB, PHASEN=	3	3	16.62	0.00
I, J, CNDB, PHASEN=	3	4	20.27	-155.41
I, J, CNDB, PHASEN=	3	5	25.01	117.93
I, J, CNDB, PHASEN=	3	6	20.27	-155.41
I, J, CNDB, PHASEN=	3	7	16.52	0.00
I, J, CNDB, PHASEN=	3	8	20.24	-155.47
I.J.CNDB.PHASEN=	4	1	20.88	-6.28
I, J, CNDB, PHASEN=	4	2	8.24	-177.97
I, J, CNDB, PHASEN=	4	3	20.27	155.41
I, J, CNDB, PHASEN=	4	4	27.08	0.00
I, J, CNDB, PHASEN=	4	5	31.71	-91.33
I, J, CNDB, PHASEN=	4	6	27.07	0.00
I, J, CNDB, PHASEN=	4	7	20.27	155.41
I, J, CNDB, PHASEN=	4	8	8.24	-177.97
I, J, CNDB, PHASEN=	5	1	20.72	0.14
I, J, CNDB, PHASEN=	5	2	20.90	6.17
I, J, CNDB, PHASEN=	5	3	25.01	-117.93
I.J.CNDB.PHASEN=	5	4	31.71	91.33
I, J, CNDB, PHASEN=	5	5	36.38	0.00
I, J, CNDB, PHASEN=	5	6	31.71	91.33
I, J, CNDB, PHASEN=	5	7	25.01	-117.93
I, J, CNDB, PHASEN=	5	8	20.90	6.17
I, J, CNDB, PHASEN=	6	1	20.88	-6.28
I, J, CNDB, PHASEN=	6	2	8.24	-177.97
I, J, CNDB, PHASEN=	6	3	20.27	155.41
I, J, CNDB, PHASEN=	6	4	27.07	0.00
I, J, CNDB, PHASEN=	6	5	31.71	-91.33
I, J, CNDB, PHASEN=	6	6	27.08	0.00
I, J, CNDB, PHASEN=	6	7	20.27	155.41
I, J, CNDB, PHASEN=	6	8	8.24	-177.97
I, J, CNDB, PHASEN=	7	1	24.98	117.94
I, J, CNDB, PHASEN=	7	2	20.24	-155.47
I, J, CNDE, PHASEM-	7	3	16.52	0.00
I, J, CNDB, PHASEN=	7	4	20.27	-155.41
I, J, CNDB, PHASEN=	7	5	25.01	117.93
I, J, CNDB, PHASEN=	7	6	20.27	-155.41
I, J, CNDB, PHASEN=	7	7	16.62	0.00
I, J, CNDB, PHASEN=	7	.3	20.24	-155.47
I, J, CNDB, PHASEN=	8	1	31.68	-91.33
I, J, CNDB, PHASEN=	8	2	27.04	0.00
I, J, CNDB, PHASEN=	8	3	20.24	155.47
I, J, CNDB, PHASEN=	8	4	8.24	177.97
I, J, CNDB, PHASEN=	8	5	20.90	-6.17

```
177.97
                                8.24
I,J,CNDB,PHASEN= 8 6
I, J, CNDB, PHASEN= 8 7
                                20.24
                                          155.47
I,J,CNDB,PHASEN= 8 8
                               27.05
                                            0.00
+>*********igenvalues**************
I.EVLDBN= 1
                 37.195
                  37.437
I,EVLDBN= 2
                     0.000
I,EVLDBN= 3
I,EVLDBN= 4
                    0.000
                    0.000
I EVLDBN= 5
                     0.000
I,EVLDBN= 6
I .EVLDBN=
           7
                     0.000
I,EVLDBN= 8
                     0.000
 CPROD1,CPROD2= 0.13835E+04-0.56843E-13 0.10000E+01 0.00000E+00
 INR NORMALIZATION PARAMETER, INRNOR= 1
                   -14.003
I= 1 WANDB=
     2 WANDB=
                    -8.939
 I=
 I= 3 WANDB=
                    -6.886
                    -8.939
 I= 4 WANDB=
 I= 5 WANDB=
                   -14.003
 I= 6 WANDB=
                    -8.939
                    -6.886
 I=
     7 WANDB=
 I= 8 WANDB=
                    -8.939
 I,WAN(I,1)= 1 0.48427E-01-0.19349E+00
I,WAN(I,1)= 2 0.35635E+00 0.26440E-01
I,WAN(I,1)= 3 0.43481E+00-0.12564E+00
I,WAN(I,1)= 4 0.35635E+00 0.26441E-01
 I.WAN(I,1)= 5 0.48427E-01-0.19349E+00
 I,WAN(I,1)= 6 0.35635E+00 0.26441E-01
I,WAN(I,1)= 7 0.43481E+00-0.12564E+00
 I,WAN(I,1)= 8 0.35635E+00 0.26441E-01
********adaptive array weights (amp.,phase)*******
                                    -75.94859
                         -7.11720
 I= 1 WANDB, WANDG=
 I= 2 WANDB, WANDG=
I= 3 WANDB, WANDG=
                          -2.05298
                                        4.24348
                                    -16.11658
                          0.00000
 I= 4 WANDB, WANDG=
                          -2.05297
                                        4.24352
                          -7.11719
 I= 5 WANDB, WANDG=
                                      ~75.94872
 I= 6 WANDB, WANDG=
                         -2.05298
     7 WANDB, WANDG=
                          0.00000
                                       -16.11658
 Ī=
 I= 8 WANDB, WANDG=
                          -2.05297
                                        4.24348
*********Cancellation***********
 CPROD1,CPROD2= 0.10001E+01-0.63768E-14 0.10000E+01 0.00000E+00
 INR NORMALIZATION PARAMETER, INRNOR= 1
                                     0.000 -31.410 DB
 INR= QUI, ADAP , CANCEL= 31.410
 CANCDB, NRAN, AVECAN= -31.40966 1 -31.40966
****file sdipjamhyperMake*****
EXECUTABLE = sdipjamhyper.out
OBJS = sdipjamhyper.o dzabgnloss.o \
        dpack2.o zabgenloss.o pack2.o fwgh.o \
        chebaf.o wwquan.o taylor.o \
        eigenv.o reordr.o ydipsubloss.o \
        plothyper.o contek.o plabel.o \
        conturek.o circsubloss.o
# use tabs for continues
IL = /usr/lib/f68881/libm.il
FLAGS = -01 -v -f68881
$(EXECUTABLE): $(OBJS)
        dislink -1F77 $(FLAGS) -o $(EXECUTABLE) $(OBJS) $(IL)
```

```
sdipjamhyper.o: /home/ajf/hyperthermia/sdipjamhyper.f
       f77 $(FLAGS) -c /home/ajf/hyperthermia/sdipjamhyper.f $(IL)
plothyper.o: /home/ajf/hyperthermia/plothyper.f
       f77 $(FLAGS) -c /home/ajf/hyperthermia/plothyper.f $(IL)
dzabgnloss.o: /home/ajf/hyperthermia/dzabgnloss.f
       f77 $(FLAGS) -c /home/ajf/hyperthermia/dzabgnloss.f $(IL)
dpack2.o: /home/ajf/monjam/dpack2.f
       f77 $(FLAGS) -c /home/ajf/monjam/dpack2.f $(IL)
contek.o: /home/ajf/monjam/contek.f
       f77 $(FLAGS) -c /home/ajf/monjam/contek.f $(IL)
plabel.o: /home/ajf/monjtr/plabel.f
       f77 $(FLAGS) -c /home/ajf/monjtr/plabel.f $(IL)
conturek.o: /home/ajf/hyperthermia/conturek.f
       f77 $(FLAGS) -c /home/ajf/hyperthermia/conturek.f $(IL)
zabgenloss.o: /home/ajf/hyperthermia/zabgenloss.f
       f77 $(FLAGS) -c /home/ajf/hyperthermia/zabgenloss.f $(IL)
pack2.o: /home/ajf/monjam/pack2.f
       f77 $(FLAGS) -c /home/ajf/monjam/pack2.f $(IL)
fwgh.o: /home/ajf/monjam/fwgh.f
       f77 $(FLAGS) -c /home/ajf/monjam/fwgh.i $(IL)
chebaf.o: /home/ajf/monjam/chebaf.f
       f77 $(FLAGS) -c /home/ajf/monjam/chebaf.f $(IL)
wwquan.o: /home/ajf/monjam/wwquan.f
       f77 $(FLAGS) -c /home/ajf/monjam/wwquan.f $(IL)
taylor.o: /home/ajf/monjam/taylor.f
       f77 $(FLAGS) -c /home/ajf/monjam/taylor.f $(IL)
eigenv.o: /home/ajf/monjam/eigenv.f
       f77 $(FLAGS) -c /home/ajf/monjam/eigenv.f $(IL)
reordr.o: /home/ajf/monjam/reordr.f
       f77 $(FLAGS) -c /home/ajf/monjam/reordr.f $(IL)
ydipsubloss.o: /home/ajf/hyperthermia/ydipsubloss.f
       f77 $(FLAGS) -c /home/ajf/hyperthermia/ydipsubloss.f $(IL)
circsubloss.o: /home/ajf/hyperthermia/circsubloss.f
       f77 $(FLAGS) -c /home/ajf/hyperthermia/circsubloss.f $(IL)
© M.I.T. LINCOLN LABORATORY 1991, ALL RIGHTS RESERVED
C***PROGRAM SDIPJAMHYPER.F --- ANALYZES FINITE ARRAYS OF DIPOLES
C***IN LOSSY DIELECTRIC OR FREE SPACE.
C...THE DIPOLES ARE ASSUMED TO BE ORIENTED PARALLEL
C***TO THE PLANE OF THE GRID FOR A PLANAR ARRAY, OR THEY CAN BE
C *** ARRANGED IN AN ANNULAR (RING) ARRAY CONFIGURATION.
C***RECEIVING CONDITIONS ARE ASSUMED.
C***DOUBLE-PRECISION VERSION
     PARAMETER (NUMCHN=8)
     PARAMETER (NUMAUX=8)
```

```
PARAMETER (NUMJAM=7)
      PARAMETER (NUMELM=8)
      PARAMETER (NUMFRQ=5)
      PARAMETER (NUMNPT=1681)
      COMPLEX PS(1825), CZ(900), VA(NUMELM), Z(NUMELM, NUMELM)
      COMPLEX VTA(NUMELM), VREFA, VRECVA(NUMELM)
      COMPLEX VCW(NUMELM), VRECVX(NUMELM), CJ, CSUMA, VXA
CC***THE ABOVE MATRICES ARE DIMENSIONED BY THE NUMBER OF ELEMENTS
      COMPLEX *16 COVNF(NUMCHN, NUMCHN), COVNFI(NUMCHN, NUMCHN)
      COMPLEX *16 COVAAJ(NUMCHN, NUMCHN), CINVCN(NUMCHN, NUMCHN)
      COMPLEX *16 VCHA(NUMCHN, NUMFRQ)
      COMPLEX *16 VMAJMA(NUMJAM, NUMFRQ)
      COMPLEX *16 VAUXJA(NUMAUX, NUMJAM, NUMFRQ)
      COMPLEX *8 VAXCWA(NUMAUX, NUMNPT)
C *** NOTE DIMENS. VCHA (NCHAN, NFREQ), VMAJMA (NJAM, NFREQ)
                                                                           DIP00200
                 VAUXJA (NAUX, NJAM, NFREQ), VAXCWA (NAUX, NPTS)
С
                                                                           DIP00210
      COMPLEX *8 VXMANA(NUMNPT), ETASP, GAMSP
      COMPLEX *16 EIGVAN(NUMCHN), EIGVEN(NUMCHN, NUMCHN)
      COMPLEX +16 WANA (NUMCHN, 1)
      COMPLEX *16 WTCTR(1,NUMCHN), WQSLC(NUMCHN,1)
      COMPLEX *16 WQNA(NUMELM,1), CMPROD(1, NUMCHN), WAN(NUMCHN,1)
      COMPLEX *16 WANDMA(NUMELM.1)
      COMPLEX +16 CWTA(NUMELM,1), ETADP, GAMDP
      DIMENSION ACALPH(NUMELM), RWTA(NUMELM), RWTB(NUMELM)
      DIMENSION VXADB(NUMNPT), VXAPH(NUMNPT)
      DIMENSION XC(NUMNPT), YC(NUMNPT), ZC(NUMNPT)
      DIMENSION THSD(1), PHSD(1), THCT(1), PHCT(1)
                                                                           DIP00320
      DIMENSION XJAMIN(NUMJAM), YJAMIN(NUMJAM), ZJAMIN(NUMJAM)
      DIMENSION XJAM(NUMJAM), YJAM(NUMJAM), ZJAM(NUMJAM)
      DIMENSION PWRJDB(NUMJAM)
                                                                           DIP00340
      DIMENSION PCHADB(NUMNPT)
      REAL *8 CNDB, PHASEN, DPDCR, WKE(144), EVLDBN
                                                                           DIP00370
      REAL *8 PREALN, PIMAGN
                                                                           DIP00380
      REAL *8 QINRDB, AINRDB, CANCDB, EXCPCB
                                                                           DIPO0390
                                                                           DIP00400
      REAL *8 DINRAA.AINRAA
      REAL *8 CNCLNA, ELSGDB, ELSGDG, AWSGDB, AWSGDG
                                                                           DIP00410
      REAL +8 SUMC, SUMAA, AVECAN, AVCAA, SAVE, SAVESR
      REAL *4 RX(100),RY(100),RZ(100),FM(41,41),DMF(41,41)
      REAL +4 AUXADB(NUMAUX), VATENA(NUMAUX)
      INTEGER IAUXA(NUMAUX), LTMP(NUMCHN), MTMP(NUMCHN)
      CHARACTER DATNAM+35, OUTNAM+35
      COMMON/A/DX,DY,NCOLX,NROWY,NEL,HZ,HL,ARAD,ZLOAD,ZCHAR
                                                                           DIP00460
      COMMON /B/ NGEN, IGEN, THETAS, PHIS, IMUT, IBLTSL, IPATRN
                                                                           DIP00470
      COMMON /D/ FGHZ, RLAMDA, IWL, IS, NSCANS
                                                                           DIP00480
      COMMON /DEG/ THDR, THDMIN
                                                                           DIP00490
      COMMON /PLT/ INEAR, IPLOTM
                                                                           DIP00500
      COMMON/XLYL/XL,YL
                                                                           DIP00510
      COMMON/NEAR/IANTX, IANTY, NPOWER, IPCONN, IPCONF, IPCUTF
                                                                           DIP00520
      COMMON/NEAR2/IPCFX, IPCFY, IPCFZ
                                                                           DIPO0530
      COMMON /P/XN,YN,ZN,RLSX,RLSY,RLSZ,NCOLXN,NROWYN,NCOLZN
                                                                           DIP00540
      COMMON /P2/EDGET, ICOMB, PUNFLX, PUNFLY
                                                                           DIP00550
      COMMON /CHEBY/ ICHEB, SLLDB
                                                                           DIPO0560
      COMMON /GROUND/ IGRNDP
                                                                           D_P00570
      COMMON /NORM/ IENORM, BIGNDB
                                                                           DIPC0580
      COMMON /NORMAL/ INRNOR
      COMMON/PCENTR/PCDXIN
                                                                           DIP00590
      COMMON /WRITE/ IWR
                                                                           DIP00600
      COMMON /SPLOSS/ ETASP, GAMSP
      COMMON /DPLOSS/ ETADP, GAMDP
```

```
COMMON /F/ FHZ, ER3, SIG3, TD3
      COMMON /CIRCLE/ ICIRC, RADIUS, HLS
      NAMELIST/DIPOLE/NCOLX, NROWY, HZIN, HLIN, ARADIN, ZLOAD, NGEN, IGEN,
                                                                             DIP00610
     1DXIN, DYIN, IMUT, NSCANS, THSD, PHSD, THSINC, IBLTSL, IPATRN, ER3, SIG3, TD3,
     2NPHCT, PHCT, NTHPT, ZCHAR, IWL, FCH2, BWFHZ, NFREQ, THDR, THDMIN, ICIRC,
     3XNIN, YNIN, ZNIN, NCOLXN, NROWYN, NCOLZN, INEAR, IPRCOM, IANGLP, CRADIN,
     41ANTX, IANTY, NPOWER, NFCOLX, NFROWY, IPCONN, IPCONF, IPCUTF, IPOL, HLSIN,
     51PCF1, 1PCFY, 1PCF2, 1TEK, 1GRNDP, 1GAIN, 1CHEB, 1WR, 1ENORM, BIGNDB,
                                                                             DIP00660
     6NXDUM, NYDUM, NPDX, EDGTDB, RLSXIN, RLSYIN, RLSZIN, ISLC, INUNIF,
     7XFOCIN, ZFOCIN, XJAMIN, YJAMIN, ZJAMIN, IATTEN, AUXADB, INRNOR,
     8NJAMS, PWRCDB, PWRJDB, NAUX, IAUXA, IAUXB, ITLTPR, ITLTDP,
     91QUAN, IRNERR, ELERDB, ELERDG, AWERDB, AWERDG, NBMOD, NBADWT, NRAN, SLLDB DIPO0700
C ** NOTE; IF IWL=0 (INCHES), IWL=2 (METERS)
      WRITE(6,2959)
2959 FORMAT(1X, 'ENTER INPUT DATA FILE NAME (typ. sdipjamhyper.data)')
      READ(5,*) DATNAM
      OPEN(4,FILE=DATNAM,FORM='FORMATTED')
      WRITE(6.3959)
3959 FORMAT(1X, 'ENTER OUTPUT DATA FILE NAME (typ sdipjamhyper.output)')
      READ(5,*) OUTNAM
      OPEN(8,FILE=OUTNAM,FORM='FORMATTED')
      CALL GETCP2(CPU1)
      PI=3.141592654
                                                                             DIP00710
      DCR=PI/180.
                                                                             DTP00720
      DPDCR=DCR
                                                                             DIP00730
                                                                             DIP00740
      CJ=(0..1.)
      CINMTR=0.0254
      IGRNDP=1
                                                                             DIP00750
      IGAIN=0
                                                                             DIP00760
      ZLOAD=0.0
                                                                             DIP00770
      ICHEB=0
                                                                             DIP00780
      IWR=0
                                                                             DIP00790
      IENORM=1
                                                                             DIP00800
      EXCPCB=-1.0D0
                                                                             DIP00810
      TOUAN=0
                                                                             DIP00820
      NRAN=1
                                                                             DIP00830
      TILTPR=0.0
                                                                             DIP00840
      TILTDP=0.0
                                                                             DIP00850
      IATTEN=0
                                                                             DIP00860
      IPOL=1
                                                                             DIP00870
      ER3=1.0
      SIG3=0.0
      TD3=-1.0
      ICIRC=0
      CRADIN=0.0
      HLSIN=0.0
      ISLC=1
      INUNIF=0
      INRNOR=0
      READ(4,DIPOLE)
                                                                             DIP00880
      IF(NFREQ.EQ 1) BWFHZ=0.0
      IF(ISLC.EQ.0) WRITE(6,81100)
81100 FORMAT(1X, 'FULLY ADAPTIVE ARRAY')
      IF(ISLC.EQ.O.AND.ICHEB.EQ.O) INUNIF=1
      IF(ICIRC.EQ.1) WRITE(6,7000)
7000 FORMAT(1X, 'RING ARRAY GEOMETRY')
      FHZ=FCH2
      IXZ=0
```

```
TYY=0
      IF(NROWYN.EQ.1) IXZ=1
      IF(NCOLZN.EQ.1) IXY=1
                                                                         D1P00890
      IF(IPOL.EQ.2) IPRCOM=0
      IF(IPRCOM.EQ.O) WRITE(6,7898)
                                                                         DIP00900
                                                                         DIP00910
7898 FORMAT(11,'NO PROBE COMPENSATION')
                                                                         DIP00920
      WRITE(6,9276)IGRNDP
                                                                         DIP00930
9276 FORMAT(1X'GROUND PLANE PARAMETER, IGRNDP=', 14)
                                                                         DIPO0940
      DO 1615 IX-1, NAUX
CC
      IF(IATTEN.EQ.O) AUXADB(IX)=0.0
                                                                         DTP00950
CC
                                                                         DIP00960
      WRITE(6,2318)IX,AUXADB(IX)
C2318 FORMAT(1X,'IX,AUXADB(IX)=',I4,2X,F12.2)
                                                                         DIP00970
                                                                         DIP00980
     VATENA(IX)=10.++(AUXADB(IX)/20.)
                                                                         DIP00990
C1615 CONTINUE
                                                                         DIP01000
      IF(NJAMS.GT.1) WRITE(8,2009)XJAHIN(1),YJAMIN(1),ZJAHIN(1)
2009 FORMAT(1X,'XJAMIN(1),YJAMIN(1),ZJAMIN(1)=',3F12.3)
                                                                         DIP01010
                                                                         DIP01020
      WRITE(6,6987) BWFHZ,NFREQ
6987 FORMAT(1X,'BWFHZ,NFREQ=',E12.5,2X,I5)
                                                                         DIP01030
                                                                         DIP01040
      IF(ITEK.EQ.O)CALL COMPRS
                                                                         DIP01050
      IF(ITEK.EQ.1)CALL TEKALL(4014,480,0,1,0)
                                                                         DIP01060
      IF(INEAR EQ.O) CALL PRNTDA
CC
                                                                         DIP01070
      IF(IMUT.EQ.0) ZLOAD=1.0
                                                                         DIP01080
      IF(ITLTPR.EQ.1) TILTPR*45.
                                                                         DIP01090
      IF(ITLTDP.EQ.1) TILTDP=45.
      CTPR=COS(TILTPR+DCR)
                                                                         DIP01100
                                                                         DIP01110
      STPR=SIN(TILTPR*DCR)
                                                                         DIP01120
      CTDP=COS(TILTDP+DCR)
                                                                         DIP01130
      STDP=SIN(TILTDP+DCR)
                                                                         DIP01140
      FGHZ=FCHZ/1.0E9
                                                                         DIP01150
      NEL=NCOLX+NROWY
      RNEL=NEL
                                                                         DIP01160
      XLIN=DXIN=(NCOLX-1)
                                                                         DIP01170
      YLIN=DYIN+(NROV:-1)
                                                                         DIP01180
      NACOLX=NCOLX-NPDX-2=NXDUM
                                                                         DIP01190
      NAROWY-NROWY-2*NYDUM
      NAEL=NACOLX+NAROWY
                                                                         DIP01200
      PCDXIN=NPDX+DXIN
                                                                         DIP01220
      NAUXP1=NAUX+1
                                                                         DIP01230
      NAUXP2=NAUX+2
      IF(ISLC.EQ.1) NMAX=NAUXP1
      IF(ISLC.EQ.O) NMAX=NEL
                                                                         DIP01250
      NMAXP1=NMAX+1
                                                                         DIP01260
      ELSGDB=ELERDB+SQRT(3.)
                                                                         DIP01270
      ELSGDG=ELERDG+SQRT(3.)
                                                                         DIP01280
      AUSGDB=AWERDB+SQRT(3.)
                                                                         DIP01290
      AWSGDG=AWERDG+SQRT(3.)
                                                                         DIP01300
      INITRN=1
      IF(ISLC.EQ.O) THEN
        NAUX=NEL
        NAUXP1=NAUX+1
        NAUXP2=NAUX+2
        DO 79130 I=1.NEL
         IAUXA(I)=I
        WRITE(6,76767)I, IAUXA(I)
       FORMAT(1X,'I, IAUXA(I)=',215)
76767
79130 CONTINUE
      ENDIF
      DO 1615 IX=1,NAUX
      IF(IATTEN.EQ.O; AUXADB(IX)=0.0
```

```
WRITE(6,2318)IX, AUXADB(IX)
2318 FORMAT(1X,'IX,AUXADB(IX)=',I4,2X,F12.2)
      VATENA(IX)=10.**(AUXADB(IX)/20.)
1615 CONTINUE
                                                                        DIP01310
    CONTINUE
30
                                                                        DIP01320
     NR=1
                                                                        DIP01330
      IF(NGEN.EQ.1) NSCANS=0
                                                                        DIP01340
     ICC=NEL
                                                                        DIP01350
     ICC1=ICC
                                                                        DIP01360
      IF(IBLTSL.EQ.1) ICC1=NROWY
                                                                        DIP01370
     IF(IWL.EQ.1)G0 TO 50
CC
                                                                        DIP01380
        WRITE(8,40)FGHZ
                                                                        DIP01390
       FORMAT(/,1X,'FGHZ=',F15.7)
40
C****COMPUTE FREE SPACE LAMBDA******
C*****ALL UNITS HAVE TO BE IN METERS
     RLAMDA=2.997925E10/FCHZ/2.54
C***PARAMETER CONVERSIONS TO PROPER UNITS
                                                                         DIP01420
55 CONTINUE
      IF(IWL.EQ.2) GO TO 50
      DX=DXIN+CINMTR
      DY=DYIN+CINMTR
      HL=HLIN+CINMTR
      HLS=HLSIN+CINMTR
      ARAD=ARADIN+CINMTR
      HZ=HZIN+CINMTR
      RLSX=RLSXIN+CINMTR
      RLSY=RLSYIN+CINMTR
      RLSZ=RLSZIN+CINMTR
      IN=ININ+CINMTR
      YN=YNIN=CINMTR
      ZN=ZNIN+CINMTR
      XFOC=XFOCIN+CINMTR
      ZFOC=ZFOCIN+CINMTR
      XL=XLIN+CINMTR
      YL=YLIN+CINMTR
      PCDX=PCDXIN+CINMTR
      RADIUS=CRADIN+CINMTR
                                                                         DIP01580
      CONTINUE
50
                                                                         DIP01590
      WRITE(8,60)DX,DY,HL,ARAD
                                                                         DIP01600
      FORMAT(1X,'DX,DY,HL,ARAD=',2X,4F14.5)
 60
                                                                         DIP01610
                                                                         DIP01620
 C***COMPUTE CALIBRATION CONSTANTS (PHASE CNLY) TO
 C ** * MAXIMIZE GAIN (FOCUS ANTENNA) TO NEAR FIELD RANGE
                                                                         DIP01630
                                                                         DIP01640
 C***PHASE CENTER 'A' VOLTAGE EXCITATION
                                                                         DIP01650
                                                                         DIP01660
       XOA=-XL/2.+NXDUM*DX+(NACOLX-1)/2.*DX
                                                                         DIP01700
       YP=0.0
       IF(ICIRC.EQ.0) WRITE(6,24690)
 24690 FORMAT(1X, 'CALLING NFDPX2')
       TF(ICIRC.EQ.0)
      2CALL NFDPX2(CTPR,STPR,CTDP,STDP,XFOC,YP,ZFOC,XOA,O,VA,VREFA)
       IF(ICIRC.EQ.1)
      2CALL NFDPC2(XFOC, YP, ZFOC, XOA, O, VA, VREFA)
                                                                          DIP01720
 C***SAVE INCIDENT VOLTAGES
                                                                          DIP01730
       DO 65 IV=1, NEL
                                                                          DIP01740
       VTA(IV)=VA(IV)
                                                                          DIP01750
       WRITE(6,3757)IV,VTA(IV)
       WRITE(8,3757)IV,VTA(IV)
```

```
DIP01760
3757 FORMAT(1X,'IV, VTA=', I4,2X,2E12.4)
                                                                         DIP01770
      CONTINUE
                                                                          DIP01780
      WRITE(6,9876)VTA(2)
9876 FORMAT(1X,'VTA(2)=',2E12.4)
                                                                          DIP01790
                                                                          DIP01800
      MB=NCOLX
                                                                          DIP01810
      IDMB=NROWY
                                                                          DIP01820
      IDM=NEL
                                                                          DIP01830
      ICC=NEL
                                                                          DIP01840
      IBLT=IBLTSL
                                                                          DIP01850
      IDM1=IDM3
                                                                          DIP01860
      IF(IBLT.EQ.O) IDM1=NEL
                                                                          DIP01890
      I2=1
8889 CONTINUE
                                                                          DIP02350
      CALL ZMATRX(CTDP,STDP,IBLTSL,ICC1,ICC,CZ,Z)
C *** SOLVE SYSTEM OF EQUATIONS FOR THE UNKNOWN CURRENTS
                                                                          DIP02370
      IF(IBLT.EQ.1) GO TO 240
                                                                          DIP02380
      TSYM=0
                                                                          DIP02390
      I12=1
                                                                          DIP02400
      WRITE(6,9876)VTA(2)
                                                                          DIP02410
      WRITE(6.6110)
                                                                          DIP02420
6110 FORMAT(1X, 'CALL CROUT')
      CALL CROUT(Z, VA, ICC1, ICC, ISYM, IWR, I12, NEL)
                                                                          DIP02430
                                                                          DIP02440
      I12=2
      WRITE(6,9876)VTA(2)
                                                                          DIP02450
                                                                          DIP02460
      GO TO 255
                                                                          DIP02470
240 IENTRY=4
                                                                          DIP02480
      CALL BLTSOL(CZ, VA, PS, NCOLX, IDMB, IENTRY)
                                                                          DIP02490
      IENTRY=3
255 WRITE(8,270)
                                                                          DIP02500
270 FORMAT(1X, 'CURRENTS')
                                                                          DIP02510
      IF(NEL.LT.40) CALL CNORM(VA,NEL)
                                                                          DIP02520
                                                                          DIP02530
      WRITE(6,280)
280 FORMAT(1X, 'AFTER CURRENTS SOLUTION')
                                                                          DIP02540
                                                                          DIP02550
C+++VA ARE CURRENTS (AMPERES) NOW
                                                                          DIP02560
C+
                                                                          DIP02570
C***COMPUTE RECEIVED VOLTAGES
                                                                          DIP02580
      DO 285 IC=1.NEL
                                                                          DIP02590
      WRITE(6,7531)ZLOAD
                                                                          DIP02600
7531 FORMAT(1X, 'ZLOAD=',F12.5)
      WRITE(6,1134)IC, VTA(IC)
                                                                          DIP02610
                                                                          DIP02620
      VRECVA(IC)=VA(IC)+ZLOAD
      WRITE(6,1174)IC, VRECVA(IC)
                                                                          DIP02630
                                                                          DIP02640
1134 FORMAT(1X, 'IC, VTA=', I4, 2X, 2E12.4)
1174 FORMAT(1X, 'IC, VRECVA=', I4, 2X, 2E12.4)
                                                                          DIP02650
                                                                          DIP02660
      WRITE(6,1135)IC, VA(IC)
                                                                          DIPO2670
1135 FORMAT(1X.'IC.VA=', I4,2X,2E12.4)
                                                                          DIP02680
      VRADB=20. *ALOG10(CABS(VRECVA(IC)))
      VRAPH=ATAN2(AIMAG(VRECVA(IC)), REAL(VRECVA(IC)))/DCR
                                                                          DIP02690
                                                                          DIP02700
C***COMPUTE CALIBRATION CONSTANTS (PHASE ONLY)
                                                                          DIP02710
      ACALPH(IC) = - VRAPH
                                                                          DIP02720
285
     CONTINUE
                                                                          DIP02730
      DO 7777 IC=1,NEL
                                                                          DIP02740
      WRITE(6,6667)IC, ACALPH(IC)
                                                                          DIP02750
      WRITE(8,6667)IC, ACALPH(IC)
                                                                          DIP02760
6667 FORMAT(1X,'IC, ACALPH=', I4, 2X, F10.2,' DEGS')
                                                                          DIP02770
7777 CONTINUE
C***COMPUTE NEAR FIELD PATTERN OF FOCUSED ARRAY
                                                                          DIP02780
                                                                          DIP02790
C***COMPUTE BEAMFORMER WEIGHTS (I.E. TAPER)
                                                                          DIP02800
```

```
IF(INUNIF.EQ.0)
                                                                           DIP02810
     2CALL VRCVWT (NACOLX, NAROWY, NXDUM, NYDUM, EDGTDB, RWTA, RWTB)
     DO 1199 IC=1.NET.
                                                                          DIP02820
      IF(INUNIF.EQ.1.AND.ISLC.EQ.0) RWTA(IC)=1./SQRT(RNEL)
     WRITE(6,5111)IC,RWTA(IC)
    WRITE(8,5111)IC,RWTA(IC)
5111 FORMAT(1X, 'IC, RWTA=', I4, 2X, E12.5)
1199 CONTINUE
                                                                          DIP02850
     DO 5333 KC=1,NEL
                                                                          DIP02860
     CWTA(KC,1)=RWTA(KC)+CEXP(CJ+ACALPH(KC)+DCR)
                                                                          DIP02870
5333 CONTINUE
                                                                          DIP02880
                                                                          DIP02890
     IF(IQUAN.ED.1)
     2CALL ADQUAN(NEL, CWTA, NBMOD, IRNERR, INITRN, ELSGDB, ELSGDG)
                                                                          DIP02900
     BIGWDB=-299.0
     DO 91020 I=1, NEL
     WRITE(6,61910) I, CWTA(I,1)
     WRITE(8,61910) I, CWTA(I,1)
61910 FORMAT(1X, 'I, CWTA(I,1)=', I4,2X,2E12.5)
     CWTADB=20.*DLDG10(CDABS(CWTA(I,1)))
      IF(CWTADB.GT.BIGWDB) BIGWDB=CWTADB
91020 CONTINUE
     DO 91120 I=1.NEL
      CWTADB=20.*DLOG10(CDABS(CWTA(I,1)))-BIGWDB
      CWTADG=DATAN2(DIMAG(CWTA(I,1)),DREAL(CWTA(I,1)))/DCR
     WRITE(6,35990)I,CWTADB,CWTADG
     WRITE(8,35990)I,CWTADB,CWTADG
35990 FORMAT(1X, 'I=', I4, 2X, 'CWTADB, CWTADG=', 1X, 2F14.5)
91120 CONTINUE
     INITRN=2
                                                                          DIP02910
C***PERFORM NEAR FIELD SCAN WITH CW RADIATING DIPOLE
                                                                          DIP02920
      IF(INEAR.EQ.O) GO TO 390
                                                                          DIP02930
     WRITE(6,350)
                                                                          DIP02940
      WRITE(6,360)XNIN, YNIN, ZNIN, RLSXIN, RLSXIN, RLSZIN, NCOLXN, NROWYN,
                                                                          DIP02950
                                                                          DTP02960
    2NCOLZN
     WRITE(6.320)IWL
                                                                          DIP02970
320 FORMAT(1X,'IWL=',14)
                                                                          DIP02980
                                                                          DIP02990
     WRITE(6.330)
   FORMAT(1X,'CHANGE NEAR FIELD SCAN PARAMETERS?, ICHANG=1')
                                                                          DIP03000
                                                                          DIP03010
     READ(5.*.END=370)ICHANG
      IF(ICHANG.EQ.0) GO TO 370
                                                                          DIP03020
340 WRITE(6,350)
                                                                          DIP03030
    FORMAT(11, 'XN, YN, ZN, RLSXIN, RLSYIN, TISZIN, NCOLIN, NROWYN, NCOLIN=') DIPO3040
350
     READ(5,*)XNIN, YNIN, ZNIN, RLSXIN, RLSXIN, RLSZIN, NCOLXN, NROWYN, NCOLZN DIPO3050
     WRITE(6,360) XNIN, YNIN, ZNIN, RLSXIN, RLSYIN, RLSZIN, NCOLXN, DWYN,
                                                                          DIP03060
     2NCOLZN
                                                                          DIP03070
                                                                          DIP03080
360 FORMAT(1X,6F10.3,2X,3I5)
     IF(IWL.EQ.2) GO TO 370
                                                                          DIP03090
     IN=ININ+CINMTR
      YN=YNIN+CINMTR
     ZN=ZNIN+CINMTR
     RLSX=RLSXIN+CINMTR
      RLSY=RLSYIN+CINMTR
     RLSZ=RLSZIN*CINMTR
370 CONTINUE
                                                                          DIP03160
C**ALL DIMENSIONS IN METERS
     NELN=NCOLXN+NROWYN
                                                                          DIP03180
                                                                          DIP03190
      NPTSN=NELN+NCGLZN
C***SET DEFAULT VALUES FOR DXN.DYN.DZN
                                                                          DIP03200
     DXN=0.0
```

```
DTP03210
      DYN=0.1
      DZN=0.0
                                                                          DIP03220
                                                                         DTP03230
      IF(NCOLIN.GT.1)DIN=RLSI/(NCOLIN-1)
      IF (NROWYN.GT.1) DYN=RLSY/(NROWYN-1)
                                                                         DIP03240
      IF(NCOLZN.GT.1)DZN=RLSZ/(NCOLZN-1)
                                                                         DIP03250
                                                                         DIP03260
     TC=0
      BIGNDB=-299.0
                                                                         DIP03270
                                                                         DIP03280
      DO 3000 ICOLZN=1,NCOLZN
      ZPOS=ZN+DZN+(ICOLZN-1)
                                                                          DIP03290
      DO 3000 IROWYN=1 NROWYN
                                                                         DIP03300
                                                                         DIP03310
      Y=YN+DYN+(IROWYN-1)
      DO 3000 ICOLXN=1,NCOLXN
                                                                         DIP03320
      CALL GETCP2(CPU2)
      CPUSUB=CPU2-CPU1
      WRITE(6,7319)CPUSUB
7319 FORMAT(1X, 'CPU SUBTOTAL=',F14.2)
                                                                          DIP03330
      IC=IC+1
                                                                          DIP03340
      I=IN+DIN+(ICOLIN-1)
                                                                          DIP03350
      IC(IC)=X
      YC(IC)=Y
                                                                          DIP03360
      ZC(IC)=ZPOS
                                                                          DIPO3370
      WRITE(6,6969) IC, XC(IC), YC(IC), ZC(IC)
                                                                         DIP03410
6969 FORMAT(1X,'IC,XC,YC,ZC=',I4,2X,3F12.3)
                                                                         DIP03420
      IF(ICIRC.EQ.0)
     2CALL NFDPX2(CTPR,STPR,CTDP,STDP,X,Y,ZPOS,0.0,0,VCW,VREFCW)
                                                                          DIP03430
     IF(ICIRC.ED.1)
     2CALL NFDPC2(X,Y,ZPOS,0.0,0,VCW,VREFCW)
     IF(IMUT.EQ.0) GO TO 4255
                                                                          DIP03440
C***SOLVE EACH SYSTEM OF EQUATIONS FOR THE UNKNOWN CURRENTS
                                                                          DIP03450
      IF(IBLT.EQ.1) GO TO 4040
                                                                          DIP03460
                                                                          DIP03470
      I12=2
      WRITE(6,6110)
                                                                          DIP03480
     CALL CROUT(Z, VCW, ICC1, ICC, ISYM, IWR, I12, NEL)
                                                                          DIP03490
     GD TO 4255
                                                                          DIP03500
                                                                         DTP03510
4040 IENTRY=3
      CALL BLTSOL(CZ, VCW, PS, NCOLX, IDMB, IENTRY)
                                                                          DIP03520
4255 CONTINUE
                                                                          DIP03530
C***COMPUTE RECEIVED VOLTAGES FOR PRESENT SCAN
                                                                          DIP03540
     DO 3285 IIIC=1,NEL
                                                                          DIP03550
      VRECVX(IIIC)=VCW(IIIC)+ZLOAD
                                                                          DIP03560
3285 CONTINUE
                                                                          DIPO3570
C***STORE AUX. CHANNEL VOLTAGES
                                                                          DIP03580
     DO 1681 IAX=1,NAUX
                                                                          DIP03590
C ** MODIFICATION TO INCLUDE AUX. ATTEN.
                                                                          DIP03600
     VAXCWA(IAX,IC)=VRECVX(IAUXA(IAX)) *VATENA(IAX)
                                                                          DIP03610
      WRITE(6,3231) IAX, IC, XC(IC), VAXCWA(IAX, IC)
                                                                          DIP03620
3231 FORMAT(1X, 'IAX, IC, XC, VAXCWA(IAX, IC)=', 214, 2X, F12.3, 2X, 2E12.5)
                                                                          DIP03630
1681 CONTINUE
                                                                          DIP03640
C+++PERFORM BEAM FORMATION
                                                                          DIP03650
      CSUMA=(0.,0.)
                                                                          DIP03660
     DD 5444 KC=1.NEL
                                                                          DIP03670
      CSUMA=CSUMA+VRECVX(KC)+CWTA(KC,1)
                                                                          DIP03680
5444 CONTINUE
                                                                          DIP03690
      VXA=CSUMA
                                                                          DIP03700
      IF(CABS(VXA).EQ.O.)VXA=(1.E-10,0.)
                                                                          DIP03710
      VIADB(IC)=20. *ALOG10(CABS(VXA))
                                                                          DIP03720
      IF(VXADB(IC).GT.BIC"DB) BIGNDB=VXADB(IC)
                                                                          DIP03730
      VXAPH(IC)=ATAN2(AIMAG(VXA), REAL(VXA))/DCR
                                                                          DIP03740
      VIMANA(IC)=10.**(VXADB(IC)/20.)*CEXP(CJ*VXAPH(IC)*DCR)
                                                                          DIP03750
     WRITE(6,4457)IC, VXMANA(IC)
CC
                                                                          DIP03760
```

```
DIP03760
     WRITE(6,4457)IC, VIMANA(IC)
CC
                                                                          DIP03770
4457 FORMAT(1X, 'IC, VXMANA=', I4, 2X, 2E12.5)
                                                                          DIP03780
      WRITE(6,6429)IC, IC(IC), VIADB(IC)
                                                                          DIP03790
6429 FORMAT(1X,'IC,X,VXADB=',I4,2X,F10.2,2X,F12.2)
                                                                         DIP03800
3000 CONTINUE
                                                                        DIP03810
     IF(IENORM.EQ.0) GO TO 2500
                                                                          DIP03820
      WRITE(6,3765)BIGNDB
                                                                          DIP03830
3765 FORMAT(1X, 'BIGNDB=', F12.2)
                                                                          DIP03840
      DO 3020 IC=1,NPTSN
                                                                          DIP03850
     WRITE(6,3343)IC, VXADB(IC)
CC
                                                                          DIP03860
      VIADB(IC)=VIADB(IC)-BIGNDB
                                                                          DIP03870
     WRITE(6,3343)IC, VIADB(IC)
CC
                                                                          DIP03880
3020 CONTINUE
                                                                          DIP03890
2500 WRITE(6,3030)
3030 FORMAT(11, 'WANT TO PLOT NEAR FIELD CUTS, IPLOTN=1')
                                                                          DIP03900
                                                                          DIP03910
      READ(5,*)IPLOTN
C***NEXT LINE ADDED TO AVOID RUN TIME ERROR
                                                                          DIP03920
                                                                          DIP03930
      IF(IPCONN.EQ.1) IPLOTN=0
      IF(IPCONN.EQ.O) THEN
        WRITE(30,18889)
        WRITE(30,18888)XN,NCOLXN,DXN,ZN,NCOLZN,DZN
        WRITE(30,4547)
4547 FORMAT(1X,'IZ,IX,VXADB(IC)')
      IPP=0
      DO 7788 IZ=1,NCOLZN
      DO 7788 IY=1,NROWYN
      DO 7788 IX=1,NCOLXN
      IPP=IPP+1
     DO 7788 IPP=1.NPTSN
                                                                          DIP03950
CC
                                                                          DIP03960
      WRITE(6,3343) IPP, VXADB(IPP)
                                                                          DIP03970
3343 FORMAT(1X, 'IPP, VXADB=', I4, 2X, F12.2)
      IF(IPCONN.EQ.O) WRITE(30,*) IZ,IX,VXADB(IPP)
                                                                          DIP03980
7788 CONTINUE
3040 IF((NCOLIN.GT.1.OR.NROWYN.GT.1).AND.IPLOTN.EQ.1)
                                                                          DIP03990
     1CALL PLOTTR(NCOLXN, NROWYN, NCOLZN, XC, YC, ZC, NPTSN, VXADB, VXAPH,
                                                                          DIP04000
                                                                          DTP04010
     2VXADB, VXAPH)
      IF (NELN.EQ.1.AND.IPLOTN.EQ 1)
                                                                          DIP04020
     1CALL PLOTAX (NPTSN, ZC, VXADB, VXAPH, VXADB, VXAPH)
                                                                          DIP04030
                                                                          DIP04040
      IF(IPCONN.EQ.O) GO TO 3939
                                                                          DIPO4050
C***THIS SECTION FOR CONTOUR PLOTS
                                                                          DIP04060
      NCLXN2=NCOLXN+2
                                                                          DIPO4070
      NRUYN2=NROUYN+2
      NCLZN2=NCOLZN+2
                                                                          DIP04080
      DO 3777 IX=1,NCOLXN
      RX(IX)=(XN+DXN+(IX-1))/CINMTR
                                                                          DIP04100
3777 CONTINUE
                                                                          DIP04110
      DO 3008 IY=1, NROWYN
      RY(IY)=(YN+DYN+(IY-1))/CINMTR
                                                                          DIP04130
3008 CONTINUE
      DO 3009 IZ=1,NCOLZN
      RZ(IZ)=(ZN+DZN+(IZ-1))/CINMTR
3009 CONTINUE
                                                                          DIP04140
      TC=0
      IF(IXZ.EQ.1) WRITE(30,18889)
18889 FORMAT(1X,'XN,NCOLXN,DXN,ZN,NCOLZN,DZN=')
      IF(IXZ.EQ.1) WRITE(30,18888) XN, NCOLXN, DXN, ZN, NCOLZN, DZN
18888 FORMAT(1X,E1..5,I5,2X,E14.5,2X,E14.5,I5,2X,E14.5)
       IF(IXZ.EQ.1) WRITE(30,4546)
```

```
4546 FORMAT(1X.'IZ.IX.FM(IZ.IX)')
      DO 3022 IZ=1,NCOLZN
                                                                         DTP04150
      DO 3022 IY=1,NROWYN
      DO 3022 IX=1.NCOLXN
                                                                         DIP04160
      IC=IC+1
                                                                         DIP04170
     IF(IXY.EQ.1) FM(IY,IX)=VXADB(IC)
      IF(IXZ.EQ.1) FM(IZ,IX)=VXADB(IC)
      IF(IX2.EQ.1.AND.IPCONN.EQ.1) WRITE(30,*) IZ,IX,FM(IZ,IX)
                                                                         DIPO4190
3022 CONTINUE
      IF(IXY.EQ.1)
     2CALL PLCONT (NCOLIN, NROWYN, NCLIN2, NRWYN2, RX, RY, FM, DMF, -50.,
                                                                         DIPO4200
     310.,5,1)
                                                                         DIPO4210
     IF(IXZ.Ed.1)
     2CALL PLCONT (NCOLXN, NCOLZN, NCLXN2, NCLZN2, RX, RZ, FM, DMF, ~50.,
     310.,5,1)
3939 CONTINUE
                                                                         DIP04220
     WRITE(25,*)NCOLAN
                                                                         DTP04470
CC
      DO 1767 IDDD=1,NCOLXN
                                                                         DIP04480
     WRITE(25,*)VXADB(IDDD)
CC
1767 CONTINUE
                                                                         DIPO4500
      WRITE(6,3050)
                                                                         DTP04510
3050 FORMAT(1X 'PLOT NEAR FIELD AGAIN?, IPLA=1')
                                                                         DIP04520
      READ(5,*)IPLA
                                                                         DIP04530
      IF(IPLA.EQ.1) GO TO 3040
                                                                         DIPO4540
390
     CONTINUE
                                                                         DIPO4550
      ISTOP=0
                                                                         DIP04560
      IF(ISTOP.EQ.1) GO TO 9999
                                                                         DIP04570
C***CALL PRNTDA (PRINT PARAMETERS)
                                                                         DIP04580
C370 CALL PRINTDA
                                                                         DIP04590
      IF(IPATRN.EQ.O) GO TO 440
                                                                         DIP04630
410 WRITE(6.420)
                                                                         DIP04640
420 FORMAT(1X, 'ISYLMBL FOR PLOTTING, LT. O THEN NOT USED')
                                                                         DIP04650
     READ(5,*)ISYMBL
                                                                         DIP04660
      WRITE(6,430)
                                                                         DIP04700
430 FORMAT(1X,'WANT TO PLOT PATTERNS AGAIN?, IPFNA=1')
                                                                         DIP04710
     READ(5,*)IPFNA
                                                                         DIP04720
      IF(IPFNA.EQ.1) GO TO 410
                                                                         DIP04730
440 CONTINUE
                                                                         DIPO4740
C***THIS SECTION FOR COVARTANCE MATRIX COMPUTATION
                                                                         DIP04750
     IF(NJAMS.EQ.0) GD TO 9999
      DO 7999 ICH=1,NMAX
                                                                         DIP04770
     DO 7999 JCH=1.NMAX
                                                                         DIP04780
      COVNF(ICH, JCH) = DCMPLX(0.0D0,0.0D0)
                                                                         DIPO4790
7999 CONTINUE
                                                                         DIP04800
8888 IF(NJAMS.EQ.O) GO TO 4444
                                                                         DIPO4810
C***THIS SECTION FOR JAMMER COVARIANCE MATRIX
                                                                         DIP04820
     FMINHZ=FCHZ-BWFHZ/2.
                                                                         DIP04830
      DELFHZ=0.0
     IF(NFREQ.GT.1) DELFHZ=BWFHZ/(NFREQ-1)
      DO 600 IFR=1,NFREQ
                                                                         DIP04850
      FHZ=FMINHZ+DELFHZ*(IFR-1)
                                                                         DIP04860
      FGHZ=FHZ/1.0E9
                                                                         DIP04870
C ** * COMPUTE FREE SPACE WAVELENGTH AT EACH FREQUENCY
     RLAMDA=2.997925E10/FHZ/2.54
C+++NOTE: GAMMA= ALPHA +J BETA
          AND RLAMDA=2 PI/ BETA (REF. HAYT PG. 334)
C+++
C***THUS NEED TO COMPUTE GAMMA, AND ETA FOR EACH FREQ.
     IF(IWL.EQ.2) GO TO 8789
      DX=DXIN+CINMTR
```

```
DY=DYIN+CINMTR
      HL=HLIN+CINMTR
      ARAD=ARADIN+CINMTR
      HZ=HZIN+CINMTR
8789 CONTINUE
      CALL ZMATRX(CTDP,STDP,IBLTSL,ICC1,ICC,CZ,Z)
                                                                           DIP04940
C***COMPUTE ELEMENT INDUCED VOLTAGES DUE TO JAMMER SOURCES
                                                                          DTP04950
      CALL VJAMMR(NJAMS, NEL, PWRJDB, ICC1, ICC, PS, CZ, Z, CWTA
                                                                           DIP04960
     1,NB,IDMB,NAUX,IAUXA,ZLOAD,XJAMIN,YJAMIN,ZJAMIN,
                                                                          DIP04970
     2IFR, NFREQ, CTPR, STPR, CTDP, STDP, VMAJMA, VAUXJA)
                                                                           DIPO4980
      WRITE(6,1234)IFR, VMAJMA(1,IFR), VAUXJA(1,1,IFR)
                                                                           DIP04990
      WRITE(8,1234)IFR, VMAJMA(1, IFR), VAUXJA(1,1, IFR)
                                                                           DIP05000
1234 FORMAT(1X, 'IFR, VMAJMA, VAUXA=', I4, 2X, 4E11.4)
                                                                           DIP05010
                                                                           DIP05020
600 CONTINUE
C***FORM RECEIVED VOLTAGE MATRIX
                                                                           DIP05030
C***VRECVM(MAIN A, AUX A1, AUX A2,... AUX AN :)
                                                                           DIPO5050
      DO 9001 IJAM=1,NJAMS
      DO 9002 IFR=1,NFREQ
                                                                           DIP05060
      IF(ISLC.EQ.1) VCHA(1,IFR)=VMAJMA(IJAM,IFR)
      WRITE(6,6789)IJAM, IFR, VMAJMA(IJAM, IFR)
                                                                           DIP05080
6789 FORMAT(1X, 'IJAM, IFR, VMAJMA(IJAM, IFR) = ', 2I4, 2X, 2E12.5)
                                                                           DIP05090
      DO 9003 IA=1,NAUX
                                                                           DIPO5100
      IAP1=IA+1
                                                                           DIP05110
      IF(ISLC.EQ.1) VCHA(IAP1,IFR)=VAUXJA(IA,IJAM,IFR)+VATENA(IA)
      IF(ISLC.EQ.O) VCHA(IA, IFR)=VAUXJA(IA, IJAM, IFR)+VATENA(IA)
      WRITE(6,8876) IJAM, IFR, IA, VAUXJA(IA, IJAM, IFR)
      WRITE(8,8876)IJAM, IFR, IA, VAUXJA(IA, IJAM, IFR)
8876 FORMAT(1X, 'IJAM, IFR, IA, VAUXJA=', 314, 2X, 2E12.4)
                                                                           DIP05150
9003 CONTINUE
                                                                           DIP05160
9002 CONTINUE
                                                                           DIP05170
      DO 5555 KKK=1,4
                                                                           DIP05180
     DO 5555 LLL=1,NFREQ
                                                                           DTP05190
CCC WRITE(6,6655)KKK,LLL,VCHA(KKK,LLL)
                                                                           DIP05200
5555 CONTINUE
                                                                           DIP05210
      WRITE(6,5533)
                                                                           DIP05220
5533 FORMAT(1X,'NOW COMPUTE JAMMER COVARIANCE MATRIX')
                                                                          DIP05230
C***COMPUTE COVARIANCE MATRIX FOR JTH JAMMER SOURCE
                                                                           DIP05240
      IF(NFREQ.GT.1) CALL COVSWC(VCHA, VCHA, NMAX, NFREQ, BWFHZ, COVAAJ)
      DO 9005 ICH=1.NMAX
      DO 9005 JCH=1,NMAX
     COVNF(ICH, JCH) = COVNF(ICH, JCH) + COVAAJ(ICH, JCH)
                                                                           DIP05280
     WRITE(6,2299) ICH, JCH, COVNF (ICH, JCH)
                                                                           DIP05290
9005 CONTINUE
                                                                           DIP05300
9011 CONTINUE
                                                                           DIP05310
C *** ADD RECEIVER NOISE TO DIAGONAL ELEMENTS
                                                                           DIPO5320
4444 DO 8006 ICH=1.NMAX
                                                                           DIP05330
      COVNF(ICH, ICH) = COVNF(ICH, ICH) +1.0D0
                                                                           DIP05340
8006 CONTINUE
                                                                           DIP05350
      DO 2727 I=1,NMAX
                                                                           DIP05360
      DO 2727 J=1,NMAX
                                                                          DIP05370
      IF(CDABS(COVNF(I,J)).EQ.O.ODO) GO TO 2727
                                                                          DIP05380
      CNDB=10.*DLOG10(CDABS(COVNF(I,J)))
                                                                           DIP05390
      PREALN=DREAL(COVNF(I,J))
                                                                           DIP05400
      PIMAGN=DIMAG(COVNF(I,J))
                                                                           DIP05410
      PHASEN=DATAN2(PIMAGN, PREALN)/DCR
                                                                           DIP05420
      WRITE(8,4411)I, J, CNDB, PHASEN
                                                                           DIP05430
      IF(I.EQ.1)WRITE(6,4411)I,J,CNDB,PHASEN
                                                                           DIP05440
4411 FORMAT(1X,'I,J,CNDB,PHASEN=',2I4,2X,2F12.2)
                                                                           DIP05450
2727 CONTINUE
                                                                           DIP05460
```

C===C	OMPUTE COVARIANCE MATRIX INVERSE	DIP05470	
	DO 8007 ICH=1,NMAX	DIP05480	
	DO 8007 JCH=1,NMAX	DIP05490	
	COVNFI(ICH, JCH)=COVNF(ICH, JCH)	DTP05500	
8007		DIP05510	
	CALL DCHINV(COVNFI, LTMP, MTMP, NMAX, NMAX)	DIP05520	
C***C	HECK MATRIX INVERSION ACCURACY (CINVERSE+C=I)	DIP05530	
	CALL CHMULT(COVNFI, COVNF, NMAX, NMAX, NMAX, CINVCN)	DIP05540	
	DO 8008 ICH=1,NMAX	DIP05550	
	)O 8008 JCH=1,NMAX	DIP055 30	
	IF(ICH.EQ.1)WRITE(6,8009)ICH, JCH, CINVCN(ICH, JCH)	DIP05570	
8009	FORMAT(1X,'ICH, JCH, CINVCN=', 214, 2X, 2E12.5)	DIP05580	
8008	CONTINUE	DIP05590	
C+++C	OMPUTE EIGENVALUES (ALSO EIGENVECTORS AND PERFORMANCE INDEX)	DIP05600	
	IJOB=2	DIP05610	
	CALL EIGCC(COVNF, NMAX, NMAX, IJOB, EIGVAN, EIGVEN, NMAX, WKE, IER)	DIP05620	
	WRITE(6,155)IER	DIP05630	
155	FORMAT(1X,'AFTER COMPUTE EIGENVALUES, IER=',15)	DIP05640	
	DO 200 I=1,NMAX	<b>DIP056</b> 50	
	WRITE(6,300)I,EIGVAN(I)	DIP05660	
300	FORMAT(1X,'I,EIGVAN=',2X,I4,2E10.3)	DIP05670	
200	CONTINUE	DIP05680	
	DO 205 I=1,NMAX	DIP05690	
	EVLDBN=10.*DLOG10(CDABS(EIGVAN(I)))	<b>DIP0571</b> 0	
	WRITE(8,207)I,EVLDBN	DIP05720	
	WRITE(6,207)I,EVLDBN	DIP05730	
207	FORMAT(1X,'I,EVLDBN=',I4,2X,2F12.3)	DIP05740	
205	CONTINUE	DIP05750	
	CALL PRNTDA	DIP05760	
C++F	ILL-IN SIDELOBE CANCELLER QUIESCENT WEIGHTS	DIP05770	
	DO 8985 I=1,NMAX	DIP05780	
	WQSLC(I,1)=DCMPLX(0.0D0,0.0D0)	DIP05790	
8985	CONTINUE	DIP05800	
	WQSLC(1,1)=DCMPLX(1.0D0,0.0D0)	DIP05810	
	DO 7769 I=1,NMAX	DIP05820	
	WRITE(6,1212)I, WQSLC(I,1)	DIP05830	
	WRITE(8,1212)I, WQSLC(I,1)	DIP05840	
	FORMAT(1X,'I,WQSLC=',I4,2X,2E12.5)	DIP05850	
	CONTINUE	DIP05860	
C***C	OMPUTE QUIESCENT INR	DIP05870	
	IF(ISLC.EQ.1) CALL INRTIO(WQSLC,COVNF,NMAX,WTCTR,CMPROD,QINRDB)		
	IF(ISLC.EQ.O) CALL INRTIO(CWTA, COVNF, NMAX, WTCTR, CMPROD, QINRDB)	DIDOFGOO	
CC	QINRAA=10.*DLOG10(CDABS(COVNF(1,1))) OMPUTE AVERAGE CANCELLAFION	DIP05890	
(===(		DIBOE010	
	SUMC=0.0 SUMAA=0.0	DIP05910	
	DO 1829 IR=1,NRAN	DIP05920 DIP05930	
C7	ERO-OUT ADAPTIVE WEIGHTS INITIALLY	DIP05930 DIP05940	
C2	DO 57 I=1.NMAX	DIP05950	
	WAN(I,1)=DCMPLX(0.0D0,0.0D0)	DIP05960	
57	CONTINUE	DIP05900	
	OMPUTE ADAPTIVE ARRAY WEIGHTS	DIP05980	
IF(ISLC.EQ.1) CALL CMMULT(COVNFI, WQSLC, NMAX, NMAX, 1, WAN)			
IF(ISLC.EQ.O) CALL CMMULT(COVNFI, CWTA, NHAX, NHAX, 1, WAN)			
C ***QUANTIZE AND RANDOMIZE ADAPTIVE WEIGHT SETTINGS DIF			
IF(IQUAN.EQ.1.AND.NBADWT.LT.20)			
2CALL ADQUAN(NMAX, WAN, NBADWT, IRNERR, INITRN, AWSGDB, AWSGDG)			
	ORMALIZE FULLY ADAPTIVE WEIGHTS	DIP06020	
,	IF(ISLC.EQ.0) THEN		

```
SAVE=0.0D0
        DO 33345 I=1,NEL
        SAVE=SAVE+CDABS(WAN(I,1))**2
33345
        CONTINUE
        SAVESR=DSQRT(SAVE)
        DO 44456 I=1.NEL
        WAN(I,1)=WAN(I,1)/SAVESR
        WANDB=20. *DLOG10(CDABS(WAN(I,1)))
        IF(IR.EQ.1) WRITE(6,3599)I, WANDB
        IF(IR.EQ.1) WRITE(8,3599)I, WANDB
3599
        FORMAT(1X,'I=', I4,2X,'WANDB=',1X,F12.3)
44456
       CONTINUE
      ENDIF
C****TO PRINT NORMALIZED WEIGHTS
                                                                          DIP06030
      BIGWDB=-299.0
                                                                          DIPO6040
      DO 9102 I=1,NMAX
                                                                          DIP06050
      IF(IR.EQ.1) WRITE(6.6191)I.WAN(I.1)
                                                                          DIP06100
      IF(IR.EQ.1) WRITE(8,6191)I,WAN(I,1)
                                                                          DIPO6110
6191 FORMAT(1X, 'I, WAN(I,1) = ', I4,2X,2E12.5)
      WANDBM=20. +DLOG10(CDABS(WAN(I,1)))
                                                                          DIP06060
      IF(WANDBM.GT.BIGWDB) BIGWDB=WANDBM
                                                                          DIP06070
9102 CONTINUE
                                                                          DIP06080
      DO 9112 I=1,NMAX
                                                                          DIP06090
      WANDB=20. *DLOG10(CDABS(WAN(I,1)))-BIGWDB
                                                                          DIP06130
      WANDG=DATAN2(DIMAG(WAN(I,1)), DREAL(WAN(I,1)))/DCR
      IF(IR.EQ.1) WRITE(6,35991)I, WANDB, WANDG
      IF(IR.EQ.1) WRITE(8,35991)I, WANDB, WANDG
35991 FORMAT(1X,'I=',I4,2X,'WANDB,WANDG=',1X,2F14.5)
9112 CONTINUE
                                                                          DIP06170
C***COMPUTE ADAPTED INR
                                                                          DIP06180
      CALL INRTIO(WAN, COVNF, NMAX, WTCTR, CMPROD, AINRDB)
                                                                          DIP06190
C***COMPUTE CANCELLATION
                                                                          DIP06200
      CANCER=AINREP-GINRER
                                                                          DIP06210
      WRITE(6,3007) GINRDB, AINRDB, CANCDB
                                                                          DIP06220
      WRITE(8,3007)QINRDB,AINRDB,CANCDB
                                                                          DIP06230
3007 FORMAT(1X,'INR= QUI,ADAP ,CANCEL=',2X,3F10.3,2X,' DB')
                                                                          DIP06240
      SUMC=SUMC+CANCDB
                                                                          DIP06250
C***COMPUTE ADAPTED INR CH. A
                                                                          DIP06260
      DO 2255 I=1,NMAX
                                                                          DIP06270
      WANDMA(I,1)=WAN(I,1)
                                                                          DIP06280
      IF(I.GT.NMAX) WANDMA(I,1)=(0.0D0,0.0D0)
                                                                          DIP06290
2255 CONTINUE
                                                                          DIP06300
      CALL INRTIO (WANDMA, COVNF, NMAX, WTCTR, CMPROD, AINRAA)
                                                                          DIP06310
C ** * CANCELLATION CH. A
                                                                          DIP06320
      CNCLNA=AINRAA-QINRAA
                                                                          DIP06330
      WRITE(6,3738)QINRAA,AINRAA,CNCLNA
                                                                          DIP06340
      WRITE(8,3739)QINRAA,AINRAA,CNCLNA
                                                                          DIP06350
3738 FORMAT(1X, 'SIDELOBE CANCELLER CH. A INR= QUI, ADAP, CAN=',2X,3F10.3)
      SUMAA=SUMAA+CNCLNA
                                                                          DIP06370
1829 CONTINUE
                                                                          DIP06380
      AVECAN=SUMC/NRAN
                                                                          DIP06390
      WRITE(6,4456)CANCDB, NRAN, AVECAN
                                                                          DIP06400
4456 FORMAT(1X, 'CANCDB, NRAN, AVECAN=', F12.5, 2X, I5, 2X, F12.5)
                                                                          DIP06410
      WRITE(8,4456)CANCDB, NRAN, AVECAN
                                                                          DIPO6420
      AVECAA=SUMAA/NRAN
                                                                          DIP06430
      WRITE(6,2220)AVECAA
                                                                          DIPO6440
      WRITE(8,2220)AVECAA
                                                                          DIP06450
2220 FORMAT(1X, 'AVE. CANEL, CH. A, =',2X,F12.5)
                                                                          DIP06460
C***SECTION TO COMPUTE ADAPTIVE ARRAY RADIATION PATTERNS
                                                                          DIP06470
```

```
CC
      IF(INEAR.EQ.O.OR.IANGLP.EQ.O) GO TO 9990
                                                                         DTP06480
      IF(INEAR.EQ.0) GG TO 9990
                                                                         DIP06490
      BIGADB=-299.
                                                                         DIP06500
                                                                         DIP06510
      IC=0
      DO 8919 IZ=1.NCOLZN
      DO 8919 IY=1.NROWYN
                                                                         DTP06520
      DO 8919 IX=1,NCOLXN
                                                                         DIP06530
                                                                         DTP06540
      IC=IC+1
      CSUMA=(0.,0.)
                                                                         DIP06550
      IF(ISLC.EQ.1) CSUMA=CSUMA+DCONJG(WAN(1,1))+VIMANA(IC)
      DO 7921 IAX=1.NAUX
                                                                         DIP06570
      IAXP1=IAX+1
                                                                         DIP06580
      IF(ISLC.EQ.1) CSUMA=CSUMA+DCONJG(WAN(IAXP1,1)) *VAXCWA(IAX.IC)
      IF(ISLC.EQ.O) CSUMA=CSUMA+DCONJG(WAN(IAX,1))=VAXCWA(IAX,IC)
7921 CONTINUE
                                                                         DTP06610
      IF(CABS(CSUMA).EQ.O.) CSUMA=(1.E-10,0.)
                                                                         DIP06620
      PCHADB(IC)=20. *ALOG10(CABS(CSUMA))
                                                                         DIP06630
      IF(PCHADB(IC).GT.BIGADB) BIGADB=PCHADB(IC)
                                                                         DIP06650
      WRITE(6,4999)IC,PCHADB(IC)
                                                                         DIP06660
8919 CONTINUE
                                                                         DIP06670
C***NORMALIZE ADAPTIVE PATTERNS
                                                                         DIP06680
      IF(IIZ.EQ.1) WRITE(31,18889)
      IF(IXZ.EQ.1) WRITE(31,18888)XN,NCOLXN,DXN,ZN,NCOLZN,DZN
      IF(IXZ.EQ.1) WRITE(31,4546)
      IC=0
                                                                         DIP06690
      DO 3459 IZ=1,NCOLZN
      DO 3459 IY=1,NROWYN
                                                                         DIPO6700
      DO 3459 IX=1,NCOLXN
                                                                         DIP06710
      IC=IC+1
                                                                         DIP06720
      PCHADB(IC)=PCHADB(IC)-BIGADB
                                                                         DIP06730
      IF(IPCONN.EQ.0) WRITE(31,*)IZ,IX,PCHADB(IC)
      IF(IPCONN.EQ.O) GO TO 3459
                                                                         DIP06760
      IF(IXY.EQ.1) FM(IY,IX)=PCHADB(IC)
      IF(IXZ.EQ.1) FM(IZ,IX)=PCHADB(IC)
      IF(IXZ.EQ.1.AND.IPCONN.EQ.1) WRITE(31,*) IZ,IX,FM(IZ,IX)
3459 CONTINUE
                                                                         DIP06780
      WRITE(25,*)NCOLXN
                                                                         DIP06800
      DO 1879 IDDD=1.NCOLXN
                                                                         DIP06810
      WRITE(25,*)PCHADB(IDDD)
1879 CONTINUE
                                                                         DTP06830
      IF ((NCOLXN.GT.1.OR.NROWYN.GT.1).AND.IPLOTN.EQ.1.AND.IANGLP.EQ.0)
                                                                         DIP06840
     2CALL PLOTTR(NCOLXN, NROWYN, NCOLZN, XC, YC, ZC, NPTSN, PCHADB, PCHADB,
                                                                         DTP06850
                                                                         DIP06860
     3PCHADB, PCHADB)
      IF(IPCONN.EQ.1.AND.IXY.EQ.1)
     2CALL PLCONT (NCOLXN, NROWYN, NCLXN2, NRWYN2, RX, RY,
     3FM,DMF,-50.,10.,5,1)
                                                                         DIP06880
      IF(IPCONN.EQ.1.AND.IXZ.EQ.1)
     2CALL PLCONT(NCOLXN, NCOLZN, NCLXN2, NCLZN2, RX, RZ,
    3FM,DMF,-50.,10.,5,1)
9990 CONTINUE
                                                                         DIP06890
9999 CONTINUE
                                                                         DTP06900
      CALL DONEPL
                                                                         DIP06910
      CALL GETCP2(CPUL)
      CPUTOT=CPUL-CPU1
      WRITE(6,2006)CPUTGT
2006 FORMAT(1X, 'TOTAL CPU TIME=',F15.2)
                                                                         DIP06920
                                                                         DIP06930
      FMD
C***SUBROUTINE TO COMPUTE IMPEDANCE MATRIX
                                                                         DIP07720
```

```
DIP07730
      SUBROUTINE ZMATRX (CTDP, STDP, IBLTSL, ICC1, ICC, CZ, Z)
      COMPLEX Z(ICC1,ICC),CZ(1)
                                                                          DIP07740
                                                                          DIP07750
      COMMON/A/DX, DY, NCOLX, NROWY, NEL, HZ, HL, ARAD, ZLOAD, ZCHAR
                                                                           DIP07760
      COMMON /GROUND/ IGRNDP
      COMMON /CIRCLE/ ICIRC, RADIUS, HLS
                                                                           DIP07770
      NB=NCOLX
                                                                          DIP07780
      IDMB=NROWY
C***COMPUTE MUTUAL IMPEDANCES Z(1,1),Z(1,2),Z(1,3),...,Z(1,NEL).
                                                                          DIP07790
                                                                           DIP07800
CC
     CALL RGDZMN(IGRNDP, ICC1, ICC, Z)
                                                                           DIP07810
      CALL RGDZAB(CTDP,STDP,IGRNDP,ICC1,ICC,Z)
      IF(ICIRC.EQ.O) CALL RGDZA2(CTDP, STDP, IGRNDP, ICC1, ICC, Z)
      IF(ICIRC.EQ.1) CALL CADZA2(CTDP,STDP,IGRNDP,ICC1,ICC,Z)
      WRITE(8,10) NEL, NCOLX, NRUWY
                                                                           DIP07830
                                                                           DIP07840
      FORMAT(1X, 'NEL, NCOLX, NROWY=', 315)
                                                                           DIP07850
      ICOUNT=0
                                                                           DIP07860
      DO 20 I=1,NCOLX
                                                                           DIP07870
      DO 20 J=1,NROWY
                                                                           DTP07880
      ICOUNT=ICOUNT+1
                                                                           DIP07890
      WRITE(6,80) I,J,Z(1,ICOUNT)
      FORMAT(1X,'I,J,Z(1,ICOUNT)=',2I4,2X,2E12.5)
                                                                           DIP07900
80
                                                                           DIP07910
20
      CONTINUE
                                                                           DIP07920
      IDM=NEI.
                                                                           DIP07930
      ICC=NEL
C***FILL THE IMPEDANCE MATRIX
                                                                           DIP07940
                                                                           DIP07950
      IF(NCOLX.LE.1) GO TO 70
                                                                           DIP07960
      IBLT=IBLTSL
                                                                           DIP07970
      IDM1=IDMB
                                                                           DIP07980
      IF(IBLT.EQ.O) IDM1=NEL
                                                                           DTP07990
      IF(NROWY.GT.1) GO TO 40
                                                                           DIP08000
C***FILL TOEPLITZ MATRIX
      DO 30 I=2, NEL
                                                                           DIP08010
                                                                           DIP08020
      DO 30 J=I, NEL
      K=1+J-I
                                                                           DIP08030
                                                                           DIPC8040
      Z(I,J)=Z(1,K)
                                                                           DIP08050
      CONTINUE
30
                                                                           DIP08060
      GO TO 50
                                                                           DIP08070
      CALL BTOEPL(IBLT, NB, IDMB, IDM1, IDM, Z)
40
                                                                           DIP08080
      CONTINUE
                                                                           DIP08090
      IF(IBLT.EQ.0) GO TO 70
                                                                           DTP08100
      DO 60 I=1,IDMB
                                                                           DIP08110
      DO 60 J=1,IDM
                                                                           DIP08120
        IC=(J-1)+IDMB+I
                                                                           DIP08130
        CZ(IC)=Z(I,J)
CC
      WRITE(8,7878)I,J,Z(I,J)
                                                                           DIPO8140
                                                                           DIP08150
7878 FORMAT(1X, 'I, J, Z(I, J)=', 2I4, 2X, 2E12.5)
                                                                           DIPO8160
      CONTINUE
60
                                                                           DIP08170
70
      CONTINUE
                                                                           DIP08180
      RETURN
                                                                           DIP08190
C***SUBROUTINE TO COMPUTE RECEIVE BEAMFORMER WEIGHTS
                                                                           DIP09150
      SUBROUTINE VRCVWT(NACOLX, NAROWY, NXDUM, NYDUM, EDGTDB, WA, WB)
                                                                           DIP09160
                                                                           DIP09170
      DIMENSION WA(1), WB(1)
                                                                           DIP09180
      DIMENSION WT(180)
      COMMON/A/DX,DY,NCOLX,NROWY,NEL,HZ,HL,ARAD,ZLOAD,ZCHAR
                                                                           DIP09200
      COMMON /CHEBY/ ICHEB, SLLDB
                                                                           DIP09210
                                                                           DIP09220
      PI=3.141592654
      TFI=2.*PI
                                                                           DIP09230
                                                                           DIP09240
      DCR=PI/180.
                                                                           DIP09250
      SLLDB=40
CC
```

```
DIP09260
     TAP=10. ** (EDGTDB/20.)
                                                                          DTP09270
      WRITE(6,*)TAP
                                                                          DIP09280
      AMP=(1.-TAP)/2.
      NPDX=NCQLX-NACQLX-2=NXDUM
                                                                          DIP09290
                                                                          DIP09300
      AXL=DX+(NACOLX-1)
                                                                          DIP09310
      AYL=DY+(NAROWY-1)
                                                                          DIP09320
      IF(TAP.NE.1.0)FX=AXL/2.*PI/ACOS(TAP)
                                                                          DIP09330
      IF(TAP.NE.1.0)FY=AYL/2.*PI/ACOS(TAP)
                                                                          DIP09340
      WRITE(6.20)
      FORMAT(1X,'BEFORE CALL CHEBWT')
                                                                          DIP09350
                                                                          DIP09360
      WRITE(6,30)SLLDB
                                                                          DIP09370
      FORMAT(1X, 'SLLDB=', 2X, F10.2)
30
      IF(ICHEB.EQ.1)CALL CHEBWT(NACOLX,SLLDB,WT,RLOSS)
                                                                          DTP09380
                                                                          DIP09390
      WRITE(6,40)
                                                                          DIP09400
      FORMAT(1X,'AFTER CHEBWT')
40
                                                                          DIP09410
      DO 1000 IC=1, NEL
                                                                          DIP09420
      WA(IC)=0.0
                                                                          DIP09430
      WB(IC)=0.0
                                                                          DIP09440
1000 CONTINUE
                                                                          DTP09450
C***COMPUTE EFFECTIVE DIPOLE CENTER COORDS. FOR BOTH PHASE CENTERS
                                                                          DIP09460
      XO=-AXL/2.
                                                                          DTP09470
      YO=-AYL/2.
                                                                          DIP09480
      IC=0
                                                                          DIP09490
      DO 80 I=1, NACOLX
                                                                          DIP09500
      DO 80 J=1, NAROWY
                                                                          DIP09510
      IC=IC+1
                                                                          DIP09520
      X=X0+DX+(I-1)
                                                                          DIP09530
      Y=Y0+DY=(J-1)
                                                                          DIP09540
      TAPERX=1.0
                                                                          DIP09550
      TAPERY=1.0
                                                                          DTP09560
      IF(ICHEB.EQ.1) GO TO 70
                                                                          DIP09570
      IF(TAP.NE.1.0)TAPERX=COS(PI*X/FX)
      IF(TAP.NE.1.0.AND.FY.NE.0.0)TAPERY=COS(PI+Y/FY)
                                                                          DTP09580
                                                                          DIP09590
60
      TAPER=TAPERX+TAPERY
                                                                          DIP09600
      WT(IC)=TAPER
                                                                          DIP09610
70
      CONTINUE
                                                                          DIP09620
80
      CONTINUE
                                                                          DIP09630
C***TRANSFORM FROM SUB-APERTURES TO FULL ARRAY
                                                                           DIP09640
      IBGNA=NXDUM+NROWY+NYDUM
                                                                           DIP09650
      IBGNB=IBGNA+NPDX+NROWY
                                                                           DIF09660
      IC=0
                                                                           DIP09670
      DO 1010 IX=1,NACOLX
                                                                           DIP09680
      DO 1010 IY=1,NAROWY
                                                                           DTP09690
      IC=IC+1
                                                                           DIP09700
      IA=IBGNA+(IX-1)*NROWY+IY
                                                                           DIP09710
      IB=IBGNB+(IX-1) +NROWY+IY
                                                                           DIP09720
      WA(IA)=WT(IC)
                                                                           DIP09730
      WB(IB)=WI(IC)
                                                                           DIP09740
1010 CONTINUE
                                                                           DIP09750
      RETURN
                                                                           DIP09760
      END
C *** SUBROUTINE TO COVARIANCE MATRIX BASED ON NUMERICAL INTEGRATION
                                                                           DIP09770
C *** IN THE FREQUENCY DOMAIN ACCORDING TO SIMPSON'S RULE
                                                                           DTP09780
                                                                           DIP09790
      SUBROUTINE COVSWC(VA, VB, NCHAN, NFREQ, BWFHZ, COVAB)
       COMPLEX *16 VA(NCHAN, NFREQ), VB(NCHAN, NFREQ)
                                                                           DIP09800
                                                                           DIP09810
      COMPLEX *16 COVAB(NCHAN, NCHAN), CSUM, DCABF
                                                                           DIP09820
      REAL +8 SWC(101)
                                                                           DIP09830
      DELTAF=BWFHZ/(NFREQ-1)
                                                                           DIP09840
       CALL SIMUC(NFREQ, SWC)
```

```
DO 10 ICH=1,NCHAN
                                                                           DIP09850
                                                                           DIP09860
      DO 10 JCH=1,NCHAN
      CSUM=(0.CD0,0.0D0)
                                                                           DIP09870
      DO 20 IFR=1,NFREQ
                                                                           DIP09880
      DCABF=VA(ICH, IFR) *DCONJG(VB(JCH, IFR))
                                                                           DIP09890
      CSUM=CSUM+DCABF+SWC(IFR)
                                                                           DIP09900
CC
      WRITE(6,6767) IFR, CSUM
                                                                           DTP09910
6767 FORMAT(1X,'IFR,CSUM=',14,2X,2E12.5)
                                                                           DIP09920
20
      CONTINUE
                                                                           DIP09930
      COVAB(ICH, JCH) = DELTAF /2. *CSUM
                                                                           DIF 09940
C***NEW LINE TO NORMALIZE COVAB
                                                                           DIP09950
      COVAB(ICH, JCH) = COVAB(ICH, JCH) / BWFHZ
                                                                           DTP09960
      WRITE(6,4567) ICH, JCH, COVAB(ICH, JCH)
                                                                           DIP09970
4567 FORMAT(1X,'ICH,JCH,COVAB(ICH,JCH)=',2I4,2X,2E12.5)
                                                                           DTP09980
10
      CONTINUE
                                                                           DIP09990
                                                                           DIP10000
      RETURN
      END
                                                                           DIP10010
C***SUBROUTINE TO GENERATE SIMPSON'S 1/3 RULE WEIGHTING COEF.
                                                                           DIP10020
C---INTGERAL F(X)DX=(DELTAX/3.)+(F(1)+4+F(2)+2+F(3)+4+F(4)+ ..+F(00D)) DTP10030
C*** THE SERIES 1 4 2 4 2 4 ..... ARE SIMPSON'S COEF.
                                                                           DIP10040
      SUBROUTINE SIMWC(NCOEF, SWC)
                                                                           DTP10050
                                                                           DIP10060
      REAL +8 SWC(NCOEF)
      DO 10 N=1.NCDEF
                                                                           DIP10070
      XNN=FLOAT(N)
                                                                           DIP10080
      NN=N/2
                                                                           DIP10090
      TT=INN/2
                                                                           DTP10100
      DIF=TT-FLOAT(NN)
                                                                           DIP10110
      NC=2
                                                                           DTP10120
      IF(DIF.EQ.O.) NC=4
                                                                           DIP10130
      IF(N.EQ.1.OR.N.EQ.NCOEF) NC=1
                                                                           DIP10140
      SWC(N)=NC
                                                                           DIP10150
10
      CONTINUE
                                                                           DIP10160
                                                                           DIP10170
      RETURN
                                                                           DIP10180
C *** SUBROUTINE TO COMPUTE RECEIVED VOLTAGES DUE TO JAMMER SOURCES
                                                                           DTP10190
      SUBROUTINE VJAMMR(NJM, NEL, PWRJDB, ICC1, ICC, PS, CZ, Z, CWTA
                                                                           DTP10200
     1, NB, IDMB, NAUX, IAUXA, ZLOAD, XJAMIN, YJAMIN, ZJAMIN,
                                                                           DIP10210
     21FR, NFR, CTPR, STPR, CTDP, STDP, VMAINA, VAUXA)
                                                                           DIP10220
      COMPLEX *16 CWTA(NEL,1)
                                                                           DIP10230
      COMPLEX PS(1),CZ(1),Z(ICC1,ICC)
                                                                           DTP10240
      COMPLEX VJM(180), VREFJM
                                                                           DIP10250
                                                                           DIP10260
      COMPLEX CJ.CSUMA
      COMPLEX +16 VMAINA(NJM,NFR)
                                                                           DIP10270
      COMPLEX *16 VAUXA(NAUX, NJM, NFR)
                                                                           DIP10280
      DIMENSION PWRJDB(1), PWRJ(10)
                                                                           DIP10290
      DIMENSION XJAMIN(1), YJAMIN(1), ZJAMIN(1)
                                                                           DIP10300
                                                                           DTP10310
      INTEGER IAUXA(1)
      COMMON/PCFNTR/PCDXIN
                                                                           DIP10320
      COMMON /B/ NGEN, IGEN, THETAS, PHIS, IMUT, IBLTSL, IPATRN
                                                                           DIP10330
      COMMON /D/ FGHZ, RLAMDA, IWL, IS, NSCANS
                                                                           DIP10340
      COMMON /WRITE/ IWR
                                                                           DIP10350
      COMMON /CIRCLE/ ICIRC, RADIUS, HLS
      PI=3.141592654
                                                                           DIP10360
                                                                           DTP10370
      DCR=PI/180.
                                                                           DIP10380
      CJ=(0.,1.)
      CINMTR=0.0254
      ISYM=0
                                                                           DIP10390
C *** CONVERT DB TO POWER (RELATIVE TO NOISE)
                                                                           DIP10400
                                                                           DIP10410
      DO 10 II=1,NJM
```

	PWRJ(II)=10.**(PWRJDB(II)/10.)	DIP10420
	WRITE(6,66)II,PWRJ(II)	DIP10430
	WRITE(8,66)II,PWRJ(II)	
	WRITE(0,00)11,FWD(11)	DIP10440
<b>6</b> 6	FORMAT(11, 'II, PWRJ IN POWER=', I4, 21, E12.5)	DIP10450
10	CONTINUE	DIP10460
	IF(IFR.GT.1) GO TO 55	
	WRITE(6,20)IWL	DIP10470
	DO 5544 IJAM=1,NJM	DIP10480
	WRITE(6 40)	DIP10490
	WRITE(6,50)XJAMIN(IJAM), YJAMIN(IJAM), ZJAMIN(IJAM)	DIP10500
5544	CONTINUE	DIP10510
	FORMAT(1X,'IWL*',I4)	DIP10520
20		DIP10530
	WRITE(6,30) FORMAT(1X, 'CHANGE NEAR FIELD JAMMER POSITIONS (INCHES)?, ICH=1')	
30		DIP10550
	READ(5,*)ICH	DIP10560
	IF(ICH.EQ.O) GD TO 55	
	DO 8887 IJAM=1,NJM	DIP10570
	WRITE(6,40)	DIP10580
40	FORMAT(1X, 'XJAMIN, YJAMIN ZJAMIN =')	DIP10590
	READ(5,*,END=55)XJAMIN(IJAM),YJAMIN(IJAM),ZJAMIN(IJAM)	DIP10600
	WRITE(6,50)XJAMIN(IJAM), YJAMIN(IJAM), ZJAMIN(IJAM)	DIP10610
50	FORMAT(11,5F10.3,21,2I5)	DIP10620
	CONTINUE	DIP10630
8887		DIP10640
55	CONTINUE	
	PCDX=PCDXIN+CINMTR	DIP10660
	WRITE(6,4757)IFR	DIP10670
4757	FORMAT(1X, 'FREQ. INDEX, IFR=', 14)	
60	CONTINUE	DIP10680
C**AL	L DIMENSIONS IN METERS	
	DO 180 IPHACN=1,1	DIP10700
	WRITE(6,2223)IPHACN	DIP10710
2223	FORMAT(1X,'IPHACN=',I4)	DIP10720
	IREFDP=-PCDX/2.+(IPHACN-1)*PCDX	DIP10730
C+++D	ERFORM JAMMER SOURCE SCAN	DIP10740
C+++F	DO 180 IJAM=1,NJM	DIP10750
	I=XJAMIN(IJAM)*CINMTR+PCDX*(IPHACN-1)	
	Y=YJAMIN(IJAM)+CINMTR	
	ZPOS=ZJAMIN(IJAM)+CINMTR	DIP10790
	WRITE(6,6688)I,Y,ZPOS,IREFDP	D11 10/30
6688	FORMAT(1X,'X,Y,ZPOS, XREFDP (METERS)=',4F12.4)	
	IF(ICIRC.EQ.O)	
	2CALL NFDPX2(CTPR,STPR,CTDP,STDP,X,Y,ZPOS,XREFDP,1,VJM,VREFJM)	
	IF(ICIRC.EQ.1)	
	2CALL NFDPC2(X,Y,ZPOS,XREFDP,1,VJM,VREFJM)	
Cessi	NORMALIZE INCIDENT JAMMER POWER	DIP10820
<b>0.</b>	DO 70 INORM=1, NEL	DIP10830
	VJM(INORM)=VJM(INORM)/VREFJM+SQRT(PWRJ(IJAM))	DIP10840
	WRITE(6.77) INORM, VJM(INORM)	DIP10850
	FORMAT(1X,'INORM, VJM(INORM) VOLTAGE*', 14,2X,2E12.5)	DIP10860
77		DIP10870
70	CONTINUE	DIP10880
	IF(IMUT.EQ.0) GO TO 90	DIP10890
C+++5	SOLVE EACH SYSTEM OF EQUATIONS FOR THE UNKNOWN CURRENTS	-
	IF(IBLTSL.EQ.1) GO TO 80	DIP10900
	I12=1	DIP10910
	IF(IJAM.GT.1.OR.IPHACN.GT.1) I12=2	DIP10920
	WRITE(6,6110)	DIP10930
6110	THE PARTY OF THE P	DIP10940
0110	CALL CROUT(Z,VJM,ICC1,ICC,ISYM,IWR,I12,NEL)	DIP10950
	GO TO 90	DIP10960
	90 10 30	

80	IENTRY=4	DIP10970
00	IF(IJAM.GT.1.OR.IPHACN.GT.1) IENTRY=3	DIP10970
	CALL BLTSOL(CZ,VJM,PS,NB,IDMB,IENTRY)	
	IF(NEL.LT.40) CALL CHORM(VJM,NEL)	DIP10990
00		DIP11000
90	CONTINUE COMPUTE RECEIVED VOLTAGES FOR PRESENT SOURCE POSITION	DIP11010
((		DIP11020
	DO 100 IEL=1, NEL	DIP11030
CC	VJM(IEL)=VJM(IEL)+ZLOAD	DIP11040
CC	WRITE(6,4456) IEL, VJM(IEL)	DIP11050
100	FORMAT(1X,'IEL,VJM(IEL) RECEIVED VOLT.=',I4,2X,2E12.5) CONTINUE	DIP11060
	ERFORM BEAM FORMATION FOR MAIN A	DIP11070
	HASE CENTER A	DIP11080
(r		DIP11090
	CSUMA=(0.0D0,0.0D0)	DIP11100
	DO 110 KC=1,NEL	DIP11110
110	CSUMA=CSUMA+VJM(KC)+CWTA(KC,1)	DIP11120
110	CONTINUE	DIP11130
	VMAINA(IJAM,IFR)=CSUMA	DIP11140
2222	WRITE(6,222)IJAM, IFR, VMAINA(IJAM, IFR)	DIP11150
	FORMAT(1X,'IJAM,IFR,VMAINA= (AFTER B.F.)',214,2X,2E12.5)	DIP11160
C+++C	OMPUTE AUXILIARY CHANNEL VOLTAGES	DIP11170
	DO 7000 IAUX=1,NAUX	DIP11180
7000	VAUXA(IAUX,IJAM,IFR)=VJM(IAUXA(IAUX)) CONTINUE	DIP11190 DIP11200
180	CONTINUE	DIP11200 DIP11210
160	DO 9000 IJAM=1.NJAMS	DIP11210 DIP11220
ccc	WRITE(6,4433)IJAM, IFR, VMAINA(IJAM, IFR)	DIP11230
4433	FORMAT(1X, 'VJAMMR: IJAM, IFR, VMAINA=', 214, 2X, 2E12.5)	DIP11230
1100	DD 9001 IAUX=1.NAUX	DIP11240
	WRITE(6,3333)IAUX, IJAM, IFR, VAUXA(IAUX, IJAM, IFR)	DIP11260
3333	FORMAT(1X, 'VJAMMR:IAUX,IJAM,IFR,VAUXA=',314,2X,2E12.5)	DIP11270
9001	CONTINUE	DIP11280
	CONTINUE	DIP11290
	RETURN	DIP11300
	END	DIP11310
C***S	UBROUTINE TO COMPUTE INTERFERENCE TO NOISE RATIO	DIP11320
	SUBROUTINE INRTIO(WT,COV,NEL,WTCTR,CMPROD,DBINR)	DIP11330
	IMPLICIT REAL *8 (A-H, D-Z)	DIP11340
	COMPLEX *16 WT(NEL,1),COV(NEL,NEL),WTCTR(1,NEL)	DIP11350
	COMPLEX *16 CMPROD(1,NEL),CPROD1(1,1),CPROD2(1,1),CINR	DIP11360
	COMMON /NORMAL/ INRNOR	
	WRITE(12.2222)	DIP11370
2222	FORMAT(1X,'INSIDE INRTIO SUBROUTINE')	DIP11380
	DO 110 I=1, NEL	DIP11390
CC	WRITE(12,445)I,WT(I,1)	DIP11400
445	FORMAT(1X,'I,WT(I,1)=',2X,I4,2X,2E12.5)	DIP11410
110	CONTINUE	DIP11420
	CALL CONJTR(WT, NEL, 1, WTCTR)	DIP11430
	DO 111 I=1,NEL	DIP11440
CC	WRITE(12,666)I,WT(1,1),WTCTR(1,I)	DIP11450
666	FORMAT(1X,'I,WT,WTCTR=',2X,I4,2X,4E12.5)	DIP11460
111	CONTINUE	DIP11470
	CALL CHMULT(WTCTR,COV,1,NEL,NEL,CMPROD)	DIP11480
	CALL COMULT(CMPROD, WT, 1, NEL, 1, CPROD1)	DIP11490
	CALL CHMULT(WTCTR, wT,1,NEL,1,CPROD2)	DIP11500
	WRITE(6,333)CPROD1(1,1),CPROD2(1,1)	
	WRITE(8,333)CPROD1(1,1),CPROD2(1,1)	
333	FORMAT(1X,'CPROD1,CPROD2=',4E12.5)	DIP11520
C+++I	NR NORMALIZED	DIP11530

```
CC
      INRNOR=0
      WRITE(6,7739)INRNOR
                                                                         DIP11550
                                                                         DIP11560
      WRITE(8,7739)INRNOR
7739 FORMAT(1X, 'INR NORMALIZATION PARAMETER, INRNOR=', 14)
                                                                         DIP11570
                                                                         DIP11580
     IF(INRNOR.EQ.1) CINR=CPROD1(1,1)/CPROD2(1,1)
C+++INR NOT NORMALIZED FOR INRNOR=0
                                                                         DIP11600
      IF(INRNOR.EQ.0) CINR=CPROD1(1,1)
      DBINR=10. *DLOG10(CDABS(CINR))
                                                                         DIP11610
                                                                         DIP11620
      RETURN
                                                                         DIP11630
      END
      SUBROUTINE CHMULT(A,B,L,M,N,C)
                                                                         DIP11640
                                                                         DIP11650
      COMPLEX *16 A(L,M),B(M,N),C(L,N)
                                                                         DIP11660
      DO 20 I=1.L
                                                                         DIP11670
      DO 20 J=1,N
                                                                         DIP11680
      C(I,J)=DCMPLX(0.0D0,0.0D0)
                                                                         DIP11690
      DO 20 K=1,M
      C(I,J)=C(I,J)+(A(I,K)*B(K,J))
                                                                         DIP11700
                                                                         DIP11710
    CONTINUE
                                                                         DIP11720
      RETURN
                                                                         DIP11730
C***SUBROUTINE TO COMPUTE CONJUGATE TRANSPOSE OF A MATRIX
                                                                         DIP11740
      SUBROUTINE CONJTR(A,L,M,ACTR)
                                                                         DIP11750
                                                                         DIP11760
      COMPLEX *16 A(L,M), ACTR(M,L)
                                                                         DIP11770
      DO 10 I=1,M
                                                                         DIP11780
      DO 10 J=1,L
      ACTR(1,J)=DCONJG(A(J,I))
                                                                         DIP11790
                                                                         DIP11800
      WRITE(12,7766)I,J,A(J,I),ACTR(I,J)
7766 FORMAT(1X,'I,J,A(J I),ACTR(I,J)=',2X,2I4,2X,4E12.5)
                                                                         DIP11810
                                                                         DIP11820
      CONTINUE
      RETURN
                                                                         DIP11830
                                                                         DIP11840
      END
      SUBROUTINE DCMINV(A,L,M,IDM,NEQ)
                                                                         DIP11850
                                                                         DIP11860
      COMPLEX *16 A(IDM,IDM),BIGA,HOLD,DET
      INTEGER L(IDM),M(IDM)
                                                                         DIP11870
                                                                         DIP11880
      N=NEQ
                                                                         DIP11890
      DET=DCMPLX(1.0D0,0.0D0)
      DO 80 K=1,N
                                                                         DIT11900
                                                                         DIP11910
      L(K)≈K
                                                                         DIP11920
      M(K)=K
                                                                         DIP11930
      BIGA=A(K,K)
                                                                         DIP11940
      DO 20 J=K,N
                                                                         DIP11950
      DO 20 I=K,N
      IF(CDABS(BIGA)-CDABS(A(I,J)))15,19,19
                                                                         DIP11960
10
                                                                         DIP11970
      BIGA=A(I,J)
      L(K)=I
                                                                         DIP11980
                                                                         DIP11990
      M(K) = J
                                                                         DIP12000
      CONTINUE
19
                                                                         DIP12010
      CONTINUE
20
                                                                         DIP12020
      J=L(K)
      IF(J-K)35,35,25
                                                                         DIP12030
                                                                         DIP12040
25
      CONTINUE
      DO 30 I=1,N
                                                                         DIP12050
                                                                         DIP12060
      HOLD = -A(K,I)
      A(K,I) = A(J,I)
                                                                         DIP12070
                                                                         DIP12080
30
      A(J,I)=HOLD
      I=M(K)
                                                                         DIP12090
35
                                                                         DIP12100
      IF(I-K) 45,45,38
      CONTINUE
                                                                         DIP12110
38
                                                                         DIP12120
      Đũ 10 J=1,N
```

	HOLD=-A(J,K)	DIP12130
	$\mathbf{A}(\mathbf{J},\mathbf{K}) = \mathbf{A}(\mathbf{J},\mathbf{I})$	DIP12140
40	A(J,I)=HOLD	DIP12150
45	CONTINUE	DIP12160
	DO 55 I=1,N	DIP12170
	IF(I-K)50,55,50	DIP12180
50	A(I,K)*A(I,K)/(-BIGA)	DIP12190
55	CONTINUE	DIP12200
	DO 65 I=1,N	DIP12210
	DO 65 J=1,N	DIP12220
	IF(I-K)60,64,60	DIP12230
60	IF(J-K)62,64,62	DIP12240
62	A(I,J)=A(I,K)+A(K,J)+A(I,J)	DIP12250
64	CONTINUE	DIP12260
65	CONTINUE	DIP1227C
	DO 75 J=1,N	DIP12280
	IF(J-K)70,75,70	DIP12290
70	A(K,J)=A(K,J)/BIGA	DIP12300
75	CONTINUE	DIP12310
	DET*DET*BIGA	DIP12320
	A(K,K)=1.ODO/BIGA	DIP12330
80	CONTINUE	DIP12340
	K=N	DIP12350
100	K=K-1	DIP12360
	IF(K)150,150,105	DIP12370
105	I=L(K)	DIP12380
	IF(I-K)120,120,108	DIP12390
108	CONTINUE	DIP12400
	DO 110 J=1,N	DIP12410
	HOLD=A(J,K)	DIP12420
	$\mathbf{A}(\mathbf{J},\mathbf{K}) = \mathbf{A}(\mathbf{J},\mathbf{I})$	DIP12430
110	A(J,I)=HOLD	DIP12440
120	J≖M(K)	DIP12450
	IF(J-K)100,100,125	DIP12460
125	CONTINUE	DIP12470
	DO 130 I=1,N	DIP12480
	HOLD=A(K,I)	DIP12490
	$\mathbf{A}(\mathbf{K},\mathbf{I}) = -\mathbf{A}(\mathbf{J},\mathbf{I})$	DIP12500
130	A(J,I)=HOLD	DIP12510
	GO TO 100	DIP12520
150	RETURN	DIP12530
	END	DIP12540
	UBROUTINE CNORM	DIP12550
	OMPUTES A NORMALIZED COMPLEX COLUMN VECTOR SUCH THAT THE	DIP12560
C***M	AXIMUM ELEMENT HAS UNITY MAGNITUDE.	DIP12570
С		DIP12580
	RINTS MAGNITUDE AND PHASE OF NORMALIZED VECTOR	DIP12590
С		DIP12600
	SUBROUTINE CNORM(V,N)	DIP12610
C***	V IS THE INPUT COMPLEX COLUMN VECTOR	DIP12620
C***	N IS THE LENGTH OF Y	DIP12630
	COMPLEX V(1),SS	DIP12640
	CNOR=0.0	DIP12650
	DO 10 K=1,N	DIP12660
	SA=CABS(V(K))	DIP12670
	IF(SA.GT CNOR) CNOR=SA	DIP12680
10	CONTINUE	DIP12690
	IF(CNOR.LE.O.) CNOR=1.0	DIP12700
	DO 30 K=1,N	DIP12710

```
SS=V(K)
                                                                         DIP12720
      SA=CABS(SS)
                                                                         DIP12730
      SNOR=SA/CNOR
                                                                         DIP12740
      PHR=0.
                                                                         DIP12750
      IF(SA.GT.O.) PHR=ATAN2(AIMAG(SS), REAL(SS))
                                                                         DIP12760
      PH=57.29578*PHR
                                                                         DIP12770
      WRITE(8,20) K, SNOR, SA, PH
                                                                         DIP12780
                                                                         DIP12790
20
     FORMAT(11,15,F10.6,31,E15.3,F10.0)
30
      CONTINUE
                                                                         DIP12800
      RETURN
                                                                         DIP12810
      END
                                                                         DIP12820
      SUBROUTINE GETCP2(RCPU)
      REAL ETIME, TARRAY(2), RCPU
      TIME=ETIME(TARRAY)
      RCPU=TARRAY(1)
      RETURN
      END
C***SUBROUTINE TO COMPUTE INDUCED VOLTAGE BETWEEN PROBE AND
                                                                         YD100010
C***DIPOLE ARRAY ELEMENTS
C***RING ARRAY
      SUBROUTINE NFDPC2(XP,YP,ZP,XREF,IREF,V,VREF)
C***DISTANCES ARE IN METERS
      COMPLEX V(1), VD, VR1, VREF
      COMPLEX *16 DZABG
                                                                         YD100060
      COMMON/A/DX,DY,NCOLX,NROWY,NEL,HZ,HL,ARAD,ZLOAD,ZCHAR
                                                                         YD100070
      COMMON /CIRCLE/ ICIRC, RADIUS, HLS
C***NOTE: AA,BB,... FOR PROBE
                                11,22,... FOR DIPOLE
                                                                         YDI00120
      PI=3.1415926535
      DCR=PI/180.
C***PROBE DIMENSIONS
      XAA=XP
                                                                         YD100260
      IBB=IP
                                                                         YD100270
      XCC=XP
                                                                         YD100280
      YAA=YF-HLS
      YBB=YP
                                                                         YD100300
      YCC=YP+HLS
785 ZAA=ZP
      ZBB*ZP
                                                                         YD100330
      ZCC=ZP
      IC=0
                                                                         YD100430
C****CIRCULAR DIPOLE ARRAY ELEMENTS
     DELPHI=360./NEL
      DO 10 IX=1, NEL
      PHI=DELPHI*(IX-1)
      ID=RADIUS+COS(PHI+DCR)
      ZD=RADIUS*SIN(PHI*DCR)
      DIST=SQRT((XP-XD)**2+(ZP-ZD)**2)
      IF(DIST.LT.ARAD) XD=XD+ARAD
      X11=XD
      122=1D
                                                                         YD100490
      X33=XD
      Y11=YP-HL
      Y22=YP
      Y33=YP+HL
      Z11=ZD
      Z22=ZD
      Z33=ZD
9446 IC=IC+1
                                                                         YD100650
    WRITE(6,2233)IX,IY,XD,YD
                                                                         YD100660
```

```
2233 FORMAT(1X,'IX,IY,XD,YD=',2X,2I4,2X,2F12.5)
                                                                           YD100670
      WRITE(6,7854)XAA,XBB,XCC,YAA,YBB,YCC,ZAA,ZBB,ZCC
7854 FORMAT(1X, 'XYZABC=',9E12.5)
CC
      WRITE(6,7855)X11,X22,X33,Y11,Y22,Y33,Z11,Z22,Z33
7855 FORMAT(1X,'XYZ123=',9E12.5)
      CALL DSZABG(XAA, XBB, XCC, YAA, YBB, YCC, ZAA, ZBB, ZCC, X11, X22, X33,
                                                                           YD100680
     2Y11,Y22,Y33,Z11,Z22,Z33,DZABG)
                                                                           YD100690
      VD=DZABG
                                                                           YD100700
CC
      WRITE(6,33)IC,VD
      FORMAT(1X,'IC, VD=',2X, I5, 2E12.5)
      V(IC)=VD
      VAMPDB=20. *ALOG10(CABS(V(IC)))
                                                                           YD100780
      VPHASE=ATAN2(AIMAG(V(IC)), REAL(V(IC))) + 180./3.141592654
                                                                           YD100790
      WRITE(6,4455)IC, VAMPDB, VPHASE
                                                                           YD100800
4455 FORMAT(1X,'IC, VAMDB, VPHASE*', 2X, I4, 2X, 2F12.2)
                                                                           YD100810
                                                                           YD100820
10
      CONTINUE
C *** COMPUTE REFERENCE VOLTAGE (FICTITIOUS ELEMENT AT IREF, Y=0)
                                                                           YD100830
C***SET IREF=RADIUS, THIS MAKES ELEMENT 1 THE REFERENCE
      IREF=RADIUS
C***SET IREF=ARAD, (REF. CLOSE TO THE ORIGIN)
      IREF=ARAD
      X11=XREF
      I22≈IREF
                                                                           YD100870
      X33=XREF
      Y11=-HL
      Y22=0.0
                                                                           YD100980
      Y33=HL
      Z11=0.0
      Z22=0.0
      Z33=0.0
9267 CALL DSZABG(XAA, XBB, XCC, YAA, YBB, YCC, ZAA, ZBB, ZCC, X11, X22, X33,
                                                                           YDI01000
     2Y11, Y22, Y33, Z11, Z22, Z33, DZABG)
                                                                           YDI01010
                                                                           YDI01020
      VR1=DZABG
      VREF=VR1
      RETURN
                                                                           YD101080
99
                                                                           YDI01090
C***SUBROUTINE TO COMPUTE MUT. IMPED. BETWEEN STRAIGHT DIPOLES
C***ARRANGED IN A RING ARRAY (CIRCLE)
      SUBROUTINE CADZA2(CTDP,STDP,IGRNDP,ICC1,ICC,Z)
      COMPLEX ZMA, ZABG, Z(ICC1, ICC)
      COMMON /A/ DX, DY, NCOLX, NROWY, NEL, HZ, HL, ARAD, ZLOAD, ZCHAR
      COMMON /CIRCLE/ ICIRC, RADIUS, HLS
C***ALL DIMENSIONS IN METERS
      PI=3.1415926535
      DCR=PI/180.
C***FIXED POSITION FOR ELEMENT 1
      XAA=RADIUS
      IBB=RADIUS
      XCC=RADIUS
      YAA=-HL
      YBB=0.0
      YCC=HL
      ZAA=0.0
      ZBB=0.0
      ZCC=0.0
     Y11=YAA
      Y22=YBB
      Y33=YCC
```

```
C***VARIABLE POSITION FOR RING ARRAY ELEMENTS
      DELPHI=360./NEL
      IC=0
      DO 20 I=1,NEL
      IC=IC+1
      PHI=DELPHI+(I-1)+DCR
      ID=RADIUS*COS(PHI)
      ZD=RADIUS*SIN(PHI)
      X11=XD
      122=ID
      X33=XD
      IF(I.EQ.1) X11=X11+ARAD
      IF(I.EQ.1) X22=X22+ARAD
      IF(I.EQ.1) X33=X33+ARAD
      Z11=ZD
      722=7D
      Z33=ZD
CC
     WRITE(6,87)IC
87
     FORMAT(1X,'IC=',I5)
CC
      WRITE(6,88)XAA,XBB,XCC,YAA,YBB,YCC,ZAA,ZBB,ZCC
     FORMAT(1X,'XYZABC=',9E12.4)
88
CC
      WRITE(6,89) X11, X22, X33, Y11, Y22, Y33, Z11, Z22, Z33
     FORMAT(1X,'XYZ123=',9E12.4)
89
      ZMA=ZABG(XA:, XBB, XCC, YAA, YBB, YCC, ZAA, ZBB, ZCC, X11, X22, X33,
     2Y11, Y22, Y33, Z11, Z22, Z33)
      Z(1,IC)=ZMA
      IF(IC.EQ.1) Z(1,1)=Z(1,1)+ZLOAD
      CONTINUE
20
      DO 40 I=1,NEL
      IF(I.LT.9)WRITE(6,30)I,Z(1,I)
      WRITE(8,30)I,Z(1,I)
      FORMAT(1X,'Z(1,',I4,')=',2E12.5)
30
     CONTINUE
40
      RETURN
      END
ZAB00010
C***PROGRAM TO CALCULATE MUTUAL IMPEDANCE BETWEEN TWO DIPOLES
C***WITH ARBITRARY LENGTH AND ORIENTATION. A PIECEWISE-
                                                                       ZAB00020
C***SINUSOIDAL CURRENT DISTRIBUTION IS ASSUMED.
                                                                        ZAB00030
      COMPLEX FUNCTION ZABG(X1,X2,X3,Y1,Y2,Y3,Z1,Z2,Z3,XA,XB,XC,YA,YB,YCZAB00040
     2,ZA,ZB,ZC)
                                                                        ZAB00050
      COMPLEX P11,P12,P21,P22,Q11,Q12,Q21,Q22,R11,R12,R21,R22
                                                                       ZAB00060
      COMPLEX S11,S12,S21,S22,JCOM,GAM,CGDS,SGDS,SGDT,ETA,EP3
                                                                       ZAB00070
      COMPLEX EGDS, EGDT
      COMMON /F/ FHZ.ER3.SIG3.TD3
C***ALL DIMENSIONS IN METERS
      PI=3.141592654
      TPI=2.*PI
      B=TPI
                                                                        ZAB00090
      JCOM=(0.,1.)
      E0=8.854E-12
      UO=1.2566E 3
      OMEGA=TPI+FHZ
                                                                       WIR01670
      IF(SIG3.LT.0.)EP3=ER3*E0*CMPLX(1.,-TD3)
                                                                       WIR01720
      IF(TD3.LT.0.)EP3=CMPLX(ER3+E0,-SIG3/OMEGA)
                                                                        WIR01730
                                                                       WIR01740
      ETA=CSQRT(UO/EP3)
      GAM=OMEGA+CSQRT(-U0+EP3)
                                                                       WIR01750
```

```
IF(CABS(GAM+AM).GT.0.06) WRITE(6,7923)AM
7923 FORMAT(11, 'CABS(GAM*AM) IS GREATER THAN 0.06, AM=', E14.5)
      TNT=0
                                                                            ZAR00130
      XBA=XB-IA
                                                                            ZAB00140
      YBA=YB-YA
                                                                            ZAB00150
      784=78-74
                                                                            ZAB00160
      X21=X2-X1
                                                                            ZAB00170
      Y21=Y2-Y1
                                                                            Z4R00180
      Z21=Z2-Z1
                                                                            ZAB00190
      DS=SQRT(XBA+XBA+YBA+YBA+ZBA+ZBA)
                                                                            ZAB00200
      DT=SQRT(X2.*X21+Y21*Y21+Z21*Z21)
                                                                            ZAB00210
CC
      DSK=B+DS
                                                                            ZAB00220
                                                                            ZAB00230
CC
     DTK=B+DT
     CGDS=CMPLX(COS(DSK),0.0)
                                                                            ZAB00240
CC
CC
                                                                            74500250
     SGDS=CMPLX(0.0,SIN(DSK))
     SGDT=CMF (0.0.SIN(DTK))
                                                                            ZAB00260
C***FOR LOSSY MEDIUM THE NEXT LINES ARE APPROPRIATE
      EGDS=CEXP(GAM+DS)
      EGDT=CEXP(GAM+DT)
      CGDS=(EGDS+1./EGDS)/2.
      SGDS=(EGDS-1./EGDT)/2.
      SGDT=(EGDT-1./EGDT)/2.
     WRITE(6,1345)XA,XB,YA,YB,ZA,ZB
1345 FORMAT(1X,'XYZAB=',6E14.5)
      WRITE(6,1346)X1,X2,Y1,Y2,Z1,Z2
1346 FORMAT(1X,'XYZ12=',6E14.5)
      CALL GGS(XA,YA,ZA,XB,YB,ZB,X1,Y1,Z1,X2,Y2,Z2,AM.
                                                                            ZAB00270
     2DS, CGDS, SGDS, DT, SGDT, INT, ETA, GAM, P11, P12, P21, P22)
                                                                            ZAB00280
      CALL GGS(XA, YA, ZA, XE, YB, ZB, X2, Y2, Z2, X3, Y3, Z3, AM,
                                                                            Z4B00290
     2DS, CGDS, SGDS, DT, SGDT, INT, ETA, GAM, Q11, Q12, Q21, Q22)
                                                                            ZAB00300
      CALL GGS(XB, YB, ZB, XC, YC, ZC, X1, Y1, Z1, X2, Y2, Z2, AM,
                                                                            ZAB00310
     2DS, CGDS, SGDS, DT, SGDT, INT, ETA, GAM, R11, R12, R21, R22)
                                                                            ZAB00320
      CALL GGS(XB,YB,ZB,XC,YC,ZC,X2,Y2,Z2,X3,Y3,Z3,AM,
                                                                            ZAB00330
     2DS, CGDS, SGDS, DT, SGDT, INT, ETA, GAM, S11, S12, S21, S22)
                                                                            ZAB00340
      ZABG=P22+Q21+R12+S11
                                                                            ZAB00350
CC
      WRITE(6.7898)ZABG
7898 FORMAT(1X, 'EXITING ZABG WITH ZABG=', 2E14.5)
      RITURN
                                                                            ZAB00360
      END
                                                                            ZAB00370
*********file dzabgnloss.f******
C****DOUBLE PRECISION VERSION
C***PROGRAM TO CALCULATE MUTUAL IMPEDANCE BETWEEN TWO DIPOLES
                                                                            DZA00020
                                                                            DZA00030
C***WITH ARBITRARY LENGTH AND ORIENTATION. A PIECEWISE-
C***SINUSOIDAL CURRENT DISTRIBUTION IS ASSUMED.
                                                                            DZA00040
      SUBROUTINE DSZABG(SX1,SX2,SX3,SY1,SY2,SY3,SZ1,SZ2,SZ3,
                                                                            DZ400050
     2SXA, SXB, SXC, SYA, SYB, SYU, SZA, SZB, SZC, DZABG)
                                                                            DZA00060
      IMPLICIT REAL+8 (A-11,0-Z)
                                                                            DZA00070
      COMPLEX*16 P11,P12,P21,P22,Q11,Q12,Q21,Q22,R11,R12,R21,R22
                                                                            DZ#00080
      COMPLEX*16 S11,S12,S21,S22,JCOM,GAM,CGDS,SGDS,SGDT,ETA
                                                                            DZA00090
      COMPLEX * 16 DZABG, EP3, EGDS, EGDT
      REAL*4 SX1, SX2, SX3, SY1, SY2, SY3, SZ1, SZ2, SZ3
                                                                            DZA00110
      REAL *4 SXA, SXB, SXC, SYA, SYB, SYC, SZA, SZB, SZC
                                                                            DZ400120
      REAL+4 FHZ, ER3, SIG3, TD3
      COMMON /F/ FHZ, ER3, SIG3, TD3
      JCOM=(0.D0.1.D0)
                                                                            DZ400140
      PI=3.1415926535898D0
      TPI=2.0D0*PI
```

```
B=TPI
      E0=8.854D-12
      HO=1.2566D-6
                                                                          WIR01670
      OMEGA=TPI+FHZ
                                                                          WIR01680
CC
      WRITE(6,2843)OMEGA
                                                                          W1801690
2843 FORMAT(1X.'OMEGA=',E12.5)
      IF(SIG3.LT.0.ODO)EP3=ER3*EO*DCMPLX(1.ODO,-TD3)
                                                                          WIR01720
                                                                          WIR01730
      IF(TD3.L7.0.OD0)EP3=DCMPLX(ER3=E0,-SIG3/OMEGA)
      WRITE(6,7755)ER3,E0,EP3
7755 FORMAT(1X, 'ER3, EO, EP3=', 4E12.5)
                                                                          WIR01740
      ETA=CDSQRT(UO/EP3)
      GAM=DMEGA+CDSQRT(-U0+EP3)
                                                                          WIR01750
C**COMPUTE GAMMA BY EQUATION IN HAYT PAGE 333
      GAM=JCOM+OMEGA+CDSQRT(UO+EP3)+CDSQRT(1.DO-JCOM+SIG3/(JMEGA+EP3))
CC
      WRITE(6,8888)GAM
CC
8888 FORMAT(1X,'HYAT GAM=',2E12.5)
      AM=0.0001D0
                                                                          DZA00180
      INT=0
                                                                          DZA00190
      X1=SX1
                                                                          DZA00200
      12=SX2
                                                                          DZA00210
      X3=SX3
                                                                          DZ#00220
      Y1=SY1
                                                                          DZA00230
      Y2=SY2
                                                                          DZ400240
      Y3=SY3
                                                                          DZA00250
      Z1=SZ1
                                                                          DZA00260
      72=572
                                                                          DZA00270
      Z3=SZ3
                                                                          D7A00280
      IA-SIA
                                                                          DZ400290
      IB=SIB
                                                                          DZA00300
      IC=SIC
                                                                          DZA00310
      YA-SYA
                                                                          DZA00320
      YB=SYB
                                                                          DZA00330
      VC=SYC
                                                                          DZ400340
      ZA=SZA
                                                                           DZA00350
      ZB=SZB
                                                                          DZA00360
      ZC=SZC
                                                                          DZA00370
      IBA=XB-XA
                                                                          DZ400380
      YBA=YB-YA
                                                                           DZA00390
      ZBA=ZB-ZA
                                                                           DZA00400
      X21=X2-X1
                                                                           DZA00410
      Y21=Y2-Y1
                                                                           DZA00420
      Z21=Z2~Z1
                                                                           DZ400430
      DS*DSQRT(XBA*XBA+YBA*YBA+ZBA*ZBA)
                                                                           DZA00440
      DT=DSQRT(X21+X21+Y21+Y21+Z21+Z21)
                                                                           DZA00450
      DSK=B+DS
CC
                                                                           DZA00460
      DTK=B+DT
      CGDS=DCMPLX(DCOS(DSK),0.0D0)
                                                                           DZ400470
CC
                                                                           DZA00480
      SGDS=DCMPLX(0.0D0,DSIN(DSK))
CC
                                                                           DZA00490
      SGDT=DCMPLX(0.0D0,DSIN(DTK))
CC
      EGDS=CDEXP(GAM*DS)
      EGDT=CDEXP(GAM+DT)
      CGDS=(EGDS+1.0D0/EGDS)/2.0D0
      SGDS=(EGDS-1.0D0/EGDS)/2.0D0
      SCDT=(EGDT-1.0D0/EGDT)/2.0D0
                                                                           DZA00500
      CALL DGGS(XA, YA, ZA, XB, YB, ZB, X1, Y1, Z1, X2, Y2, Z2, AM,
     2DS,CGDS,SGDS,DT,SGDT,INT,ETA,GAM,P11,P12,P21,P22)
                                                                           DZA00510
                                                                           DZA00520
      CALL DGGS(XA,YA,ZA,XB,YB,ZB,X2,Y2,Z2,X3,Y3,Z3,AM,
     2DS,CGDS,SGDS,DT,SGDT,INT,ETA,GAM,Q11,Q12,Q21,Q22)
                                                                           DZA00530
```

	CALL DGGS(XB,YB,ZB,XC,YC,ZC,X1,Y1,Z1,X2,Y2,Z2,AM,	DZA00540
	2DS, CGDS, SGDS, DT, SGDT, INT, ETA, GAM, R11, R12, R21, R22)	DZA00550
	CALL DGGS(IB,YB,ZB,XC,YC ZC,X2,Y2,Z2,X3,Y3,Z3,AM,	DZAOC560
	2DS, CGPS, SGDS, DT, SGDT, INT, ETA, GAM, S11, S12, S21, S22)	DZA00570
	DZABG=P22+Q21+R12+511	DZA00580
CC	WRITE(6,8899)DZABG	
8899	FORMAT(1X, 'EXITING DZABG, DZA'3G=',2E14.5)	
	RETURN	DZA00590
	END	DZA00600

#### REFERENCES

- 1. C.A. Perez and L.W. Brady, *Principles and Practice of Radiation Oncology*, J.B. Lippincott Co., Philadelphia, 1987.
- 2. A.W. Guy, "History of biological effects and medical applications of microwave energy," *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-32, No. 9, pp. 1182-1200, September 1984.
- 3. J. Overgaard, "The effect of local hyperthermia alone, and in combination with radiation on solid tumors," in Cancer Therapy by Hyperthermia and Radiation, Proceedings of the 2<sup>nd</sup> International Symposium. C. Streffer (Editor), 2-4 June 1977, pp. 49-61, Urban & Schwarzenberg, Inc., Baltimore, 1978.
- 4. J. Overgaard, "Clinical hyperthermia, an update," Proceedings of the 8th International Congress of Radiation Research, Vol. 2, pp. 942-947, July 1987.
- 5. "Special issue on phased arrays for hyperthermia treatment of cancer," *IEEE Transactions* on Microwave Theory and Techniques, Vol. MTT-34, No. 5, May 1986.
- 6. "Special issue on hyperthermia and cancer therapy," *IEEE Transactions on Biomedical Engineering*, Vol. BME-31, No. 1, January 1984.
- 7. V. Sathiaseelan, M.F. Iskander, G.C.W. Howard, and N.M. Bleehen. "Theoretical analysis and clinical demonstration of the effect of power control using the annular phased-array hyperthermia system." *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-34, No. 5, pp. 514-519. May 1986.
- 8. V. Sathiaseelan. "Potential for patient-specific optimization of deep heating patterns through manipulation of amplitude and phase," *Strahlentherapie Onkologie*, Vol. 165, No. 10, pp. 743-745, 1989.
- 9. D. Sullivan. "Three-dimensional computer simulation in deep regional hyperthermia using the finite-difference time-domain method," *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-38, No. 2, pp. 204-211, Feb. 1990.
- 10. B.S. Trembly, A.H. Wilson, M.J. Sullivan, A.D. Stein, T.Z. Wong, and J.W. Strohbehn, "Control of the SAR pattern within an interstitial microwave array through variation of antenna driving phase," *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-34, No. 5, pp. 568-571, May 1986.
- 11. G. Sato. C. Shibata. S. Sekimukai. H. Wakabayashi, K. Mitsuka, and K. Giga, "Phase-controlled circular array heating equipment for deep-seated tumors: preliminary experiments," *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-34, No. 5, pp. 520-525, May 1986.

- 12. P.A. Cudd, A.P. Anderson, M.S. Hawley, and J. Conway, "Phased-array design considerations for deep hyperthermia through layered tissue," *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-34, No. 5, pp. 526-531, May 1986.
- N. Morita, T. Hamasaki, and N. Kumagai, "An optimal excitation method in multiapplicator systems for forming a hot zone inside the human body," *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-34, No. 5, pp. 532-538, May 1986.
- C. De Wagter, "Optimization of simulated two-dimensional temperature distributions induced by multiple electromagnetic applicators," *IEEE Transactions on Microwave Theory and Tech*niques, Vol. MTT-34, No. 5, pp. 589-596, May 1986.
- 15. M. Knudsen and U. Hartmann, "Optimal temperature control with phased-array hyperthermia system," *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-34, No. 5, pp. 597-603, May 1986.
- C.F. Babbs, V.A. Vaguine, and J.T. Jones, "A predictive-adaptive, multipoint feedback controller for local heat therapy of solid tumors," *IEEE Transactions on Microwave Theory and Techniques*, Vol. MTT-34, No. 5, pp. 604-611, May 1986.
- 17. R.B. Roemer, "Physical and engineering aspects of hyperthermia," *Proceedings of the 8th International Congress of Radiation Research*, Vol. 2, pp. 948-953, July 1987.
- R.B. Roemer, K. Hynynen, C. Johnson, R. Kress, "Feedback control and optimization of hyperthermia heating patterns: present status and future needs," *IEEE Eighth Annual Con*ference of the Engineering in Medicine and Biology Society, pp. 1496-1499, November 1986.
- A. Boag and Y. Leviatan, "Optimal excitation of multiapplicator systems for deep regional hyperthermia." *IEEE Transactions on Biomedical Engineering*, Vol. BME-37, No. 10, pp. 927-995, October 1990.
- J.T. Loane III and S.W. Lee, "Gain optimization of a near-field focusing array for hyperthermia applications," *IEEE Transactions on Microwave Theory and Techniques*, Vol. 37, No. 10, pp. 1629-1635, October 1989.
- 21. P.F. Turner, T. Schaefermeyer, and T. Saxton, "Future trends in heating technology of deep-seated tumors," Recent Results in Cancer Research, Vol. 107, pp. 249-262, 1988.
- 22. P.F. Turner, A. Tumeh, T. Schaefermeyer, "BSD-2000 approach for deep local and regional hyperthermia: physics and technology," *Strahlentherapie Onkologie*, Vol. 165, No. 10, pp. 738-741, 1989.
- 23. R.T. Compton, Jr., Adaptive Antennas, Concepts and Performance, Prentice Hall, New Jersey, 1988.

- 24. A. J. Fenn. "Theory and analysis of near-field adaptive nulling." 1986 IEEE Antennas and Propagation Society International Symposium Digest, Vol. 2, IEEE, New York, pp. 579-582, 1986.
- 25. A. J. Fenn, "Theory and analysis of near-field adaptive nulling," 1986 Asilomar Conf. on Signals, Systems and Computers, Computer Society Press of the IEEE, Washington, D.C., pp. 105-109, 1986.
- 26. A. J. Fenn, "Theoretical near-field clutter and interference cancellation for an adaptive phased-array antenna." 1987 IEEE Antennas and Propagation Society International Symposium Digest, Vol. 1, IEEE, New York, pp. 46-49, 1987.
- 27. A.J. Fenn. "Moment-method analysis of near-field adaptive nulling." *IEEE Sixth International Conference on Antennas and Propagation, ICAP 89*, 4-7 April 1989, pp. 295-301.
- 28. A.J. Fenn. "Evaluation of adaptive phased-array antenna far-field nulling performance in the near-field region." *IEEE Transactions on Antennas and Propagation*, Vol. 38, No. 2, pp. 173-185. February 1990.
- A.J. Fenn, "Near-field testing of adaptive radar systems," Proc. 12th Annual Antenna Measurement Techniques Association Meeting and Symposium. 8-11 October 1990, pp. 13 9-13 14.
- A.J. Fenn, H.M. Aumann, F.G. Willwerth, and J.R. Johnson. "Focused near-field adaptive nulling: experimental investigation." 1990 IEEE Antennas and Propagation Society International Symposium Digest, Vol. 1, 7-11 May 1990, pp. 186-189.
- 31. M.I. Skolnik. *Introduction to Radar Systems*, Second Edition, McGraw-Hill Book Company, 1980.
- 32. J.W. Strohbehn and R.B. Roemer, "A survey of computer simulations of hyperthermia treatments," *IEEE Transactions on Biomedical Engineering*, Vol. BME-31, No. 1, pp. 136-149, January 1984.
- 33. K.B. Ocheltree and L.A. Frizzell, "Determination of power deposition patterns for localized hyperthermia: a transient analysis," *International Journal of Hyperthermia*, Vol. 4, No. 3, pp. 281-296, 1988.
- 34. C.K. Chou, G.W. Chen, A.W. Guy, and K.H. Luk, "Formulas for preparing phantom muscle tissue at various radio frequencies," *Bioelectromagnetics*, Vol. 5, No. 4, pp. 435-441, 1984.
- 35. Y. Zhang. W.T. Joines, J.R. Oleson, "The calculated and measured temperature distribution of a phased interstitial antenna array," *IEEE Transactions on Microwave Theory and Techniques*, Vol. 38, No. 1, pp. 69-77, January 1990.

- J.H. Richmond, "Radiation and scattering by thin-wire structures in the complex frequency domain." Ohio State University, ElectroScience Laboratory, Technical Report 2902-10, July 1973.
- J.H. Richmond, "Computer program for thin-wire structures in a homogeneous conducting medium," Ohio State University, ElectroScience Laboratory, Technical Report 2902-12, August 1973.
- 38. J.H. Richmond, "Radiation and scattering by thin-wire structures in a homogeneous conducting medium (computer program description)," *IEEE Trans. Antennas Propagat.*, Vol. AP-22, no. 2, p. 365, March 1974.
- 39. A.J. Fenn. "Focused near-field nulling for adaptive electromagnetic hyperthermia applications," submitted to 1991 Progress in Electromagnetics Research Symposium, 1-5 July 1991.
- 40. A.J. Fenn, "Adaptive hyperthermia for improved thermal dose distribution," submitted to 9th International Congress of Radiation Research, 7-12 July 1991.
- 41. J.R. Johnson, A.J. Fenn, H.M. Aumann, and F.G. Willwerth, "An experimental adaptive nulling receiver utilizing the sample matrix inversion algorithm with channel equalization," *IEEE Transactions on Microwave Theory and Techniques*, Vol. 39, No. 5, pp. 798-808 May 1991.
- 42. G. Strang, Linear Algebra and Its Applications, New York: Academic, 1976.
- 43. J.T. Mayhan, "Some techniques for evaluating the bandwidth characteristics of adaptive nulling systems," *IEEE Trans. Antennas Propagat.*, Vol. AP-27, No. 3, pp. 363-373, May, 1979.
- 44. A.J. Fenn, "Maximizing jammer effectiveness for evaluating the performance of adaptive nulling array antennas," *IEEE Trans. Antennas Propagat.*, Vol. AP-33, No. 10, pp. 1131-1142, Oct. 1985.
- 45. W.L. Stutzman and G.A. Thiele, Antenna Theory and Design, New York: Wiley, 1981.
- 46. I.J. Gupta and A.A. Ksienski, "Effect of mutual coupling on the performance of adaptive arrays," *IEEE Trans. Antennas Propagat.*, Vol. AP-31, No. 5, pp. 785-791
- 47. A.J. Fenn. "Moment-method analysis of near-field adaptive nulling," MIT Lincoln Laboratory, Technical Report 842, 7 April 1989.
- 48. J.D. Kraus, Antennas, 2nd Edition, McGraw-Hill Book Company, New York, 1988.
- 49. D.T. Paris and F.K. Hurd, Basic Electromagnetic Theory, Chapter 7, McGraw-Hill Book Company, New York, 1969.

- 50. S.A. Schelkunoff and H.T. Friis, Antennas: Theory and Practice, Wiley, New York, 1952.
- 51. "Program manual for the simplified space payload thermal analyzer (Version 3.0/VAX)," Arthur D. Little, Inc., Cambridge, Massachusetts, Prepared for National Aeronautics and Space Administration, Contract NAS5-27606.
- 52. J.T. Bartoszek, B. Huckins, M. Coyle, "A simplified shuttle payload thermal analyzer (SSPTA) program," AIAA 14th Thermophysics Conference, June 4-6, 1979
- 53. "Thermal Network Modeling Handbook," TRW Systems Group, prepared for National Aeronautics and Space Administration, Contract NAS 9-10435, January 1972.
- D. Sullivan. "Mathematical methods for treatment planning in deep regional hyperthermia,"
   *IEEE Transactions on Microwave Theory and Techniques*, Vol. 39, No. 5, pp. 864-872, May 1991.
- 55. R.C. Johnson and H. Jasik, Antenna Engineering Handbook, Second Edition, Chapter 22, McGraw-Hill Book Company, New York, 1984.

#### REPORT DOCUMENTATION PAGE

Form Approved OMB No. 0704-0188

Pubic reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching evision data sources, gathering and maintaining the data needed and complisting and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, so Washington Headquarters Services. Directorates for information Operations and Reports, 1215 Jefferson Davis Highway, Surfe 1204, Arlington, VA 22202-4302, and to the Office of Management and Budget, Paperwork Reduction Project (0704-0188) Washington, DC 20503

AGENCY USE ONLY (Leave blank)	2. REPORT DATE 3 July 1991		3. REPORT TYPE AND DATES COVERED Technical Report	
4. TITLE AND SUBTITLE			5. FUNDING NUMBERS	
Application of Adaptive Nulling to E Dose Distribution in Cancer Therapy		for Improved Thermal	G. Disconding	
6. AUTHOR(S)			C — F19628-90-0002 PE — 63250F	
			PR — 227	
Alan J. Fenn		-		
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES)			8. PERFORMING ORGANIZATION REPORT NUMBER	
Lincoln Laboratory, MIT			<b>7</b> 7 0.0	
P.O. Box 73 Lexington, MA 02173-9108			TR-917	
9. SPONSORING/MONITORING AGEN	CY NAME(S) AND ADDRESS(E	S)	10. SPONSORING/MONITORING AGENCY REPORT NUMBER	
U.S. Air Force Space Systems Division			ESD TR OLOTO	
PO Box 92960			ESD-TR-91-079	
Los Angeles, CA 90009-2960				
11.SUPPLEMENTARY NOTES				
None				
12a. DISTRIBUTION/AVAILABILITY ST	ATEMENT		12b. DISTRIBUTION CODE	
Approved for public release; distribution is unlimited.				
13. ABSTRACT (Maximum 200 words)				
hyperthermia threatment of car With the proposed design conce target body and simultaneously ture regions (hot spots) occur. algorithm would be used to rapic Analysis of an annular phas merit of combining adaptive nul	ncer. A system design concep pt, it may be possible to max minimize or reduce the electr In a clinical situation, either tly form the adaptive null (or led-array antenna embedded lling with conventional near-f	t for implementing adap imize the applied electri ic field at target position a gradient search algor nulls) prior to any signif in an infinite homogene ield focusing used in hyp	se distribution in electromagnetic otive hyperthermia is introduced. ic field at a tumor position in the as where undesired high temperatithm or sample matrix inversion icant tissue heating.	

15. NUMBER OF PAGES 14. SUBJECT TERMS electromagnetic hyperthermia phased arrays 119 adaptive nulling cancer therapy 16. PRICE CODE method-of-moments adaptive arrays 17. SECURITY CLASSIFICATION 18. SECURITY CLASSIFICATION 19. SECURITY CLASSIFICATION 20. LIMITATION OF OF REPORT OF THIS PAGE OF ABSTRACT **ABSTRACT** SAR Unclassified Unclassified Unclassified

theory and software used to compute the moment-method received voltage at a short-dipole probe due to a transmitting dipole array are documented. Computer simulations show that adaptive nulling can prevent undesired high-tempera-

ture regions from occurring while simultaneously heating a deep-seated tumor site.