AD-A259 908

(12)

ERDEC-CR-002

DETECTION OF PICOGRAM QUANTITIES OF BOTULINUS TOXIN-B USING THE LIGHT ADDRESSABLE POTENTIOMETRIC SENSOR



James P. Chambers

UNIVERSITY OF TEXAS AT SAN ANTONIO San Antonio, TX 78285

James J. Valdes

RESEARCH AND TECHNOLOGY DIRECTORATE

December 1992

Approved for public release; distribution is unlimited.

93-02054

U.S. ARMY CHEMICAL AND BIOLOGICAL DEFENSE AGENCY



Disclaimer

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorizing documents.

REPORT DOCUMENTATION PAGE

Form Approved
OMB No 0704-0188

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson (Nazis Highway, Suite 1204, Arlington, VAI 22202-4302, and to the Office of Management and Budget, Paperwork Reduction Project (0704-0188), Washington, DC 20503

1. AGENCY USE ONLY (Leave blank)	2. REPORT DATE 1992 December	3. REPORT TYPE AND Final, 91 Sep			
4. TITLE AND SUBTITLE	L	5	. FUNDING NUMBERS		
Detection of Picogram Qua	antities of Botulinu	s Toxin-B	C-DAAA15-C-0008		
Using the Light Addressable Potentiometric Sensor			51		
6. AUTHOR(S)					
Chambers, James P. (University of Texas at San					
Antonio), and Valdes, Jam	nes J. (ERDEC)				
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES)			PERFORMING ORGANIZATION REPORT NUMBER		
University of Texas at San Antonio, San Antonio, TX 78285					
9. SPONSORING/MONITORING AGENCY	NAME(S) AND ADDRESS(ES)	1	0. SPONSORING / MONITORING		
DIR, ERDEC, * ATTN: SCBRD-RT, APG, MD 21010-5423			AGENCY REPORT NUMBER ERDEC-CR-002		
11. SUPPLEMENTARY NOTES					
*At the time this work wa	as completed, ERDEC	was known as CRD	EC.		
12a. DISTRIBUTION / AVAILABILITY STAT	EMENT	1	26. DISTRIBUTION CODE		
Approved for public releatunlimited.	ase; distribution is	3			
13. ABSTRACT (Maximum 200 words)		 			
Detection of picogram qua purified biotin- and fluo					
report. All assays were					
Sensor. Results indicate	e a near linear dose	response curve	for increasing the		
	amounts of BoTX (10-1,000 pg). The presence of 1-30 ng of toxin resulted in a				
•	dose response curve that yielded high output signals (>4,000 μ V/s) at saturation				
(approximately 25 ng). V					
indicated very low output signals (15% control), which indicated minimal biosensor					
nonspecificity.		•			
,					
Į					
1					
<u> </u>					
14. SUBJECT TERMS			15. NUMBER OF PAGES		
Receptors Toxins		13			
Biosensors			16. PRICE CODE		
	SECURITY CLASSIFICATION OF THIS PAGE	19. SECURITY CLASSIFICA OF ABSTRACT	ATION 20. LIMITATION OF ABSTRACT		

UNCLASSIFIED

 \mathbf{u}

UNCLASSIFIED

Blank

PREFACE

The work described in this report was authorized under Contract No. DAAA15-C-0008. This work was started in September 1991 and completed in January 1992.

The use of trade names or manufacturers' names in this report does not constitute an official endorsement of any commercial products. This report may not be cited for purposes of advertisement.

Reproduction of this document in whole or in part is prohibited except with permission of the Director, Edgewood Research, Development and Engineering Center (ERDEC), ATTN: SCBRD-RTD, Aberdeen Proving Ground, MD 21010-5423. However, the Defense Technical Information Center and the National Technical Information Service are authorized to reproduce the document for U.S. Government purposes.

This report has been approved for release to the public.

<u>Acknowledgments</u>

The authors thank Dr. Eppie Rael (University of Texas at El Paso, El Paso, TX) for supplying the immunoreagents; James Colston and Dr. Pramod Kumas (University of Texas at San Antonio, San Antonio, TX) for technical support; and Marcy Goforth (Research and Technology Directorate, ERDEC) for preparing the final manuscript.

Anna	sion Pe	
	GRALI	2
DTIC	TAB	
Unanc	ounced	ā
Justi	fication	0
;		
By	ance of the same o	
	ibution	/
Aval	1361116	
•	ivall a	10/06
Dist	Speci	al l
10 - 11	1	
IH I		2004
11.		· E

Blank

14.4 C

k ia Selici 4

CONTENTS

		Page
1.	INTRODUCTION	7
2.	MATERIALS AND METHODS	7
2.1	Reagents	7
2.2	Instrumentation	7
2.3	Purification and Labelling of Immunoglobulin (IgG)	8
2.4	Assay Procedure	8
3.	RESULTS	9
4.	DISCUSSION	9
	LITERATURE CITED	12

Blank

DETECTION OF PICOGRAM QUANTITIES OF BOTULINUS TOXIN-B USING THE LIGHT ADDRESSABLE POTENTIOMETRIC SENSOR

1. INTRODUCTION

The extremely potent exotoxins produce by <u>Clostridium botulinum</u> [also known as botulinus toxin (BoTX)] cause fatal neurotoxicity in both animals and humans. Seven serologically distinct BoTXs were isolated from this organism and each share similar structural characteristics. A 145,000-150,000 M_r protoxin is enzymatically cleaved into the active toxin, which consists of two peptides of approximately 50,000 and 100,000 M_r linked by a disulfide bridge. The active toxin binds specifically to gangliosides and prevents presynaptic release of acetylcholine by an unknown mechanism. $_2$

Sensitive and rapid detection of BoTX is invaluable for testing contamination of foodstuffs. In this report, we describe an immunoligand assay that used a Light Addressable Potentiometric Sensor (LAPS), purchased from Molecular Devices Corporation (Menlo Park, CA), to detect picogram (pg) quantities of BoTX. The LAPS is a recently developed, sensitive and reliable technology utilizing nitrocellulose filters onto which fluorescein-labelled immunocomplexes are immobilized and labelled with an antifluorescein urease conjugated antibody. 3,4,5 Quantitation of toxin is accomplished by inserting immobilized, urease-labelled immunocomplex into a reader containing a solution of urea. Enzymatic breakdown of urea produces a potentiometric shift proportional to the amount of urease, hence, toxin, which is detected rapidly with high sensitivity and reproducibility. Biotin-and fluorescein-labelled anti-BoTX polyclonal antibodies were used in this study to detect picogram quantities of toxin on the LAPS.

2. MATERIALS AND METHODS

2.1 Reagents.

Goat polyclonal antiserum against BoTX-B was obtained from Biodesign International (Kenne Bunkport, ME). Antibodies were biotinylated and fluoresceinated using a labelling kit that contained the N-hydroxysuccinimide esters of DNP-biotin and carboxyfluorescein. Labelling reagents were reconstituted with dimethylformamide. Samples were concentrated using Centricon-30 Microconcentrators (Amicon Division, W.R. Grace and Company, Beverly, MA). All general laboratory reagents were of the highest quality. Venoms from the Heloderma horridum, Naja nigricollis, Agkistrodon acutus, Apis mellifera, Leiurus quinquestriatus, Bitis nasicornus, Naja atra, Crotalus, viridis helleri, and Bufo americus were purchased from Sigma Chemical Company (St. Louis, MO).

2.2 <u>Instrumentation</u>.

All assays were performed on the LAPS, and the data were collected and stored digitally on a microcomputer.

2.3 <u>Purification and Labelling of Immunoglobulin (IqG)</u>.

IgG was purified from serum by passage over a Protein A column, following the procedures described by the manufacturer (Promab, JRH Biosciences, Lenexa, KS). Antibodies were concentrated using a Centricon-30 microconcentrator (20,000 x g for 20 min) and washed two times with phosphate buffered saline (PBS) until a volume of 1 mL remained. The final concentration of purified IgG was 1.36 mg/mL as determined by the Bradford method⁶ using bovine serum albumin as a standard.

Antibody preparations were separated into two aliquots. One aliquot was biotinylated [DNP biotin labelling reagent for 2 hr at room temperature using a molar coupling ration (MCR) of 30:1]. The remaining aliquot was fluoresceinated (fluorescein labelling reagent for 2 hr at room temperature using an MCR of 30:1). We determined (data not shown) that a n MCR of 30:1 provides adequate and efficient labelling of IgG antibodies. Unreacted haptens were removed by dialysis following six, 1-L changes of PBS. Dialyzed samples were stored at 4 °C.

2.4 <u>Assay Procedure</u>.

For each test, BoTX and venoms were dissolved in a $100-\mu L$ assay buffer (Molecular Devices Corporation). Labelled antibodies were then diluted in the assay buffer to a final concentration of 2 ng/ μL . Prior to assaying, biotinylated and fluoresceinated antibodies were combined in a 1:1 ratio in that each 100 μL contained 100 ng of each labelled antibody. We determined this ratio and concentration of antibodies to be optimal for sensitive detection of BoTX (data not shown).

Primary immunocomplexes were formed by adding 100 μ l antibody solution to the reaction mixture containing toxin. After gently vortexing, assay mixtures were incubated at 37 °C for 45 min. After incubation, 1 mL of capture reagent (Molecular Devices Corporation) containing streptavidin was added to each reaction mixture followed by a brief mixing. Each assay mixture was then transferred immediately to separate test wells, and vacuum filtered through a biotinylated nitrocellulose filter mounted on a test stick (Molecular Devices Corporation). Immunocomplexes were captured (immobilized) by formation of a nitrocellulose biotin-streptavidin-biotin conjugate. Unbound material was removed by adding a 2-mL wash buffer (Molecular Devices Corporation) to each test well and followed by vacuum filtration.

Urease conjugated, antifluorescein, antibody solution (Molecular Devices Corporation) was added to each test well and filtered. A final wash with a 2-mL wash buffer was performed, and the nitrocellulose test sticks were removed and then placed in a beaker containing wash buffer until read. Each nitrocellulose stick was inserted into the LAPS reader, and the respective immunoconjugates were monitored via hydrolysis of urea.

3. RESULTS

The effect of increasing amounts of BoTX on output signal is shown in Figure 1. A linear response (r = 0.98) over the lowest range of toxin from 100 to 1,000 pg was observed (shown in A, Figure 1). A lower detection limit of 250 pg (78.7 \pm 1.5 $\mu V/s$) is attainable under these conditions. A dose response curve of increasing amounts of BoTX (1-30 ng) that indicates saturation to occur at approximately 25 ng with a corresponding mean output signal of 4219.4 \pm 301.4 $\mu V/s$ is shown in B, Figure 1. Further, the effect of increasing BoTX over an extended concentration range (10 $^{-10}$ - 10 $^{-6}$ g) is shown in C, Figure 1. Figure 2 shows the effect of various venoms on the output signal, and the results indicate very low, assay nonspecificity (<15% of control).

4. DISCUSSION

Active BoTX-B consists of heavy and light chains (M_r approximately 100,00 and 50,000, respectively) linked by a disulfide that are bridged together with several <u>Botulinum</u>-associated protein that yield a complex a with a molecular weight of nearly 500,000 datlons. Toxins produced from nonproteolytic cultures require exogenous protease treatment for protoxin activation. However, in this study, BoTX-B was not treated with exogenous protease. Previous studies in our laboratory indicated the toxin complex to be extensively nicked (5 pieces), which upon Western blotting, revealed approximately 70% of the immunoreactivity to reside in 2 peptides ($M_r = 37,480$ and 50,470).

Detection of picogram quantities of BoTX was achieved using labelled polyclonal IgG antibodies. Data in Figure 1 indicate the lowest range of detection under these assay conditions to be between 100 and 250 pg. This amount of toxin approaches the minimal lethal dose in the mouse. 10,11 Although the output signals observed when detecting picogram quantities of BoTX are rather low (<200 $\mu V/s$), they are very reproducible, which are indicated by the fact that the standard error bars are not visibly distinct from the data points shown in A, Figure 1.

Output signal produced by toxin in the 1-30 ng range (B, Figure 1) indicated saturation to occur in approximately 25 ng toxin. BoTX in excess of 30 ng results in a sharp decrease in the signal shown in C, Figure 1, is most likely due to dilution of antibody by excess toxin. It is essential that antibody concentration be in relatively large excess when compared to that of toxin for reliable detection. Detection of toxin in excess of 30 ng is possible, but larger amounts of labelled antibody is required, which results in a reduction of sensitivity.

It was important to determine the extent of biosensor nonspecificity. To accomplish this, the effect of nine different challenge venoms on the signal output was determined. No attempt was made to isolate and purify individual proteins from the venoms for the purpose of exposing the assay to the maximum amount of toxic venom component(s). When compared to an




Figure 1. Detection of BoTX




Figure 2. Nonspecific Output Signal Produced by Nine Difference Venoms

LITERATURE CITED

- 1. DasGupta, B.R., and Sugiyama, H., <u>Biochem. Biophys. Res.</u>
 <u>Commun.</u> Vol. 48, p 108 (1972).
- 2. Burgen, A.S.V., Dickens, F., and Eatman, L.J., <u>J. Physiol.</u> <u>London</u> Vol. 109, p 10 (1949).
 - 3. Briggs, J., Nature, Vol. 329, p 565, 1987.
- 4. Hafeman, D., Parce, J., and McConnell, H., <u>Science</u>, Vol. 240, p 1182, 1988.
 - 5. Parce, J., et al., Science, Vol. 240, p 1182, 1988.
 - 6. Bradford, M.A., Anal. Biochem. Vol. 72, p 248 (1976).
- 7. Wagman, J., and Bateman, J.B., Arch. Biochem. Biophys. Vol. 31, p 424 (1951).
 - 8. DasGupta, B.R., J. <u>J. Physiol. Paris</u> Vol. 84, p 220 (1990).
- 9. Kumar, P., et al., "Detection of Botulinus Toxin Using an Evanescent Wave Immunosenor," <u>Recent Advances in Toxicology Research</u>, Vol. III, pp 361-374, 1992.
- 10. Hanig, J.P., and Lamanna, C., <u>J. Theor. Biol.</u> Vol. 77, p 107 (1979).
- 11. Schantz, E.J., and Johnson, E.A., <u>Microbiol.</u> Vol. 56, p 80 (1992).